# VIASURE


## Real Time PCR Detection Kits

by CerTest  
BIOTEC

### ***Herpes virus 1, Herpes virus 2 & Treponema pallidum***

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|                                                                                                                                     |            |
|-------------------------------------------------------------------------------------------------------------------------------------|------------|
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile   | VS-HHT106L |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile  | VS-HHT106H |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile  | VS-HHT112L |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile | VS-HHT112H |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 96-well plate, low profile       | VS-HHT113L |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 96-well plate, high profile      | VS-HHT113H |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®   | VS-HHT136  |
| VIASURE <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®  | VS-HHT172  |

## ENGLISH

### 1. Intended use

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* in clinical specimens from patients with signs and symptoms of genital ulcer. This test is intended to be used as an aid in the diagnosis of *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific for *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum*.

### 2. Summary and Explanation

*Treponema pallidum* and herpes simplex virus types 1 or 2 (HSV-1/ 2) are the most common sexually transmitted pathogens that cause ulcers of the genital, anal and oropharyngeal region. Genital ulcers are a significant risk factor for transmission and acquisition of human immunodeficiency virus (HIV). Identification of ulcer-causing pathogens is a prerequisite for an effective therapy and reducing the risk of HIV infection.

Genital herpes can be caused by two very similar viruses, herpes simplex virus HSV-1 or HSV-2, being the infections caused by first one more often. Herpes simplex virus (HSV)-1 and -2 are large, double-stranded DNA viruses that cause lifelong persistent infections characterized by periods of quiescence and recurrent disease. These two HSV types cannot be distinguished clinically. In fact, they share a high degree of genetic homology, but they also have specific regions with small nucleotide variations which may allow discrimination. HSV-1 and HSV-2 infection occurs via inoculation of virus particles into susceptible mucosal surfaces. Afterwards, these neurotropic viruses can become latent in the local sensory ganglion, periodically reactivating to cause symptomatic lesions, or undergo asymptomatic viral release, with the potential for disease transmission and infection. The same treatment is used for both HSV-1 and HSV-2 infections, the location of the lesions and the chronicity of the infection (primary or recurrent) determine dosage and frequency.

The spirochete *Treponema pallidum* (Tp), the etiologic agent of syphilis, causes a multistage sexually transmitted infection (STI). Pathogenic treponemes cause venereal syphilis, yaws, endemic syphilis, and pinta—multistage, infections that, although similar, can be differentiated based on clinical, epidemiologic, and geographic criteria. Only venereal syphilis is transmitted by sexual activity. The pathogenic treponemes are uncultivable, slow-growing microorganisms with identical flat-wave morphologies. They poorly tolerate desiccation, elevated temperature, and ambient oxygen tension, traits that explain why efficient transmission requires close personal contact.

Several types of laboratory test are available for diagnosis of *Treponema pallidum* and herpes simplex virus types 1 or 2. Serology, dark field microscopy and PCR are the gold standards to diagnose these pathogens. However, Real Time PCR assays have been shown to be a more sensitive and specific diagnosis tool for their detection.

### 3. Principle of the procedure

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit is designed for the diagnosis of *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* in clinical samples. After DNA isolation, the identification of *Herpes virus 1*, *Herpes virus 2* and/or *Treponema pallidum* is performed by the amplification of a conserved region of the *US4* gene for *Herpes virus 1*, *US6* gene for *Herpes virus 2* and 16S rRNA gene for *Treponema pallidum*, using specific primers and a fluorescent-labelled probe.

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in an stabilized format. *Herpes virus 1* DNA targets are amplified and detected in the FAM channel, *Herpes virus 2* DNA targets are amplified and detected in the ROX channel, *Treponema pallidum* DNA targets are amplified and detected in the Cy5 channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

### 4. Reagents provided

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

| Reagent/Material                                                                              | Description                                                                                            | Colour      | Amount              |
|-----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------------|---------------------|
| <i>Herpes virus 1</i> , <i>Herpes virus 2</i> & <i>Treponema pallidum</i><br>8-well strips    | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White       | 6/12 x 8-well strip |
| Rehydration Buffer                                                                            | Solution to reconstitute the stabilized product                                                        | Blue        | 1 vial x 1.8 mL     |
| <i>Herpes virus 1</i> , <i>Herpes virus 2</i> & <i>Treponema pallidum</i><br>Positive Control | Non-infectious synthetic lyophilized DNA                                                               | Red         | 1 vial              |
| Negative control                                                                              | Non template control                                                                                   | Violet      | 1 vial x 1 mL       |
| Water RNase/DNase free                                                                        | RNase/DNase free water                                                                                 | White       | 1 vial x 1 mL       |
| Tear-off 8-cap strips                                                                         | Optical caps for sealing wells during thermal cycling                                                  | Transparent | 6/12 X 8-cap strip  |

Table 1. Reagents and materials provided in VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit with Ref. VS-HHT106L, VS-HHT106H, VS-HHT112L and VS-HHT112H.

| Reagent/Material                                                                | Description                                                                                            | Color       | Amount           |
|---------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------------|------------------|
| <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> 96-well plate    | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White       | 1 plate          |
| Rehydration Buffer                                                              | Solution to reconstitute the stabilized product                                                        | Blue        | 1 vial x 1.8 mL  |
| <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Positive Control | Non-infectious synthetic lyophilized DNA                                                               | Red         | 1 vial           |
| Negative control                                                                | Non template control                                                                                   | Violet      | 1 vial x 1 mL    |
| Water RNase/DNase free                                                          | RNase/DNase free water                                                                                 | White       | 1 vial x 1 mL    |
| Tear-off 8-cap strips                                                           | Optical caps for sealing plate during thermal cycling                                                  | Transparent | 12 X 8-cap strip |

Table 2. Reagents and materials provided in VIASURE *Herpes virus 1, Herpes virus 2 & Treponema pallidum* Real Time PCR Detection Kit with Ref VS-HHT113L and VS-HHT113H.

| Reagent/Material                                                                | Description                                                                                            | Colour      | Amount              |
|---------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------------|---------------------|
| <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> 4-well strips    | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | Transparent | 9/18 x 4-well strip |
| Rehydration Buffer                                                              | Solution to reconstitute the stabilized product                                                        | Blue        | 1 vial x 1.8 mL     |
| <i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i> Positive Control | Non-infectious synthetic lyophilized DNA                                                               | Red         | 1 vial              |
| Negative control                                                                | Non template control                                                                                   | Violet      | 1 vial x 1 mL       |
| Water RNase/DNase free                                                          | RNase/DNase free water                                                                                 | White       | 1 vial x 1 mL       |
| 4-cap strips                                                                    | Optical caps for sealing wells during thermal cycling                                                  | Transparent | 9/18 X 4-cap strip  |

Table 3. Reagents and materials provided in VIASURE *Herpes virus 1, Herpes virus 2 & Treponema pallidum* Real Time PCR Detection Kit with Ref. VS-HHT136 and VS-HHT172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Herpes virus 1, Herpes virus 2 & Treponema pallidum* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-HHT113L, VS-HHT113H, VS-HHT136 and VS-HHT172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-HHT136 and VS-HHT172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.

- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

## 8. Test procedure

### 8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from clinical samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- Invisorb® Spin Universal Kit (Strattec).

### 8.2. Lyophilized positive control

*Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

| Cycles | Step                                   | Time   | Temperature |
|--------|----------------------------------------|--------|-------------|
| 1      | Initial denaturation                   | 2 min  | 95°C        |
| 45     | Denaturation                           | 10 seg | 95°C        |
|        | Annealing/Extension (Data collection*) | 50 seg | 60°C        |

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (*Herpes virus 1*), ROX (*Herpes virus 2*), Cy5 (*Treponema pallidum*) and HEX, JOE or VIC channels (IC). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

| Herpes virus 1 (FAM) | Herpes virus 2 (ROX) | <i>Treponema pallidum</i> (Cy5) | Internal control (HEX) | Negative Control | Positive Control | Interpretation                                                                     |
|----------------------|----------------------|---------------------------------|------------------------|------------------|------------------|------------------------------------------------------------------------------------|
| +                    | +                    | +                               | +/-                    | -                | +                | Herpes virus 1, Herpes virus 2 and <i>Treponema pallidum</i> Positive              |
| -                    | -                    | -                               | +                      | -                | +                | Herpes virus 1, Herpes virus 2 and <i>Treponema pallidum</i> Negative              |
| +                    | -                    | -                               | +/-                    | -                | +                | Herpes virus 1 Positive, Herpes virus 2 and <i>Treponema pallidum</i> Negative     |
| +                    | +                    | -                               | +/-                    | -                | +                | Herpes virus 1 and Herpes virus 2 Positive, and <i>Treponema pallidum</i> Negative |
| +                    | -                    | +                               | +/-                    | -                | +                | Herpes virus 1 and <i>Treponema pallidum</i> Positive, and Herpes virus 2 Negative |
| -                    | +                    | -                               | +/-                    | -                | +                | Herpes virus 2 Positive, Herpes virus 1 and <i>Treponema pallidum</i> Negative     |
| -                    | +                    | +                               | +/-                    | -                | +                | Herpes virus 2 and <i>Treponema pallidum</i> Positive, Herpes virus 1 Negative     |
| -                    | -                    | +                               | +/-                    | -                | +                | <i>Treponema pallidum</i> Positive, Herpes virus 1 and Herpes virus 2 Negative     |
| -                    | -                    | -                               | -                      | -                | +                | Experiment fail                                                                    |
| +                    | +                    | +                               | +                      | +                | -                | Experiment fail                                                                    |

Table 5. Sample interpretation

+: Amplification curve

-: No amplification curve


A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from samples in transport medium, endocervical swabs, urine, plasma, cell lysates and synthetic CSF.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum*, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit was tested using 8 endocervical swabs from symptomatic patients. The comparison was carried out using the results obtained from the routine analysis performed in the hospital with the molecular detection method "CLART® ENTHERPEX" (Genómica). Two samples were HSV-1 positive and 6 were detected as HSV-2 positive. The concordance of both kits was 100%.

The clinical performance of VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit was tested using 89 samples from QCMD, UK NEQAS and INSTAND EQA programs. These samples included urine, plasma, cell lysates, synthetic CSF, genital, anal or oropharyngeal ulcer-causing pathogens. These results were compared with those obtained with a molecular detection method (FTD Genital ulcer (Fast Track Diagnostics)) and/or with the corresponding reports. All samples were detected correctly. 23/81 were HSV-1 positive, 21/81 were HSV-2 positive and 15/81 were *Treponema pallidum* positive.

In conclusion, the results show a high sensitivity and specificity to detect *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* using VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit.

### 12.2. Analytical sensitivity

VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *Herpes virus 1*, *Herpes virus 2* and *Treponema pallidum* (Figure 2, 3 and 4).

Figure 2. Dilution series of Herpes virus 1 ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).


Figure 3. Dilution series of Herpes virus 2 ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).



Figure 4. Dilution series of *Treponema pallidum* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

### 12.3. Analytical specificity

The specificity of the *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* assay was confirmed by testing panel consisting of different microorganisms representing the most common enteric and genitourinary pathogens or flora present in the gastrointestinal tract or urogenital system. No cross-reactivity was detected between any of the following microorganisms tested.

| Cross-reactivity testing           |   |                                           |   |                                                  |   |
|------------------------------------|---|-------------------------------------------|---|--------------------------------------------------|---|
| <i>Candida albicans</i>            | - | <i>Enterococcus faecium</i>               | - | <i>Proteus mirabilis</i>                         | - |
| <i>Candida parapsilosis</i>        | - | <i>Escherichia coli</i> 0.1285; O18:H7:K1 | - | <i>Pseudomonas aeruginosa</i>                    | - |
| <i>Candida tropicalis</i>          | - | <i>Gardnerella vaginalis</i>              | - | <i>Serratia marcescens</i>                       | - |
| <i>Candida glabrata</i>            | - | <i>Haemophilus influenzae</i>             | - | <i>Staphylococcus aureus</i>                     | - |
| <i>Candida krusei</i>              | - | <i>Haemophilus ducrey</i>                 | - | <i>Stenotrophomonas maltophilia</i>              | - |
| <i>Candida dubliniensis</i>        | - | <i>Klebsiella oxytoca</i>                 | - | <i>Streptococcus pneumoniae</i>                  | - |
| <i>Aspergillus fumigatus</i>       | - | <i>Klebsiella pneumoniae</i>              | - | <i>Streptococcus agalactiae</i>                  | - |
| <i>Bacteroides fragilis</i>        | - | <i>Listeria ivanovii</i>                  | - | <i>Trichomonas vaginalis</i>                     | - |
| <i>Acinetobacter baumannii</i>     | - | <i>Listeria monocytogenes</i>             | - | <i>Ureaplasma parvum</i>                         | - |
| <i>Chlamydia trachomatis</i> (LGV) | - | <i>Listeria innocua</i>                   | - | <i>Ureaplasma urealyticum</i>                    | - |
| <i>Chlamydia trachomatis</i> (SW)  | - | <i>Mycoplasma genitalium</i>              | - | Cytomegalovirus                                  | - |
| <i>Enterobacter cloacae</i>        | - | <i>Mycoplasma hominis</i>                 | - | Hepatitis A                                      | - |
| <i>Enterobacter aerogenes</i>      | - | <i>Neisseria gonorrhoeae</i>              | - | Human papillomavirus genotypes<br>HP16 and HPV18 | - |
| <i>Enterococcus faecalis</i>       | - | <i>Neisseria meningitidis</i>             | - |                                                  |   |

Table 6. Sexual transmitted pathogens and microbial flora used in this study.

### 12.4. Analytical reactivity

The reactivity of VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit for *Herpes virus 1* was evaluated against HSV 1 MacIntyre strain showing positive result.

The reactivity of VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit for *Herpes virus 2* was evaluated against HSV 2 MS strain showing positive result.

The reactivity of VIASURE *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Real Time PCR Detection Kit for *Treponema pallidum* was evaluated against *Treponema pallidum* showing positive result.

## ANNEX 1

## COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

| Table A.1 LOW PROFILE BLOCK THERMOCYCLERS |                                                             |
|-------------------------------------------|-------------------------------------------------------------|
| Manufacturer                              | Model                                                       |
| Agilent Technologies                      | AriaMx/AriaDx Real-Time PCR System                          |
| Applied Biosystems                        | 7500 Fast Real-Time PCR System <sup>(1)</sup>               |
| Applied Biosystems                        | 7500 Fast Dx Real-Time PCR System <sup>(1)</sup>            |
| Applied Biosystems                        | QuantStudio™ 12K Flex 96-well Fast                          |
| Applied Biosystems                        | QuantStudio™ 6 Flex 96-well Fast                            |
| Applied Biosystems                        | QuantStudio™ 7 Flex 96-well Fast                            |
| Applied Biosystems                        | QuantStudio™ 3 Fast Real-Time PCR System <sup>(5)</sup>     |
| Applied Biosystems                        | QuantStudio™ 5 Fast Real-Time PCR System                    |
| Applied Biosystems                        | StepOne Plus™ Real-Time PCR System <sup>(5)</sup>           |
| Applied Biosystems                        | StepOne™ Real-Time PCR System <sup>(5)</sup>                |
| Applied Biosystems                        | ViiA™ 7 Fast Real-Time PCR System                           |
| BIONEER                                   | Exicycler™ 96                                               |
| Bio-Rad                                   | CFX96™ / CFX96™ IVD Real-Time PCR Detection System          |
| Bio-Rad                                   | Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup> |
| Cepheid                                   | SmartCycler® <sup>(3)</sup>                                 |
| Qiagen                                    | Rotor-Gene® Q <sup>(3)</sup>                                |
| Roche                                     | LightCycler® 480 Real-Time PCR System <sup>(4)</sup>        |
| Roche                                     | LightCycler® 96 Real-Time PCR System <sup>(4)</sup>         |
| Roche                                     | Cobas z480 Analyzer <sup>(4)</sup>                          |

| Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS |                                                                        |
|--------------------------------------------|------------------------------------------------------------------------|
| Manufacturer                               | Model                                                                  |
| Abbott                                     | Abbott m2000 RealTime System                                           |
| Applied Biosystems                         | 7300 Real-Time PCR System <sup>(5)</sup>                               |
| Applied Biosystems                         | 7500 Real-Time PCR System                                              |
| Applied Biosystems                         | 7900 HT Real-Time PCR System <sup>(5)</sup>                            |
| Applied Biosystems                         | ABI PRISM 7000 <sup>(6)</sup>                                          |
| Applied Biosystems                         | ABI PRISM 7700 <sup>(5)</sup>                                          |
| Applied Biosystems                         | QuantStudio™ 12K Flex 96-well                                          |
| Applied Biosystems                         | QuantStudio™ 6 Flex 96-well                                            |
| Applied Biosystems                         | QuantStudio™ 7 Flex 96-well                                            |
| Applied Biosystems                         | QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>                     |
| Applied Biosystems                         | QuantStudio™ 5 Real-Time PCR System                                    |
| Applied Biosystems                         | ViiA™ 7 Real-Time PCR System                                           |
| Analytik Jena Biometra                     | TOptical                                                               |
| Analytik Jena Biometra                     | qTOWER 2.0                                                             |
| BIONEER                                    | Exicycler™ 96                                                          |
| Bio-Rad                                    | CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System |
| Bio-Rad                                    | iCycler iQ™ Real-Time PCR Detection System                             |
| Bio-Rad                                    | iCycler iQ™5 Real-Time PCR Detection System                            |
| Bio-Rad                                    | MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>                    |
| Bio-Rad                                    | MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>                   |
| Cepheid                                    | SmartCycler® <sup>(3)</sup>                                            |
| DNA-Technology                             | DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>              |
| DNA-Technology                             | DTlite Real-Time PCR System <sup>(2)</sup>                             |
| Eppendorf                                  | Mastercycler™ep realplex                                               |
| Qiagen                                     | Rotor-Gene® Q <sup>(3)</sup>                                           |
| Stratagene / Agilent Technologies          | Mx3000P™ Real Time PCR System                                          |
| Stratagene / Agilent Technologies          | Mx3005P™ Real Time PCR System                                          |
| VIASURE                                    | VIASURE 48 Real Time PCR System <sup>(2)</sup>                         |
| VIASURE                                    | VIASURE 96 Real Time PCR System <sup>(2)</sup>                         |

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

## DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

| REAL-TIME PCR THERMOCYCLER    | VIASURE CHANNEL | DETECTION CHANNEL | OBSERVATIONS                                                                                                                                                                                                                                                                                                                                                    |
|-------------------------------|-----------------|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Bio-Rad CFX96™                | FAM             | FAM               | Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.                                                                                      |
|                               | HEX             | HEX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| ABI 7500 Applied Biosystems   | FAM             | FAM               | Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.         |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Roche Lightcycler®480II       | FAM             | 465/510           | Colour Compensation is required                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | 533/580           |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | 533/610           |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | 618/660           |                                                                                                                                                                                                                                                                                                                                                                 |
| Smartcycler® Cepheid          | FAM             | Channel 1         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | Channel 2         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Channel 3         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Channel 4         |                                                                                                                                                                                                                                                                                                                                                                 |
| Abbott m2000rt                | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Mx3000P™ Mx 3005P™ Stratagene | FAM             | FAM               | Passive reference option for ROX must be none                                                                                                                                                                                                                                                                                                                   |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| AriaMx Agilent                | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | HEX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Rotor-Gene®Q Qiagen           | FAM             | Green             | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".                                                                                                   |
|                               | HEX             | Yellow            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Orange            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Red               |                                                                                                                                                                                                                                                                                                                                                                 |
| Mic Real Time PCR Cycler bms  | FAM             | Green             | In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10) |
|                               | HEX             | Yellow            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Orange            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Red               |                                                                                                                                                                                                                                                                                                                                                                 |
| Exicycler™ 96 BIONEER         | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | JOE               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

**OPTICAL MEASUREMENT EXPOSURE SETTING**

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cyclers (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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## 13. Bibliography/Bibliografía

1. J. D. Radolf *et al.* *Treponema pallidum*, the syphilis spirochete: making a living as a stealth pathogen. *Nature Reviews Microbiology* 2016; 14(12): 744–759.
2. M. Lieveld *et al.* A high resolution melting (HRM) technology-based assay for cost-efficient clinical detection and genotyping of herpes simplex virus (HSV)-1 and HSV-2. *Journal of Virological Methods* 2017; Volume 248, 181-186.
3. M. Costa-Silva *et al.* Cross-sectional study of *Treponema pallidum* PCR in diagnosis of primary and secondary syphilis. *International Journal of Dermatology* 2017; doi: 10.1111/ijd.13823.

4. M.A. Minaya *et al.* Molecular evolution of herpes simplex virus 2 complete genomes: Comparison between primary and recurrent infections. *American Society for Microbiology* 2017; doi:10.1128/JVI.00942-17.
5. C.O. Onyango *et al.* Evaluation of a TaqMan Array Card for Detection of Central Nervous System Infections. *Journal of Clinical Microbiology* 2017 Jul;55(7): 2035-2044.

## 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: March 2019

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.cer-test.es](http://www.cer-test.es)

VIASURE online

F-362 rev01

**VIASURE**

CE

IVD

Real Time PCR Detection Kits

**CerTest**  
BIOTEC