VIASURE

Real Time PCR Detection Kit

Group A Streptococcus

CE IVD

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Content

1.	Intended use	6
2.	Summary and Explanation	6
3.	Principle of the procedure	6
4.	Reagents provided	7
5.	Reagents and equipment to be supplied by the user	7
6.	Transport and storage conditions.....	8
7.	Precautions for users	8
8.	Test procedure	9
9.	Result interpretation	10
10.	Limitations of the test	13
11.	Quality control.....	14
12.	Performance characteristics.....	14
	ANNEX 1	17
A1.1	Principle of the procedure.....	17
A1.2	Reagents provided.....	17
A1.3	Test procedure	18
	ANNEX 2	20
A2.1	Principle of the procedure.....	20
A2.2	Reagents provided.....	20
A2.3	Test procedure	20
	ANNEX 3	23
A3.1	Principle of the procedure.....	23
A3.2	Reagents provided.....	23
A3.3	Test procedure	25
	ANNEX 4	27
A4.1	Principle of the procedure.....	27
A4.2	Reagents provided.....	27
A4.3	Test procedure	27

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ENGLISH

1. Intended use

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit is designed for the specific detection of DNA of the *Streptococcus pyogenes* (group A (beta-hemolytic) *Streptococcus* (Strep-A)) in pharyngeal and otic exudate swabs from patients with signs and symptoms of respiratory infection. This test is intended to be used as an aid in the diagnosis of *Streptococcus pyogenes* infection in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Streptococcus pyogenes*.

2. Summary and Explanation

Group A *Streptococcus*, also known as group A beta-hemolytic streptococcus (GAS), is a gram-positive bacterium that grows on chains, producing white to grayish colonies surrounded by a halo of hemolysis when grown in blood agar medium. It differs from other beta-hemolytic streptococci by possessing, in its cell wall, a group A specific polysaccharide antigen (Lancefield antigen).

It is spread mainly by contact, for example, through sneezing or kissing. Some people can carry the bacteria without showing signs of illness.

Group A *Streptococcus* is usually responsible for upper respiratory tract infections and skin infections. In some cases, infection with these microorganisms can go beyond acute episodes and trigger serious post-infectious syndromes such as acute rheumatic fever and glomerulonephritis. Diseases caused by this microorganism are usually classified as suppurative diseases (pharyngotonsillitis, scarlet fever, impetigo, erysipelas, toxic shock syndrome and necrotizing fasciitis) and non-suppurative (rheumatic fever and acute glomerulonephritis).

In order to avoid these complications, a timely and adequate diagnosis and treatment is essential to avoid the risk of the appearance of suppurative sequelae and, especially, non-suppurative sequelae.

The reference test is the culture of tonsillar exudate, whose sensitivity and specificity are high but has limitations such as the delay in the result, and which does not allow us to distinguish between acute infection and carrier status. That is why, since the 1980s, rapid immunological techniques have been developed that are simple to use and low cost, which allow streptococcal antigen to be detected in a few minutes, and currently molecular diagnostic methods have also been developed as they provide greater sensitivity and specificity.

3. Principle of the procedure

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit is designed for the diagnosis of *Streptococcus pyogenes* in clinical samples. After DNA isolation, the identification of *Streptococcus pyogenes* is performed by the amplification of a conserved region of the *speB* gene, using specific primers and a fluorescent-labelled probe.

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent

signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format.

4. Reagents provided

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Group A *Streptococcus* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.

- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *Group A Streptococcus* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at – 20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is indented for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-GAS113L, VS-GAS113H, VS-GAS136 and VS-GAS172, VS-GAS113LE, VS-GAS113HE, VS-GAS136E and VS-GAS172E). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-GAS136 and VS-GAS172, VS-GAS136E and VS-GAS172E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.

- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE Group A *Streptococcus* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), "VIASURE Real Time PCR Detection Kits" does not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because it does not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for Open and Rotor-Gene format with extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE Group A *Streptococcus* Real Time PCR Detection Kit has been validated tested in pharyngeal swabs in Stuart-Amies media (Deltalab) and otic exudate collected from swabs in Amies liquid transport medium (Transystem™, COPAN). Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at 2 to 8°C for up to 72 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 72 hours), It is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 72 hours or frozen at -20°C or ideally at -70°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids

Otic exudates can be collected from standard swab with transport medium or with sterile wide-mouth, hermetically sealed bottle (samples obtained by aspiration). It is recommended to process the sample as soon as possible. The samples can be stored at room temperature for up to 2 hours. For long term store, it is recommended to keep the samples at -20°C or lower.

8.2. DNA extraction

Prepare the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from respiratory samples and otic exudates, you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co).
- Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (ROCHE)
- Maxwell® RSC Viral Total Nucleic Acid Purification Kit
- Invisorb® Spin Universal Kit (Invitex.).

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Streptococcus pyogenes* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>Streptococcus pyogenes</i> (FAM) ¹	Internal Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≤40	Valid

Table 1. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40 in control wells (PC and NC)).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

<i>Streptococcus pyogenes</i> (FAM)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
≤40	≤40 ≥40 or no signal ¹	≥40 or no signal	≤40	<i>Streptococcus pyogenes</i> Positive
≥40 or no signal	≤35 ²	≥40 or no signal	≤40	<i>Streptococcus pyogenes</i> Negative
≥40 or no signal	≥35 ² or no signal	≥40 or no signal	≥40 or no signal	Invalid
≤40	≤35 ²	≤40	≤40	Invalid

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curves.

1 The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.


2 In the case of Group A *Streptococcus* target genes negative, IC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each RT-qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

9.2. References with extraction control (references in Annex 3 and 4)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Streptococcus pyogenes* in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>Streptococcus pyogenes</i> (FAM) ¹	Extraction Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≤40	Valid

Table 3. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The Extraction Control (EC) should show an amplification signal (Ct ≤40 in control wells (PC and NC)). Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following criteria for read and analyze the results:

<i>Streptococcus pyogenes</i> (FAM)	External control (HEX)	Negative Control	Positive Control	Interpretation
≤40	≤40 ≥40 or no signal ¹	≥40 or no signal	≤40	<i>Streptococcus pyogenes</i> Positive
≥40 or no signal	≤35 ²	≥40 or no signal	≤40	<i>Streptococcus pyogenes</i> Negative
≥40 or no signal	≥35 or no signal ₂	≥40 or no signal	≥40 or no signal	Invalid
≤40	≤35 ²	≤40	≤40	Invalid


Table 4. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

2 In the case of *Streptococcus pyogenes* target genes negative, EC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each RT-qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from pharyngeal and otic exudate collected from swabs.
- The quality of the test depends on the quality of the sample; properly extracted nucleic acid from clinical samples must be extracted.

- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Group A Streptococcus*, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Group A Streptococcus* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNAse free (white vial). Each procedure must take place in established order and in separate laboratory areas.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the bacterial DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown *Streptococcus pyogenes* variant.
 - A bacterial load in the specimen below the limit of detection for the assay.
 - The presence of RT-qPCR inhibitors or other types of interfering substances. The impacts of antibiotics used to prevent the infection or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable bacteria and does not imply that these bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets bacterial sequences.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *Group A Streptococcus* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Group A Streptococcus* Real Time PCR Detection Kit was tested using clinical samples (pharyngeal and otic exudate collected from swabs) from patients with suspected respiratory infection. In order to determine the clinical diagnostic accuracy, different multicenter evaluations have been conducted through collaboration with national entities. A summary of the sites, sample type and workflow is included in the following table. The results were as follows:

	Site	Sample type	Workflow	Target
1	Hospital Clínico Universitario Lozano Blesa (Zaragoza, Spain)	Pharyngeal swabs	Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep / MagLEAD® 12gC nucleic acid extractor using MagDEA® Dx reagents / Maxwell® RSC Viral Total Nucleic Acid Purification Kit and Invisorb® Spin Universal Kit + AriaMx Real-time PCR System/Applied BioSystem 7500 RT PCR System	<i>Streptococcus pyogenes</i>
2	Hospital Clínico Universitario Lozano Blesa (Zaragoza, Spain)	Otic exudate swabs	Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep / MagLEAD® 12gC nucleic acid extractor using MagDEA® Dx reagents / Maxwell® RSC Viral Total Nucleic Acid Purification Kit and Invisorb® Spin Universal Kit + AriaMx Real-time PCR System/Applied BioSystem 7500 RT PCR System	<i>Streptococcus pyogenes</i>

Table 5. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE Group A *Streptococcus* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following tables:


Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	<i>Streptococcus pyogenes</i> Real-TM Quant kit molecular assay.	<i>Streptococcus pyogenes</i>	95	288	1	0	1 (0.96-1)	0.99 (0.98-1)	1 (0.96-1)	0.99 (0.98-1)
2	<i>Streptococcus pyogenes</i> Real-TM Quant kit molecular assay.	<i>Streptococcus pyogenes</i>	15	7	0	0	1 (0.78-1)	1 (0.59-1)	1 (0.78-1)	1 (0.59-1)

Table 6. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE Group A *Streptococcus* Real Time PCR Detection Kit.

Results show high agreement to detect *Streptococcus pyogenes* using VIASURE Group A *Streptococcus* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for *Streptococcus pyogenes* (Figure 3).

Figure 3. Dilution series of *Streptococcus pyogenes* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

12.3. Analytical specificity

The specificity of the *Streptococcus pyogenes* assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	<i>Moraxella catarrhalis</i>	-	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-
<i>Bordetella parapertussis</i>	-	<i>Streptococcus pneumoniae</i> Z022	-	Influenza A/Anhui/1/2013 (H7N9) virus	-
<i>Bordetella holmesii</i>	-	<i>Streptococcus salivarius</i>	-	Influenza B/Brisbane/60/2008 virus	-
<i>Bordetella bronchiseptica</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Influenza B/Florida/04/06 virus	-
<i>Chlamydia caviae</i>	-	<i>Staphylococcus epidermidis</i>	-	Influenza B/Phuket/3073/2013 virus	-
<i>Chlamydia psittaci</i> genotypes A and C	-	Influenza A/New Caledonia/20/99 (H1N1) virus	-	Human Adenovirus Types 1-5, 8, 15, 31, 40 and 41	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/California/7/2009 (H1N1) virus	-	Human Bocavirus	-
<i>Haemophilus influenzae</i> MinnA	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	Human coronavirus 229E, OC43 and NL63	-
<i>Legionella bozemanii</i>	-	Influenza A/Singapore/GP1908/2015 IVR-180 (H1N1)pdm09 virus	-	MERS Coronavirus	-
<i>Legionella dumoffii</i>	-	Influenza A/Victoria/210/2009 (H3N2) virus	-	Human metapneumovirus A and B	-
<i>Legionella longbeache</i>	-	Influenza A/Thüringen/5/2017 (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Legionella micdadei</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Respiratory Syncytial virus (RSV) A and B	-
<i>Legionella pneumophila</i>	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Human rhinovirus type C	-
<i>Mycoplasma pneumoniae</i>	-				

Table 7. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE Group A *Streptococcus* Real Time PCR Detection Kit was evaluated against *Streptococcus pyogenes* Rosenbach 1884, showing positive results.

ANNEX 1

OPEN FORMAT AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-GAS106L
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-GAS106H
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-GAS112L
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-GAS112H
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-GAS113L
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-GAS113H
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-GAS136
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-GAS172
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-GAS101L
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-GAS101H
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS-GAS101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition.

Target	Channel	Gene
<i>Streptococcus pyogenes</i>	FAM	<i>speB</i> gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.ceratest.es.

A1.2 Reagents provided

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.ceratest.es).

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS101L, VS-GAS101H, VS-GAS106L, VS-GAS106H, VS-GAS112L and VS-GAS112H.

Reagent/Material	Description	Color	Amount
Group A <i>Streptococcus</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref VS-GAS113L and VS-GAS113H.

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS101, VS-GAS136 and VS-GAS172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Group A *Streptococcus* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized Group A *Streptococcus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Group A Streptococcus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Streptococcus pyogenes*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es). In the Applied Biosystems 7500 Fast Real-Time PCR System, StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-GAS196T

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Streptococcus pyogenes</i>	FAM	speB gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A2. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.ceratest.es.

A2.2 Reagents provided

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

Group A *Streptococcus* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized Group A *Streptococcus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Group A Streptococcus* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Group A Streptococcus* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Group A Streptococcus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.cer-test.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*Streptococcus pyogenes*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (to

check the most common detection channels consult website www.certest.es). In the Applied Biosystems 7500 Fast Real-Time PCR System, StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 3

OPEN FORMAT AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-GAS106LE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-GAS106HE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-GAS112LE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-GAS112HE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-GAS113LE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-GAS113HE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-GAS136E
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-GAS172E
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-GAS101LE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-GAS101HE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS-GAS101E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Streptococcus pyogenes</i>	FAM	speB gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A3.2 Reagents provided

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3, A3.4 and A3.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A3.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A3.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es)

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A3. 3. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS101LE, VS-GAS101HE, VS-GAS106LE, VS-GAS106HE, VS-GAS112LE and VS-GAS112HE.

Reagent/Material	Description	Color	Amount
Group A <i>Streptococcus</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A3. 4. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref VS-GAS113LE and VS-GAS113HE.

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A3. 5. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS101E, VS-GAS136E and VS-GAS172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Group A Streptococcus* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

Group A Streptococcus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Group A Streptococcus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Group A Streptococcus* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A3. 6. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*Streptococcus pyogenes*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es). In the Applied Biosystems 7500 Fast Real-Time PCR System, StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Group A <i>Streptococcus</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-GAS196TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Streptococcus pyogenes</i>	FAM	speB gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.ceratest.es.

A4.2 Reagents provided

VIASURE Group A *Streptococcus* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
Group A <i>Streptococcus</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Group A <i>Streptococcus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE Group A *Streptococcus* Real Time PCR Detection Kit with Ref. VS-GAS196TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Group A Streptococcus* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

Group A Streptococcus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Group A Streptococcus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Group A Streptococcus* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Group A Streptococcus* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Group A Streptococcus* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*Streptococcus pyogenes*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es). In the Applied Biosystems 7500 Fast Real-Time PCR System, StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

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Bibliography/Bibliografía

1. Dunne EM *et al.* Detection of group a streptococcal pharyngitis by quantitative PCR. *BMC Infectious Diseases*. 2013 Jul 11;13:312. doi: 10.1186/1471-2334-13-312.
2. Felsenstein S *et al.* Molecular and Clinical Diagnosis of Group A Streptococcal Pharyngitis in Children. *Journal of Clinical Microbiology*. 2014 Nov;52(11):3884-9. doi: 10.1128/JCM.01489-14.
3. Church DL *et al.* Evaluation of Simplexa Group A Strep Direct Kit Compared to Hologic Group A Streptococcal Direct Assay for Detection of Group A Streptococcus in Throat Swabs. *Journal of Clinical Microbiology*. 2018 Feb 22;56(3). pii: e01666-17. doi: 10.1128/JCM.01666-17.

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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