VIASURE

Real Time PCR Detection Kit

CTX, TEM, SHV & mcr

CE IVD

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EN For download IFUS from other languages, please enter in certest.es/viasure/labeling. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use.....	6
2.	Summary and Explanation	6
3.	Principle of the procedure	7
4.	Reagents provided	7
5.	Reagents and equipment to be supplied by the user	7
6.	Transport and storage conditions.....	8
7.	Precautions for users	8
8.	Test procedure	9
8.1.	Specimen collection, transport and storage	9
8.2.	DNA extraction.....	10
9.	Result interpretation	10
10.	Limitations of the test	12
11.	Quality control.....	13
12.	Performance characteristics.....	13
12.1.	Clinical sensitivity and specificity.....	13
12.2.	Analytical sensitivity	15
12.3.	Analytical specificity	17
12.4.	Analytical reactivity	19
	ANNEX 1	20
A1.1	Principle of the procedure.....	20
A1.2	Reagents provided.....	20
A1.3	Test procedure	21
A1.3.1	Lyophilized positive control.....	21
A1.3.2	PCR protocol.....	22
	ANNEX 2	23
A2.1	Principle of the procedure.....	23
A2.2	Reagents provided.....	23
A2.3	Test procedure	23
A2.3.1	Lyophilized positive control.....	23
A2.3.2	Lyophilized reaction mix tube.....	24
A2.3.3	PCR protocol.....	24

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ENGLISH

1. Intended use

VIASURE *CTX*, *TEM*, *SHV* & *mcr* Real Time PCR Detection Kit is a real time PCR assay designed for the qualitative detection and differentiation of family-specific β -lactamases *CTX-M*, *TEM* and *SHV*; and *mcr-1* colistin genes in gram-negative bacteria associated to resistances in blood culture and swab samples such as biopsy, gluteal swabs, ulcer swabs, wound swabs, conjunctive swab, tongue swab, joint fluid swab, oral swabs, hair swabs, perineal swab, urethral swab, abscess swabs from patients with suspect of bacterial infection and/or multi-resistant infection or colonization and swabs for epidemiological control such as nasal, pharyngeal, rectal swabs, axillary swabs, inguinal swabs, epithelial swabs, colostomy swabs, stoma swab, BAS, BAL and sputum specimens by their healthcare professional (HCP). This test is intended to be used as an aid in the diagnosis of family-specific *CTX-M*, *TEM* and *SHV*; and *mcr-1* genes in gram-negative bacteria in clinical samples in both, surveillance and timely identification of antibiotic-resistant infections and colonization, in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens amplified using real time PCR and detected using fluorescent reporter dye probes specific for *CTX-M*, *TEM*, *SHV* and *mcr-1* genes resistances.

2. Summary and Explanation

β -Lactams are the most widely used class of antibiotics. Since the discovery of benzylpenicillin in the 1920s, new penicillin derivatives and related β -lactam classes of cephalosporins, cephamycins, monobactams, and carbapenems have been discovered. Each new class of β -lactam has been developed either to increase the spectrum of activity to include additional bacterial species or to address specific resistance mechanisms that have arisen in the targeted bacterial population. Resistance to β -lactams is primarily because of bacterially produced β -lactamase enzymes that hydrolyze the β -lactam ring, thereby inactivating the drug. The newest effort to circumvent resistance is the development of novel broad-spectrum β -lactamase inhibitors that work against many problematic β -lactamases, including cephalosporinases and serine-based carbapenemases, which severely limit therapeutic options.

On the other hand, only two polymyxins, polymyxin E (colistin) and polymyxin B, are currently commercially available. Colistin has re-emerged as a last-hope treatment in the mid-1990s against multidrug-resistant Gram-negative pathogens due to the development of extensively drug-resistant Gram-negative bacteria. Unfortunately, rapid global resistance towards colistin has emerged following its resurgence. Different mechanisms of colistin resistance have been characterized, including intrinsic, mutational, and transferable mechanisms.

Multidrug-resistant pathogens are a serious problem not only making treatment difficult but also worsening the prognosis of infected patients.

The detection of the common ESBL genes such as *TEM*, *SHV* and *CTX-M* by molecular methods in the ESBL producing bacteria and their patterns of antimicrobial resistance can provide useful information about their epidemiology and can aid a rational antimicrobial therapy.

3. Principle of the procedure

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit is designed for the diagnosis family-specific CTX-M, TEM and SHV genes in gram-negative bacteria associated to resistances in clinical samples. After DNA isolation, the detection of β -lactamases (CTX-M, TEM, SHV) and colistin (*mcr-1*) resistance genes is performed by the amplification of a conserved region of the *Bl_{CTX-M1}*, *Bl_{CTX-M9}*, (*Cluster A or Bl_{CTX-M-A}*), *Bl_{CTX-M-2}*, *Bl_{CTX-M-8}*, *Bl_{CTX-M25}* (*Cluster B or Bl_{CTX-M-B}*), *Bl_{TEM}*, *Bl_{SHV}* genes for β -lactamases resistances and *mcr-1* gene for colistin resistances, using specific primers and fluorescent-labelled probes.

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer, and polymerase) in a stabilized format.

4. Reagents provided

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format products, Annex 2 for tube format products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tube format (Annex 2 and 4).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: -30°C to -10°C and/or \leq -70°C.
- Centrifuge for 1.5mL tubes and PCR well strips or 96-well plate (if available).
- Vortex.
- Micropipetes (0.5-20 μ L, 20-200 μ L).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit has been validated on the following equipments: 7500 Fast Real-Time PCR System (Applied Biosystems), CFX96™ Real-Time PCR Detection System (Bio-Rad), DTprime Real-time Detection Thermal Cycler (DNA-Technology), DTLite Real-time Detection Thermal Cycler (DNA-

Technology), Rotor-Gene® Q (Qiagen), Cobas z 480 Analyzer (Roche Molecular Diagnosis), VIASURE V-Lab96 Cycler (CerTest Biotec S.L.) and AriaMx Real-Time PCR System (Agilent Technologies).

To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the CTX, TEM, SHV & mcr Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at – 20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-BLC113L, VS-BLC113H, VS-BLC136 and VS-BLC172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.

- For references VS-BLC101, VS-BLC136 and VS-BLC172 (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit and any additional reagents or equipment required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink, or smoke in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions, and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format products Test Procedure, Annex 2 for Tube format products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit has been tested in blood culture and swab samples (gluteal, ulcer, wound, conjunctiva, tongue, joint fluid, oral, hair, perineal, urethral, abscess, axillary, inguinal, epithelial and colostomy, nasal, pharyngeal, rectal, stoma swabs, biopsy, BAS, BAL and sputum). Swab

specimens have been collected in AMIES transport medium (DeltaLab). Other types of samples must be validated by the user.

Collection, storage, and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at Room Temperature (RT) for up to 2 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 2 hours). It is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored frozen at -20°C or ideally at -80°C (if medium allows it) for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For details, refer to the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94), and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from clinical samples, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- Invisorb® Spin Universal Kit (Invitex).

9. Result interpretation

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for CTX, TEM, SHV and mcr in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

(ROX) ¹ Controls	CTX resistance gene (<i>Bla</i> _{CTX-M-A} and <i>Bla</i> _{CTX-M-B}) (FAM) ¹	<i>mcr-1</i> resistance gene (Cy5) ¹	SHV resistance gene (<i>Bla</i> _{SHV}) (ROX) ¹	TEM resistance gene (<i>Bla</i> _{TEM}) (HEX) ¹	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	>40 or no signal	>40 or no signal	>40 or no signal	Valid

Table 1. Expected Performance of Controls. Ct values. no signal = no amplification curves.

¹ In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read, and analyse the results:

CTX resistance gene (<i>Bla</i> _{CTX-M-A} and <i>Bla</i> _{CTX-M-B}) (FAM)	<i>mcr-1</i> resistance gene (Cy5)	SHV resistance gene (<i>Bla</i> _{SHV}) (ROX)	TEM resistance gene (<i>Bla</i> _{TEM}) (HEX)	Interpretation for patient's individual samples	
≤37	>40 or no signal	>35 or no signal	>30 or no signal	Valid	CTX resistance gene DNA Detected
≤37	≤40	>35 or no signal	>30 or no signal	Valid	CTX and <i>mcr-1</i> resistance gene DNA Detected
≤37	≤40	≤35	>30 or no signal	Valid	CTX, <i>mcr-1</i> and SHV resistance gene DNA Detected
≤37	≤40	≤35	≤30	Valid	CTX, <i>mcr-1</i>, SHV and TEM resistance gene DNA Detected
≤37	>40 or no signal	≤35	>30 or no signal	Valid	CTX and SHV resistance gene DNA Detected
≤37	>40 or no signal	>35 or no signal	≤30	Valid	CTX and TEM resistance gene DNA Detected
>37 or no signal	≤40	>35 or no signal	>30 or no signal	Valid	<i>mcr-1</i> resistance gene DNA Detected
>37 or no signal	≤40	≤35	>30 or no signal	Valid	<i>mcr-1</i> and SHV resistance gene DNA Detected
>37 or no signal	≤40	>35 or no signal	≤30	Valid	<i>mcr-1</i> and TEM resistance gene DNA Detected
>37 or no signal	≤40	≤35	≤30	Valid	<i>mcr-1</i>, SHV and TEM resistance gene DNA Detected
>37 or no signal	>40 or no signal	≤35	>30 or no signal	Valid	SHV resistance gene DNA Detected
>37 or no signal	>40 or no signal	≤35	≤30	Valid	SHV and TEM resistance gene DNA Detected
≤37	>40 or no signal	≤35	≤30	Valid	CTX, SHV and TEM resistance gene DNA Detected


>37 or no signal	>40 or no signal	>35 or no signal	≤30	Valid	TEM resistance gene DNA Detected
>37 or no signal	>40 or no signal	>35 or no signal	>30 or no signal	Valid	Targets not Detected

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curves.

In case of a continued ambiguous result, it is recommended to review the instructions for use; the extraction process used by the user; to verify the correct performance of each qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- repeat qPCR with the same isolated DNA sample, or
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from blood cultures and swab samples such as gluteal, ulcer, wound, conjunctiva, tongue, joint fluid, oral, hair, perineal, urethral, abscess, axillary, inguinal, epithelial and colostomy, nasal, pharyngeal, rectal, stoma swabs, biopsy, BAS, BAL and sputum.
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by β -lactamases (*CTX-M*, *TEM*, *SHV*) and colistin (*mcr-1*) resistances genes, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).

- Degradation of the bacterial DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown by β -lactamases (*CTX-M*, *TEM*, *SHV*) and colistin (*mcr-1*) resistances genes mutations.
 - A bacterial load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of antibiotics used to prevent the infection or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable bacteria and does not imply that these bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of resistance targets bacteria.
 - Negative results do not preclude by β -lactamases (*CTX-M*, *TEM*, *SHV*) and colistin (*mcr-1*) resistances genes bacteria infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak bacterial levels during infections caused by β -lactamases (*CTX-M*, *TEM*, *SHV*) and colistin (*mcr-1*) resistances bacterias have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the resistance genes.
 - If diagnostic tests for β -lactamases (*CTX-M*, *TEM*, *SHV*) and colistin (*mcr-1*) genes resistances are negative and the patient's clinical presentation and epidemiological information suggest that bacteria resistances infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
 - Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *CTX*, *TEM*, *SHV* & *mcr* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *CTX*, *TEM*, *SHV* & *mcr* Real Time PCR Detection Kit was tested using clinical samples from patients with suspected of bacterial resistance infection and for epidemiological control. In order to determine the clinical diagnostic accuracy, different multicenter evaluations have been conducted through collaboration with national entities. A summary of the sites, sample type and workflow are included in the following table. The results were as follows:

Site	Sample type	Workflow	Target
1 Hospital Royo Villanova (Zaragoza, Spain) and CerTest facilities.	a Blood culture	MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.) + DTprime Real-time Detection Thermal Cycler (DNA-Technology)	CTX
			TEM
			SHV
			<i>mcr-1</i>
	b Swab samples*	MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.) + DTprime Real-time Detection Thermal Cycler (DNA-Technology)	CTX
			TEM
			SHV
			<i>mcr-1</i>

Table 3. Site, sample type, workflow, and target.

* Swabs samples as gluteal, ulcer, wound, conjunctiva, tongue, joint fluid, oral, hair, perineal, urethral, abscess, axillary, inguinal, epithelial and colostomy, nasal, pharyngeal, rectal, stoma swabs, biopsy, BAS, BAL and sputum.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE CTX, TEM, SHV & *mcr* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1a	Culture	CTX	17	401	0	0	1 (0.80-1)	1 (0.99-1)	1 (0.81-1)	1 (0.99-1)
		TEM	54	364	0	0	1 (0.93-1)	1 (0.99-1)	1 (0.93-1)	1 (0.99-1)
		SHV	26	392	0	0	1 (0.86-1)	1 (0.99-1)	1 (0.87-1)	1 (0.99-1)
	ARM-D® Kit, TEM/SHV/GES Detected Variants (STRECK) and ARM-D® Kit, β-Lactamase Detected Variants (STRECK)	CTX	17	68	0	0	1 (0.80-1)	1 (0.94-1)	1 (0.80-1)	1 (0.94-1)
		TEM	54	31	0	0	1 (0.93-1)	1 (0.88-1)	1 (0.93-1)	1 (0.88-1)
		SHV	26	59	0	0	1 (0.86-1)	1 (0.93-1)	1 (0.86-1)	1 (0.93-1)
1b	Culture	CTX	97	548	0	0	1 (0.96-1)	1 (0.99-1)	1 (0.96-1)	1 (0.99-1)
		TEM	181	464	0	0	1 (0.98-1)	1 (0.99-1)	1 (0.97-1)	1 (0.99-1)
		SHV	126	519	0	0	1 (0.97-1)	1 (0.99-1)	1 (0.97-1)	1 (0.99-1)
	ARM-D® Kit, TEM/SHV/GES Detected Variants (STRECK) and ARM-D® Kit, β-Lactamase Detected Variants (STRECK)	CTX	97	548	0	0	1 (0.96-1)	1 (0.99-1)	1 (0.96-1)	1 (0.99-1)
		TEM	181	464	0	0	1 (0.98-1)	1 (0.99-1)	1 (0.97-1)	1 (0.99-1)
		SHV	126	519	0	0	1 (0.97-1)	1 (0.99-1)	1 (0.97-1)	1 (0.99-1)

Table 4. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE CTX, TEM, SHV & *mcr* Real Time PCR Detection Kit.

Results show high agreement to detect β-lactamases (CTX-M, TEM, SHV) and colistin (*mcr-1*) resistances genes using VIASURE CTX, TEM, SHV & *mcr* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit has a detection limit of 0.02 CFU per reaction for β -lactamase resistance gene CTX (*BlaCTX-M1*, *BlaCTX-M9*, (Cluster A or *BlaCTX-M-A*) and *BlaCTX-M2*, *BlaCTX-M8*, *BlaCTX-M25* (Cluster B or *BlaCTX-M-B*)), 0.08 CFU per reaction for β -lactamase resistance gene type TEM (*BlaTEM* gene), 8.75 copies per reaction for β -lactamase resistance gene type SHV (*Bla_{SHV}* gene), and 0.02 CFU per reaction for colistin resistance gene *mcr-1*, with a positive rate of 95%.

Figure 3. Dilution series of CTX-M-1 and CTX-M-9 gene (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).


Figure 4. Dilution series of CTX-M-2, CTX-M-8 and CTX-M-25 gene (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).


Figure 5. Dilution series of *SHV* gene (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel ROX).


Figure 6. Dilution series of *mcr-1* genes (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel Cy5).



Figure 7. Dilution series of *TEM* genes (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel HEX).

12.3. Analytical specificity

The specificity of VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit was confirmed by testing a panel consisting of different microorganisms associated to gastrointestinal diseases, microorganisms transmitted by blood and microorganisms associated with resistances. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing (wet testing)					
Human Adenovirus type 1, 3 4, 5 8, 13, 31, 40 and 41	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>enteritidis</i>	-	Japanese encephalitis	-
<i>Aeromonas caviae</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhimurium</i> . Serotype 4,5,12:i:1,2	-	Japanese encephalitis virus strain Nakayama	-
<i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>gallinarum</i>	-	<i>Plasmodium falciparum</i> strain 3D7	-
<i>Arcobacter butzleri</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>pullorum</i>	-	Rift Valley Fever Virus AR21229	-
Astrovirus genotype I-VIII	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi A</i>	-	Rift Valley Fever Virus MP12	-
<i>Bacteroides fragilis</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi B</i>	-	St Louis Encephalitis Virus	-
<i>Campylobacter coli</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhi</i>	-	<i>Trypanosoma cruzi</i>	-
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	-	Sapovirus	-	<i>Usutu Virus</i>	-
<i>Campylobacter hyointestinalis</i>	-	<i>Serratia liquefaciens</i>	-	West Nile virus strain NY99	-
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-	<i>Shigella dysenteriae</i> serotype 1	-	West Nile virus strain Heja	-
<i>Campylobacter lari</i>	-	<i>Shigella flexneri</i>	-	West Nile virus strain Ug 37	-
<i>Campylobacter upsaliensis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Yellow Fever Virus strain 17D	-
<i>Candida albicans</i>	-	<i>Vibrio parahaemolyticus</i> serotype O1:K1	-	Yellow Fever Virus strain French Neurotropic	-
<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3	-	Zika virus strain FB-GWUH-2016	-
<i>Clostridium difficile</i>	-	<i>Yersinia enterocolitica</i> O:9	-	Zika Virus African strain	-
<i>Clostridium difficile</i> O:27	-	<i>Anaplasma marginale</i>	-	Zika Virus Asian strain PF13/251013-18	-
<i>Clostridium perfringens</i>	-	<i>Bartonella henselae</i> strain Houston-1	-	Zika Virus French Polynesian strain 11474/16	-
Coxsackievirus A24, A9 and B3	-	<i>Borrelia hermsii</i>	-	Zika Virus French Polynesian strain 11468/16	-
<i>Cryptosporidium parvum/hominis</i>	-	<i>Borrelia lusitanae</i>	-	<i>Citrobacter brakii</i> carbapenemase positive (VIM-1)	-
<i>Dientamoba fragilis</i> / <i>Blastocystis hominis</i>	-	<i>Borrelia valaisiana</i>	-	<i>Enterobacter cloacae</i> ESBL/Carbapenemase positive (TEM-1 (non ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL), NDM-1)	+/-

<i>Echovirus 30</i>	-	<i>Borrelia azfeli</i> strain P-Ko/1984	-	<i>Enterococcus avium</i> vancomycin A resistant	-
<i>Entamoeba dispar</i> / <i>Blastocystis hominis</i>	-	<i>Borrelia bavariensis</i>	-	<i>Enterococcus faecalis</i> vancomycin A resistant	-
<i>Entamoeba histolytica</i>	-	<i>Borrelia bisetti</i>	-	<i>Enterococcus faecalis</i> vancomycin B resistant (equivalent to ATCC51299)	-
<i>Enterococcus faecalis</i>	-	<i>Borrelia burgdorferi</i> strain IRS	-	<i>Enterococcus faecium</i> vancomycin A resistant LMG16165	-
<i>Enterovirus 68</i> and <i>71</i>	-	<i>Borrelia burgdorferi sensu stricto</i> strain B31	-	<i>Enterococcus faecium</i> vancomycin A resistant IOWA1	-
Enterohemorrhagic <i>Escherichia coli</i> serotype O157:H7	-	<i>Borrelia garinii</i>	-	<i>Enterococcus faecium</i> S + <i>Enterococcus gallinarum</i> vancomycin C resistant (MI12043391 + LMG16289)	-
Enteroinvasive <i>Escherichia coli</i>	-	<i>Borrelia japonica</i>	-	<i>Enterococcus gallinarum</i> vancomycin B and C resistant ENT20120142	-
Enteropathogenic <i>Escherichia coli</i>	-	<i>Borrelia miyamotoi</i>	-	<i>Escherichia coli</i> carbapenemase positive (OXA-244)	-
Enterotoxigenic <i>Escherichia coli</i> serotype O25:H42	-	<i>Borrelia spielmanii</i>	-	<i>Escherichia coli</i> carbapenemase positive (TEM-1 (non ESBL), IMP-1)	+/-
<i>Giardia intestinalis</i> strain WB clone C6	-	<i>Coxiella burnetii</i> strain Nine Mile Q	-	<i>Helicobacter pylori</i> claritromycin resistant A2146G	-
<i>Helicobacter cinaedi</i>	-	<i>Leptospira</i>	-	<i>Helicobacter pylori</i> claritromycin resistant A2147G	-
<i>Helicobacter heilmannii</i>	-	<i>Rickettsia conorii</i> strain Moroccan	-	<i>Klebsiella pneumoniae</i> carbapenemase positive (SHV-1 (non ESBL), KPC-3, OXA-48)	+/-
<i>Helicobacter hepaticus</i>	-	<i>Theileria annulata</i>	-	<i>Klebsiella pneumoniae</i> ESBL/carbapenemase positive (TEM-1 (non ESBL), SHV-1 (non ESBL), CTX-M-2 (ESBL), KPC-2)	+/-
<i>Helicobacter pylori</i> J99	-	Tick-Borne encephalitis virus strain Neudorf1	-	<i>Klebsiella pneumoniae</i> (TEM, SHV, CTX-M-1, NDM)	+/-
Human Rotavirus A	-	<i>Treponema phagedenis</i>	-	<i>Staphylococcus aureus</i> meticillin resistant N315	-
<i>Klebsiella oxytoca</i>	-	Chikungunya virus S27 Petersfield	-	<i>Staphylococcus aureus</i> meticillin resistant ST398	-
<i>Listeria monocytogenes</i> serovar 1/2c	-	Chikungunya Virus Martinique	-	<i>Staphylococcus aureus</i> meticillin resistant mecC	-
Norovirus GII (II.4 RNA)	-	Chikungunya Virus F24	-	<i>Escherichia coli</i> (Encodes CTX-M-15 ESBL. Plasmid Pek499)	+/-
<i>Proteus vulgaris</i>	-	Chikungunya Virus WHO IS (R91064)	-	<i>Escherichia coli</i> (colistin resistant. MCR-1 positive)	+/-

<i>Pseudomonas aeruginosa</i>	-	Dengue virus type 1 strain Hawaii A	-	Dengue virus type 3 strain H87	-
<i>Salmonella bongori</i> serovar 66:z41	-	Dengue virus type 2 strain New Guinea C	-	Dengue virus type 4 strain H241	-

Table 5. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit for CTX-M β -lactamase resistance gene was evaluated against the following strains (as templates), showing positive results:

- *Enterobacter cloacae* ESBL/Carbapenemase positive (SHV-12 (ESBL), **CTX-M-9** (ESBL), OXA-48).
- *Enterobacter cloacae* ESBL/Carbapenemase positive (TEM-1 (non ESBL), SHV-12 (ESBL), **CTX-M-15** (ESBL), NDM-1).
- *Klebsiella pneumoniae* ESBL/carbapenemase positive (TEM-1 (non ESBL), SHV-1 (non ESBL), **CTX-M-2** (ESBL), KPC-2).
- *Escherichia coli* (encodes **CTX-M-15** ESBL. Plasmid pEK499) (NCTC 13400).
- *Klebsiella pneumoniae* (TEM, SHV, **CTX-M-1**, NDM).

The reactivity of VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit for TEM β -lactamase resistance gene was evaluated against the following strains (as templates), showing positive results:

- *Enterobacter cloacae* ESBL/Carbapenemase positive (**TEM-1** (non ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL), NDM-1).
- *Escherichia coli* carbapenemase positive (**TEM-1** (non ESBL), IMP-1).
- *Klebsiella pneumoniae* ESBL/carbapenemase positive (**TEM-1** (non ESBL), SHV-1 (non ESBL), CTX-M-2 (ESBL), KPC-2).
- *Escherichia coli* (encodes CTX-M-15 ESBL. **Plasmid pEK499**) (NCTC 13400).

The reactivity of VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit for SHV β -lactamase resistance gene was evaluated against the following strains (as templates), showing positive results:

- *Enterobacter cloacae* ESBL/Carbapenemase positive (**SHV-12** (ESBL), CTX-M-9 (ESBL), OXA-48).
- *Enterobacter cloacae* ESBL/Carbapenemase positive (TEM-1 (non ESBL), **SHV-12** (ESBL), CTX-M-15 (ESBL), NDM-1).
- *Klebsiella pneumoniae* carbapenemase positive (**SHV-1** (non ESBL), KPC-3, OXA-48).
- *Klebsiella pneumoniae* ESBL/carbapenemase positive (TEM-1 (non ESBL), **SHV-1** (non ESBL), CTX-M-2 (ESBL), KPC-2).

The reactivity of VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit for Colistin resistance was evaluated against the following strain (as template), showing positive results:

- *Escherichia coli* colistin resistant (MCR-1 positive) (NCTC 13846).

ANNEX 1

OPEN FORMAT AND ROTOR-GENE FORMAT

Annex for the following references:

PRODUCT	REFERENCE
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS- BLC101L
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS- BLC101H
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS- BLC106L
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS- BLC106H
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS- BLC112L
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS- BLC112H
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 96-well plate, low profile	VS- BLC113L
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 96-well plate, high profile	VS- BLC113H
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS- BLC101
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS- BLC136
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS- BLC172

Table A1 1. References

A1.1 Principle of the procedure

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format.

Target	Channel	Gene
CTX resistance gene	FAM	<i>Bla_{CTX-M1}</i> , <i>Bla_{CTX-M9}</i> , (Cluster A or <i>Bla_{CTX-M-A}</i>) and <i>Bla_{CTX-M2}</i> , <i>Bla_{CTX-M8}</i> , <i>Bla_{CTX-M25}</i> (Cluster B or <i>Bla_{CTX-M-B}</i>)
<i>mcr-1</i> resistance gene	Cy5	<i>mcr-1</i> gene
SHV resistance gene	ROX	<i>Bla_{SHV}</i> gene
TEM resistance gene	HEX, VIC or JOE*	<i>Bla_{TEM}</i> gene

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.cerTEST.es.

A1.2 Reagents provided

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.cerTEST.es).

Reagent/Material	Description	Colour	Amount
CTX, TEM, SHV & mcr 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CTX, TEM, SHV & mcr Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit with Ref. VS-BLC101L, VS-BLC101H, VS-BLC106L, VS-BLC106H, VS-BLC112L and VS-BLC112H.

Reagent/Material	Description	Colour	Amount
CTX, TEM, SHV & mcr 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers a in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CTX, TEM, SHV & mcr Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit with Ref VS-BLC113L and VS-BLC113H.

Reagent/Material	Description	Colour	Amount
CTX, TEM, SHV & mcr 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CTX, TEM, SHV & mcr Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	2/9/18 X 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit with Ref. VS-BLC101, VS-BLC136 and VS-BLC172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

CTX, TEM, SHV & mcr Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized CTX, TEM, SHV & mcr Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *CTX*, *TEM*, *SHV* & *mcr* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*CTX* gene resistance), ROX (*SHV* gene resistance), Cy5 (*mcr-1* gene resistance) and HEX, JOE or VIC (*TEM* resistance gene). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT

Annex for the following references:

PRODUCT	REFERENCE
VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS- BLC196T

Table A2. 1.References.

A2.1 Principle of the procedure

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format.

Target	Channel	Gene
CTX resistance gene	FAM	<i>Bla_{CTX-M1}</i> , <i>Bla_{CTX-M9}</i> , (Cluster A or <i>Bla_{CTX-M-A}</i>) and <i>Bla_{CTX-M2}</i> , <i>Bla_{CTX-M8}</i> , <i>Bla_{CTX-M25}</i> (Cluster B or <i>Bla_{CTX-M-B}</i>)
<i>mcr-1</i> resistance gene	Cy5	<i>mcr-1</i>
SHV resistance gene	ROX	<i>Bla_{SHV}</i> gene
TEM resistance gene	HEX, VIC or JOE*	<i>Bla_{TEM}</i> gene

Table A2. 2.Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.ceratest.es.

A2.2 Reagents provided

VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
CTX, TEM, SHV & mcr Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CTX, TEM, SHV & mcr Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE CTX, TEM, SHV & mcr Real Time PCR Detection Kit with Ref. VS-BLC196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

CTX, TEM, SHV & mcr Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized CTX,

TEM, *SHV* & *mcr* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *CTX*, *TEM*, *SHV* & *mcr* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *CTX*, *TEM*, *SHV* & *mcr* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *CTX*, *TEM*, *SHV* & *mcr* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*CTX* gene resistance), ROX (*SHV* gene resistance), Cy5 (*mcr-1* gene resistance) and HEX, JOE or VIC (*TEM* resistance gene). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

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Bibliography/Bibliografía

- Centers for Disease Control and Prevention (<https://www.cdc.gov/>)
- World Health Organization (<http://www.who.int/en/>)
- Roschanski N, Fischer J, Guerra B, Roesler U. Development of a multiplex real-time PCR for the rapid detection of the predominant beta-lactamase genes CTX-M, SHV, TEM and CIT-type AmpCs in Enterobacteriaceae. PLoS One. 2014 Jul 17;9(7):e100956. doi: 10.1371/journal.pone.0100956. PMID: 25033234; PMCID: PMC4102473.
- Irrgang A, Roschanski N, Tenhagen BA, Grobbel M, Skladnikiewicz-Ziemer T, Thomas K, Roesler U, Käsbohrer A. Prevalence of mcr-1 in E. coli from Livestock and Food in Germany, 2010-2015. PLoS One. 2016 Jul 25;11(7):e0159863. doi: 10.1371/journal.pone.0159863. PMID: 27454527; PMCID: PMC4959773.
- Tolosi R. Rapid detection and quantification of plasmid-mediated colistin resistance genes (mcr-1 to mcr-5) by real-time PCR in bacterial and environmental samples. Journal of Applied Microbiology. 2020 2 June. doi:10.1111/jam.14738.
- Cantón R, González-Alba JM, Galán JC. CTX-M enzymes: origin and diffusion. Frontiers in microbiology. 2012 April. doi: 10.3389/fmicb.2012.00110.
- El-Sayed Ahmed, M. A. E. G., Zhong, L. L., Shen, C., Yang, Y., Doi, Y., & Tian, G. B. (2020). Colistin and its role in the Era of antibiotic resistance: an extended review (2000-2019). *Emerging Microbes & Infections*, 9(1), 868–885. <https://doi.org/10.1080/22221751.2020.1754133>
- Bush, K., & Bradford, P. A. (2016). β -lactams and β -lactamase inhibitors: An overview. *Cold Spring Harbor Perspectives in Medicine*, 6(8). <https://doi.org/10.1101/cshperspect.a025247>
- Occurrence of the CTX-M, SHV and the TEM Genes Among the Extended Spectrum β -Lactamase Producing Isolates of Enterobacteriaceae in a Tertiary Care Hospital of North India. *J Clin Diagn Res*. 2013 Apr; 7(4): 642–645. Published online 2013 Apr 1. doi: [10.7860/JCDR/2013/5081.2872](https://doi.org/10.7860/JCDR/2013/5081.2872)

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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