VIASURE

Real Time PCR Detection Kit

CMV q

CE

0318

IVD

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Content

1.	Intended use.....	5
2.	Summary and Explanation	5
3.	Principle of the procedure	6
4.	Reagents provided	6
5.	Reagents and equipment to be supplied by the user	6
6.	Transport and storage conditions.....	7
7.	Precautions for users	7
8.	Test procedure	9
8.1.	Specimen collection, transport and storage	9
8.2.	DNA extraction	9
8.3.	Reference workflow of VIASURE CMV q Real Time PCR Detection Kit for Conversion Factor (CF) calculation	10
8.4.	Standard curve preparation for quantitative assay	10
9.	Result interpretation	11
10.	Limitations of the test	14
11.	Quality control.....	15
12.	Performance characteristics.....	15
12.1.	Clinical sensitivity and specificity	15
12.2.	Linear range	17
12.3.	Analytical sensitivity	18
12.4.	Lower Limit of Quantification (LLOQ)	19
12.5.	Analytical specificity.....	19
12.6.	Analytical reactivity	19
	ANNEX 1	20
A1.1	Principle of the procedure.....	20
A1.2	Reagents provided.....	20
A1.3	Test procedure.....	21
A1.3.1	Conversion factor protocol	21
A1.3.2	Lyophilized quantitative standards	22
A1.3.3	PCR protocol.....	22
	ANNEX 2	24
A2.1	Principle of the procedure.....	24
A2.2	Reagents provided.....	24
A2.3	Test procedure	25
A2.3.1	Conversion factor protocol	25
A2.3.2	Lyophilized quantitative standards	25
A2.3.3	Lyophilized reaction mix tube.....	26
A2.3.4	PCR protocol.....	26

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ENGLISH

1. Intended use

VIASURE CMV q Real Time PCR Detection Kit is a real-time PCR test designed for the detection and quantitative measurement of Cytomegalovirus (CMV) DNA in human EDTA-plasma samples. This test is intended for use as an aid in the universal screening of the CMV infection, in the diagnosis and management of CMV in solid organ transplant patients and in hematopoietic stem cell transplant patients. The test can be used in these populations to assess the need to initiate antiviral treatment, and in patients receiving anti-CMV therapy, serial DNA measurements can be used to assess viral response to treatment. The results from the VIASURE CMV q Real Time PCR Detection Kit must be interpreted within the context of all relevant clinical and laboratory findings. DNA is extracted from clinical specimens, amplified using real-time PCR, and detected using fluorescent reporter dye probes specific for CMV.

2. Summary and Explanation

Cytomegalovirus (CMV) is a wide-spread double-stranded DNA virus which belongs to the viral family known as herpesviruses, *Herpesviridae*, or human herpesvirus-5 (HHV-5). Human cytomegalovirus infections may be asymptomatic in healthy people, but it can be life-threatening in an immunocompromised patient. It is possible to be transmitted transplacentally to neonates, or through breast milk of an infected and shedding mother (congenital CMV infection), as well as by intimate contact and by transplantation from (or sharing syringes with) an infected individual. Like other herpesviruses, after recovery of the initial infection, CMV often remains latent, but it can reactivate at any time, during the compromise of the immune system with immunosuppression. Of all herpesviruses, CMV harbours the largest number of genes dedicated to evading innate and adaptative immunity in the host. CMV infects between 60% to 70% of adults in industrialized countries and close to the 100% in emerging countries.

In some cases, infection in healthy people can cause mild illness that may include: fever, sore throat, fatigue or swollen glands. Occasionally, CMV can cause mononucleosis or hepatitis (liver problems). People with weakened immune system who get CMV can have more serious symptoms affecting the eyes, lungs, liver, esophagus, stomach, and intestines.

Congenital CMV infection is a leading infectious cause of deafness, learning disabilities, and intellectual disability in infants worldwide. Congenital CMV infection occurs when a baby is born with CMV infection, and it is estimated that one out of every 200 babies is born infected with CMV. CMV-positive babies can have brain, liver, spleen, lung and growth problems, but the most common long-term health problem in babies born with congenital CMV infection is hearing loss.

The diagnosis of human cytomegalovirus infection is confirmed by demonstrating the presence of infectious virus, viral antigens, or viral DNA in clinical samples (including blood and other body fluids) from infected patients. Serologic tests, highlighting ELISA test, are commonly used for diagnosis of CMV infection for people older than 12 months, but PCR-based assays are currently considered the standard diagnostic method to confirm congenital CMV infection in newborns.

3. Principle of the procedure

VIASURE CMV *q* Real Time PCR Detection Kit is designed for the detection and quantitative measurement of CMV DNA in human EDTA-plasma samples. After DNA isolation, the detection of CMV is performed by the amplification of a conserved region of the *pp65* and *gB* genes, using specific primers and fluorescent-labelled probes.

VIASURE CMV *q* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE CMV *q* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an **endogenous internal control** to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an endogenous internal control (EIC) (human *RNase P* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

In addition, this kit contains 5 lyophilized vials (CMV1Q1 to CMV1Q5) with different standardized concentrations of CMV specific DNA (CMV *q* Quantitative Standards) for the quantification of CMV viral DNA in clinical samples. These quantitative standards are calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162), ensuring the quality of the results and the comparison with results from other laboratories.

4. Reagents provided

VIASURE CMV *q* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open format with internal control products and Annex 2 for tube format with internal control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE CMV *q* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

Optional:

- 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code 09/162), or another reference material (in case the final user needs to calculate its own conversion factor when using a different workflow with VIASURE CMV q Real Time PCR Detection Kit).

VIASURE CMV q Real Time PCR Detection Kit has been validated with different workflows, although the MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad) workflow in plasma samples (hereinafter referred to as VIASURE CMV q reference workflow in this document) was the one considered for the analytical/clinical validation, and therefore, for the interpretation of results.

The different workflows validated are listed below:

	Matrix	Workflow
1 (recommended)	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)
2	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + VIASURE V-Lab96 Cyclor (CerTest Biotec S.L)
3	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + 7500 Fast Real-Time PCR System (Applied Biosystems)
4	Plasma	Maxwell® RSC 16 Viral TNA/ Maxwell® RSC (Promega) + CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
5	Plasma	QIAamp® MiniElute® Virus Spin Kit (QIAGEN) + CFX96™ Real-Time PCR Detection System (Bio-Rad)

Table 1. Sample matrix extraction-amplification methods (workflows) validated using VIASURE CMV q Real Time PCR Detection Kit.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the CMV q Quantitative Standards has been re-suspended, they are stable at room temperature during the time needed to dispense a 96-well plate, but they should be stored at 4 °C for up to 5 months, or at -20° C or lower for longer period of time. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Quantitative Standards has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the CMV q Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.

- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-CMV113L and VS-CMV113H). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- **Do not mix CMV q Quantitative Standards from different kits and / or lots.** Each quantitative standard set has been validated for each lot of master mix with which it is supplied.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE CMV q Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Use calibrated micropipettes and proper tips with VIASURE CMV q Real Time PCR Detection Kit and stick to test procedure of this instructions of use. Inaccurate volumes due to imprecise pipetting or malpractices may result in an erroneous quantification.
- Standard curve should be included in the same run of samples to be quantify and should be generated by running each of CMV q Quantitative Standard in triplicate.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.

- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.
- When using 1st WHO International Standard for human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC 09/162), please note that the stability of this standard can be affected by some use practices. Therefore, it is recommended to reconstitute the standard and prepared the aliquots maintaining the cold chain. Additionally, aliquots should be for single use only since freeze-thawing cycles affect the integrity of the reference material.

8. Test procedure

Please, see Annex 1 for Open format with internal control products Test Procedure and Annex 2 for Tube format with internal control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE CMV q Real Time PCR Detection Kit has been tested in human EDTA-plasma clinical samples. Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at room temperature for up to 2 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 2 hours), we recommend shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For details, refer to the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94).

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

VIASURE CMV q Real Time PCR Detection Kit has been validated with EDTA-plasma samples extracted with MagDEA Dx SV Kit / magLEAD® 12gC instrument (Precision System Science Co.) and amplified on CFX96™ Real-Time PCR Detection System (Bio-Rad).

However, although they are not considered as reference workflows, other extraction-amplification workflows have been validated when using EDTA-plasma samples (see table 1 in section 5).

8.3. Reference workflow of VIASURE CMV q Real Time PCR Detection Kit for Conversion Factor (CF) calculation

Before CMV quantification from clinical samples, an extraction-amplification procedure (workflow) must be decided in advance and maintained for further analysis. Because CMV concentration must be expressed in International Units (IU), and the quantitative results are obtained in copies per mL after interpolation in the linear calibration curve, a conversion factor (CF) is required to be calculated in order to relate copies/mL values with IU/mL data. **The conversion factor (CF) depends on the type of sample tested and the extraction and amplification procedure used (workflow).**

VIASURE CMV q Real Time PCR Detection Kit has been validated using EDTA-plasma samples extracted with MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) and run on CFX96™ Real-Time PCR Detection System (Bio-Rad). The conversion factor value was calculated using the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162) and established at **1.263 IU/copies**.

$$\text{IU/mL} = \text{copies/mL} \times 1.263$$

For instance, a CMV viral load determined at 10,000 copies/mL in a sample would correspond to 12,630 IU/mL.

Other workflows were validated using VIASURE CMV q Real Time PCR Detection Kit. The CF value obtained for each one is summarized in tables A1.5 and A2.4 (see Annexes 1 and 2, respectively), but please note that these CF values are applicable only for EDTA-plasma samples as sample matrix. Nevertheless, **each laboratory must establish its own conversion factor**.

If you wish to have a quick guide for how to calculate your CF value and an excel template to perform it, visit the website www.certest.es.

NOTE: It has been noticed that the **stability of WHO International Standard can be affected by some use practices**. Therefore, it is recommended to reconstitute the standard and prepare the aliquots maintaining the cold chain. Additionally, aliquots should be for single use only since freeze-thawing cycles affect the integrity of the reference material.

8.4. Standard curve preparation for quantitative assay

This kit contains 5 lyophilized vials (CMV1Q1 to CMV1Q5) with different standardized concentrations of CMV specific DNA (CMV q Quantitative Standards) for the quantification of CMV viral DNA in clinical samples. These Quantitative Standards were calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162) and can be used to generate a standard curve with which to determine the amount of CMV viral DNA in a sample in IU/mL. Lyophilized CMV q Quantitative Standard panel of VIASURE CMV q Real Time PCR Detection Kit are ready to use after being reconstituted in RNase/DNase free Water.

The CMV q Quantitative Standards have the following concentrations:

CMV q Quantitative Standard	Concentration (copies/mL)
CMV1Q1	5×10^6
CMV1Q2	5×10^5
CMV1Q3	5×10^4
CMV1Q4	5×10^3
CMV1Q5	5×10^2


Table 2. CMV q Quantitative Standards and concentrations.

To perform a quantitative assay, take CMV q Quantitative Standard set (aluminium pouches, red vials), and prepare the standard curve as follows:

- Rehydrate CMV q Quantitative Standard set. Pipette 100 μ L of RNase/DNase free Water in each lyophilized CMV q Quantitative Standard, changing the micropipette tip between vials.
- Vortex each vial for 20 seconds.

CMV q Quantitative Standards should be manipulated in an area separated from the extraction procedure, and it is recommended to have them ready to be added for when samples are to be loaded.

Figure 1. Preparation of CMV q Quantitative Standard set for quantitative assay.

9. Result interpretation

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (EIC) signal to verify the extraction and correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the CMV q Quantification Standards amplification curves as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the

fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of CMV q Quantitative Standards and negative control in each run validate the reaction by checking the absence of signal in the negative control well and the presence of signal (i.e. exponential) in CMV q Quantitative Standards wells.

For a valid diagnostic test run, the following control conditions must be met:

Controls	CMV (FAM) ¹	Endogenous Internal Control (EIC) (HEX) ²	Interpretation of Controls
CMV q Quantitative Standards (CMV1Q1-CMV1Q5)	+	-	Valid
Negative Control (NC)	-	-	Valid

Table 2. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the CMV q Quantitative Standards wells for the target channel), all results are reported as 'Invalid' and retesting is required.

2 Due to human housekeeping *RNase P* gene target is present in nucleated human cells, amplification signal of Endogenous Internal Control (EIC) is not expected in CMV q Quantitative Standards.

Assessment of clinical samples test results should be performed after the CMV q Quantitative Standards and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

A standard curve can be generated from the dilution series of CMV q Quantitative Standards (CMV1Q1-5) following the formula:

$$Ct = m \log (Q) + b$$

Where: Ct = Threshold Cycle; m = Slope; Q = Concentration; and b = Intercept.

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{(Ct-b)/m}$$

The obtained quantification result is reported in the concentration units of the calibrators set.

NOTE: Standard curve should be included in the same run of samples to be quantify and should be generated by running each of CMV q Quantitative Standard in triplicate.

The quantification results are valid if the generated standard curve reaches the following control parameter values:

- Efficiency: 80 - 120%
- R^2 : ≥ 0.98

The concentration of your "Positive sample" refers to the concentration in the eluted DNA after extraction, not to the original clinical sample. **To determine the target concentration of the original sample, consider the dilutions of the extraction procedure and PCR Set-up.**

For interpretation of individual patient sample results, use the following table, read and analyzed the results if VIASURE CMV q reference workflow (MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)) has been performed:

CMV (FAM)	Endogenous internal control (HEX)	Interpretation of patients' individual samples	
< 298.12 IU/mL	≤40 or no signal ¹	Valid	Non-quantifiable CMV positive sample
298.12 to 631.5 IU/mL	≤40 or no signal ¹	Valid	Quantifiable CMV positive sample by extrapolation ²
631.5 to 631.5 x 10 ⁴ IU/mL	≤40 or no signal ¹	Valid	Quantifiable CMV positive sample in the quantification range
> 631.5 x 10 ⁴ IU/mL	≤40 or no signal ¹	Valid	Quantifiable CMV positive sample by extrapolation ³
No signal	≤ 35 ³	Valid	Target not detected
No signal	>35 or no signal ³	No valid	Test Failure – Repeat Testing ⁴

Table 3. Interpretation of individual patient sample results. No signal = no amplification curves.

1 The Endogenous Internal Control (EIC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 CMV target genes positive below the quantification range let quantification of CMV viral load by extrapolation of the standard curve, but results could not be accurate.


3 CMV target genes positive above the quantification range let quantification of CMV viral load by extrapolation of the standard curve, but results could not be accurate. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100.

4 In the case of CMV target genes negative, EIC must show an amplification signal with Ct ≤35. The Ct value could be very variable due to the endogenous internal control is a human housekeeping gene that should be present in all human nucleated cells in the original samples. If there is an absence of signal or Ct value > 35 of the endogenous internal control, the result is considered as 'invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use; the correct performance of each workflow steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- a) repeat qPCR with the same isolated DNA sample, or
- b) re-extract and retest another aliquot of the same specimen or,
- c) obtain a new specimen and retest.

Figure 3. Standard curve from the dilution series (blue) and quantification of a positive sample of unknown concentration (orange) run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with DNA extracted EDTA-plasma samples.
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- Lower levels of target below the low limit of quantification might be detected, but results may not be accurately quantifiable.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by CMV, either the high number of DNA template copies which contains each CMV q Quantitative Standard vial, samples containing high concentrations of target DNA, or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the viral DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown CMV variants.
 - A viral load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent infection of CMV or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- Inaccurate viral DNA results (loss of viral load) may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Improper maintenance of commonly used equipment, especially micropipettes and extraction/amplification systems.
 - Use of uncalibrated micropipettes and/or improper tips.

- Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targeted viral sequences.
- Negative results do not preclude CMV infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by CMV have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other immunodepressive-related illnesses are negative and the patient's clinical presentation and epidemiological information suggest that CMV infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.
- The conversion factor (CF) depends on the type of sample tested and the extraction and amplification instruments used (workflow). Therefore, each sample type-extraction-amplification workflow requires the calculation of its own CF.

11. Quality control

VIASURE CMV q Real Time PCR Detection Kit contains CMV q Quantitative Standards and a negative control that must be included in each run to correctly interpret and quantify the results. Also, the endogenous internal control (EIC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE CMV q Real Time PCR Detection Kit was tested using plasma-EDTA clinical samples from patients who underwent to a universal screening of the CMV infection, for the analysis and management of CMV in solid organ transplant patients, or in hematopoietic stem cell transplant patients before the antiviral treatment. The study also included patients that have been received anti-CMV therapy.

In order to determine the clinical diagnostic accuracy, an evaluation has been conducted through collaboration with national and international entities. A summary of the sites, sample type and workflow are included in the following table:

	Site	Sample type	Workflow	Target
1	CerTest Biotec S.L (Zaragoza, Spain) in collaboration with Cerba Xpert (Saint Ouen L'Aumone, France)	EDTA-Plasma	MagDEA® Dx SV Kit using the magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	CMV
2	CerTest Biotec S.L (Zaragoza, Spain) in collaboration with Cerba Xpert (Saint Ouen L'Aumone, France)	EDTA-Plasma	QIAamp® MiniElute® Virus Spin Kit (QIAGEN) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	CMV

Table 4. Site, sample type, workflow, and target.

Correlation between extraction methods using VIASURE CMV q Real Time PCR Detection Kit with comparator assays were assessed as shown in the following figures:


Figure 2. Correlation between the QIAamp® MiniElute® Virus Spin Kit and MagDEA® Dx SV Kit / magLEAD® 12gC extraction method results using VIASURE CMV q Real Time PCR Detection Kit.


Figure 3. Correlation between VIASURE CMV q Real Time PCR Detection Kit and RealStar® CMV 1.0 Kit (Altona Diagnostics) results.


Figure 4. Correlation between VIASURE CMV q Real Time PCR Detection Kit and RealStar® CMV 1.0 Kit (Altona Diagnostics) results using QIAamp® MiniElute® Virus Spin Kit (QIAGEN).


Figure 5. Correlation between VIASURE CMV q Real Time PCR Detection Kit using QIAamp® MiniElute® Virus Spin Kit (QIAGEN) and cobas® CMV (Roche Diagnostics) results.


Figure 5. Correlation between VIASURE CMV q Real Time PCR Detection Kit and cobas® CMV (Roche Diagnostics) results using MagDEA® Dx SV Kit / magLEAD® 12gC and quantifiable samples.


Figure 6. Correlation between VIASURE CMV q Real Time PCR Detection Kit and cobas® CMV (Roche Diagnostics) results using MagDEA® Dx SV Kit / magLEAD® 12gC and samples over VIASURE LoD (231 IU/mL).

Lastly, it was possible to calculate the clinical sensitivity and specificity and the negative and positive predictive values (with 95% confidence intervals) using VIASURE CMV q Real Time PCR Detection Kit. Briefly, when a sigmoid curve was observed in a clinical sample when analyzed with VIASURE assay, it was considered positive, and samples with absence of curve were negative.

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	cobas® CMV (Roche Diagnostics)	CMV	120	153	0	0	1 (0.97-1)	1 (0.97-1)	1 (0.96-1)	1 (0.97-1)

Table 5. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE CMV q Real Time PCR Detection Kit, using MagDEA® Dx SV Kit / magLEAD® 12gC + CFX96™ Real-Time PCR Detection System (Bio-Rad).


Results show high agreement to detect and quantify CMV using VIASURE CMV q Real Time PCR Detection Kit.

12.2. Linear range

The linearity of VIASURE CMV q Real Time PCR Detection Kit was determined using specific and synthetic DNA belonging to CMV, and further confirmed with reference material in plasma samples.

First, a series of ten-fold dilutions containing a known concentration (ranging from 10^7 to 10^1 copies/ μ L) of specific CMV DNA (CMV1Q) was amplified, and Ct values in all dilution tested were obtained.

Figure 7. Dilution series of CMV1Q (10^7 - 10^1 copies/ μ L) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).

Afterwards, linear range was further confirmed using negative plasma samples spiked with a known concentration of CMV reference material (1st WHO IS for human Cytomegalovirus NAT (NIBSC code 09/162) and NATtrol™ Cytomegalovirus Linearity Panel (Zeptometrix Corporation, Ref: NATCMV-LIN).

Using VIASURE CMV q Real Time PCR Detection Kit, amplifications in all dilutions tested were obtained. End-point fluorescence level was maintained, and no crosstalk was observed. Assays with 1st WHO IS for human Cytomegalovirus and NATtrol™ Cytomegalovirus Linearity Panel showed statistical significance.

Figure 8. Dilution series of 1st WHO IS for human CMV template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).



Figure 9. Dilution series of NATtrol™ Cytomegalovirus Linearity Panel template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).

12.3. Analytical sensitivity

To determine the LoD of the product, a dilution series of the 1st WHO International Standard for human CMV NAT (NIBSC code 09/162) and a tentative LoD concentration of NATtrol™ Cytomegalovirus Strain AD-169 (Zeptometrix, Ref: NATCMV-0002) were set up.

A total of 71 replicates per concentration level of 1st WHO IS for human CMV reference material and 89 replicates of the NATtrol™ Cytomegalovirus Strain AD-169 reference material were extracted with the MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) and amplified in CFX96™ Real-Time PCR Detection System (Bio-Rad). Using the 1st WHO IS for human CMV NAT (NIBSC code 09/162) and Probit Model to analyse the results, VIASURE CMV q Real Time PCR Detection Kit showed a detection limit of 236.04 copies /mL, equivalent to 298.12 IU/mL for CMV, with a positive rate of $\geq 95\%$. Using NATtrol™ Cytomegalovirus Strain AD-169 material, VIASURE CMV q Real Time PCR Detection Kit showed a detection limit of 100 copies/mL, equivalent to 504 IU/mL for CMV, with a positive rate of $\geq 95\%$.

CMV concentration (copies/mL)	CMV concentration (IU/mL)	N° of replicates	N° positive replicates	Hit rate (%)
500	631.5	71	71	100
400	505.2	71	71	100
300	378.9	71	70	98.6
200	252.6	71	63	88.7

Table 6. Analytical sensitivity (limit of detection or LoD) evaluated with negative plasma samples spiked with the 1st WHO International Standard for human CMV NAT (NIBSC code 09/162) and VIASURE CMV q Real Time PCR Detection Kit, (IU) = International Units.

CMV concentration (copies/mL)	CMV concentration (IU/mL)	N° of replicates	N° positive replicates	Hit rate (%)
100	504	89	85	95.5

Table 7. Analytical sensitivity (limit of detection or LoD) evaluated with negative plasma samples spiked with the NATtrol™ Cytomegalovirus Strain AD-169 (Zeptometrix, Ref: NATCMV-0002) and VIASURE CMV q Real Time PCR Detection Kit, (IU) = International Units.

12.4. Lower Limit of Quantification (LLOQ)

The lower limit of quantification (LLOQ) is the lowest concentration value that CMV can be detected and at which the total error meets the acceptance criterium of being below 40%. The LLOQ was determined to be 500 copies/mL, equivalent to 631.5 IU/mL, for CMV.

12.5. Analytical specificity

The specificity of the VIASURE CMV q Real Time PCR Detection Kit was confirmed by testing a panel consisting of different microorganisms associated to common disorders in transplant recipients or immunocompromised individuals. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing (wet testing)					
Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Hepatitis A virus	-	<i>Plasmodium falciparum</i> strain 3D7	-
<i>Acinetobacter baumannii</i>	-	Herpes simplex virus 1 (HSV-1) strain MacIntyre	-	<i>Pneumocystis jirovecii</i> type A1 and g885652	-
<i>Aspergillus fumigatus</i>	-	Herpes simplex virus 2 (HSV-2) strain MS	-	<i>Pseudomonas aeruginosa</i>	-
<i>Bacteroides fragilis</i>	-	Human herpesvirus 6 (HHV6) strain Z29 and types A and B	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
BK virus types 1b-2 and IV	-	JC virus types 1A and 2B	-	<i>Toxoplasma gondii</i> type II	-
<i>Candida albicans</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Treponema pallidum</i>	-
<i>Enterobacter aerogenes</i> serotype Cloaca B	-	<i>Listeria monocytogenes</i> serotype 1/2b and serovar 4b	-	Varicella-Zoster Virus strain Ellen	-
<i>Enterobacter cloacae</i> serotype Cloaca A	-	<i>Listeria innocua</i> serotype 6a	-	<i>Yersinia enterocolitica</i> O:3	-
Epstein-Barr virus	-	<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	-	Parvovirus B19 genotypes 1, 2 and 3a	-
<i>Escherichia coli</i> strain 0.1285;O18:H7:K1	-				

Table 8. Reference pathogenic microorganisms used in this study.

12.6. Analytical reactivity

The reactivity of VIASURE CMV q Real Time PCR Detection Kit was evaluated against DNA from Cytomegalovirus, strains Merlin and AD-169 (as template), showing positive results.

ANNEX 1

OPEN FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE CMV q Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-CMV106L
VIASURE CMV q Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-CMV106H
VIASURE CMV q Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-CMV112L
VIASURE CMV q Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-CMV112H
VIASURE CMV q Real Time PCR Detection Kit 96-well plate, low profile	VS-CMV113L
VIASURE CMV q Real Time PCR Detection Kit 96-well plate, high profile	VS-CMV113H

Table A1 1. References

A1.1 Principle of the procedure

VIASURE CMV q Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an Endogenous Internal Control (EIC) (human *RNase P* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

In addition, this kit contains 5 lyophilized vials (CMV1Q1 to CMV1Q5) with different standardized concentrations of CMV specific DNA (CMV q Quantitative Standards) for the quantification of CMV viral DNA in clinical samples. These quantitative standards are calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162), ensuring the quality of the results and the comparison with results from other laboratories.

Target	Channel	Gene
CMV	FAM	<i>pp65</i> and <i>gB</i>
Endogenous Internal control (EIC)	HEX, VIC or JOE *	human <i>RNase P</i>

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.certest.es.

A1.2 Reagents provided

VIASURE CMV q Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3 and A1.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices.

Reagent/Material	Description	Colour	Amount
CMV q 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CMV q Quantitative Standards (CMV1Q1-CMV1Q5)	Non-infectious synthetic lyophilized DNA	Red	5 (1 vial per concentration)
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE CMV q Real Time PCR Detection Kit with Ref. VS-CMV106L, VS-CMV106H, VS-CMV112L and VS-CMV112H.

Reagent/Material	Description	Colour	Amount
CMV q 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CMV q Quantitative Standards (CMV1Q1-CMV1Q5)	Non-infectious synthetic lyophilized DNA	Red	5 (1 vial per concentration)
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE CMV q Real Time PCR Detection Kit with Ref VS-CMV113L and VS-CMV113H.

A1.3 Test procedure

A1.3.1 Conversion factor protocol

CMV q Quantitative Standards has been calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162) and can be used to generate a standard curve with which to determine the CMV viral DNA in a sample in International Units (IU) per mL after calculating the conversion factor (CF). **The conversion factor (CF) depends on the type of sample tested and the extraction and amplification instruments used (workflow).** Therefore, the combination for which the conversion factor is to be tested must be decided in advance and maintained in all corresponding analyses.

The reference workflow established to use with VIASURE CMV q Real Time PCR Detection Kit is EDTA-plasma samples extracted with **MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) and run on CFX96™ Real-Time PCR Detection System (Bio-Rad)**. If this workflow is performed, the concentration of samples in IU/mL can be calculated using the CF value **1.263 IU/copies** as shown in the following formula:

$$\text{IU/mL} = \text{copies/mL} \times 1.263$$

Other alternative workflows were validated using VIASURE CMV q Real Time PCR Detection Kit. The CF value obtaining for each one is summarized in table A1.5:

	Matrix	Workflow	Conversion Factor (CF) (IU/copies)
1 (recommended)	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	1.263
2	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + VIASURE V-Lab96 Cyclor (CerTest Biotec S.L)	1.03
3	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + 7500 Fast Real-Time PCR System (Applied Biosystems)	1.30
4	Plasma	Maxwell® RSC 16 Viral TNA/ Maxwell® RSC (Promega) + CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)	2.12
5	Plasma	QIAamp® MiniElute® Virus Spin Kit (QIAGEN) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	1.86

Table A1 5. Calculation of conversion factor (CF) with different extraction-amplification methods using plasma samples spiked with a known concentration of 1st WHO International Standard for CMV DNA (NIBSC 09/162). (IU) = International Units.

If another workflow is followed, then CF value must be calculated by the user. If you wish to have a quick guide for how to calculate your CF value and an excel template to perform it, visit the website www.certest.es.

A1.3.2 Lyophilized quantitative standards

CMV q Quantitative Standards contain high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized CMV q Quantitative Standards (red vials) by adding 100 µL of the supplied Water RNase/DNase free (white vial) into each vial and vortex thoroughly for 20 seconds.

Once the CMV q Quantitative Standards has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. At least one replicate of each CMV q Quantitative Standard and one replicate of negative control must be included in each run for each assay. However, standard curve should be generated by running each of CMV q Quantitative Standard in triplicate, as well as to be included in the same run of samples to be quantified.

Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, each reconstituted CMV q Quantitative Standards (CMV1Q1-5) (aluminium pouches, red vials), and Negative Control (violet vial) and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler.

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (CMV) and HEX, JOE or VIC (Endogenous Internal Control (EIC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE CMV q Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-CMV196T

Table A2. 1.References.

A2.1 Principle of the procedure

VIASURE CMV q Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an Endogenous Internal Control (EIC) (human *RNase P* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

In addition, this kit contains 5 lyophilized vials (CMV1Q1 to CMV1Q5) with different standardized concentrations of CMV specific DNA (CMV q Quantitative Standards) for the quantification of CMV viral DNA in clinical samples. These quantitative standards are calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162), ensuring the quality of the results and the comparison with results from other laboratories.

Target	Channel	Gene
CMV	FAM	<i>pp65</i> and <i>gB</i>
Endogenous Internal control (EIC)	HEX, VIC or JOE *	human <i>RNase P</i>

Table A2 6. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.cerfest.es.

A2.2 Reagents provided

VIASURE CMV q Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
CMV q Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
CMV q Quantitative Standards (CMV1Q1-CMV1Q5)	Non-infectious synthetic lyophilized DNA	Red	5 (1 vial per concentration)
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 2. Reagents and materials provided in VIASURE CMV q Real Time PCR Detection Kit with Ref. VS-CMV196T.

A2.3 Test procedure

A2.3.1 Conversion factor protocol

CMV q Quantitative Standards has been calibrated against the 1st WHO International Standard for Human Cytomegalovirus for Nucleic Acid Amplification Techniques (NIBSC code: 09/162) and can be used to generate a standard curve with which to determine the CMV viral DNA in a sample in International Units (IU) per mL after calculating the conversion factor (CF). **The conversion factor (CF) depends on the type of sample tested and the extraction and amplification instruments used (workflow).** Therefore, the combination for which the conversion factor is to be tested must be decided in advance and maintained in all corresponding analyses.

The reference workflow established to use with VIASURE CMV q Real Time PCR Detection Kit is EDTA-plasma samples extracted with **MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) and run on CFX96™ Real-Time PCR Detection System (Bio-Rad)**. If this workflow is performed, the concentration of samples in IU/mL can be calculated using the CF value **1.263 IU/copies** as shown in the following formula:

$$\text{IU/mL} = \text{copies/mL} \times 1.263$$

Other alternative workflows were validated using VIASURE CMV q Real Time PCR Detection Kit. The CF value obtaining for each one is summarized in table A2.4:

	Matrix	Workflow	Conversion Factor (CF) (IU/copies)
1 (recommended)	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	1.263
2	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + VIASURE V-Lab96 Cyclor (CerTest Biotec S.L)	1.03
3	Plasma	MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) + 7500 Fast Real-Time PCR System (Applied Biosystems)	1.30
4	Plasma	Maxwell® RSC 16 Viral TNA/ Maxwell® RSC (Promega) + CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)	2.12
5	Plasma	QIAamp® MiniElute® Virus Spin Kit (QIAGEN) + CFX96™ Real-Time PCR Detection System (Bio-Rad)	1.86

Table A2 4. Calculation of conversion factor (CF) with different extraction-amplification methods using plasma samples spiked with a known concentration of 1st WHO International Standard for CMV DNA (NIBSC 09/162). (IU) = International Units.

If another workflow is followed, then CF value must be calculated by the user. If you wish to have a quick guide for how to calculate your CF value and an excel template to perform it, visit the website www.certest.es.

A2.3.2 Lyophilized quantitative standards

CMV q Quantitative Standards contain high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized CMV q Quantitative Standards (red vials) by adding 100 µL of the supplied Water RNase/DNase free (white vial) into each vial and vortex thoroughly for 20 seconds.

Once the CMV q Quantitative Standards has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (at least one replicate of each CMV q Quantitative Standard and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the CMV q Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.4 PCR protocol

Determine and separate the number of required reactions including samples and controls. At least one replicate of each CMV q Quantitative Standard and one replicate of negative control must be included in each run for each assay. However, standard curve should be generated by running each of CMV q Quantitative Standard in triplicate, as well as to be included in the same run of samples to be quantified.

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated CMV q Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, each reconstituted CMV q Quantitative Standards (CMV1Q1-5) (aluminium pouches, red vials), and Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler.

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 5. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (CMV) and HEX, JOE or VIC (Endogenous Internal Control (EIC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

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Bibliography/Bibliografía

1. CMV Infection Laboratory Testing | CDC. (n.d.). Retrieved April 6, 2022, from <https://www.cdc.gov/cmvp/clinical/lab-tests.html>
2. Gupta, M., & Shorman, M. (2021). Cytomegalovirus. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK459185/>
3. Plotkin, S. A., & Boppana, S. B. (2019). Vaccination against the human cytomegalovirus. *Vaccine*, 37, 7437–7442. <https://doi.org/10.1016/j.vaccine.2018.02.089>
4. Stockdale, L., Nash, S., Nalwoga, A., Painter, H., Asiki, G., Fletcher, H., & Newton, R. (2018). Human cytomegalovirus epidemiology and relationship to tuberculosis and cardiovascular disease risk factors in a rural Ugandan cohort. *PLoS ONE*, 13(2), e0192086. <https://doi.org/10.1371/journal.pone.0192086>
5. MM6-A Quantitative Molecular Methods for Infectious Diseases; Approved Guideline CLSI.
6. EP05-A Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline CLSI — Second Edition.
7. EP17-A Protocols for determination of limits of detection and limits of quantitation; Approved Guideline CLSI
8. WHO manual for the preparation of secondary reference materials for in vitro diagnostic assays designed for infectious disease nucleic acid or antigen detection: calibration to WHO International Standards, Annex 6, TRS No 1004

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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