

# HbA1c (GHb) Calibrators

**REF** 01-04-0018

**CAL** 1 x 400 µL Level 1 Calibrator  
1 x 400 µL Level 2 Calibrator

Pour d'autres langues  
Für andere Sprachen  
Para otras lenguas  
Per le altre lingue  
Dla innych języków

Para outras línguas  
Για τις άλλες  
Λόγους  
För andra språk  
For andre språk



[www.trinitybiotech.com](http://www.trinitybiotech.com)

## INTENDED USE

The Hemoglobin A1c Calibrator is intended for use as a quality control material to monitor the precision of laboratory testing procedures for HbA1c quantitation. The Calibrators are for the calibration of quantitative HbA1c affinity assay systems. For *in vitro* diagnostic use only. **IVD**

## SUMMARY AND EXPLANATION OF THE TEST

Recognition of the usefulness of HbA1c in the control of diabetes has led to a substantial demand for its measurement by clinical laboratories.<sup>1,2,3,4,5</sup> Reliable use of the data generated necessitates good accuracy and precision in the HbA1c measurement. Factors such as pH, ionic strength, temperature, dilution and column equilibration can affect the performance of the assay. In order to assure the accuracy and precision of this assay, high quality calibrator material should be used. The Trinity Biotech HbA1c (GHb) Calibrators were developed to assist with these requirements.

Good laboratory practice provides a system of management controls to ensure the consistency and reliability of analytical results (e.g., standard operating procedures, uniform sample collection and handling practices, performance training of personnel, statistical evaluation of control results, proper storage of test kits and controls, a permanent record of control results, etc.). Use of external control samples assures the test reagents are working properly, that the Trinity Biotech Affinity HbA1c analyzer is calibrated properly, and that the operator has performed the test correctly. If the controls do not perform as expected (i.e., within given acceptance ranges), review the instructions for use to see if the test was performed correctly; repeat the test or contact Trinity Biotech Technical Service before testing patient samples.

The staff at each laboratory site will benefit by establishing a quality assurance plan, based on their institution's policies. Run quality control specimens, for example, under the following conditions:

- At regular intervals determined by the laboratory procedures
- With each new shipment of reagents and with each new lot of reagents
- Each time a calibration is performed
- To train and confirm performance acceptability for new analysts
- When results do not match the patient's clinical condition or symptoms

Good laboratory practices include a well-designed and implemented quality control process. These practices, for example, may involve:

- Proper storage and handling of reagents kits
- Careful sample collection and handling procedures
- Training of testing personnel
- Routine review of sample and control results
- Periodic quality system reviews
- Retention of quality control testing records

If the problem cannot be corrected, or the reason for an out-of-limits result cannot be determined, contact Trinity Biotech or the Trinity Biotech Distributor nearest you.

## PRINCIPAL OF THE PROCEDURE





Utilizing the principles of boronate affinity and high-performance liquid chromatography (HPLC), the assay separates and quantitates Glycated Hemoglobin in whole blood and hemolysates of whole blood. Glycated Hemoglobin is separated from non-Glycated Hemoglobin which is then quantitated with an HbA1c result generated based on that ratio.

## REAGENTS / COMPONENTS

1 vial lyophilized Trinity Biotech HbA1c (GHb) Level 1 Calibrator  
1 vial lyophilized Trinity Biotech HbA1c (GHb) Level 2 Calibrator  
1 Package insert

The Trinity Biotech Level 1 and Level 2 HbA1c (GHb) Calibrators included in the HbA1c (GHb) Calibrator Kit are prepared from non-diabetic and glycated non-diabetic whole blood. Each sample provides a clear, cherry red hemolysate containing oxyhemoglobin with low and high HbA1c values as indicated in the assay data section. Once reconstituted, each vial's content is 400 µL (microliters) of HbA1c (GHb) Calibrator material.

## STORAGE AND STABILITY

	<b>Lyophilized:</b> Lyophilized vials of Trinity Biotech HbA1c (GHb) Calibrator stored at 2°C – 8°C are stable until the expiration date on the label.
	<b>Reconstituted:</b> Once reconstituted, the Trinity Biotech HbA1c (GHb) Calibrator is stable at 2°C – 8°C for 30 days.
	<b>Diluted:</b> Diluted samples of Trinity Biotech HbA1c (GHb) Calibrator are stable for 8 hours at 20°C – 25°C.
	<b>DO NOT USE</b> after the expiration date.

## PRECAUTIONS



**CAUTION:** For *In Vitro* Diagnostic **IVD** Use ONLY

**SAFETY GLASSES, GLOVES AND LAB COAT ARE RECOMMENDED WHEN USING THE TRINITY BIOTECH HbA1c (GHb) CALIBRATORS.**

**DO NOT USE:** If diluted sample turns dark brown.



## POTENTIALLY BIOHAZARDOUS MATERIAL

Human sourced materials were used in the manufacturing of this product. This product was found to be non-reactive for hepatitis B surface antigen (HBsAg), antibodies to hepatitis C (HCV), and antibodies to Human Immunodeficiency Viruses (HIV-1 and HIV-2), when tested by FDA cleared methods. No known test method can offer assurance that products derived from human blood will not transmit disease, and material should be handled as such.

## PREPARATION PROCEDURE

### RECONSTITUTION

1. Tap bottom of vial gently to settle material in vial.
2. Open the vial.
3. Add 400µL of Trinity Biotech Diluent to the vial. Trinity Biotech Diluent **REF** 01-03-0013 (DIL, 940mL), 01-03-0056 (2DIL, 3.8L), 01-03-0059 (2DIL, 940mL), 01-03-0066 (PDQ DIL), 01-03-0101 (Premier Hb9210™ Sample Diluent, 940 mL), or 01-03-0097 (Premier Hb9210™ DIL Reagent, 3.8L). Replace cap. Allow the vial to stand for ten minutes, then rotate gently until completely dissolved.
4. For further dilution, treat the HbA1c (GHb) calibrator as whole blood.

### DILUTIONS

Following instructions below, dilute the reconstituted HbA1c (GHb) Calibrators using Trinity Biotech Diluent in the same manner required in the assay for a patient sample. The reconstituted HbA1c (GHb) Calibrators have the same hemoglobin concentration as whole blood.

System/Injector with Injection Volume (µL)	Dilution Ratio	Typical Dilution µL Calibrator : µL Diluent	Vial Type
Premier Hb9210™ (5µL)	1:150	10:1490	Untreated Test Tube
PDQ / PDQ PLUS (10µL)	1:100	20:1980	Untreated PDQ Tube
ultra <sup>2</sup> w/215 (20µL)	1:200	5:995 (or 8:1592)	Shell or Crimp Top

Prior to making each calibrator dilution, rotate the HbA1c (GHb) Calibrator vial gently and ensure material is uniformly mixed and in solution.

## TEST PROCEDURE

After reconstitution and dilution of the HbA1c (GHb) Calibrator material, it should be analyzed as part of the system calibration.

## RESULTS AND INTERPRETATION OF RESULTS

These Trinity Biotech HbA1c (GHb) Calibrator materials have been assigned values utilizing NGSP and/or IFCC HbA1c reference materials. According to these reference materials, each Trinity Biotech HbA1c system has been assigned a system-specific calibrator set point to optimize accuracy. When assayed by the Trinity Biotech HPLC Affinity method, the set points on each Trinity Biotech analyzer are as follows:

## Trinity Biotech HbA1c (GHb) Calibrator Set Points:

Trinity Biotech HbA1c (GHb) Calibrator Kit <b>LOT</b> 6550 <b>EXP 2017-04-30</b>			
	Units	Calibrator 1 <b>LOT</b> 6551	Calibrator 2 <b>LOT</b> 6552
<b>Premier Hb9210™</b>			
HbA1c(NGSP) <sup>6</sup>	%	5.2	12.2
HbA1c(IFCC) <sup>7</sup>	mMol HbA1c / Mol Hb	33	110
<b>PDQ</b>			
HbA1c(NGSP) <sup>6</sup>	%	5.2	11.8
HbA1c(IFCC) <sup>7</sup>	mMol HbA1c / Mol Hb	33	105
<b>ultra<sup>2</sup></b>			
HbA1c(NGSP) <sup>6</sup>	%	5.4	12.3
HbA1c(IFCC) <sup>7</sup>	mMol HbA1c / Mol Hb	36	111

Users of other methods, including HPLC ion-exchange, should determine their own values.

**Premier Hb9210™ Primary Reporting Method Barcodes**

Please ensure you scan the relevant barcode set below applicable to your Primary Reporting Method. Note that if you are dual-reporting on the Premier Hb9210™ you do not need to scan the barcodes for the secondary reporting method.

**For %HbA1c Values (NGSP/DCCT Referenced)**

Premier Hb9210™ HbA1c (GHb) Calibrator 1 Barcode, lot 6551  yrn*:4{n}	
Premier Hb9210™ HbA1c (GHb) Calibrator 2 Barcode, lot 6552  Zrn_:dY3	

**For mmol HbA1c/mol Hb Values (IFCC Referenced)**

Premier Hb9210™ HbA1c (GHb) Calibrator 1 Barcode, lot 6551  yrn*:4{n}	
Premier Hb9210™ HbA1c (GHb) Calibrator 2 Barcode, lot 6552  Zrn_:dY3	

**Validation of Assigned Values**

The Premier Hb9210™ NGSP & IFCC set points detailed for this HbA1c (GHb) Calibrator kit have been independently verified by the IFCC Network Laboratory of the Queen Beatrix Hospital in Winterswijk, The Netherlands (Dr. Cas Weykamp).

**LIMITATIONS**

- This product should not be used past the expiration date.
- If there is evidence of microbial contamination, brown color or excessive turbidity in the reconstituted HbA1c (GHb) Calibrator, discard the vial.
- This product is not intended for use as a standard.

	<b>Important Information</b>	Prior to running HbA1c (GHb) calibrations following column changes, please verify that the baseline is smooth and quiet prior to beginning calibration. Do not proceed if excessive noise is present. Please refer to the system Operator's Manual chapter titled 'Results and Interpretation' for additional information regarding baseline verification and chromatography verification checks.
--	------------------------------	---

**REFERENCES**

1. Trivelli, *New Engl. J. of Med.*, 284, 353, (1971).
2. Gabbay, K.H., *et al.*, *J. Clin. End & Metab.*, 44, 859, (1977).
3. Dix, D. *et al. Clin. Chem.*, 25, 877, (1979).
4. Abraham, E.D. *Diabetes*, 27, 950, (1978).
5. Goldstein, D.E. *et al. Diabetes*, 50, 70, (1982).
6. Goldstein, D.E. *et al.*, in *Methods in Diabetes Research, Vol II*, 475-504 (1986).
7. Hoelzel W, *et al.*, *Clin.Chem.* 2004; 50: 166-174

**ORDERING INFORMATION**

Catalog No.	Item	Quantity
01-04-0018	Kit, HbA1c (GHb) Calibrators	1 x 400 µL Calibrator 1 1 x 400 µL Calibrator 2

Manufactured  
**EC REP**  
Authorized Representative  
**REF**  
Product Number  
**LOT**  
Lot  
**EXP**  
Use by YYYY-MM-DD  
8°C  
2°C  
Store Lyophilized Material at 2-8°C  
25°C  
20°C  
Store Diluted Material at 20-25°C

European Conformity  
  
Consult accompanying documents  
**CAL**  
Calibrator  
**IVD**  
For In Vitro Diagnostic use  
  
Caution, consult accompanying documents  
8°C  
2°C  
Store Reconstituted Material at 2-8°C  
  
Biohazard

Rev J 04/15

Trinity Biotech  
Kansas City, MO 64132  
Tel. 1 800-325-3424  
Fax: 1 816-361-1974

**EC REP**  
Trinity Biotech plc  
Bray Co. Wicklow, Ireland  
Tel. 353 1 2769800  
Fax 353 1 2769888  
www.trinitybiotech.com



*Handwritten signature and date:*  
12/10/16  
P.H.