**CHORUS
VIT B12**

REF 86400

DIESSE

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INSTRUCTIONS FOR USE

CHORUS VIT B12

For the quantitative determination of vitamin B12 in human serum

For *in vitro* diagnostic use only

1. USE

CHORUS VIT B12 is an immunological kit for the automated quantitative determination of vitamin B12.

The test is conducted on human serum using a single-use device applied to the CHORUS TRIO instrument.

The kit, intended to monitor the level of vitamin B12, is an aid in the diagnosis of vitamin B12 deficiency.

It must only be used by professional laboratory personnel.

2. INTRODUCTION

Vitamin B12 is one of the nine water-soluble vitamins, which intervenes in man in the formation of red blood cells and myelin sheaths around the nerves. Vitamin B12 cannot be synthesized in the human body; the main sources of absorption are meat, fish, eggs, and dairy products. Absorption in the gastro-intestinal tract depends on the intrinsic factor to which the vitamin binds to maintain its stable form. The most frequent cause of severe vitamin B12 deficiency is the lack of intrinsic factor due to autoimmune atrophic gastritis. The disease is historically called "pernicious anaemia" even though many patients have mainly neurological manifestations. Other causes of vitamin B12 deficiency are malabsorption due to gastrectomy, inflammatory bowel disease or dietary deficiency. In general, the prevalence increases with age due to a worsening of absorption in the gastro-intestinal tract. Symptoms of vitamin B12 deficiency can be very ambiguous and your deficiency may manifest itself even after months or years. Diagnosis is essential because of the latent nature of this disorder and the risk of permanent neurological damage. Generally, the primary test performed to confirm the diagnosis of vitamin B12 deficiency is the measurement of the serum level of vitamin B12. In addition, other biomarkers such as folate, methylmalonic acid (MMA), homocysteine and holotranscobalamin should also be evaluated.

3. PRINCIPLE OF THE METHOD

The CHORUS VIT B12 device is ready to use for the determination of vitamin B12 in human serum in the Chorus TRIO instrument.

The test is based on a competition ELISA method.

The microplate wells are coated with streptavidin. The calibrator, positive control and samples are added to a device. The sample and biotinylated anti-B12 antibody are added to a plate and during a first incubation, an immunocomplex between any vitamin B12 present in the serum and the biotinylated anti-B12 antibody is formed, which in turn binds to the immobilized streptavidin on the solid phase.

After a short incubation, the conjugate (containing vitamin B12 marked with peroxidase) is added (this delay allows increasing the sensitivity for low-concentration samples). The conjugate competes with the vitamin B12 present in the serum under examination for the binding sites to the anti-B12 antibody. After washing to eliminate the excess of unbound reagents, the peroxidase substrate is added. The enzymatic reaction is stopped by adding stop solution, which causes the solution to turn to yellow. The colour that develops is inversely proportional to the amount of vitamin B12 in the patient sample.

The single-use devices contain all the reagents to perform the test when applied to the Chorus TRIO instrument.

The results are expressed in pg/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin that have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. Follow all the precautions normally adopted in laboratory practice when handling material of human origin.

Disposal of residues: the serum samples, calibrators and strips used must be treated as infected residues and disposed of in accordance with the provisions of current laws.

Health and safety information

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. An exposure of 30 minutes to 1% sodium hypochlorite may be necessary to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present in the spillage, dry the contaminated area before using sodium hypochlorite. All materials used to decontaminate any spillage, including gloves, must be disposed of as potentially infected waste. Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) for at least 30 minutes and use within 60 minutes.


1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed on the bottom.
3. Check the actual presence of the reagents in the device and the integrity of the device itself; do not use devices that lack any reagents and/or have foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.
The kit can only be used with CHORUS TRIO with an updated version of the software. Make sure that the software installed on the instrument is a version (Rel.) equal to or later than the one shown in the table published on the Diesse website
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode on the grip of the device so that it can correctly be read by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Any defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipaemic or icteric samples of not fully coagulated serum or samples with microbial contamination may be a source of error.
11. Do not use the device after its expiry date.
12. **Check that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packs of 6 devices each

Description:

Position 9: Space available for barcode label

Position 8: Empty

Position 7: MICROPLATE WELL

Coated with streptavidin

Position 6: Empty

Position 5: BIOTINYLATED VITAMIN B12

Content: biotinylated anti-vitamin B12 antibodies stabilized in buffer and preservative.

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine and hydrogen peroxide stabilized in buffer.

Position 3: STOP SOLUTION

Content: 0.3 M sulphuric acid solution

Position 2: VIT B12 CONJUGATE

Content: Vitamin B12 marked with peroxidase stabilized in buffer

Position 1: EMPTY WELL

to transfer the sample to

Use: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR **1 x 0.450 ml**
Content: Human serum containing vitamin B12 and preservative.
Liquid, ready to use.

CONTROL + POSITIVE CONTROL **1 x 0.650 ml**
Content: Human serum containing vitamin B12 and preservative.
Liquid, ready to use.

SAMPLE DILUENT/NEGATIVE CONTROL
SAMPLE DILUENT/NEGATIVE CONTROL **1 x 1.100 ml**
Content: Human serum not containing vitamin B12 and preservative.
Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED

- CHORUS VITAMIN EXTRACTION **REF** 86500
- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Vortex
- Micropipettes capable of accurately withdrawing volumes of 50-200 µl.
- Single-use gloves
- Sodium hypochlorite solution (5%)
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see section 9: Test validation).

The expiry date is printed on each component and on the kit label.

The reagents have limited stability after opening:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C
SAMPLE DIL/NEG CTRL	8 weeks at 2-8°C

7. SAMPLE COLLECTION AND STORAGE

The sample is serum obtained from blood collected by normal venous blood sampling and handled with the precautions required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C;

for longer storage periods, freeze at temperatures $\leq -20^{\circ}\text{C}$.

The sample can be thawed for a maximum of 3 times. Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

The extracted sample can be used within 8 hours from extraction if kept at 4°C.

It is strongly recommended to avoid any intensive sample processing.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

The CHORUS VIT B12 test (REF. 86400) can be performed simultaneously ONLY with the CHORUS PROMONITOR tests.

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Prepare the samples, calibrator, positive control and negative control (when necessary) using the product CHORUS VITAMIN EXTRACTION (REF 86500) according to the extraction protocol given in the relative instructions for use.
Pretreatment cannot be performed in the CHORUS AUTOSAMPLER instrument (REF 81210).
4. Dispense into well no. 1 of each device:

SAMPLE	100 µl/device
CALIBRATOR	100 µl/device
POSITIVE CONTROL	100 µl/device

At each batch change use one device for the calibrator.

5. Place the devices on the CHORUS TRIO instrument. Perform the calibration (if required) and the test as indicated in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the positive control has a value outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the positive control result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF RESULTS

The Chorus TRIO instrument provides the result in pg/ml calculated on the basis of a batch-dependent graph stored in the instrument.

11. LIMITATIONS

All values obtained require a careful interpretation without ever disregarding other indicators relating to the same patient. In fact, the test cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. CALIBRATION RANGE

Calibration range: 150.0-2000.0 pg/ml.

For samples >2000.0 pg/ml, repeat the test diluting the sample with Vit B12 Sample Diluent/Negative Control provided with the kit. It is advisable to make 1:5 and 1:10 dilutions of the off-scale samples; then proceed with extraction of the samples before analysis.

13. REFERENCE RANGE

120 samples (56 men, 64 women) from an apparently healthy population between the ages of 21 and 89 years were analysed. Pregnant women were excluded.

The reference range of the Chorus VIT B12 kit is 237.5 - 764.8 pg/ml.

These values should only be used as guidelines.

Each laboratory should investigate the transferability of the expected values to their patient population and, if necessary, determine its own reference range.

14. ANALYTICAL SPECIFICITY

3 samples with increasing concentration of vitamin B12 were tested to which the following interferents were added:

Bilirubin (4.5 - 45 mg/dl)
Triglycerides (250 - 1500 mg/dl)
Haemoglobin (2.5 - 10 mg/ml)
Rheumatoid factor (44 - 220 IU/ml)

The presence of the above listed interfering substances in the serum under examination does not alter the test result, with the exception of haemoglobin.

1 sample with vitamin B12 concentration in the normal range was tested to which increasing concentrations of Biotin (0.5 - 1000.0 ng/ml) were added.

The presence of Biotin in the serum under examination does not alter the test result at concentrations up to 35 ng/ml. Higher Biotin concentrations determine incorrect results.

15. COMPARISON STUDIES

In one trial, 167 samples were analysed with a Diesse kit and another commercially available kit: the correlation data are summarized in the following table:

Correlation	R	95% CI
Pearson	0.96	0.94-0.97
Spearman	0.89	0.85-0.92

Pearson and Spearman's correlation coefficients respectively fall within the very high and high correlation ranges.

16. PRECISION AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (pg/ml)	CV%	Mean (pg/ml)	CV%
1	1876.2	3.1	1681.8	13.1
2	912.8	1.8	872.9	0.8
3	540.0	2.6	480.5	4.4
4	405.2	1.6	457.1	1.2
5	308.8	4.5	349.3	6.7
6	254.0	2.2	216.1	11.3

Sample	Between batches		Between instruments	
	Mean (pg/ml)	CV%	Mean (pg/ml)	CV%
1	1758.1	3.5	1752.0	12.3
2	908.9	8.1	840.0	13.0
3	484.5	3.0	540.9	10.4
4	473.5	12.2	485.8	8.4
5	365.3	10.8	336.2	11.8
6	204.4	13.9	216.6	11.3

17. LINEARITY

The linearity of the test was evaluated using 3 samples with known concentrations of vitamin B12. Serial dilutions were made for each sample, which demonstrated the linearity of the 150.0-2000.0 pg/ml method.

18. ACCURACY

The accuracy of the test was evaluated by measuring the International Standard Vitamin B12, Serum Folate, HOLO TC (NIBSC code 03/178). Two kit batches were tested on one instrument. The mean recovery with respect to the expected standard value (480 pg/ml) is 111%.

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