CERTEST

Influenza A+B+ *RSV+Adeno Resp.*

ONE STEP *Influenza type A,*
Influenza type B, Respiratory
Syncytial Virus and Adenovirus
Respiratory
COMBO CARD TEST

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ENGLISH

INTENDED USE

CerTest *Influenza A+B+RSV+Adeno Resp.* one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Influenza* type A and type B, *Respiratory Syncytial Virus* and *Adenovirus* antigens from nasopharyngeal swab, nasopharyngeal wash or aspirate specimens.

CerTest *Influenza A+B+RSV+Adeno Resp.* combo card test is a simple and a highly sensitive screening assay to make a presumptive diagnosis of *Influenza* type A and B, *Respiratory Syncytial Virus* and/or *Adenovirus* respiratory infection. The clinical presentation of infections caused by the heterogeneous group of the respiratory viruses can be very similar. Thus, the implementation of this qualitative assay that rapidly identifies the most important viruses involved shows a great interest.

INTRODUCTION

Influenza is caused by a virus that attacks mainly the upper respiratory tract - the nose, throat and bronchi and rarely also the lungs. The infection usually lasts for about a week. It is characterized by sudden onset of high fever, myalgia, headache and severe malaise, non-productive cough, sore throat and rhinitis. Most people recover within one to two weeks without requiring any medical treatment. In the very young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, *Influenza* poses a serious risk. In these people, the infection may lead to severe complications of underlying diseases, pneumonia and death.

The currently circulating *Influenza* viruses that cause human disease are divided into two groups: A and B. *Influenza* A has 3 subtypes which are important for humans: A (H3N2), A (H1N1) and A (H5N1), of which the former is currently associated with most deaths. *Influenza* viruses are defined by 2 different protein components, known as antigens, on the surface of the virus. They are spike-like features called haemagglutinin (H) and neuraminidase (N) components.

Respiratory Syncytial Virus (RSV) is the most common cause of bronchiolitis and pneumonia among infants and children under 1 year of age. Illness begins most frequently with fever, runny nose, cough and sometimes wheezing. Severe lower respiratory tract disease may occur at any age, especially among the elderly or among those with compromised cardiac, pulmonary or immune systems. *RSV* is spread from respiratory secretions through close contact with infected persons or contact with contaminated surfaces or objects.

Adenoviruses most commonly cause respiratory illness; however, depending on the infecting serotype, they may also cause various other illnesses, such as gastroenteritis, conjunctivitis, cystitis and rash illness. Symptoms of respiratory illness caused by *Adenovirus* infection range from the common cold syndrome to pneumonia, croup and bronchitis. Patients with compromised immune systems are especially susceptible to severe complications of *Adenovirus* infection. Adenovirus is transmitted by direct contact, fecal-oral transmission and occasionally waterborne transmission. Some types are capable of establishing persistent asymptomatic infections in tonsils, adenoids and intestines of infected hosts and shedding can occur for months or years.

TEST PRINCIPLE

CerTest *Influenza A+B+RSV+Adeno Resp.* is based on the principle of a qualitative immunochromatographic assay for the determination of *Influenza* type A in the strip A, *Influenza* type B in the strip B, *Respiratory Syncytial Virus* in the strip C and *Adenovirus* in the strip D from nasopharyngeal swab, nasopharyngeal wash or aspirate specimens.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza* type A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-type A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza* type B and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-type B) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip C consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *RSV* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*RSV*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip D consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Adenovirus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Adenovirus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

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If the sample is *Influenza A* positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-type A monoclonal antibodies-red polystyrene microspheres) in the strip A and if the sample is *Influenza B* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-type B monoclonal antibodies-red polystyrene microspheres) in the strip B, which were previously pre-dried on the absorbent pad. If the sample is *RSV* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*RSV* monoclonal antibodies-red polystyrene microspheres) in the strip C and if the sample is *Adenovirus* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*Adenovirus* monoclonal antibodies-red polystyrene microspheres) in the strip D, which were previously pre-dried on the absorbent pad. The mixtures then move upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-type A antibodies present on the membrane of strip A (test line), the anti-type B antibodies present on the membrane of strip B (test line), the anti-*RSV* antibodies present on the membrane of strip C (test line) and the anti-*Adenovirus* antibodies present on the membrane of strip D (test line) capture the coloured conjugate and a red line will be visible in the strips. These bands are used to interpret the result.

If the sample is negative, there is no *Influenza* type A or type B, *RSV* and *Adenovirus* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Influenza* type A, the anti-*Influenza* type B, the anti-*RSV* and the anti-*Adenovirus* antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the four strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on the four membranes will capture control green-conjugate complex and the four control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Influenza A+B+RSV+Adeno Resp.* combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.
6. Take out the cap of the sample collection tube with Reagent B (1).
7. Introduce the swab into the sample collection tube (2) and mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution. Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
8. Close the tube with the reagent and sample. Shake the tube in order to assure good sample dispersion. Mix 60s (3).

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

Specimen collection

(1) Take out the cap of sample collection tube

(2) Introduce the swab into the tube, rotating 1 minute and extract the liquid

(3) Mix 60s

Nasopharyngeal aspirate method (suction apparatus, sterile suction catheter):

For adult:

1. Place the irrigator up to the nose.
2. Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side. Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
3. Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

For children:

1. Use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head.
2. Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
3. Repeat for the other nostril and transfer the fluid into the same specimen container.

Sample procedure for both samples (adult and children):


4. Add 6 drops (1) of the nasal wash or aspirate samples with a pipette and 9 drops (2) of Reagent B in a testing tube. Mixer with vortex at least 1 minute (3) to homogenize.

Specimen collection

(1) Add 6 drops of nasal wash

(2) Add 9 drops of Reagent B

(3) Mix the solution with vortex 1 minute

MATERIALS

MATERIALS PROVIDED

- CerTest *Influenza A+B+RSV+Adeno Resp.* combo card tests
- Sample collection tubes with diluent (Reagent B)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer
- Vortex


TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

- Procedure A using nasopharyngeal swab samples:

1. Proceed to shake the sample collection tube in order to assure good sample dispersion.
2. Remove the CerTest *Influenza A+B+RSV+Adeno Resp.* combo card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and, using each time the same tube, dispense 4 drops in the circular window marked with the letter A (5), 4 drops in the circular window marked with the letter B (6), 4 drops in the circular window marked with the letter C (7) and 4 drops in the circular window marked with the letter D (8).
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of reagent B until seeing the liquid running through the reaction zone.

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




- Procedure B using nasopharyngeal wash or aspirate samples:


1. Mix the solution with the vortex at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (3).
2. Remove the CerTest *Influenza A+B+RSV+Adeno Resp.* combo card test from its sealed bag just before using it.
3. Dispense, each time from the same testing tube, 4 drops in the circular window marked with letter A (4), 4 drops in the circular window marked with the letter B (5), 4 drops in the circular window marked with the letter C (6) and 4 drops in the circular window marked with the letter D (7).
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of reagent B until seeing the liquid running through the reaction zone.


Influenza A Strip A-procedure

Add 4 drops in A-circular window


Influenza B Strip B-procedure

Add 4 drops in B-circular window


RSV Strip C-procedure

Add 4 drops in C-circular window


Adenovirus Strip D-procedure

Add 4 drops in D-circular window


INTERPRETATION OF RESULTS (please refer to the illustration below)

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













	A (Influenza A)	B (Influenza B)	C (RSV)	D (Adeno)	Interpretation of the results
1.	-	-	-	-	No infection caused by <i>Influenza</i> type A and type B, <i>Respiratory Syncytial Virus</i> and <i>Adenovirus</i> .
	GREEN	GREEN	GREEN	GREEN	
2.	+	+	+	+	Simultaneous infection caused by <i>Influenza</i> type A, <i>Influenza</i> type B, <i>Respiratory Syncytial Virus</i> and <i>Adenovirus</i> .
	GREEN-RED	GREEN-RED	GREEN-RED	GREEN-RED	
3.	+	+	+	-	Simultaneous infection caused by <i>Influenza</i> type A, <i>Influenza</i> type B and <i>Respiratory Syncytial Virus</i> .
	GREEN-RED	GREEN-RED	GREEN-RED	GREEN	
4.	+	+	-	+	Simultaneous infection caused by <i>Influenza</i> type A, <i>Influenza</i> type B and <i>Adenovirus</i> .
	GREEN-RED	GREEN-RED	GREEN	GREEN-RED	
5.	+	-	+	+	Simultaneous infection caused by <i>Influenza</i> type A, <i>Respiratory Syncytial Virus</i> and <i>Adenovirus</i> .
	GREEN-RED	GREEN	GREEN-RED	GREEN-RED	
6.	-	+	+	+	Simultaneous infection caused by <i>Influenza</i> type B, <i>Respiratory Syncytial Virus</i> and <i>Adenovirus</i> .
	GREEN	GREEN-RED	GREEN-RED	GREEN-RED	
7.	+	+	-	-	Simultaneous infection caused by <i>Influenza</i> type A and type B.
	GREEN-RED	GREEN-RED	GREEN	GREEN	
8.	-	-	+	+	Simultaneous infection caused by <i>Respiratory Syncytial Virus</i> and <i>Adenovirus</i> .
	GREEN	GREEN	GREEN-RED	GREEN-RED	
9.	-	+	-	+	Simultaneous infection caused by <i>Influenza</i> type B and <i>Adenovirus</i> .
	GREEN	GREEN-RED	GREEN	GREEN-RED	
10.	-	+	+	-	Simultaneous infection caused by <i>Influenza</i> type B and <i>Respiratory Syncytial Virus</i> .
	GREEN	GREEN-RED	GREEN-RED	GREEN	
11.	+	-	+	-	Simultaneous infection caused by <i>Influenza</i> type A and <i>Respiratory Syncytial Virus</i> .
	GREEN-RED	GREEN	GREEN-RED	GREEN	
12.	+	-	-	+	Simultaneous infection caused by <i>Influenza</i> type A and <i>Adenovirus</i> .
	GREEN-RED	GREEN	GREEN	GREEN-RED	
13.	+	-	-	-	Infection caused by <i>Influenza</i> type A.
	GREEN-RED	GREEN	GREEN	GREEN	
14.	-	+	-	-	Infection caused by <i>Influenza</i> type B.
	GREEN	GREEN-RED	GREEN	GREEN	
15.	-	-	+	-	Infection caused by <i>RSV</i> .
	GREEN	GREEN	GREEN-RED	GREEN	
16.	-	-	-	+	Infection caused by <i>Adenovirus</i> .
	GREEN	GREEN	GREEN	GREEN-RED	
17.	Any other result				Invalid result: either A, B, C or D, we recommend repeating the assay using the same sample with another test.

Whatever result detailed in the above-mentioned table (2 to 16) should be followed up with additional confirmatory diagnostic procedures.

Whatever results detailed between 2 and 6 are not very common. Dual simultaneous respiratory infections are more frequent than triple or quadruple respiratory virus infection.

INVALID: Total absence of control coloured lines (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

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NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest *Influenza A+B+RSV+Adeno Resp.* combo card test should be used only with nasopharyngeal swab, nasopharyngeal wash and aspirate samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of *Influenza* type A, *Influenza* type B, *RSV* and/or *Adenovirus* respiratory infection. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal swab, nasopharyngeal wash and aspirate samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture and PCR.

EXPECTED VALUES

Multiple viral infections are frequent in hospitalized children with respiratory tract disease (17.4%). Multiple viral infections are linked to higher fever, longer hospital stays and more frequent use of antibiotics than in the case of infants with single *RSV* infections.

Respiratory Syncytial Virus (RSV) is the leading cause of emergency visits and hospitalization for acute lower respiratory tract infections (LRTI) in infants and young children worldwide.

In adults and older, healthy children, the symptoms of *Respiratory Syncytial Virus* are mild and typically mimic the common cold. Self-care measures are usually all that's needed to relieve any discomfort.

Infection with *Respiratory Syncytial Virus* can be severe in some cases, especially in premature babies and infants with underlying health conditions. *RSV* can also become serious in older adults, adults with heart and lung diseases, or anyone with a very weak immune system (immunocompromised).

The exact prevalence and incidence of adenoviral infections are unknown, because most cases are seen by general practitioners and optometrists. *Adenovirus* is a very common infection, estimated to be responsible for between 2% and 5% of all respiratory infections. In winter, infection with type 4 or 7 causes recognisable illness in military recruits, with about 25% requiring hospitalisation for fever and lower respiratory tract disease.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

The evaluation, with nasopharyngeal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Influenza A+B+RSV+Adeno Resp.*, CerTest) and others commercial available immunochromatographic test (BinaxNOW® *Influenza A&B* (Alere), BinaxNOW® *RSV* (Alere) and *Adenovirus Respi*, (CorisBioConcept)) and a commercial available Immunofluorescence test (PathoDx® *Adenovirus*, Remel). The results were as follows:

IC test: CerTest Influenza A+B+RSV+Adeno Resp. (Influenza A and B)		IC test: BinaxNOW® Influenza A&B		
		+	-	Total
	+	5	0	5
	-	0	6	6
	Total	5	6	11

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The results showed a high sensitivity and specificity to detect *Influenza* type A and type B, *Respiratory Syncytial Virus* and *Adenovirus* using CerTest Influenza A+B+RSV+Adeno Resp.

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Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Influenza A+B+RSV+Adeno Resp.*; no cross reactivity against organisms that cause other respiratory infections:

- *Influenza type A strip no cross reactivity with:*

<i>Influenza type B</i>	<i>Respiratory Syncytial Virus</i>	<i>Adenovirus</i>
-------------------------	------------------------------------	-------------------
- *Influenza type B strip no cross reaction with:*

<i>Influenza type A</i>	<i>Respiratory Syncytial Virus</i>	<i>Adenovirus</i>
-------------------------	------------------------------------	-------------------
- *RSV strip no cross reaction with:*

<i>Influenza type A</i>	<i>Influenza type B</i>	<i>Adenovirus</i>
-------------------------	-------------------------	-------------------
- *Adenovirus strip no cross reaction with:*

<i>Influenza type A</i>	<i>Influenza type B</i>	<i>Respiratory Syncytial Virus</i>
-------------------------	-------------------------	------------------------------------

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- DANIEL E. NOYOLA, et al. « Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens". *Journal of Clinical Microbiology*, Mar. 2000, p. 1161–1165 Vol. 38, No. 3
- ADRIANA WEINBERG AND MIRANDA L. WALKER. "Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza Virus A and B". *Clinical And Diagnostic Laboratory Immunology*, Mar. 2005, p. 367–370 Vol. 12, No. 3
- ANDREEA C. CAZACU, GAIL J. DEMMLER et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. *Journal of Clinical Microbiology*, Aug. 2004, p. 3661–3664 Vol. 42, No. 8.
- KENNETH E. IRMEN AND JAMES J. KELLEHER. "Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital". *Clinical and Diagnostic Laboratory Immunology*, May 2000, p. 396–403 Vol. 7, No. 3.
- MARCELA ECHAVARRIA, JOSE L. SANCHEZ, et al. "Rapid Detection of *Adenovirus* in Throat Swab Specimens by PCR during Respiratory Disease Outbreaks among Military Recruits", *Journal of Clinical Microbiology*, Feb. 2003, Vol. 41, No. 2, p. 810–812
- MARILYN J. AUGUST' AND ANN L. WARFORD; "Evaluation of a Commercial Monoclonal Antibody for Detection of *Adenovirus* Antigen" *Journal of Clinical Microbiology*, Nov. 1987, Vol. 25, No. 11; p. 2233-2235
- ANN R. FALSEY and EDWARD E. WALSH. "Respiratory Syncytial Infection in Adults", *Clinical Microbiology Reviews*, July 2000, Vol. 13, No.3, p. 371-384.
- DIANE C. HALSTEAD, SANDRA TODD, and GALE FRITCH; "Evaluation of Five Methods for Respiratory Syncytial Virus Detection" *Journal of Clinical Microbiology*, May 1990, Vol. 28, No. 5; p. 1021-1025.
- LINA B, VALETTE M, FORAY S, et al. "Surveillance of community-acquired viral infections due to respiratory viruses in Rhone-Alpes (France) during winter 1994 to 1995", *J Clin Microbiol* 1996;34:3007–3011
- DREWS AL, ATMAR RL, GLEZEN WP, BAXTER BD, PIEDRA PA, GREENBERG SB. "Dual respiratory virus infections", *Clin Infect Dis* 1997;25:1421–1429.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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