CHORUS

RUBELLA IgM

REF 81031 (36 tests)

Prodotto da/Manufactured by/Fabricado por:
DIESSE Diagnostica Senese
 Via delle Rose 10
 53035 Monteriggioni (Siena) - Italy

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Chantler and C.J. Evans: Selection and performance of monoclonal and polyclonal antibodies in an IgM antibody capture enzyme immunoassay for rubella. J. Immunol. Meth. 87: 109 (1986).
3. R.S. Tedder and J.L. Yao: The production of monoclonal antibodies to rubella haemagglutinin and their use in antibody-capture assays for rubella-specific IgM. J. Hyg. Camb. 88: 335 (1982).
4. B. Forghani, C. Myoraku, N. Schmidt: Use of monoclonal antibodies to human immunoglobulin M in "capture" assays for measles and Rubella immunoglobulin M. J. Clin. Microbiol. 18, 652 (1983).
5. T.Zhang , C.A. Mauracher , L.A. Mitchell , : Detedtion of Rubella Virus-Specific Immunoglobulin G (IgG) , IgM , and IgA Antibodies by Immunoblot Assays. J.Clin.Microbiol. Apr.1992.p.824-830
6. P.Partanen, H. Seppanen , Jukka Suni , A. Vaheri : Selective Reactivity of antibodies to Human Immunoglobulins G , M and A with Rubella Virus Proteins. J.clin.microbial.,May 1985 , p.800-802

INSTRUCTIONS FOR USE

CHORUS RUBELLA IgM

REF 81031

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO RUBELLA VIRUS IN HUMAN SERUM USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENTS

2. INTRODUCTION

Rubella virus causes a very common infectious disease in man, accompanied by a typical diffused rash. The course of the illness is generally short and without complications. However, it is harmful for the fetus if contracted during pregnancy, because the virus is teratogenic and may cause very serious malformations especially during the first month of pregnancy. It is therefore essential to determine the patient's immune status before the start of pregnancy, if possible, in order to perform vaccination or, if this is impossible, to monitor the course of the pregnancy to observe whether seroconversion occurs. Several methods are employed for the assay of specific immunoglobulins (1). ELISA is one of these, and has the advantage of being completely automated. The serological diagnosis of Rubella may be performed with the assay of specific IgM because these antibodies, together with the clinical symptoms, are of considerable diagnostic relevance.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The antigen, composed of purified and inactivated Rubella Virus, is bound to the solid phase. Through incubation with human serum diluted in a diluent blocking the IgG, the specific IgM are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgM monoclonal antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.


The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL coated with Rubella Virus purified Antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: DILUENT FOR THE SAMPLES

Contents: Proteic solution containing anti human IgG serum, phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgM monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT [REF] 83607
- Chorus / Chorus Trio Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to

clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to be used to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the test sample and that of the Cut-Off) which can be used as a quantitative measure, as it is proportional to the amount of specific IgM present.

The results can be interpreted as follows:

NEGATIVE: when the result is < 0.9

DOUBTFUL: when the result is between 0.9 and 1.1

POSITIVE: when the ratio is > 1.1

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

11. LIMITATIONS OF THE PROCEDURE

Specimens collected earlier than 5 days after onset of symptoms may not contain sufficient IgM antibody to be positive in the test.

Since maternal IgM antibody does not cross the placenta, virus-specific IgM antibody in the new-born can be assumed to be a result of congenital infection. Because of the possibility of contamination of cord blood with maternal IgM, it is prudent to confirm positive viral IgM antibody results on cord blood samples by testing a follow-up specimen from the infant. Rubella IgM antibody responses recently have been reported in a small proportion of patients with infectious mononucleosis; some of these showed narrow specificity for certain strains of Rubella virus. Polyclonal stimulation of B memory cells by EBV was considered in some cases to be a possible mechanism for heterotypic IgM antibody responses.

The results must always be interpreted together with the clinical picture and data from other diagnostic investigations. Sera from patients in an early or late stage of the disease could give a repeatedly negative result close to the Cut-off value. In such cases, a confirmation of the result is recommended, by testing a new sample drawn from the same patient at 7-10 days interval from the first sample.

12. ANALYTICAL SPECIFICITY

The presence in the samples of potentially interfering substances such as Rheumatoid factor, Bilirubin, Triglycerides and haemolytic samples does not influence the result. No influence on the result was verified in the test of sera samples from women during pregnancy.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

431 samples were tested in the experimentation both with the Diesse kit and with another commercial reference kit.

The samples that gave positive results or were in disagreement, were confirmed with another test.

The results are summarized in the following table.

Table 1.

		Reference		Total
		+	-	
Diesse	+	74	4	78
	-	2	351	353
Total		76	355	431

The kit shows a specificity of 98.9 % (CI_{95%}: 97.1 – 99.6) and a sensitivity of 97.4 % (CI_{95%}: 90.9 – 99.3).

14. PRECISION

Table 2. Within-run precision

Sample	INDEX	SD	CV%
Positive control	1.7	0.16	9
RVM 1	0.5	0.08	15
RVM 2	0.9	0.12	13
RVM 3	2.2	0.14	6

Table 3. Between-run precision

Sample	Mean Index					Mean	SD	CV%
	Run 1	Run 2	Run 3	Run 4	Run 5			
RVM 3	1.6	1.6	1.5	1.8	1.5	1.6	0.12	8
RVM 2	0.7	0.8	0.8	0.7	0.8	0.8	0.05	6
Serum 32	0.9	1.0	1.0	1.0	0.8	0.9	0.09	10
Serum 2208	1.1	0.8	1.0	1.1	1.1	1.0	0.13	13
Pos. con.	1.4	1.4	1.4	1.4	1.4	1.4	0.00	-

15. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Chantler and C.J. Evans: Selection and performance of monoclonal and polyclonal antibodies in an IgM antibody capture enzyme immunoassay for rubella. J. Immunol. Meth. 87: 109 (1986).
3. R.S. Tedder and J.L. Yao: The production of monoclonal antibodies to rubella haemagglutinin and their use in antibody-capture assays for rubella-specific IgM. J. Hyg. Camb. 88: 335 (1982).
4. B. Forghani, C. Myoraku, N. Schmidt: Use of monoclonal antibodies to human immunoglobulin M in "capture" assays for measles and Rubella immunoglobulin M. J. Clin. Microbiol. 18, 652 (1983).
5. T.Zhang , C.A. Mauracher , L.A. Mitchell ,: Detedtion of Rubella Virus-Specific Immunoglobulin G (IgG) , IgM , and IgA Antibodies by Immunoblot Assays. J.Clin.Microbiol. Apr.1992.p.824-830
6. P.Partanen, H. Seppanen , Jukka Suni , A. Vaheeri : Selective Reactivity of antibodies to Human Immunoglobulins G , M and A with Rubella Virus Proteins. J.clin.microbial.,May 1985 , p.800-802.

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Chantler and C.J. Evans: Selection and performance of monoclonal and polyclonal antibodies in an IgM antibody capture enzyme immunoassay for rubella. J. Immunol. Meth. 87: 109 (1986).
3. R.S. Tedder and J.L. Yao: The production of monoclonal antibodies to rubella haemagglutinin and their use in antibody-capture assays for rubella-specific IgM. J. Hyg. Camb. 88: 335 (1982).
4. B. Forghani, C. Myoraku, N. Schmidt: Use of monoclonal antibodies to human immunoglobulin M in "capture" assays for measles and Rubella immunoglobulin M. J. Clin. Microbiol. 18, 652 (1983).
5. T.Zhang , C.A. Mauracher , L.A. Mitchell , : Detedtion of Rubella Virus-Specific Immunoglobulin G (IgG), IgM , and IgA Antibodies by Immunoblot Assays. J.Clin.Microbiol. Apr.1992.p.824-830.
6. P.Partanen, H. Seppanen , Jukka Suni , A. Vaheri : Selective Reactivity of antibodies to Human Immunoglobulins G , M and A with Rubella Virus Proteins. J.clin.microbial.,May 1985 , p.800-802.

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Chantler and C.J. Evans: Selection and performance of monoclonal and polyclonal antibodies in an IgM antibody capture enzyme immunoassay for rubella. J. Immunol. Meth. 87: 109 (1986).
3. R.S. Tedder and J.L. Yao: The production of monoclonal antibodies to rubella haemagglutinin and their use in antibody-capture assays for rubella-specific IgM. J. Hyg. Camb. 88: 335 (1982).
4. B. Forghani, C. Myoraku, N. Schmidt: Use of monoclonal antibodies to human immunoglobulin M in "capture" assays for measles and Rubella immunoglobulin M. J. Clin. Microbiol. 18, 652 (1983).
5. T.Zhang , C.A. Mauracher , L.A. Mitchell , : Detedtion of Rubella Virus-Specific Immunoglobulin G (IgG) , IgM , and IgA Antibodies by Immunoblot Assays. J.Clin.Microbiol. Apr.1992.p.824-830.
6. P.Partanen, H. Seppanen , Jukka Suni , A. Vaheri : Selective Reactivity of antibodies to Human Immunoglobulins G , M and A with Rubella Virus Proteins. J.clin.microbial.,May 1985 , p.800-802.

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusqu'e Ημερομηνία λήξης Prazo de validade
	EN ES IT	Do not reuse No reutilizar Non riutilizzare	FR GR PT	Ne pas réutiliser Μην κάνετε επαναληπτική χρήση Não reutilizar
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbicante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

CE
0123

Diesse Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) – Italy
Tel. 0577-587111