

Product PN 10497511, 10497500, 10497501, 10497502, 10497503, 10497603,

10497600, 10497601, 10497602, 10497504, 10497506, 10497507,

10497508, 10497509, 10497510

Mod. 984

Rev. 02

Description Microbiological Monitors and Analytical Funnels

Microbiological Monitors and Analytical Funnels



PRODUCT DESCRIPTION

Microbiological Monitors (100 ml) with 47 mm or 56 mm Nitrocellulose Membrane--sterile, disposable, single-use

Produc	t Code		
47 mm	56 mm	Description	Quantity
10497511	10497603	Monitor, Nitrocellulose, 0.2 um, white/black grid	50/pk
10497500	10497600	Monitor, Nitrocellulose, 0.45 um, white/black grid	50/pk
10497501	N/A	Monitor, Nitrocellulose, 0.45 um, white/black grid,	50/pk
		individually packed	
10497502	10497601	Monitor, Nitrocellulose 0.45 um, black/white grid	50/pk
10497503	10497602	Monitor, Nitrocellulose 0.8 um, black/white grid	50/pk

Analytical Funnel (100 ml) with removable Nitrocellulose Membrane--sterile, disposable, single-use

Product	Description	Quantity
Code		
10497504	Funnel, Nitrocellulose, white/black grid, sterile, 0.45 um	50/pk
10497506	Funnel, Nitrocellulose, white/black grid, sterile, 0.45 um, individually packaged	50/pk
10497507	Funnel, Nitrocellulose, white/black grid, sterile, 0.2 um	50/pk
10497508	Funnel, Nitrocellulose, black/white grid, sterile, 0.45 um	50/pk
10497509	Funnel, Nitrocellulose, black/white grid, sterile, 0.45 um, individually packaged	50/pk
10497510	Funnel, Nitrocellulose, white/black grid, sterile, 0.2 um, individually packaged	50/pk

MANUFACTURER NAME

GVS North America 63 Community Drive Sanford, Me 04073

Phone: +1.866.736.1250

eMail: <u>CustomerCareGVSLS@gvs.com</u> - Website: <u>www.gvs.com</u>

INTENDED USE / APPLICATION

Microbiological analysis of:

- Potable water
- Beer and wine
- Waste water
- Dairy products
- Soft drinks



Product PN 10497511, 10497500, 10497501, 10497502, 10497503, 10497603,

10497600, 10497601, 10497602, 10497504, 10497506, 10497507,

10497508, 10497509, 10497510

Mod. 984 Rev. 02

Description Microbiological Monitors and Analytical Funnels

Description Mic	crobiological Monitors and Analytical Funnels
PRODUCT CHARACTERISTIC	 Rapid testing with no need to sterilize funnels or filter base between samples, testing time can be reduced by up to 70% No flaming required: Combined filtration unit minimizes the risk of cross-contamination All-in-one system: Filtration unit easily converts to a Petri dish, which can be labeled and incubated for culturing Reduced contamination: Single use materials virtually eliminate cross-contamination between funnel and membrane Reproducible results: All-in-one filtration unit reduces the chance of external error Easy handling: Ready-to-use, pre-sterilized monitors are simple to use
INSTRUCTIONS WARNINGS	Store at room temperature and avoid direct sunlight. Nitrocellulose membrane is flammable. Keep away from heat and open flame.
PRODUCT SHELF LIFE	3 years when stored properly
STERILIZATION	Product is sold sterilized by gamma irradiation at 25- 40 kgy.
PACKAGING AND LABELLING	Primary: 50 units and IFU are packaged in a labeled box; 2 layers of 25 units. Box is shrinkwrapped. GVS North America Sanford, ME 04073 - USA www.gvs.com +1 866-736-1250 Microbiological Monitors Nitrocellulose, 47mm, 0.2 μm, white/black grid REF 10497511 LOT SAMPLE QTY 100 EA Jan/2018 Jan/2021 Secondary (for shipping): Another barcode label is applied to outside of shipping box, with reference to the cumulative quantity.
CERTIFICATE OF COMPLIANCE	With each shipment, GVS NA will provide the CoC with item number, lot number and expiration date. The Quality management system is in compliance with ISO 9001.



Product PN 10497511, 10497500, 10497501, 10497502, 10497503, 10497603,

10497600, 10497601, 10497602, 10497504, 10497506, 10497507,

10497508, 10497509, 10497510

Mod. 984

Rev. 02

Description Microbiological Monitors and Analytical Funnels

VISUAL REQUIREMENTS

Visual acceptance requirements apply when inspected under below conditions:

Magnification: none

Light type: normal room lighting

	Acceptance Requirement	AQL	Sampling Plan
1	Scratches or abrasion that may affect function	None	
2	Scratches that interfere with a clear view through the lid or		
	graduation areas of the funnel	None	
3	Debris or contaminants	None	4000/
4	Excess flash or flash that may affect function, such as at the	None	100%
	base vacuum port		
5	Proper placement of media	All	
6	No torn or miss cut media	All	
7	Grid lines: alignment, visible	All	

PERFORMANCE REQUIREMENTS

Acceptance Requirement	AQL	Sampling Plan
Leak Test	No leaks after 10 minutes	12 at start-up & twice per shift

This material specification describes the properties of product above indicated.

This document contains general requirements, material description, drawing references, defect specification, biological material requirements.

REVISIONS AND APPROVALS:

DATE	REV.	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
9/1/18	0	Initial release	Debra English/Manager Biochemistry Applications	Kevin Wrigley/Director, Quality Here M. W.



Product PN	10497511, 10497500, 10497501, 10497502, 10497503, 10497603, 10497600, 10497601, 10497602, 10497504, 10497506, 10497507, 10497508, 10497509, 10497510	Mod. 984
Description	Microbiological Monitors and Analytical Funnels	TREV. UZ

CustomerApproval:

Company name	
Approved by: (Name, Function) (Signature)	
Date(Company stamp)	

Please send back this document signed for approval. If we will not receive this specification signed , we consider the first order placed as implicit approval.