CHORUS

**Herpes Simplex 1
IgG Recombinant**

REF 81023

DIESSE

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

Chorus Herpes simplex 1 IgG Recombinant is an immunoassay kit for automated qualitative detection of IgG class antibodies against Herpes simplex virus (Type 1). The test is performed in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments. The kit is intended to detect the exposure to Herpes simplex virus (Type 1) infection as an aid to the relative diagnosis. It must be used by professional laboratory users only.

2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 is a highly contagious infection, that is common and endemic throughout the world. Even if HSV-1 is mainly transmitted by oral-to-oral contact, HSV-1 can also be transmitted to the genital area through oral-genital contact.

HSV-1 causes chiefly oral-labial and oral-facial lesions, but a proportion of HSV-1 infections are genital, sexually transmitted herpes (infections in the genital or anal area).

Oral herpes infection is mostly asymptomatic. Main symptoms of oral herpes include painful blisters or open ulcers in or around the mouth. Genital sexually transmitted herpes caused by HSV-1 can be asymptomatic or can have mild symptoms that go unrecognized. When symptoms occur, genital herpes is characterized by one or more genital or anal blisters or ulcers.

Exposure to HSV can be determined indirectly by assaying for specific serum antibodies of the viruses through ELISA tests, which compared to other tests, can provide a high level of specificity, sensitivity and reliability.

Due to the significant cross-reactivity existing among viruses HSV-1 and HSV-2 is virtually impossible to reliably differentiate HSV-1 infections from those HSV-2 with these tests. Some studies have shown that tests for IgG that are based on the type-specific proteins gG-1 (for HSV-1) and gG-2 (for HSV-2) are the most accurate for discriminating between infection with HSV-1 or HSV-2.

Since the product Chorus Herpes simplex 1 IgG Recombinant is based on the presence of the recombinant protein gG-1 (used for the coating of the MT Plate), the device is used for the detection of anti-HSV-1 IgG antibodies as an aid in the HSV-1 diagnosis.

3. PRINCIPLE OF THE METHOD

The Chorus Herpes Simplex 1 IgG Recombinant device is ready to use for the detection of IgG antibodies against Herpes Simplex Virus (type 1) in the Chorus/Chorus TRIO instruments. The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The recombinant antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially

biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with recombinant HSV1 antigen, composed of recombinant G1 glycoprotein specific for Herpes Simplex Virus Type 1 (maximum concentration of 2 µg/ml).

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02%, indicator to reveal the presence of the serum and cells extracts

Positon 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgG labelled with horseradish peroxidase (maximum concentration of 1 µg/ml), in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the sample patient has to be added.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 ml

Contents: 100% human serum containing antibodies anti-Herpes Simplex Virus (type 1) IgG, Proclin 300 0.2%, Gentamicin 0.05% and Ethylene Glycol 23%.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: 100% human serum containing antibodies anti-Herpes Simplex Virus (type 1) IgG and and Proclin 300 0.2%, Gentamicin 0.05% and Ethylene Glycol 23%.

Liquid, ready for use.

The device is a qualitative test: Calibrator and Positive control cannot provide traceability to International Standard Units. Confidence in measurements is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced starting from positive human sample diluted to reach a theoretical concentration of antibodies. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference product.

MATERIALS REQUIRED BUT NOT PROVIDED

- Chorus Instruments
- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS SAMPLE DILUENT/NEGATIVE CONTROL **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated

and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum from whole blood collected by venipuncture and handled following standard laboratory techniques.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C (for at least 38 months) and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. At each change of batch, perform the calibration.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	50 µl/device
POSITIVE CONTROL	50 µl/device

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum as reported in the Instrument Operating Manual to check the validity of the obtained results. The acceptability range is reported on the enclosed Quality control certificate.

If the result of the positive control is out of the acceptability range, the calibration must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1

NEGATIVE: when the result is < 0.9

DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

The product must be used by professional laboratory users only.

The test is not applicable to sample types different from serum. All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 0.1 and 0.6 Index.

Among the affected population the expected values, which have been determined by examining 80 pre-characterized positive sera, were between 1.2 and 8.4 Index.

13. ANALYTICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (44 – 220 UI/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the serum sample of the interfering substances described above except for Rheumatoid Factor does not affect the test result.

14. CROSS-REACTIONS

56 samples positive to Rubella, Epstein-Barr Virus, Varicella, Cytomegalovirus, HSV 2 and Toxoplasma were tested.

No significant cross-reactions were found.

15. METHOD COMPARISON

In an experimentation 153 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	80	0	80
	-	6	67	73
	Total	86	67	153

Percent Positive Agreement (~Diagnostic Sensitivity):

93.0% CI_{95%}: 85.4-96.7

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 94.6-99.9

Positive Predicted Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predicted Value (NPV): 91.8% CI_{95%}: 87.4-96.1

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean Index	CV%	Mean Index	CV%
1	0.1	40.0*	0.1	-
2	0.3	-	0.3	-
3	0.5	-	0.5	12.0
4	0.8	3.8	0.8	10.0
5	1.3	5.4	1.4	7.9
6	1.8	4.4	1.8	8.3
7	3.2	4.7	3.0	7.3
8	4.5	5.1	4.2	7.4
9	6.4	3.6	6.6	10.6
10	6.9	4.5	7.2	12.5

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero

Sample	Precision between batches		Precision between instruments	
	Mean Index	CV%	Mean Index	CV%
1	0.1	-	0.1	-
2	0.3	-	0.3	-
3	0.5	-	0.5	-
4	0.8	7.5	0.8	-

5	1.5	4.0	1.5	8.0
6	1.8	5.6	1.8	-
7	3.2	4.7	3.2	4.7
8	4.5	1.3	4.5	1.3
9	7.4	9.5	7.4	7.6
10	8.2	5.4	8.2	4.4

17. REFERENCES

1. (<https://www.who.int/news-room/fact-sheets/detail/herpes-simplex-virus#hsv1>)
2. Anzivino E. et All. Herpes simplex virus infection in pregnancy and in neonate: status of art of epidemiology, diagnosis, therapy and prevention. Virology Journal 2009, 6:40.
3. Singh A. et All. Canadian STI best practice laboratory guidelines. The laboratory diagnosis of Herpes simplex virus infections. Can J Infect Dis Med Microbiol Vol 16 No 2 March/April 2005
4. (<https://www.who.int/news-room/fact-sheets/detail/herpes-simplex-virus#hsv1>).

18. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

19. SUMMARY OF SAFETY AND PERFORMANCE

This document, available on Eudamed database, is part of Technical Documentation.

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