CHORUS

VARICELLA IgM

REF 81079

REF 81079/12

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(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO VARICELLA ZOSTER VIRUS IN HUMAN SERUM USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENTS

2. INTRODUCTION

Varicella and Herpes Zoster are two clinical manifestations resulting from infection by Varicella Zoster Virus (VZV). Varicella, or chicken pox, is a highly contagious disease which is generally a consequence of primary infection by VZV and normally affects children. Infection caused by VZV during pregnancy can cause disease or malformation of the fetus; if it occurs at the end of pregnancy it may lead to the death of the neonate.

Herpes Zoster is a disease which chiefly affects adults and appears to be caused by a reactivation of the virus, which can remain latent for long periods in the spinal sensorial ganglia following primary infection. The infection causes painful cutaneous eruptions along the path of the nerves affected.

Serological methods are usually adopted to determine the immune state of subjects at risk (chiefly immunodepressed patients) and in pre- and post-natal diagnosis of infected subjects.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The antigen, composed of inactivated and purified Varicella Zoster virus, is bound to the solid phase. Through incubation with human serum diluted in a diluent blocking the IgG, the specific IgM are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgM monoclonal antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81079).


The kit is sufficient for 12 tests (REF 81079/12).

DD DEVICES

6 packages each containing 6 devices (REF 81079).

2 packages each containing 6 devices (REF 81079/12).

Use: equilibrate a package at room temperature, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL coated with Varicella Zoster purified Antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: DILUENT FOR THE SAMPLES

Contents: Proteic solution containing anti human IgG serum, phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum

Position 2: CONJUGATE

Contents: anti-human IgM monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minutes exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to be used to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the test sample and that of the Cut-Off) which can be used as a quantitative measure, as it is proportional to the amount of specific IgM present.

The test serum can be interpreted as follows:

NEGATIVE: when the result is < 0.9

DOUBTFUL/EQUIVOCAL: for all the results between 0.9 and 1.1

POSITIVE: when the result is > 1.1

In the case of a doubtful/equivocal result, repeat the test. If the test remains doubtful/equivocal, collect a new sample after 1-2 weeks.

11. LIMITATIONS OF THE PROCEDURE

The test results should be used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

12. ANALYTICAL SPECIFICITY

The presence of potentially interfering substances such as Rheumatoid Factor, Bilirubin, Triglycerides and haemolytic samples does not interfere with the test result.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

168 samples were analyzed during the experimentation both with the DIESSE kit and with another commercial reference kit.

The results are summarized in the following table:

Table 1.

		REFERENCE	
		+	-
DIESSE	+	33	6
	-	1	128
	tot	34	134

The kit has a sensitivity of 97% and a specificity of 96%.

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