CHORUS

MYCOPLASMA PNEUMONIAE IgG

REF 81034 (36 tests)

Prodotto da/Manufactured by/Fabricado por/
 Произведено от/Proizvodi / Tootja / Pagamino /
 Wyprodukowany przez / Proizvaja:
DIESSE Diagnostica Senese
Strada dei Laghi 39
53035 Monteriggioni (SIENA) – Italy

	Capitolo Section Kapitola Kapitel Κεφάλαιο Capítulo Chapitre Capítulo Capitol Peatükk Rozdział Odsek Skyrius
Modifiche introdotte nella revisione corrente Changes introduced in the current revision Změny zavedené v aktuální revizi In der aktuellen Revision eingeführte Änderungen Αλλαγές που εισάγονται στην τρέχουσα αναθεώρηση Cambios introducidos en la revisión actual Changements introduits dans la révision actuelle Alterações introduzidas na revisão atual Modificări introduse în actuala revizuire Praeguses versioonis tehtud muudatused Zmiany wprowadzone w aktualnej wersji Spremembe, uvedene v trenutni reviziji Pakeitimai, pateikti dabartinėje redakcijoje	1 – 2 – 3 – 4 – 6 – 8 – 10

INDICE / INDEX / ÌNDICE / ИНДЕКС / INDEKS / INDEKS / SPIS TREŚCI / KAZALO

1. UTILIZZAZIONE / INTENDED USE / INDICACIONES / UTILIZAÇÃO / ИСПОЛЗВАЙТЕ / KORIŠTENJE / KASUTA / NAUDOJIMAS / ZASTOSOWANIE / UPORABA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / RESUMEN Y EXPLICACIÓN DEL TEST / INTRODUÇÃO / ВЪВЕДЕНИЕ / UVOD / SISSEJUHATUS / ĮVADAS / WPROWADZENIE / UVOD
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / PRINCIPIO DEL MÉTODO / PRINCÍPIO DO MÉTODO / ПРИНЦИП НА МЕТОДА / NAČELO METODE / MEETODI PÕHIMÕTE / METODO PRINCIPAS / ZASADA METODY / NAČELO METODE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / COMPOSIÇÃO DO KIT E PREPARAÇÃO DOS REAGENTES / СЪСТАВ НА КОМПЛЕКТА И ПРИГОТВЯНЕ НА РЕАГЕНТ / SASTAV KOMPLETA I PRIPREMA REAGENSA / KOMPLEKTI KOOSTIS JA REAGENTIDE VALMISTAMINE / RINKINIO SUDĖTIS IR REAGENTŲ PARUOŠIMAS / SKŁAD ZESTAWU I PRZYGOTOWANIE ODCZYNNIKÓW / SESTAVA KOMPLETA IN PRIPRAVA REAGENTOV
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALIDADES DE CONSERVAÇÃO E ESTABILIDADE DOS REAGENTES / СЪХРАНЕНИЕ И СТАБИЛНОСТ НА РЕАКТИВИТЕ / NAČIN ČUVANJA I STABILNOST REAGENSA / REAGENTIDE SÄILITAMINE JA STABIILSUS / REAGENTŲ LAIKYMAS IR STABILUMAS / PRZECHOWYWANIE I STABILNOŚĆ ODCZYNNIKÓW / NAČIN SHRANJEVANJA IN STABILNOST REAGENTOV
6. PRECAUZIONI / PRECAUTIONS / PRECAUCIONES DE USO / PRECAUÇÕES / ПРЕДУПРЕЖДЕНИЯ / MJERE OPREZA / ETTEVAATUST / ATSARGUMO PRIEMONĖS / ŚRODKI OSTROŻNOŚCI / PREVIDNOSTNI UKREPI
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / TIPO DE MUESTRA Y CONSERVACION / TIPO DE AMOSTRA E CONSERVAÇÃO / ВИД НА ПРОБАТА И СЪХРАНЕНИЕ / VRSTA UZORAKA I ČUVANJE / PROOVIDE TÛÜP JA SÄILITAMINE / MĖGINIO TIPAS IR LAIKYMAS / RODZAJ PRÓBK I JEJ PRZECHOWYWANIE / VRSTA VZORCA IN SHRANJEVANJE
8. PROCEDIMENTO / TEST PROCEDURE / PROCEDIMENTO / MODO DE PROCEDER / ПРОЦЕДУРА / POSTUPAK / PROTSEDUUR / PROCEDŪRA / PROCEDURA / POSTOPEK
9. VALIDAZIONE DEL TEST / TEST VALIDATION / VALIDACIÓN DEL TEST / VALIDAÇÃO DO TESTE / ВАЛИДИРАНЕ НА ТЕСТА / VALIDACIJA TESTA / TESTI VALIDEERIMINE / TYRIMO PATVIRTINIMAS / WALIDACJA BADANIA / VALIDACIJA TESTA
10. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETACIÓN DE LOS RESULTADOS / INTERPRETAÇÃO DOS RESULTADOS / ТЪЛКУВАНЕ НА ТЕСТА / TUMAČENJE TESTA / TESTI TÖLGENDAMINE / TYRIMO AIŠKINIMAS / INTERPRETACJA BADANIA / RAZLAGA TESTA
11. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / LIMITACIONES / LIMITES DO TESTE / ОГРАНИЧЕНИЯ НА ПРОЦЕДУРАТА / OGRANIČENJA POSTUPKA / PROTSEDUURI PIIRANGUD / PROCEDŪROS APRIBOJIMAI / OGRANICZENIA PROCEDURY / OMEJITVE POSTOPKA
12. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ESPECIFICIDAD ANALÍTICA / ESPECIFICIDADE ANALÍTICA / АНАЛИТИЧНА СПЕЦИФИЧНОСТ / ANALITIČKA

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

INSTRUCTIONS FOR USE

CHORUS MYCOPLASMA PNEUMONIAE IgG

REF 81034

(English)

1. INTENDED USE

SEMIQUANTITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO MYCOPLASMA IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS/CHORUS TRIO INSTRUMENTS

2. INTRODUCTION

Mycoplasma pneumoniae is the most common etiological agent causing pneumonia acquired in community environments, particularly in the age range of 5 to 30 years; it can cause epidemics which develop slowly as the incubation period varies from 10 to 14 days, and the infection involves subjects living in close contact or in segregated groups (schools, military barracks, family groups). Pneumonia caused by Mycoplasma is also known as atypical primary pneumonia or Eaton's pneumonia.

Mycoplasma pneumoniae attacks and destroys ciliated epithelial cells of the respiratory tract mucosa. Microscopic examination evidences interstitial pneumonia, bronchitis and bronchiolitis.

Due to the recurrence of symptoms caused by different etiological agents, additional diagnostic measures such as serological tests are necessary for the diagnosis of an acute pneumonia infection.

The immune response, following a Mycoplasma pneumoniae infection, depends on the type of infection: the IgM-class antibodies are more frequently found in case of primary infections as demonstrated by their presence in younger patients. The IgM are low or not found in older patients with higher probability of reinfection. On the other hand the IgA represent the highest sensitivity marker in case of an infection in progress, a reinfection or a recent infection.

Specific IgA-class antibodies appear early during the infection and reach peak levels during the first 4 weeks, then declining quite rapidly, earlier than the IgG; diagnosis of an acute infection can often be made on the basis of a single sample test. The IgG appear after the IgM and IgA and they generally reach the peak level during the 5th week. A high titer together with a significant increase between two samples taken approximately 2 weeks apart, confirms that an infection is in progress.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The sonicated, solubilized Mycoplasma antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation with the conjugate (anti-human IgG monoclonal antibodies conjugated with horse radish peroxidase) is performed. The conjugate which has not been bound is eliminated and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.


The results are expressed in AU/ml calculated in reference to an internal standard curve.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: equilibrate a package at room temperature, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with *Mycoplasma pneumoniae* antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE.

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgG labeled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT [REF] 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation

DEVICES	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

6. PRECAUTIONS**ONLY FOR DIAGNOSTIC USE IN VITRO.****Attention:**

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS/CHORUS TRIO instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:

- a) The conjugate contains phenol
- b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the Chorus/Chorus TRIO instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument gives a semiquantitative result in AU/mL, calculated on the basis of a lot-dependent graph stored in the instrument.

The degree of immunity of the test serum can be interpreted as follows:

IMMUNE: when the IgG concentration in the sample is >18 AU/mL

NOT IMMUNE: if the concentration is < 12 AU/mL

DOUBTFUL/EQUIVOCAL: when the concentration is in the range 12-18 AU/mL

Titer	Interpretation of the results		CHORUS Result
> 40 AU/mL	High titer	Titer compatible with a recent or current <i>M.pneumoniae</i> infection	P
≤ 40 AU/mL >18 AU/mL	Moderate titer	The interpretation is only possible with a new collection of sample after 8-15 days.	P
≤ 18 AU/mL ≥12 AU/mL	Low titer		D/E
< 12 AU/mL	Not significant titer		N

When the sign > appears on the instrument the AU/mL are over 100 AU/mL; in this case dilute the sample in Sample Diluent/Negative Control [REF] 83607 (e.g. 1:3) and repeat the test. The result must be multiplied by the dilution used. If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS OF THE PROCEDURE

A negative result does not exclude a *Mycoplasma* infection: in case of clinical suspect collect a new sample after 15 days.

Serum samples collected during the acute phase of the infection can give a positive result to IgM antibodies only. It is suggested to determine the *Mycoplasma* IgM level using the *Mycoplasma* IgM kit (REF 81035).

Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures.

12. ANALYTICAL SPECIFICITY

The presence in the serum sample of potentially interfering substances such as Reactive Protein C, bilirubin, rheumatoid factor, haemolysis and triglycerides does not interfere with the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

61 samples were tested, 42 of which were positive and 19 negative with the Chorus kit in comparison with another commercial method. The discordant samples were tested using Western Blotting considered as reference. Doubtful samples were considered positive. The Diesse kit gave 2 false positive results.

The results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	40	2	42
	-	0	19	19
	Total	40	21	61

Sensitivity 100.0 % CI_{95%}: 91.2 to 100.0

Specificity 90.5 % CI_{95%}: 71.1 to 97.3

14. PRECISION**Table 1. Within-run precision**

	No. of replicates	Mean AU/ml	S. D.	CV%
Lot. no. 475	9	61.0	2.80	4.6
Lot. no. 459	9	45.3	4.90	10.8
Lot. no. 504	9	37.4	4.84	12.9

Table 2. Between-run and between-lots precision

Sample	Mean AU/ml			Mean	S. D.	CV%
	Lot 475	Lot 459	Lot 504			
MYG 1	<10	<10	<10	<10	0.00	0.0
MYG 2	15.3	14.7	14	14.7	0.65	4.4
MYG 3	41.7	39.6	40.7	40.7	1.05	2.6
MYG 4	65.7	78.3	73.7	72.6	6.38	8.8
MYG 5	>100	>100	>100	>100	0.00	0.0

15. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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16. POUŽITÁ LITERATURA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAFÍA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten: Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAPHIE

1. G.B. Wisdom : Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt : Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean : Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, str.1295-1299, Prosinac 1988.

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15. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Zv.369, str.1295-1299, December 1988.

	EN ES IT CS HR BG LT SL	Date of manufacture Fecha de fabricación Data di fabbricazione Datum výroby Datum proizvodnje Дата на производство Pagaminimo data Datum izdelave	FR EL PT RO DE ET PL	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico Data fabricatiei Herstellungsdatum Valmistamise kuupäev Data produkcji
	EN ES IT CS HR BG LT SL	Use By Fecha de caducidad Utilizzare entro Použitelné do Upotrijebiti do Използване в рамките на Sunaudoti per Uporabiti do	FR EL PT RO DE ET PL	Utiliser jusqu'à Ημερομηνία λήξης Prazo de validade A se folosi pana la Verwendbar bis Kasutage siseselt Data minimalnej trwałości
	EN ES IT CS HR BG LT SL	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso Pozor, čtěte přiložené dokumenty Pozor, vidi upute za uporabu Внимание, вижте инструкциите за употреба Dėmesio, žr. naudojimo instrukciją Pozor, glejite navodila za uporabo	FR EL PT RO DE ET PL	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída Atentie, consultați documentele insotitoare Achtung, Begleitdokumente beachten Tähelepanu, vt kasutusjuhendit Uwaga, patrz instrukcja obsługi
	EN ES IT CS HR BG LT SL	Manufacturer Fabricante Fabbricante Výrobce Proizvođač Производител Gamintojas Proizvajalec	FR EL PT RO DE ET PL	Fabricant Κατασκευαστής Fabricante Productator Hersteller Tootja Producent
	EN ES IT CS HR BG LT SL	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Obsah stačí na <n> testů Sadržaj dovoljan za „tot.“ Testova Достатъчно съдържание за "n" есета Turinio pakanka „n“ tyrimų Dovolj vsebine za "n" poskusov	FR EL PT RO DE ET PL	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios Continut sufficient pt <n> teste Inhalt ausreichend für <n> Prüfungen Piisav sisu "n" essee jaoks Zawiera wystarczającą ilość do „n” próbek
	EN ES IT CS HR BG LT SL	Temperature limitation Límite de temperatura Limiti di temperatura Teplotní omezení Temperaturne granice Температурни граници Temperatūros ribos Temperaturne omejitve	FR EL PT RO DE ET PL	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura Limita da temperatura Temperaturbegrenzung Temperatuuri piirangud Wartości graniczne temperatury
	EN ES IT CS HR BG LT SL	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso Čtěte návod k použití Pogledati upute za uporabu Вижте инструкциите за употреба Žr. naudojimo instrukcijas Glejite navodila za uporabo	FR EL PT RO DE ET PL	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização Pentru utilizare consultați instrucțiunile Gebrauchsanweisung beachten Vt kasutusjuhendit Patrz instrukcja obsługi
	EN ES IT CS HR BG LT SL	Catalogue number Número de catálogo Numero di catalogo Katalogové číslo Kataloški broj Каталожен номер Katalogo numeris Kataloška številka	FR EL PT RO DE ET PL	Référence du catalogue Αριθμός καταλόγου Referência de catálogo Numar de catalog Bestellnummer Kataloogi number Numer katalogowy

IVD	EN ES IT CS HR BG LT SL	<i>In vitro</i> Diagnostic Medical Device Producto sanitario para diagnóstico <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Lékařské vybavení pro diagnostiku <i>in vitro</i> Medicinski proizvod za dijagnostiku <i>in vitro</i> Диагностично медицинско изделие <i>in vitro</i> In vitro diagnostikos medicinos priemonė Diagnostični medicinski pripomoček <i>in vitro</i>	FR EL PT RO DE ET PL	Dispositif médical de diagnostic <i>in vitro</i> <i>In vitro</i> Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico <i>in vitro</i> Dizpositiv medical pentru diagnosticare <i>in vitro</i> <i>In-Vitro</i> -Diagnostikum In vitro diagnostiline meditsiiniseade Wyrób medyczny do diagnostyki <i>in vitro</i>
LOT	EN ES IT CS HR BG LT SL	Batch code Código de lote Codice del lotto Kód šarže Šifra serije Код на партидата Partijos kodas Koda serije	FR EL PT RO DE ET PL	Code du lot Αριθμός Παρτίδας Código do lote Lot Chargenbezeichnung Partii kood Kod partii
CE	IT EN CS DE HR BG LT SL	Marcatura CE di conformità CE marking of conformity Oznakowanie zgodności CE CE-Konformität Skenneichnung CE oznaka skladnosti Маркировка за съответствие CE EB atitikties ženkla Oznaka o skladnosti CE	ES EL FR PT RO ET PL	Marcado CE de conformidad Σημανση συμμορφωσης CE Marquage de conformité CE Marcação CE de conformidade Marcajul de conformitate CE CE-vastavusmārgis Oznaczenie zgodności CE