CHORUS

MYCOPLASMA PNEUMONIAE IgM

REF 81035 (36 tests)

Prodotto da/Manufactured by/Fabricado por/ Proizveo /
 Произведено от / Tootja / Pagaminta / Wyprodukowany
 przez :

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

INSTRUCTIONS FOR USE

CHORUS MYCOPLASMA IgM

REF 81035

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO MYCOPLASMA IN HUMAN SERUM USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS/CHORUS TRIO INSTRUMENTS

2. INTRODUCTION

Mycoplasma pneumoniae is the most common etiological agent causing pneumonia acquired in community environments, particularly in the age range of 5 to 30 years; it can cause epidemics which develop slowly as the incubation period varies from 10 to 14 days, and the infection involves subjects living in close contact or in segregated groups (schools, military barracks, family groups). Pneumonia caused by *Mycoplasma* is also known as atypical primary pneumonia or Eaton's pneumonia.

Mycoplasma pneumoniae attacks and destroys ciliated epithelial cells of the respiratory tract mucosa. Microscopic examination evidences interstitial pneumonia, bronchitis and bronchiolitis.

Due to the recurrence of symptoms caused by different etiological agents, additional diagnostic measures such as serological tests are necessary for the diagnosis of an acute pneumonia infection.

The immune response, following a *Mycoplasma pneumoniae* infection, depends on the type of infection: the IgM-class antibodies are more frequently found in case of primary infections as demonstrated by their presence in younger patients. The IgM are low or not found in older patients with higher probability of reinfection. On the other hand the IgA represent the highest sensitivity marker in case of an infection in progress, a reinfection or a recent infection.

Specific IgA-class antibodies appear early during the infection and reach peak levels during the first 4 weeks, then declining quite rapidly, earlier than the IgG; diagnosis of an acute infection can often be made on the basis of a single sample test. The IgG appear after the IgM and IgA and they generally reach the peak level during the 5° week. A high titer together with a significant increase between two samples taken approximately 2 weeks apart, confirms that an infection is in progress.

The antigen used for the serological determination of specific IgM is derived from a *M. pneumoniae* membrane extract containing significant amounts of P1 cytoadhesin which is the immunodominant antigen corresponding to a transmembrane protein which is chiefly responsible for the cytoadherence of *M. pneumoniae* to the respiratory epithelium in the host during the early phase of the infection.

3. PRINCIPLE OF THE METHOD

The test is based on the indirect ELISA (Enzyme Linked Immunosorbent Assay) principle.

The antigen, composed of an extract of inactivated *Mycoplasma pneumoniae* is bound to the solid phase and is incubated with the human serum sample diluted in a diluent which blocks the IgG antibodies; the specific IgM are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate constituted of anti-human IgM monoclonal antibodies conjugated to peroxidase.

The conjugate which has not reacted is eliminated and the peroxidase substrate is added. The blue colour which develops is proportional to the concentration of specific antibodies present in the test serum.


The disposable devices contain all the reagents necessary to perform the test when applied on the Chorus/Chorus TRIO instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: equilibrate a package at room temperature, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

**Description:**

Position 8: Space for application of bar code label

Position 7: NOT USED

Position 6: MICROPLATE WELL coated with *Mycoplasma pneumoniae* antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: DILUENT FOR THE SAMPLES

Contents: Proteic solution containing anti human IgG antibody, phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-IgM monoclonal antibodies labeled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum.

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation

DEVICES 8 weeks at 2/8°C

CALIBRATOR 8 weeks at 2/8°C

POSITIVE CONTROL 8 weeks at 2/8°C

6. PRECAUTIONS**ONLY FOR DIAGNOSTIC USE IN VITRO.****Attention:**

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS/CHORUS TRIO instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to be used to room temperature (18-30°C) before use; use within 60 min.

1. **Check that the substrate contained in well 4 is colourless. Discard devices which do not have this characteristic.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).
 The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1

NEGATIVE: when the result is < 0.9

DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS OF THE PROCEDURE

For a definitive diagnosis, the test results should be used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

Samples collected during the early phase of the infection may not have developed sufficient antibodies to be detectable; in the case of doubt, collect a new sample after 2-3 weeks and repeat the test.

It is advisable to centrifuge or filter lipemic or cloudy samples before testing.

12. ANALYTICAL SPECIFICITY

The presence of potentially interfering substances such as

- Chlamydia trachomatis (n=4)
- Chlamydia pneumoniae (n=15)
- Epstein Barr virus (n=6)
- Cytomegalovirus (n=10)
- Rheumatoid factor (n=7)
- Haemolysis(n=5)

does not interfere with the test result.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

73 samples were tested. Doubtful results obtained with both kits were considered positive. The Diesse kit gave 4 false positive and 1 false negative result.

The results are summarized in the following table:

Table 1.

		Reference		Total
		+	-	
Diesse	+	18	4	22
	-	1	50	51
Total		19	54	73

Sensitivity 94.7 % $CI_{95\%}$: 75.4 to 99.1

Specificity 92.6 % $CI_{95\%}$: 82.4 to 97.1

14. PRECISION

Table 2. Within-run precision

	No. of replicates	Mean Index	SD	CV%
Lot. no. 505	9	2.5	0.29	11.6
Lot. no. 506	9	2.1	0.27	12.9
Lot. no. 509	8	1.9	0.23	12.1

Table 3. Between-run and between-lots precision

Sample	Mean Index			Mean	SD	CV%
	Lot 505	Lot 506	Lot 509			
MYM 1	0.5	0.5	0.4	0.5	0.06	12.0
MYM 2	1.1	1.0	1.1	1.1	0.06	5.5
MYM 3	2.9	2.1	2.4	2.5	0.40	16.0

15. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAFÍA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAPHIE

1. G.B. Wisdom : Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt : Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean : Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. BIBLIOGRAFIA

- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
- 3.S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
- 4.J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
- 5.E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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15. **BIBLIOGRAFIJA**

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, str.1295-1299, Prosinac 1988.

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