# Pacific Hemostasis® APTT-XL

## I. Intended Use

Pacific Hemostasis APTT-XL is intended for use in performing the activated partial thromboplastin time (APTT) test, and for APTT-based factor assays using ellagic acid activator.

## II. Summary and Principles

The APTT is useful as a screening tool, and as a quantitative test for the intrinsic coagulation factors. It is a simple and versatile test which is sensitive to deficiencies of all plasma clotting factors except Factor VII. However, it is mainly used to detect deficiencies in Factors VIII, IX, XI, XII, and Prekallikrein.

The APTT is also commonly used to monitor heparin therapy since APTT prolongation is directly proportional to increasing amounts of heparin.<sup>12</sup>

The APTT test is performed by adding reagent containing a plasma activator and phospholipid to the test specimen. This mixture is incubated for 3 minutes at 37°C for optimum activation. Calcium chloride is added and clot formation is timed.

## III. Reagent

For *in vitro* diagnostic use.

**Composition:** 0.003% ellagic acid, 0.005% BSA, 0.30% phenol, 2.6% buffers, salts and stabilizers  
Store unopened vials at 2-8°C. **Do not freeze.** Opened vials are stable for 30 days when stored at 2-8°C.<sup>8</sup>

A yellow sediment may form after prolonged storage. Mix gently before use. Erratic values, quality control values outside established ranges, or product color variations could indicate deterioration. However, poor performance could also be due to other factors within the test system.

## IV. Specimen Collection

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22-24°C. For more details on specimen collection and storage, see NCCLS Document H21-A4.<sup>3</sup>

- Do not delay mixing the blood with anticoagulant.
- Avoid foaming the specimen.
- Use only plastic or siliconized borosilicate glass containers.
- Turbid, icteric, lipemic, or hemolyzed specimens may generate erroneous results.
- Freezing and thawing plasma that contains residual cells will generate damaged cell membranes that can affect results.
- Acute inflammatory reactions can shorten APTT results because of elevated fibrinogen.
- Plasma samples with hematocrits outside the range of 20-55% may be improperly anticoagulated and should be adjusted appropriately.

## V. Test Procedure

**Materials Provided:** APTT-XL Reagent, liquid, 10 x 10 mL or 10 x 4 mL

**Materials Required, But Not Provided:**

Pacific Hemostasis Calcium Chloride solution (0.02M)

Stopwatch or timer

Precision pipette: 0.1 mL

Normal and abnormal controls such as Pacific Hemostasis Coagulation Control Plasmas, Level 1, 2, and 3

APTT-XL is suitable for use with manual, mechanical, photo-optical, or other means of clot detection. Follow manufacturer's recommendations for proper use of instrumentation. For manual assays:

- A. Prewarm Calcium Chloride (0.02M) to 37°C.
- B. Add 0.1 mL test plasma to cuvette and prewarm to 37°C.
- C. Add 0.1 mL APTT-XL to the test plasma. Mix.
- D. Incubate the plasma-reagent mixture at 37°C for 3-5 minutes (activation time). For consistent results, test all plasmas with the same activation time.
- E. Forcibly add 0.1 mL prewarmed Calcium Chloride and time clot formation.

## VI. Quality Control

Normal and abnormal plasmas such as Pacific Hemostasis Coagulation Control Level 1, 2, and 3 should be tested in conjunction with patient plasmas. Level 1 is a normal plasma, and Levels 2 and 3 are adjusted to mimic moderately and severely deficient plasmas, respectively. A normal control and at least one abnormal control should be run at the initiation of testing each day and at least once each shift, or with each group of assays. Controls should also be tested with each reagent change or major instrument adjustment. Each laboratory should establish a control range to represent the allowable variation in day to day performance for each control.

## VII. Results

Report clotting times for each plasma to the nearest 0.1 second. A Normal Reference Range can also be reported for comparison. Do not report patient values relative to commercial control plasma clotting times. Controls are intended only for quality assurance of the test system.

## VIII. Limitations

The biochemistry of coagulation involves a series of reactions that are influenced by many pre-test conditions. These variables must be controlled to obtain reproducible results.<sup>3</sup>

### Technique

- Plasma pH will increase if exposed to air. Store samples stoppered.
- APTT-XL was designed to work at 37°C ± 0.5°C. Frequently check the temperature of all heating elements.
- All labware must be clean and free of trace amounts of detergents.
- Always follow instrument manufacturer's instructions for proper maintenance.

### Interfering Substances

- Sodium oxalate, EDTA, and heparin are not suitable anticoagulants.
- Oral contraceptives, estrogen, pregnancy, coumarin type drugs, heparin, asparaginase, and naloxone have been reported to influence APTT results.<sup>4</sup>

## IX. Expected Values

When APTT-XL was evaluated on a normal population, the following results were obtained:<sup>9</sup>

mean	± 2SD range
29.9	24.0 - 35.2
29.8	24.2 - 36.3

These values should only be used as a guideline. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent.<sup>5</sup> A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

## ENGLISH

## X. Performance Characteristics

### Heparin Sensitivity

The anticoagulant action of heparin depends on many factors, including an adequate level of Antithrombin-III, platelet activation and subsequent Platelet Factor-4 release during specimen preparation, *in vivo* presence of other medications, rate of heparin metabolism, mode of heparin administration, and delayed specimen handling. While recognizing these variables, the laboratory can determine the relative sensitivity of a given reagent to heparin by adding known amounts of heparin to pooled normal plasma and performing an APTT. For example, the following results were obtained on a photo-optical instrument with one lot of APTT-XL reagent:<sup>10</sup>

Heparin conc. (units/mL)	APTT (seconds)
0.0	28.8
0.1	38.3
0.2	50.1
0.3	63.1
0.4	80.9
0.5	98.0

Each laboratory should establish its own heparin sensitivity curve using the same heparin source used for therapy in that institution. Variations can result from different brands of heparin, tissue origin, and salt forms.<sup>12,5</sup>

### Factor Sensitivity

An APTT reagent with adequate sensitivity should demonstrate a prolonged clotting time in samples having ≤ 30-40% factor activity.<sup>12</sup> APTT-XL was evaluated on mildly and severely deficient plasmas with the following results:<sup>11</sup>

Factor	% activity	APTT (seconds)
VIII	<1%	82.0
VIII	20%	44.8
IX	<1%	83.5
IX	20%	40.9
XI	<1%	134.2
XI	20%	47.8
XII	<1%	>200
XII	20%	36.2
Prekallikrein	<1%	69.5

Furthermore, the sensitivity of APTT-XL to Factor VIII has been determined as follows:<sup>12</sup>

% Factor VIII	APTT (seconds)
100%	32.5
70%	34.0
50%	36.9
40%	38.9
30%	40.8
20%	44.4
10%	50.6
5%	56.1
1%	68.1
<1%	83.6

These values should only be used as guidelines. Each laboratory should establish sensitivity to individual factors using instruments, reagents, and techniques used in their laboratory.

### XI. References

1. Brandt, J.T., Triplett, D.A.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated Partial Thromboplastin Time. *Amer J Clin Path* 76:530, 1981.
2. Thompson, J.M.: The Control of Heparin Therapy by the Activated Partial Thromboplastin Time. Sensitivity of Various Thromboplastins to Heparin. *Standardization of Coagulation Assays: An Overview*. Edited by D.A. Triplett, College of American Pathologists, Skokie, Ill. 1982, pp 195.
3. NCCLS: *Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
4. Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. *Clin Chem* 18:1041, 1972.
5. Banez, E.L., Triplett, D.A., Koepke, J.: Laboratory Monitoring of Heparin Therapy. The Effect of Different Salts of Heparin on the Activated Partial Thromboplastin Time. *Amer J Clin Path* 74:569, 1980.
6. Wujastyk, J., Triplett, D.A.: Selecting Instrumentation and Reagents for the Coagulation Laboratory. *Pathologist* 37:398, 1983.
7. Christensen, R.L., Triplett, D.A.: Factor Assay (VIII and IX) Results in the College of American Pathologists Survey Program (1980-1982). *Amer J Clin Path* 80 (Suppl): 633, 1983.
- 8-12. Daten in 510(K)-Datei.

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#### XI. Bibliografia

1. Brandt, J.T., Triplett, D.A.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated Partial Thromboplastin Time. *Amer J Clin Path* 76:530, 1981.
2. Thompson, J.M.: The Control of Heparin Therapy by the Activated Partial Thromboplastin Time. Sensitivity of Various Thromboplastins to Heparin. *Standardization of Coagulation Assays: An Overview*. Edited by D.A. Triplett, College of American Pathologists, Skokie, Ill. 1982, pp 195.
3. NCCLS: *Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
4. Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. *Clin Chem* 18:1041, 1972.
5. Banez, E.L., Triplett, D.A., Koepke, J.: Laboratory Monitoring of Heparin Therapy. The Effect of Different Salts of Heparin on the Activated Partial Thromboplastin Time. *Amer J Clin Path* 74:569, 1980.
6. Wujastyk, J., Triplett, D.A.: Selecting Instrumentation and Reagents for the Coagulation Laboratory. *Pathologist* 37:398, 1983.
7. Christensen, R.L., Triplett, D.A.: Factor Assay (VIII and IX) Results in the College of American Pathologists Survey Program (1980–1982). *Amer J Clin Path* 80 (Suppl): 633, 1983.
- 8-12. Informazioni ottenute dal file 510(k).

Symbols Key	
	Manufacturer Hersteller Fabricant Fabricante Fabbricante
	IVD <i>In Vitro Diagnostic Medical Device</i> Medizinprodukt für die <i>in-vitro-Diagnostik</i> Matériel médical pour utilisation diagnostique <i>in vitro</i> Dispositivo medico per diagnosi <i>in vitro</i>
	LOT Lot Number Chargenummer Numéro de lot Número de lote Numero di lotto
	Use By Verfallsdatum Utiliser jusque Fecha de caducidad Da utilizzare entro
	Temperature Limitation Temperatur einschränkungen Limite de température Límite de temperatura Limiti di temperatura
	CE CE Mark CE-Markierung Marquage CE Marca CE Marchio CE
	REF Catalogue Number Katalognummer Référence catalogue Número de catálogo Numero di catalogo
	i Consult Instructions for Use Bedienungsanleitung lesen Consulter le manuel d'utilisation Consultar las instrucciones de uso Consultare le istruzioni per l'uso
	EC REP Authorized Representative in the European Community Autorisierte Vertretung in der Europäische Gemeinschaft Représentant agréé pour La Communauté européenne Representante autorizado en la Comunidad Europea Rappresentante autorizzato nella Comunità europea

||
||
||