CerTest Rapid Test

CerTest
BIOTEC

CERTEST SARS-CoV-2

ONE STEP SARS-CoV-2 Ag CARD TEST

CE IVD

The test should be used by professional trained in IVD devices

These instructions for use apply to the following references:

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SC820001P; SC820001PC; SC820001PCC

Initial part of the Instructions for use is common for all the references/

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Specific parts for each of the references are described in the corresponding annexes (see tables below)/

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
ENGLISH

2. Intended use	7
3. Summary and Explanation	7
4. Test principle	7
5. Storage and stability	8
6. Precautions for users	8
7. Specimen collection and preparation	9
8. Recommendations	10
9. Limitations of the test	10
10. Expected values	11
11. Performance characteristics	11
A. Analytical sensitivity (detection limit)	11
B. Clinical sensitivity and specificity	11
C. Hook Effect	12
D. Cross reactivity	13
E. Interferences	14
F. Repeatability and Reproducibility	14
ANNEX 1	15
I. Materials	15
II. Test procedures (samples)	15
III. Interpretation of results	17
IV. Quality Control	17
ANNEX 2	18
I. Materials	18
II. Test procedures (samples)	18
III. Interpretation of results	20
IV. Quality Control	21
A. External control. CerTest SARS-CoV-2 Positive Control	21
A. 1 Positive control procedure	21
A.1.1 Interpretation positive control procedure	22
ANNEX 3	23
I. Materials	23
II. Test procedures (samples)	23

III.	Interpretation of results	25
IV.	Quality Control	26
	A. External positive control. CerTest SARS-CoV-2 Positive Control.....	26
	A. 1 Positive control procedure	26
	A.1.1 Interpretation positive control procedure.....	27
	B. External negative control.....	28
	B. 1 Negative control procedure. CerTest Negative Control	28
	B.1.1. Interpretation negative control procedure.....	29

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ENGLISH

2. Intended use

CerTest SARS-CoV-2 one step card test is a coloured chromatographic immunoassay for the qualitative detection of nucleoprotein (or N protein) antigen of SARS-CoV-2 from nasopharyngeal swab or nasal samples from patients suspected of COVID-19 infection.

CerTest SARS-CoV-2 card test offers a simple and a highly sensitive screening assay as an aid in the diagnosis of SARS-CoV-2 infection. This product is for professional *in vitro* diagnosis use only.

3. Summary and Explanation

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus.

The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Coronaviruses (CoVs) are RNA, enveloped viruses, with a large genome (29.9 kb), amongst which β -CoVs and α -CoVs can infect mammals. Viruses belonging to the Coronaviridae family have the most abundant membrane protein (M) among other proteins, i.e. spike glycoprotein (S), nucleocapsid protein (N) and an envelope protein (E). Spike glycoprotein (S) is one of the targets of T cell response in the immune system. The S protein also promotes the binding of the virus envelope to the ACE2 receptor and the entry of the virus into the cell target. ACE2 receptors are present in cells of the arteries, veins, smooth muscles, small intestine, alveoli of the lungs, hair follicles, cardiac myofibroblasts, skin, brain, and kidney, thus SARS-CoV-2 could potentially infect these tissues.

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of the viral infection include a sore throat, aches, pains, and shortness of breath. In some cases, runny nose, nausea, and diarrhea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks. Due to the high infectious rate of SARS-CoV-2, detection of asymptomatic positive patients is probably one of the key points to controlling the outbreak.

4. Test principle

CerTest SARS-CoV-2 test is based on the principle of a qualitative immunochromatographic assay for the qualitative detection of SARS-CoV-2 from nasopharyngeal swab samples or nasal swab samples from patients suspected of COVID-19 infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse

monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no SARS-CoV-2 antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-SARS-CoV-2 antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

5. Storage and stability

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

6. Precautions for users

- **For professional *in vitro* diagnostic use.**
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous (especially samples from patients suspected of SARS-CoV-2 infection) and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- Use proper infection control practices when collecting and handling specimens. These practices should include, but are not limited to, personal protective equipment (PPE), such as a laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Follow local or national regulations regarding collecting and handling specimens.
- The tests (swabs, vials, devices) should be discarded in a proper biohazard container after testing. These containers should be discarded in accordance with local or national laws or regulations.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the **CerTest SARS-CoV-2** test. Do not use any other commercial kit component.

- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal or nasal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- All positive results should be processed following local or national laws and regulations.


7. Specimen collection and preparation

Samples should be processed as soon as possible after collection. Follow proper infection control practices. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for a maximum of 8 hours prior to testing. If the samples (only for nasopharyngeal samples) are preserved in validated transport media (VTM, UTM o Saline Buffer) could be preserved on it until 6 hours at room temperature or in the refrigerator (2-8°C).

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the sterile swab from its packing.
 2. Collect specimen with a sterile swab from one nostril.
 3. Insert the swab into the nostril to the nasopharynx, rotating against the wall (4 times) to ensure swab contains cells as well as mucus. Remove the swab from the nostril carefully.
 4. Consider repeating procedure using other nostril only if the protocol of the professional taking the sample required do it.
 5. Process the swab (sample) as soon as possible after collecting the specimen.
- Follow the test procedure (samples).

Nasopharyngeal Specimen collection

Direct nasopharyngeal swab is the preferred protocol.

Nasopharyngeal samples previously extracted by transport media:


Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTM, UTM or Saline Buffer.

1. Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0 mL. Using higher volumes will affect sensitivity of the system.
 2. Dilute the extracted sample 1:1 in the provided sample diluent (DIL Respiratory, 0.5mL/vial). Do not use the provided swabs for the dilution of these kind of samples in the sample diluent.
- Follow immediately with the point 3 of test procedure (samples).

Nasal swab method:

The device can be also used with samples from nasal swabs.

1. Remove the sterile swab from its packing.
 2. Collect specimen with a sterile swab from one nostril. Recommended blow the nose once using a tissue before collection of the specimen.
 3. Insert the swab into the nostril to the nasopharynx (approx. 2cm), rotating against the nasal wall several times (5 times) to ensure swab contains cells as well as mucus. Remove the swab from the nostril carefully and repeat the procedure with the other nostril.
 4. Once the process for both nostrils is finished, process the swab (sample) as soon as possible.
- Follow the test procedure (samples).

Nasal Specimen collection

8. Recommendations

Recommendations World Health Organization: Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 September 2020.

- 1) To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.
- 2) Where possible, all positive samples giving positive results should be transported to laboratories with NAAT (Nucleic Acid Amplification Test) capability for confirmatory testing.
- 3) This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT- confirmed COVID-19 outbreaks.
- 4) To monitor trends in disease incidence in communities.
- 5) For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
- 6) A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.

9. Limitations of the test

- 1) The test must be carried out within 2 hours after opening the sealed bag.
- 2) **CerTest SARS-CoV-2** test should be used only with nasopharyngeal swab samples or nasal swabs samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal or nasal specimens must be obtained.
- 3) The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- 4) Positive results determine the presence of SARS-CoV-2 antigens. A diagnosis of SARS-CoV-2 infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.

- 5) Positive results do not rule out co-infections with other pathogens.
- 6) Only VTM, UTM and Saline buffer transport media have been validated with **CerTest SARS-CoV-2** test. When using transport media, the test sensitivity can be reduced due to a higher dilution of sample. It is not recommended for Nasal samples.
- 7) A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal or nasal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique.
- 8) **CerTest SARS-CoV-2** has been used with samples from asymptomatic patients, detecting those with high viral loads.
- 9) To obtain accurate results, do not use bloody nasopharyngeal or nasal samples.
- 10) Very viscous nasopharyngeal/nasal samples may cause non-specific reaction in the test.

10. Expected values

On January 30, 2020, the World Health Organization declared COVID-19 a public health emergency of international concern. Since the onset of the pandemic, SARS-CoV-2 has affected millions of people, causing thousands of victims worldwide.

Recent studies have shown that the primary viral infection found in symptomatic patients was identical to that seen in asymptomatic patients, suggesting that COVID-19 patients with asymptomatic infections are capable of transmitting SARS-CoV-2. During the incubation period, asymptomatic patients are also likely to transmit the virus to healthy people.

Although SARS-CoV-2 infection can affect people of any age, children and the elderly with asthma are in a higher risk category for the disease. Those with underlying health problems such as cancer, chronic respiratory diseases, cardiovascular diseases, and diabetes are also reported to have a higher chance of developing serious illnesses.

11. Performance characteristics

A. Analytical sensitivity (detection limit)

Detection limit value of **CerTest SARS-CoV-2** is 1.25×10^2 PFU/mL nCoV-2019 D614G(S).

B. Clinical sensitivity and specificity

A multi-center evaluation with **nasopharyngeal samples** was performed comparing the results obtained by an immunochromatographic test (**CerTest SARS-CoV-2**, CerTest) vs qPCR technique (nasopharyngeal samples). The results were as follows:

		qPCR technique (nasopharyngeal samples)		
IC test: CerTest SARS-CoV-2 (nasopharyngeal samples)		+	-	Total
	+	93	1	94
	-	7	455	462
	Total	100	456	556

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity	93.0%	86.1–97.1%
Specificity	99.8%	98.8–100.0%
PPV	98.9%	94.2–100.0%
NPV	98.5%	96.9–99.4%

This multi-center evaluation (nasopharyngeal samples) with positive samples with Ct< 28 (*) showed the following results: sensitivity 95.1% (95% confidence interval: 88.0-98.7%) and specificity: 99.8% (95% confidence interval: 98.8-100.0%).

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays* (11 September 2020) from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

The results showed a high sensitivity and specificity to detect SARS-CoV-2 using **CerTest SARS-CoV-2**.

An evaluation using **nasal samples** was performed comparing the results obtained by an immunochromatographic test (**CerTest SARS-CoV-2**, Certest) using nasal samples vs qPCR technique (using nasopharyngeal samples). The results were as follows:

		qPCR technique (nasopharyngeal samples)		
IC test: CerTest SARS-CoV-2 (nasal samples)		+	-	Total
	+	129	7	136
	-	27	827	854
	Total	156	834	990

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity	82.7%	75.8 - 88.3%
Specificity	99.2%	98.3 - 99.7%
PPV	94.9%	89.7 - 97.9%
NPV	96.8%	95.4 - 97.9%

This evaluation (nasal samples for the evaluation with the rapid test) with positive samples with Ct< 28 showed the following results: sensitivity 96.9% (95% confidence interval: 91.2-99.4%) and specificity: 99.2% (95% confidence interval: 98.3-99.7%).

The results showed a high sensitivity and specificity to detect SARS-CoV-2 using **CerTest SARS-CoV-2**.

C. Hook Effect

CerTest SARS-CoV-2 test does not show full inhibitory hook effect at concentration as high as 10⁵xLoD.

D. Cross reactivity

An evaluation was performed to determine the cross reactivity of **CerTest SARS-CoV-2**; no cross reactivity against organism, pathogens, substances that could cause infections:

Virus	Adenovirus	Coronavirus Strain 229E	Coronavirus Strain OC43	Influenza B	Norovirus GII	Rhinovirus
	Astrovirus	Coronavirus Strain HKU1	Enterovirus	Metapneumovirus human (hMPV)	Parainfluenza virus	Rotavirus
	Bocavirus	Coronavirus Strain NL63	Influenza A	Norovirus GI	Respiratory Syncytial Virus	
Bacteria	<i>Bordetella pertussis</i>	<i>Clostridium difficile</i> Toxin B	<i>Listeria Monocytogenes</i>	<i>Salmonella enteritidis</i>	<i>Shigella dysenteriae</i>	<i>Streptococcus pneumococcal</i>
	<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> O:157	<i>Mycobacterium tuberculosis</i>	<i>Salmonella paratyphi</i>	<i>Shigella flexneri</i>	<i>Streptococcus pneumoniae</i>
	<i>Chlamydia pneumoniae</i>	<i>Haemophilus influenzae</i>	<i>Mycoplasma pneumoniae</i>	<i>Salmonella typhi</i>	<i>Shigella sonnei</i>	<i>Streptococcus pyogenes</i>
	<i>Clostridium difficile</i> antigen GDH	<i>Helicobacter pylori</i>	<i>Mycoplasma tuberculosis</i>	<i>Salmonella typhimurium</i>	<i>Staphylococcus aureus</i>	<i>Yersinia enterocolitica</i> O:3
	<i>Clostridium difficile</i> Toxin A	<i>Legionella pneumophila</i>	<i>Pneumocystis jirovecii</i>	<i>Shigella boydii</i>	<i>Staphylococcus epidermis</i>	<i>Yersinia enterocolitica</i> O:9
Parasites	<i>Cryptosporidium Parvum</i>	<i>Giardia lamblia</i>	<i>Entamoeba histolytica</i>			
Others	Human calprotectin	Human haemoglobin	Human lactoferrin	Human transferrin	Pig haemoglobin	Bovine haemoglobin
	Pooled human nasal wash-representative of normal respiratory microbial flora					

CerTest SARS-CoV-2 showed some cross reaction with SARS and MERS.

E. Interferences

An evaluation was performed to determine the possible interferences of **CerTest SARS-CoV-2**; no interferences against the substances tested were detected:

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (Fortasec)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (Hibor)	350.0 IU/mL	No interference
	Oseltamivir	3.0·10 ⁻² mg/mL	No interference	Almagato (Almax)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfomicin (Monurol)	3.0·10 ⁻³ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (Fluimucil)	3.0·10 ⁻² mg/mL	No interference
	Codeine (Toseina)	0.2 mg/mL	No interference	Dexketoprofen trometamol (Enantyum)	0.3 mg/mL	No interference
	Benzocaine (Angileptol)	3.0·10 ⁻² mg/mL	No interference	Levofloxacin	3.0·10 ⁻² mg/mL	No interference
	Cloperastine (Flutox)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine (Iniston mucolítico)	3.0·10 ⁻² mg/mL	No interference	Rifampicin (Rifaldin)	0.3 mg/mL	No interference
	Loratadine	0.3 mg/mL	No interference	Phenoxymethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (Polaramine)	0.3 mg/mL	No interference	Ambroxol hydrochloride (Mucosan)	0.3 mg/mL	No interference
	Ebastine (Ebastel)	3.0 mg/mL	No interference	Macrogol 3350 (Movicol)	3.0 mg/mL	No interference
	Acetyl Salicylic (Adiro)	0.3 mg/mL	No interference	Lysine Carbocysteinate (Pectox)	3.0·10 ⁻² mg/mL	No interference
	Ibuprofen (Epidifen)	0.3 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (Dolocatil)	5.0 mg/mL	No interference	Lorazepam	3.0·10 ⁻³ mg/mL	No interference
	Metamizole (Nolotil)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	2.0·10 ⁻³ mg/mL	No interference	Biotine	100.0 µg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.51 mg/mL	No interference
	CVS Nasal Spray (Cromolyn)	4 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
	Afrin (Oxymetazoline)	0.05 mg/mL	No interference	Mupirocin	0.025 mg/mL	No interference
	CVS Nasal Drops (Phenylephrine)	10 mg/mL	No interference	Fluticasone Propionate	0.05 mg/mL	No interference
	ZICAM	0.1 mg/mL	No interference	Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL	No interference
	Homeopathic	DIL 1/10	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	5.0 mg/mL	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	50.0 mg/mL	No interference
Others	HAMAs (Human anti-mouse antibodies)	10 µg/mL	No interference			

	Potencial interfering substance	Concentration (Ct)	Results	Potencial interfering substance	Concentration (Ct)	Results
Microbial interferences	STREP A	23.7/22.7/23.35	No interference	Influenza A	23.47/22.05/22.87	No interference
	ADENO	24.1/22.8/19.2	No interference	Influenza B	32.59/35.61/25.36	No interference
	RSV	29.28/27.78/38.82	No interference			

F. Repeatability and Reproducibility

Following the study of repeatability and reproducibility of **CerTest SARS-CoV-2**, performed using different internal samples, negative and positive, no differences have been observed within the evaluations.

ANNEX 1

KIT 20 CARD SARS-CoV-2. Catalogue reference: SC820001P

Annex for the following reference:

PRODUCT	REFERENCE	UDI-DI
CerTest SARS-CoV-2 Card (Ag 20 tests kit with swabs + 20 vials)	SC820001P	08435440215158

Table A 2. Reference for kit 20 card SARS-CoV-2.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest SARS-CoV-2 tests (card test) - 20x DIL Respiratory vials (0.5mL/vial). Sample extraction vials with diluent (sample and control diluent) - 20x Sterile swabs - 1x Laboratory rack - 1x IFU: Instructions for use (including Quick guide use)	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Vortex (optional)

II. Test procedures (samples)




Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

- Put the DIL Respiratory vials in the laboratory rack (identify properly samples from the different patients).
- Open the blue cap of the DIL respiratory vial (prefilled with sample diluent) and immediately put the swab into the DIL Respiratory vial. Mix and homogenize the solution by rotating the swab forcefully against the side of the DIL Respiratory vial to allow the solution to mix for 15 seconds. Best results are obtained when it is mixed to extract as much liquid as possible from the swab avoiding splashes and aerosols. Rotate the swab against the side of the DIL Respiratory vial as the swab is withdrawn to extract as much liquid as possible. Discard the swab. Close the extraction vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. The sample can be shaken using a vortex (this is optional) (1).
- Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.
- Be sure that use a separate DIL Respiratory vial and test for each sample. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120µL) from the DIL Respiratory

vial into the circular window marked with the letter S (3). After adding the drops into the sample window close the white-translucent cap (4).

5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

The Sample diluent buffer (DIL Respiratory) has SARS-CoV-2 inactivation capability: 99.68% after 1 minute and 99.98% after 10 minutes.


(1) Put the swab into the DIL Respiratory vial, rotating 15 seconds and extract the liquid. Close the vial with the blue cap. Please be sure the caps are properly closed. The sample can be shaken (vortex).

(2) Open the white-translucent cap

(3) Add 3 full drops in the circular window (S)


(4) Close the white-translucent cap.

Wait 10 minutes

III. Interpretation of results

(please refer to the illustration below)

	SARS-CoV-2		Interpretation of the results
1.	NEGATIVE		There is no SARS-CoV-2 antigen presence in the sample.
	GREEN		
2.	POSITIVE		There is SARS-CoV-2 antigen presence in the sample. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN-RED		
3.	INVALID		Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result		

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control.

INTERNAL CONTROL: control line

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

ANNEX 2

KIT 20 CARD SARS-CoV-2. Catalogue reference: SC820001PC

Annex for the following reference:

PRODUCT	REFERENCE	UDI-DI
CerTest SARS-CoV-2 Card (Ag 20 tests kit with swabs + 20 vials + PC)	SC820001PC	08435440214014

Table A 2. Reference for KIT 20 CARD SARS-CoV-2 with positive control

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest SARS-CoV-2 tests (card test) - 20x DIL Respiratory vials (0.5mL/vial). Sample extraction vials with diluent (sample and control diluent) - 20x Sterile swabs - 1x CerTest SARS-CoV-2 Positive Control - 1x Laboratory rack - 1x IFU: Instructions for use (including Quick guide use)	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Vortex (optional)




II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Put the DIL Respiratory vials in the laboratory rack (identify properly samples from the different patients).
2. Open the blue cap of the DIL respiratory vial (prefilled with sample diluent) and immediately put the swab into the DIL Respiratory vial. Mix and homogenize the solution by rotating the swab forcefully against the side of the DIL Respiratory vial to allow the solution to mix for 15 seconds. Best results are obtained when it is mixed to extract as much liquid as possible from the swab avoiding splashes and aerosols. Rotate the swab against the side of the DIL Respiratory vial as the swab is withdrawn to extract as much liquid as possible. Discard the swab. Close the extraction vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. The sample can be shaken using a vortex (this is optional) (1).
3. Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.

4. Be sure that use a separate DIL Respiratory vial and test for each sample. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120 μ L) from the DIL Respiratory vial into the circular window marked with the letter S (3). After adding the drops into the sample window close the white-translucent cap (4).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

The Sample diluent buffer (DIL Respiratory) has SARS-CoV-2 inactivation capability: 99.68% after 1 minute and 99.98% after 10 minutes.


(1) Put the swab into the DIL Respiratory vial, rotating 15 seconds and extract the liquid. Close the vial with the blue cap. Please be sure the caps are properly closed. The sample can be shaken (vortex).

(2) Open the white-translucent cap

(3) Add 3 full drops in the circular window (S)


Wait 10 minutes

(4) Close the white-translucent cap.

III. Interpretation of results

(please refer to the illustration below)

	SARS-CoV-2		Interpretation of the results
1.	NEGATIVE		There is no SARS-CoV-2 antigen presence in the sample.
	GREEN		
2.	POSITIVE		There is SARS-CoV-2 antigen presence in the sample. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN-RED		
3.	INVALID		Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result		

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

The **CerTest SARS-CoV-2** test presents two controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL, POSITIVE CONTROL SWAB:

External control: CerTest SARS-CoV-2 Positive Control (positive control swab).

A. External control. CerTest SARS-CoV-2 Positive Control


The **CerTest SARS-CoV-2 Positive Control** is an external quality control for **CerTest SARS-CoV-2** card test. Use this control swab to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C. CerTest SARS-CoV-2 Positive Control is stable through the expiration date printed on the sealed pouch. The positive control swab must remain in the sealed pouch until use. Do not freeze.

A. 1 Positive control procedure

Allow tests and positive control swab to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Open the blue cap of DIL Respiratory vial (prefilled with sample diluent) and immediately put CerTest SARS-CoV-2 Positive Control (swab) into the DIL Respiratory vial.
2. Mix the solution by rotating the positive control swab forcefully against the side of the DIL Respiratory vial at least 15 seconds. Extract as much liquid as possible from CerTest SARS-CoV-2 Positive Control rotating the swab against the wall of the DIL Respiratory vial as the swab is withdrawn. Discard CerTest SARS-CoV-2 Positive Control. Close the DIL Respiratory vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. CerTest SARS-CoV-2 Positive Control extraction can be shaken using a vortex (this is optional) (1).
3. Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.
4. Use a separate DIL Respiratory vial and test for each control. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120µL) from the DIL Respiratory vial into the circular window marked with the letter S (3). After adding the drops in the sample window, close the white-translucent cap.
5. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes. Any results read after more than 10 minutes will be considered invalid.

(1)


Introduce **CerTest SARS-CoV-2 Positive Control** into the DIL Respiratory vial, rotating 15 seconds and extract the liquid. Close the blue cap (be sure both caps are properly closed). The CerTest SARS-CoV-2 Positive Control extraction can be shaken (vortex optional)

(2)

Open the white-Translucent cap.


(3) Add 3 full drops in the circular window

Wait 10 minutes

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure

SARS-CoV-2 Positive Control		Interpretation of the results	
1.	POSITIVE		Control + SARS-CoV-2 POSITIVE: A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of SARS-CoV-2 antigens in the positive control swab.
	GREEN-RED		
2.	INVALID		INVALID: Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab. If the situation still persists, discontinue using the test kit and contact your local distributor.
	Any other result		

ANNEX 3

KIT 20 CARD SARS-CoV-2. Catalogue reference: SC820001PCC

Annex for the following reference:

PRODUCT	REFERENCE	UDI-DI
CerTest SARS-CoV-2 Card (Ag 20 tests kit with swabs + 20 vials + PC + NC)	SC82001PCC	08435440215097

Table A 3. Reference for KIT 20 CARD SARS-CoV-2 with positive and negative control

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest SARS-CoV-2 tests (card test) - 20x DIL Respiratory vials (0.5mL/vial). Sample extraction vials with diluent (sample and control diluent) - 20x Sterile swabs - 1x CerTest SARS-CoV-2 Positive Control - 1x CerTest Negative Control - 1x Laboratory rack - 1x IFU: Instructions for use (including Quick guide use)	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Vortex (optional)




II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Put the DIL Respiratory vials in the laboratory rack (identify properly samples from the different patients).
2. Open the blue cap of the DIL respiratory vial (prefilled with sample diluent) and immediately put the swab into the DIL Respiratory vial. Mix and homogenize the solution by rotating the swab forcefully against the side of the DIL Respiratory vial to allow the solution to mix for 15 seconds. Best results are obtained when it is mixed to extract as much liquid as possible from the swab avoiding splashes and aerosols. Rotate the swab against the side of the DIL Respiratory vial as the swab is withdrawn to extract as much liquid as possible. Discard the swab. Close the extraction vial with its blue cap, please be sure

3. that both caps (blue and white-translucent caps) are properly closed. The sample can be shaken using a vortex (this is optional) (1).
4. Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.
5. Be sure that use a separate DIL Respiratory vial and test for each sample. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120 μ L) from the DIL Respiratory vial into the circular window marked with the letter S (3). After adding the drops into the sample window close the white-translucent cap (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

The Sample diluent buffer (DIL Respiratory) has SARS-CoV-2 inactivation capability: 99.68% after 1 minute and 99.98% after 10 minutes.


(1) Put the swab into the DIL Respiratory vial, rotating 15 seconds and extract the liquid. Close the vial with the blue cap. Please be sure the caps are properly closed. The sample can be shaken (vortex).

(2) Open the white-translucent cap

(3) Add 3 full drops in the circular window (S)

Wait 10 minutes

(4) Close the white-translucent cap.

III. Interpretation of results

(please refer to the illustration below)

	SARS-CoV-2		Interpretation of the results
1.	NEGATIVE		There is no SARS-CoV-2 antigen presence in the sample.
	GREEN		
2.	POSITIVE		There is SARS-CoV-2 antigen presence in the sample. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN-RED		
3.	INVALID		Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result		

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

The **CerTest SARS-CoV-2** test presents three controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROLS: POSITIVE CONTROL SWAB AND NEGATIVE CONTROL SWAB:

External positive control: CerTest SARS-CoV-2 Positive Control (positive control swab).

External negative control: CerTest Negative Control (negative control swab).

A. **External positive control. CerTest SARS-CoV-2 Positive Control**


The CerTest SARS-CoV-2 Positive Control is an external quality control for **CerTest SARS-CoV-2 card test**. Use this control swab to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C. CerTest SARS-CoV-2 Positive Control is stable through the expiration date printed on the sealed pouch. The positive control swab must remain in the sealed pouch until use. Do not freeze.

A. 1 **Positive control procedure**

Allow tests and positive control swab to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Open the blue cap of DIL Respiratory vial (prefilled with sample diluent) and immediately put CerTest SARS-CoV-2 Positive Control (swab) into the DIL Respiratory vial.
2. Mix the solution by rotating the positive control swab forcefully against the side of the DIL Respiratory vial at least 15 seconds. Extract as much liquid as possible from CerTest SARS-CoV-2 Positive Control rotating the swab against the wall of the DIL Respiratory vial as the swab is withdrawn. Discard CerTest SARS-CoV-2 Positive Control. Close the DIL Respiratory vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. CerTest SARS-CoV-2 Positive Control extraction can be shaken using a vortex (this is optional) (1).
3. Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.
4. Use a separate DIL Respiratory vial and test for each control. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120µL) from the DIL Respiratory vial into the circular window marked with the letter S (3). After adding the drops in the sample window, close the white-translucent cap.
5. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes. Any results read after more than 10 minutes will be considered invalid.

(1)


Introduce **CerTest SARS-CoV-2 Positive Control** into the DIL Respiratory vial, rotating 15 seconds and extract the liquid. Close the blue cap (be sure both caps are properly closed). The CerTest SARS-CoV-2 Positive Control extraction can be shaken (vortex optional)

(2)

Open the white-Translucent cap.

(3) Add 3 full drops in the circular window (S)

Wait 10 minutes

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure

	SARS-CoV-2 Positive Control		Interpretation of the results
1.	POSITIVE		Control + SARS-CoV-2 POSITIVE: A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of SARS-CoV-2 antigens in the positive control swab.
	GREEN-RED		
2.	INVALID		INVALID: Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test device and a new positive control swab. If the situation still persists, discontinue using the test kit and contact your local distributor.
	Any other result		

B. External negative control


The CerTest Negative Control is an external quality control for **CerTest SARS-CoV-2** card test. Use this control swab to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C. The negative control swab is stable through the expiration date printed on the sealed pouch. The negative control swab must remain in the sealed pouch until use. Do not freeze.

B. 1 Negative control procedure. **CerTest Negative Control**

Allow tests and negative control swab to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Open the blue cap of DIL Respiratory vial (prefilled with sample diluent) and immediately put CerTest Negative Control (negative control swab) into the DIL Respiratory vial.
2. Mix the solution by rotating negative control swab forcefully against the side of the DIL Respiratory vial at least 15 seconds. Extract as much liquid as possible from CerTest Negative Control rotating the swab against the wall of the DIL Respiratory vial as the swab is withdrawn. Discard CerTest Negative Control. Close the DIL Respiratory vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. The negative control extraction can be shaken using a vortex (this is optional) (1).
3. Remove the **CerTest SARS-CoV-2** card test from its sealed bag just before using it.
4. Use a separate DIL Respiratory vial and test for each control. Open the white-translucent cap of the DIL Respiratory vial (2) and dispense 3 full drops (approx. 100-120µL) from the DIL Respiratory vial into the circular window marked with the letter S (3). After adding the drops in the sample window, close the white-translucent cap.
5. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes. Any results read after more than 10 minutes will be considered invalid.


(1)

Introduce the **CerTest Negative Control** into the DIL Respiratory vial rotating 15 seconds and extract the liquid. Close the blue cap (be sure both caps are properly closed). The negative control extraction can be shaken (vortex optional)

(2)

Open the white-Translucent cap.

(3) Add 3 full drops in the circular window (S)

Wait 10 minutes

Negative controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

B.1.1. Interpretation negative control procedure

Negative Control		Interpretation of the results
1.	VALID RESULT	Control - SARS-CoV-2 NEGATIVE: Only A GREEN line (control line (C)) appear across the results window during the test performance.
	GREEN (only control band appear)	
2.	INVALID	INVALID: Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or present of both lines (RED test and GREEN control lines). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test device and a new CerTest Negative Control. If the situation still persists, discontinue using the test kit and contact your local distributor.
	Any other result	

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