VIASURE

Real Time PCR Detection Kit

HLA Celiac

CE IVD

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PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE HLA Celiac Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS- CEL106L
VIASURE HLA Celiac Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS- CEL106H
VIASURE HLA Celiac Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS- CEL112L
VIASURE HLA Celiac Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS- CEL112H
VIASURE HLA Celiac Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS- CEL136
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS- CEL101L
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS- CEL101H
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS- CEL101

Table A 1. References for Open and Rotor-Gene format with internal control products. / *Referencias para productos Open y Rotor-Gene Format con control interno.*

TUBE FORMAT WITH INTERNAL CONTROL (SEE ANNEX 2) / FORMATO TUBO CON CONTROL INTERNO (VER ANEXO 2)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE HLA Celiac Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS- CEL148T

Table A 2. References for Tube format with internal control products. / *Referencias para productos formato Tubo con control interno.*

Content

1.	Intended use.....	6
2.	Summary and Explanation	6
3.	Principle of the procedure	7
4.	Reagents provided	8
5.	Reagents and equipment to be supplied by the user	8
6.	Transport and storage conditions.....	9
7.	Precautions for users	9
8.	Test procedure	10
8.1.	Specimen collection, transport and storage	10
8.2.	DNA extraction	11
9.	Result interpretation	11
9.1.	References with internal control (references in Annex 1 and 2).....	11
10.	Limitations of the test	13
11.	Quality control.....	14
12.	Performance characteristics.....	14
12.1.	Clinical sensitivity and specificity	14
12.2.	Analytical sensitivity	15
12.3.	Analytical specificity	18
12.4.	Analytical reactivity	18
	ANNEX 1	19
A1.1	Principle of the procedure.....	19
A1.2	Reagents provided.....	19
A1.3	Test procedure	20
A1.3.1	Lyophilized positive control.....	20
A1.3.2	PCR protocol	21
	ANNEX 2	22
A2.1	Principle of the procedure.....	22
A2.2	Reagents provided.....	22
A2.3	Test procedure	23
A2.3.1	Lyophilized positive control.....	23
A2.3.2	Lyophilized reaction mix tube.....	23
A2.3.3	PCR protocol	23

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ENGLISH

1. Intended use

VIASURE *HLA Celiac* Real Time PCR Detection Kit is a real-time RT-PCR multiplexed test designed for the simultaneous qualitative detection and differentiation of the main HLA alleles (*DQA1*05*, *DQB1*03:02*, *DQB1*02*, *DQA1*02*, *DQA1*03* and other alleles *DQB1* different from *DQB1*02*) associated with celiac disease (CD) in blood samples. This test is intended to be used as an aid in the detection of a genetic predisposition for celiac disease based on genotype identified. The absence of Celiac Disease-associated *HLA-DQ2.5*, *HLA-DQ2.2* and *HLA-DQ8* risk factors is an important exclusion criterion for CD. For the *HLA-DQB1*02* allele, it can determine the status of a sample as homozygous or heterozygous. DNA is extracted from blood specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for HLA alleles.

2. Summary and Explanation

Celiac disease is an inflammatory autoimmune condition that affects the small intestine primarily. At the molecular level, celiac disease is induced by an immune reaction to gluten, a protein found in wheat, barley and rye. Celiac disease is driven by aberrant immune responses to the presence of gluten in the diet in genetically predisposed individuals. The gold standard of celiac disease diagnosis requires analysis of the mucosa obtained from an intestinal biopsy and serological detection of anti-tissue transglutaminase, deaminated anti-gliadin peptides and / or anti-endomysial antibodies. Celiac disease has a strong hereditary component.

The celiac disease risk alleles *HLA-DQA1** and *DQB1** encode the alpha and beta chains, respectively, of the heterodimer of the HLA-DQ molecule, present on the surface of antigen presenting cells. Since each person presents two alleles from each loci *HLA-DQA1* and *HLA-DQB1*, it is possible to present up to 4 HLA-DQ heterodimers or molecules. A heterodimer encoded by the *HLA-DQA1* and *HLA-DQB1* alleles of the same chromosome is said to be cis-encoded, whereas if it is made up of alleles from different chromosomes, it is said to be trans-configured. Most celiac patients have at least one copy of *HLA-DQ2* or *HLA-DQ8*, two gene variants of the human leukocyte antigen (HLA class II or major histocompatibility complex, MHC). However, these variants are common in many populations, while only 1% develop celiac disease. Therefore, having one of the variants is necessary, but not sufficient for predisposition to the disease.

Each type of HLA is capable of presenting different peptide antigens, based on the amino acid sequence present in the peptide binding region of the HLA-DQ heterodimer. In the context of celiac disease, antigen presenting cells present gliadin-derived peptides, previously deaminated by tissue transglutaminase, as antigens, activating gluten-specific CD4 T cells. One of the hallmarks of celiac disease is a flow of activated lymphocytes into the small intestine, concomitant with villous atrophy. Immune activation is what ultimately induces the release of specific anti-gliadin and anti-tissue transglutaminase antibodies, these being detected in the diagnosis by serology.

The majority (95%) of celiac patients express the *HLA-DQ2* heterodimer (*HLA-DQA1*05* and *HLA-DQB1*02*) in cis or trans. *HLA-DQ8* (*HLA-DQA1*03* and *HLA-DQB1*03:02*) is present in most of the remaining patients (5-10%). *HLA-DQ2* and *HLA-DQ8* are estimated to contribute around 35% to the risk of developing celiac disease. The strongest association and most severe form occur with *HLA-DQA1*05*, *HLA-DQB1*02* (called *HLA-DQ2.5*), which, when inherited in homozygotes, is associated with a five-fold increase in risk of celiac disease development compared to individuals heterozygous for *HLA-DQ2.5*. Considering the high prevalence of these types of HLA in western

populations (40%), HLA typing has a very poor positive predictive value for the diagnosis of celiac disease. However, the absence of *HLA-DQ2* and *HLA-DQ8* makes the diagnosis of celiac disease highly unlikely (sensitivity > 96%). HLA typing has been included as a tool to exclude celiac disease in the European Society for Gastroenterology, Hepatology and Pediatric Nutrition (ESPGHAN) celiac diagnosis guide.

Therefore, a positive result in risk alleles in the HLA typing test is not sufficient for the diagnosis of celiac disease and requires other clinical correlations, such as positive serology with or without abnormal-looking duodenal biopsy, as indicated by the European and American clinical reference guides. Furthermore, it is particularly important that serology and duodenal biopsy studies are performed when the patient continues to include gluten in their diet, triggering the associated immune response.

The current diagnostic standard is based on the "rule of four out of five", which indicates that four of the five criteria are sufficient to establish the diagnosis of celiac disease: (1) typical signs and symptoms (diarrhea and malabsorption); (2) positive antibodies; (3) *HLA-DQ2* and / or *HLA-DQ8* positive; (4) intestinal damage (e.g., villous atrophy and minor injuries); and (5) clinical response to the gluten-free diet.

CerTest Biotec has developed a real-time PCR for the typing of the main celiac disease-related HLA haplotypes in patients with suspected gene susceptibility. This test has the advantage of requiring less time than the routine method, since being based on real-time PCR, only the time of the PCR reaction is required, without adding the hybridization time required with the routine method; in addition to reducing handling in the laboratory.

NOTE: For a better comprehension of nomenclature, we will refer to HLA alleles as follows:*

- *HLA-DQA1*05 alleles: DQA1*05*
- *HLA-DQB1*03:02 allele: DQB1*03:02.*
- *HLA-DQB1*02 alleles: DQB1*02.*
- *HLA-DQA1*02 alleles: DQA1*02.*
- *HLA-DQA1*03 alleles: DQA1*03.*
- *other alleles HLA-DQB1 different from HLA-DQB1*02: no" DQB1*02.*

*Antony Nolan Research Institute. (n.d.). HLA Nomenclature @ hla.alleles.org. Retrieved April 16, 2021, from <http://hla.alleles.org/nomenclature/naming.html>

3. Principle of the procedure

VIASURE HLA Celiac Real Time PCR Detection Kit is designed for the identification and differentiation of HLA alleles (*DQA1*05*, *DQB1*03:02*, *DQB1*02*, *DQA1*02*, *DQA1*03* and other alleles *DQB1* different from *DQB1*02*) from blood samples. The isolated DNA target is amplified on specific regions for the identification of the *HLA-DQA1* loci (*DQA1*05*, *DQA1*02*, and *DQA1*03* alleles), *HLA-DQB1* loci (*DQB1*03:02* and *DQB1*02* alleles and other alleles *DQB1* different from *DQB1*02*), using specific primers and a fluorescent-labeled probe.

VIASURE HLA Celiac Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *HLA Celiac* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an **endogenous internal control** to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an **Endogenous Internal Control (IC)** (hemoglobin subunit beta gene, *HBB* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels. Each kit includes two kinds of multiplex reaction mix and each one corresponds to one different assay. The first multiple reaction mix detects *DQA1*05*, *DQB1*03:02* and/or *DQB1*02* alleles (*HLA Celiac 1*) and contains the endogenous internal control. The second multiple reaction mix detects *DQA1*02*, *DQA1*03* alleles and/or "no" *DQB1*02* alleles (other alleles *DQB1* different from *DQB1*02*) (*HLA Celiac 2*).

4. Reagents provided

VIASURE *HLA Celiac* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, and Annex 2 for tube format with internal control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *HLA Celiac* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1).

VIASURE *HLA Celiac* Real Time PCR Detection Kit has been validated on the following equipments: Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 500, ROX channel – 1000 and Cy5 channel – 1000.

- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 200, ROX channel – 500 and Cy5 channel – 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *HLA Celiac* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR, human genetics and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for reference VS-CEL136). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-CEL136 (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.

- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *HLA Celiac* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Make sure to use a well for determining *DQA1*05*, *DQB1*03:02* and/or *DQB1*02* alleles and another well for determining *DQA1*02*, *DQA1*03* alleles and/or "no" *DQB1*02* alleles (other alleles *DQB1* different from *DQB1*02*). Be careful not to mix them throughout the process.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format with internal control products Test Procedure and Annex 2 for Tube format with internal control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *HLA Celiac* Real Time PCR Detection Kit has been tested in blood (with EDTA or citrate). Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at 4°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), it is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 4°C for up to 24

hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to Yong W. (eds) Biobanking. Methods in Molecular Biology, vol 1897 (2019) Humana Press, New York, NY.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from blood samples, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- Invisorb® Spin Universal Kit (Invitex).
- QIAAsymphony SP (Qiagen).

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for HLA alleles in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	HLA Celiac 1				HLA Celiac 2			Interpretation of Controls
	DQA1*05 (FAM) ¹	DQB1*03:02 (ROX) ¹	DQB1*02 (Cy5) ¹	Endogenous Internal Control (HEX) ²	DQA1*02 (FAM) ¹	DQA1*03 (HEX) ¹	"no" DQB1*02 (ROX) ¹	
Positive Control (PC)	≤40	≤40	≤40	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≥40 or no signal	≥40 or no signal	≥40 or no signal	≥40 or no signal	≥40 or no signal	≥40 or no signal	Valid

Table 1. Expected Performance of Controls

¹ In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid', and retesting is required.

2. The Endogenous Internal Control (IC) should show an amplification signal ($C_t \leq 40$) in PC control wells.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

HLA Celiac 1			HLA Celiac 2			Result	Risk ⁽²⁾
DQA1*05 (FAM)	DQB1*03:02 (ROX)	DQB1*02 (Cy5)	DQA1*02 (FAM)	DQA1*03 (HEX)	"no" DQB1*02 (ROX)		
+	-	+	-	+/-	+/- ⁽¹⁾	DQ2.5	High
+	+	+	-	+	+	DQ2.5/DQ8	High
+	-	+	+	-	- ⁽¹⁾	DQ2.5/DQ2.2	High
-	+	-	+/-	+	+	DQ8	Moderate
-	+	+	+	+	+	DQ8/DQ2.2	Moderate
+	+	-	-	+	+	DQ8/DQA1*05	Moderate
-	-	+	+	+/-	+/- ⁽¹⁾	DQ2.2	Low
+	-	-	+/-	+/-	+	DQA1*05	Very low
-	-	-	+/-	+/-	+	Negative	Risk free

Table 2. Interpretation of individual patient sample results. "no" DQB1*02= other alleles DQB1 different from DQB1*02.

+: Amplification curve

-: No amplification curve

1 A negative result for these targets implies a major risk (homozygous for DQB1*02 allele), following the indications described in Ruiz-Ortiz E. et al., 2014.

2 Risk genotypes at the HLA-DQA1 and HLA-DQB1 loci following the indications described in Erlichster M. et al., 2018.

A sample is considered positive if the C_t value obtained is less than 40 and the endogenous internal control shows or not an amplification signal. The endogenous internal control is a human housekeeping gene that should be present in all human nucleated cells.

A sample is considered negative if the sample shows no amplification signal in the detection system but the endogenous internal control is positive (C_t less than 35). A failure in the extraction procedure and/or an inhibition of the PCR reaction can be excluded by the amplification of endogenous internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System. (HLA Celiac 1).



Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System. (HLA Celiac 2).

The result is considered invalid if there is any signal of amplification in negative control or absence of target signals in the positive well. We recommend to repeat the assay again.

In case of absence of endogenous internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters, the sigmoid shape of the curve and fluorescence intensity. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other canonical diagnostic tests. It is recommended to use VIASURE HLA Celiac Real Time PCR Detection Kit as an aid to the algorithm used in the diagnosis of celiac diseases (serological screening methods, small intestinal biopsies from the duodenum, gluten-free diet etc.).
- Results obtained from genetic diagnosis have only negative predictive value, allowing the exclusion of celiac disease in a high percentage, but positive results do not imply disease condition or positive predictive value.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from blood samples.
- The quality of the test depends on the quality of the sample; properly extracted nucleic acid from clinical samples must be extracted.
- This test is a qualitative test and does not provide quantitative values.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by HLA Celiac Positive Control, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).

- Degradation of the DNA during sample shipping/storage and/or processing.
- A target in the specimen below the limit of detection for the assay.
- The presence of qPCR inhibitors or other types of interfering substances.
- Failure to follow instructions for use and the assay procedure.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *HLA Celiac* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Endogenous Internal Control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *HLA Celiac* Real Time PCR Detection Kit was tested using clinical samples (blood samples) from patients diagnosed with celiac disease and first-degree relatives. In order to determine the clinical diagnostic accuracy, a multicenter evaluation has been conducted through collaboration with national entities. A summary of the sites, sample type and workflow is included in the following table:

	Site	Sample type	Workflow	Target
1	Clinical Biochemistry Service of 'Hospital Universitario Miguel Servet' (HUMS) (Zaragoza, Spain)	Blood samples	QIASymphony DSP DNA Kit using QIASymphony SP instrument + Applied Biosystems GeneAmp® PCR System 9700 and Bio-Rad CFX96 Touch Real-Time PCR Detection System	HLA alleles associated to celiac disease
2	Certest Biotec and Department of Biochemistry of Nordsjællands hospital, (Denmark)	Blood samples	QIASymphony DSP DNA Kit using QIASymphony SP instrument + Qiagen Rotor-Gene® Q + + Bio-Rad CFX96 Touch Real-Time PCR Detection System	HLA alleles associated to celiac disease

Table 3. Summary of the sites, sample type and workflow carried out during the clinical evaluation.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *HLA Celiac* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	CeliacStrip kit (Operon)	HLA alleles associated to celiac disease	73	6	1	0	1 (0.95 – 1)	0.85 (0.42 – 0.99)	0.98 (0.91 – 0.99)	1 (0.51 – 1)
2	EUROArray HLA-DQ2/DQ8 Direct assay (EUROIMMUN)	DQ2.5/DQ8 haplotypes	11	278	0	0	1 (0.71 – 1)	1 (0.98 – 1)	1 (0.71 – 1)	1 (0.98 – 1)
		DQ2.5/DQ2.2 haplotypes	16	273	0	0	1 (0.79 – 1)	1 (0.79 – 1)	1 (0.79 – 1)	1 (0.98 – 1)
		DQ2.5 (heterocygote) haplotype	52	237	0	0	1 (0.93 – 1)	1 (0.98 – 1)	1 (0.93 – 1)	1 (0.98 – 1)
		DQ2.5 (homocygote) haplotype	24	265	0	0	1 (0.85 – 1)	1 (0.98 – 1)	1 (0.85 – 1)	1 (0.98 – 1)
		DQ8/DQ2.2 haplotypes	3	286	0	0	1 (0.29 – 1)	1 (0.98 – 1)	1 (0.29 – 1)	1 (0.98 – 1)
		DQ8/DQA1*05 haplotypes	8	281	0	0	1 (0.63 – 1)	1 (0.98 – 1)	1 (0.63 – 1)	1 (0.98 – 1)
		DQ8 haplotype	40	248	1	0	1 (0.91 – 1)	0.99 (0.97 – 1)	1 (0.91 – 1)	0.99 (0.97 – 1)
		DQ2.2 (heterocygote) haplotype	17	272	0	0	1 (0.80 – 1)	1 (0.98 – 1)	1 (0.80 – 1)	1 (0.98 – 1)
		DQ2.2 (homocygote) haplotype	7	282	0	0	1 (0.59 – 1)	1 (0.98 – 1)	1 (0.59 – 1)	1 (0.98 – 1)
		DQA1*05 haplotype	31	258	0	0	1 (0.88 – 1)	1 (0.98 – 1)	1 (0.88 – 1)	1 (0.98 – 1)
		No risk	79	210	0	0	1 (0.95 – 1)	1 (0.98 – 1)	1 (0.95 – 1)	1 (0.98 – 1)

Table 4. True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE HLA Celiac Real Time PCR Detection Kit.

In addition to this, 4 saliva samples from individuals diagnosed with celiac disease were also analysed with VIASURE HLA Celiac Real Time PCR Detection Kit, obtaining concordant results comparing to reference methods used.

Results show high agreement to detect HLA alleles associated to celiac disease using VIASURE HLA Celiac Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE HLA Celiac Real Time PCR Detection Kit has a detection limit of 50 DNA copies per reaction for HLA alleles DQA1*05, DQB1*03:02, and DQB1*02; and 10 DNA copies per reaction for HLA alleles DQA1*02, DQA1*03 and “no” DQB1*02 (other alleles DQB1 different from DQB1*02)) with a positive rate of $\geq 95\%$.

Figure 3. Dilution series of DQA1*05 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *HLA Celiac 1*, channel FAM).


Figure 4. Dilution series of DQB1*03:02 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *HLA Celiac 1*, channel ROX).


Figure 5. Dilution series of DQB1*02 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *HLA Celiac 1*, channel Cy5).


Figure 6. Dilution series of DQA1*02 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *HLA Celiac 2*, channel FAM)


Figure 7. Dilution series of DQA1*03 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction *HLA Celiac 2*, channel HEX).



Figure 8. Dilution series of "no" DQB1*02 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction *HLA Celiac 2*, channel ROX).

12.3. Analytical specificity

The specificity of the *HLA Celiac* assay was confirmed by testing a panel consisting different microorganisms representing the most common vector-borne pathogens, and other 18S-possessing homologous eukaryotes. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Anaplasma marginale</i>	-	<i>Leptospira</i>	-	Dengue 1 virus strain Hawaii	-
<i>Bartonella henselae</i> strain Houston-1	-	<i>Plasmodium falciparum</i> 3D7	-	Dengue 2 virus strain New Guinea C	-
<i>Borrelia hermsii</i>	-	<i>Rickettsia conorii</i> . strain Moroccan	-	Dengue 3 virus strain H87	-
<i>Borrelia lusitanae</i>	-	<i>Treponema phagedenis</i>	-	Dengue 4 virus strain H241	-
<i>Borrelia valaisiana</i>	-	<i>Trypanosoma cruzi</i>	-	Usutu Virus	-
<i>Borrelia azfeli</i> strain P-Ko/1984	-	Japanese encephalitis	-	West Nile virus strain NY99	-
<i>Borrelia bavariensis</i>	-	Japanese Encephalitis virus Nakayama	-	West Nile virus Heja	-
<i>Borrelia bissetti</i>	-	Rift Valley Fever Virus AR21229	-	West Nile virus Ug37	-
<i>Borrelia burgdorferi</i> sensu stricto strain IRS	-	Rift Valley Fever Virus MP12	-	Yellow Fever virus strain 17D	-
<i>Borrelia burgdorferi</i> sensu stricto strain B31	-	St Louis Encephalitis virus	-	Yellow Fever Virus French Neurotropic	-
<i>Borrelia garinii</i>	-	Chikungunya virus WHO IS (R91064)	-	Zika Virus strain 11474/16 (French Polynesia)	-
<i>Borrelia japonica</i>	-	Chikungunya virus strain S27 Petersfield	-	Zika Virus strain 11468/16(French Polynesia)	-
<i>Borrelia miyamotoi</i>	-	Chikungunya virus Martinique isolate	-	Zika Virus (African)	-
<i>Borrelia spielmanii</i>	-	Chikungunya virus strain F24	-	Zika virus strain PF13/251013-18 (Asian)	-
<i>Coxiella burnetii</i> strain Nine Mile Q	-	Tick Borne Encephalitis Virus (TBEV) strain Neudorfl	-	Zika virus strain FB-GWUH-2016	-
<i>Theileria annulata</i>	-				

Table 5. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *HLA Celiac* Real Time PCR Detection Kit was evaluated against *HLA-DQ2.5* heterodimer (*HLA-DQA1*05:01* and *HLA-DQB1*02:01* alleles), *HLA-DQ2.2* heterodimer (*HLA-DQA1*05:01* and *HLA-DQB1*02:01* alleles), *HLA-DQ8* heterodimer (*HLA-DQA1*03:01* and *HLA-DQB1*03:02* alleles), showing positive results.

ANNEX 1

OPEN AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE HLA Celiac Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS- CEL106L
VIASURE HLA Celiac Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS- CEL106H
VIASURE HLA Celiac Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS- CEL112L
VIASURE HLA Celiac Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS- CEL112H
VIASURE HLA Celiac Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS- CEL136
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS- CEL101L
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS- CEL101H
VIASURE HLA Celiac Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS- CEL101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE HLA Celiac Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an endogenous internal control to monitor PCR inhibition.

	Target	Channel	Gene (Alleles)
HLA Celiac 1	HLA alleles associated to celiac disease	FAM	DQA1*05
	HLA alleles associated to celiac disease	ROX	DQB1*03:02
	HLA alleles associated to celiac disease	Cy5	DQB1*02
	Endogenous Internal control (IC)	HEX, VIC or JOE#	HBB gene (β-globin)
HLA Celiac 2	HLA alleles associated to celiac disease	FAM	DQA1*02
	HLA alleles associated to celiac disease	HEX, VIC or JOE#	DQA1*03
	HLA alleles associated to celiac disease	ROX	no DQB1*02

Table A1 2. Target, channel and genes.

#Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.certest.es.

A1.2 Reagents provided

VIASURE HLA Celiac Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3 and A1.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
<i>HLA Celiac</i> 1 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	1/3/6 x 8-well strip
<i>HLA Celiac</i> 2 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	White	1/3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>HLA Celiac</i> Positive Control	synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	2/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *HLA Celiac* Real Time PCR Detection Kit with Ref. VS-CEL101L, VS- CEL101H, VS- CEL106L, VS- CEL106H, VS- CEL112L and VS- CEL112H.

Reagent/Material	Description	Colour	Amount
<i>HLA Celiac</i> 1 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	Transparent	2/9 x 4-well strip
<i>HLA Celiac</i> 2 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	Transparent	2/9 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>HLA Celiac</i> Positive Control	synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	4/18 x 4-cap strip

Table A1 4. Reagents and materials provided in VIASURE *HLA Celiac* Real Time PCR Detection Kit with Ref. VS-CEL101 and VS-CEL136. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

HLA Celiac Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *HLA Celiac* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *HLA Celiac* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A1 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*DQA1*05* and *DQA1*02*), HEX, JOE or VIC (Endogenous Internal Control (IC) and *DQA1*03*), ROX (*DQB1*03:02* and "no" *DQB1*02*), and Cy5 channels (*DQB1*02*). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE HLA Celiac Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-CEL148T

Table A2. 1.References.

A2.1 Principle of the procedure

VIASURE HLA Celiac Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an endogenous internal control to discard the inhibition of the polymerase activity.

	Target	Channel	Gene (Alleles)
HLA Celiac 1	HLA alleles associated to celiac disease	FAM	DQA1*05
	HLA alleles associated to celiac disease	ROX	DQB1*03:02
	HLA alleles associated to celiac disease	Cy5	DQB1*02
	Endogenous Internal control (IC)	HEX, VIC or JOE [#]	HBB gene (β-globin)
HLA Celiac 2	HLA alleles associated to celiac disease	FAM	DQA1*02
	HLA alleles associated to celiac disease	HEX, VIC or JOE [#]	DQA1*03
	HLA alleles associated to celiac disease	ROX	no DQB1*02

Table A2. 2.Target, channel and genes.

[#]Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.cer-test.es.

A2.2 Reagents provided

VIASURE HLA Celiac Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
HLA Celiac 1 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	2 vials
HLA Celiac 2 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	White	2 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
HLA Celiac Positive Control	synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE HLA Celiac Real Time PCR Detection Kit with Ref. VS-CEL148T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

HLA Celiac Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *HLA Celiac* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *HLA Celiac* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *HLA Celiac* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *HLA Celiac* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*DQA1*05* and *DQA1*02*), HEX, JOE or VIC (Endogenous Internal Control (IC) and *DQA1*03*), ROX (*DQB1*03:02* and "no" *DQB1*02*), and Cy5 channels (*DQB1*02*). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.