VIASURE

Real Time PCR Detection Kit

Carbapenemase-producing *Enterobacteriaceae*

CE IVD

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Content

1.	Intended use	7
2.	Summary and Explanation	7
3.	Principle of the procedure	8
4.	Reagents provided	9
5.	Reagents and equipment to be supplied by the user	9
6.	Transport and storage conditions	10
7.	Precautions for users	10
8.	Test procedure.....	11
8.1.	Specimen collection, transport and storage	12
8.2.	DNA extraction	12
9.	Result interpretation	13
9.1.	References with internal control (references in Annex 1 and 2).....	13
9.2.	References with extraction control (references in Annex 3 and 4)	15
10.	Limitations of the test	17
11.	Quality control.....	18
12.	Performance characteristics.....	18
12.1.	Clinical sensitivity and specificity	18
12.2.	Analytical sensitivity.....	20
12.3.	Analytical specificity	22
12.4.	Analytical reactivity	22
	ANNEX 1	23
A1.1	Principle of the procedure	23
A1.2	Reagents provided	23
A1.3	Test procedure	24
A1.3.1	Lyophilized positive control.....	24
A1.3.2	PCR protocol	25
	ANNEX 2	26
A2.1	Principle of the procedure	26
A2.2	Reagents provided	26
A2.3	Test procedure	27
A2.3.1	Lyophilized positive control.....	27
A2.3.2	Lyophilized reaction mix tube.....	27
A2.3.3	PCR protocol	27

ANNEX 3	29
A3.1 Principle of the procedure	29
A3.2 Reagents provided	29
A3.3 Test procedure	31
A3.3.1 Lyophilized extraction control	31
A3.3.2 Lyophilized positive control.....	31
A3.3.3 PCR protocol	31
ANNEX 4	33
A4.1 Principle of the procedure	33
A4.2 Reagents provided	33
A4.3 Test procedure	34
A4.3.1 Lyophilized extraction control	34
A4.3.2 Lyophilized positive control.....	34
A4.3.3 Lyophilized reaction mix tube.....	34
A4.3.4 PCR protocol	34

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ENGLISH

1. Intended use

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit is designed for the simultaneous qualitative detection and differentiation of the main carbapenemase-encoding genes (*NDM*, *VIM*, *OXA*, *KPC* and/or *IMP*) from bacterial isolates of clinical samples and directly from rectal swabs from individuals suspected of infection by carbapenem-resistant pathogens by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of infection caused by carbapenem-resistant *Enterobacteriaceae* in combination with clinical and epidemiological risk factors. DNA is extracted from bacterial isolates and rectal swabs samples, amplified using real time PCR, and detected using fluorescent reporter dye probes specific for carbapenemase-encoding genes.

2. Summary and Explanation

Antimicrobial resistance is a major public health problem nowadays, but the resistance to carbapenem antimicrobials is becoming particularly concerning. Carbapenem antibiotics are commonly considered to be the most potent group of antimicrobial agents with proven efficacy in the treatment of patients with severe bacterial infections, including those caused by antimicrobial-resistant strains. However, the recent increase in the global detection of carbapenem-resistant *Enterobacteriaceae* (CRE) represents one of the most pressing threats to public health. This is mainly due to their implication in numerous outbreaks and Health Care Associated Infection (HAI), which imply longer lengths of stay, increased healthcare costs, and higher mortality than carbapenem-susceptible infections. The current definition of CRE recommended by CDC is an *Enterobacteriaceae* isolate that is resistant to ertapenem, imipenem, meropenem, or doripenem, according to the current Clinical and Laboratory Standards Institute (CLSI) breakpoints, or documentation that the isolate produces carbapenemase (carbapenemase-producing *Enterobacteriaceae* – CPE –). Both the terms CRE and CPE are used sometimes interchangeably, but the former refers to the resistance phenotype, whereas the latter implies the mechanism underlying the phenotype. CPE are mostly CRE that are resistant to carbapenems due to carbapenemase production, which is the most prominent mechanism underlying carbapenem resistance in Gram-negative pathogens.

Enterobacteriaceae is an ubiquitous family of Gram-negative bacteria, from which certain species are part of the normal human flora, but frequently associated to diarrheal disease and extraintestinal infections. They are responsible of a variety of community and healthcare-acquired infections, and a major driving force of resistance is the production of β -lactamases (encoded by *bla* gene). The most frequently occurring species of CRE are *Klebsiella pneumoniae* (*K. pneumoniae*) and *Escherichia coli* (*E. coli*). While some strains are innately resistant to carbapenems, others contain mobile genetic elements (for example, plasmids, transposons) which can be readily horizontal transmitted between among various bacterial species and genera. The surge in CRE is mostly driven by the emergence and spread of carbapenemases, a specific group of β -lactamases that are capable of hydrolyzing carbapenems, along with other β -lactams. There are currently five major carbapenemases, for which are all on mobile genetic elements: *Klebsiella pneumoniae* carbapenemase (KPC), New Delhi metallo- β -lactamase (NDM), the carbapenem hydrolyzing oxacillinases (OXA-48 and OXA-48-like), Verona integron-encoded metallo- β -lactamase (VIM) and metallo- β -lactamase active on imipenem (IMP).

The prevalence of CRE and types of carbapenemases is highly dependent on geography, being *K. pneumoniae* the predominant specie of *Enterobacteriaceae* across regions. Increasing trends of invasive isolates of carbapenem-resistant *K. pneumoniae* were observed from 2012-2015. According to a point prevalence survey of HAI and antimicrobial use in Europe, 18 of 28 countries reported CRE and three countries reported HAI with more than 20% of resistant isolates, with the highest percentage (39.9%) in Greece. In the United States of America (USA), 49.5% of *A. baumannii*, 19.2% of *P. aeruginosa*, 7.9% of *K. pneumoniae* and 0.6% of *E. coli* invasive isolates submitted to the National Healthcare Safety Network were resistant to carbapenems in 2014. KPC is the most common observed carbapenemase, being prevalent in countries like Greece, Italy, Brazil, China, and several others including USA and Colombia. However, while in Indian subcontinent (India, Pakistan, Bangladesh) NDM is the most common carbapenemase, OXA-48 is concentrated in European and Mediterranean countries including North Africa.

The most important source for CRE dissemination is intestinal carriers (infected or colonized people) and their prompt identification is crucial to reduce cross-transmission of CRE through direct person-to-person contact in community and healthcare settings. They primarily affect patients in acute and long-term healthcare settings, who are being treated for another condition, patients who have compromised immune systems or have invasive devices like tubes going into their body. Once these bacteria have spread outside the gut, they could cause serious infections, such as pneumonia, bacteremia, urinary tract infections, wound infections, surgical site infections and meningitis.

The phenotypic methods for identification of CPE, like modified Hodge test, have been considered 'reference' methods, but they may be complex, occasionally inconclusive and time-consuming owing to the wide range of carbapenem minimum inhibitory concentration. Therefore, molecular methods for rapid detection of carbapenemase-producing *Enterobacteriaceae* (CPE) represent an advantageous tool for faster identification of CRE and allow precise differentiation of the carbapenemase type.

NOTE: For a better comprehension of nomenclature, we will refer to carbapenemase-encoding genes as follows:

- *bla_{NDM}*: NDM.
- *bla_{VIM}*: VIM.
- *bla_{OXA-48 like}*: OXA.
- *bla_{KPC}*: KPC.
- *bla_{IMP}*: IMP (IMP-1 cluster + IMP-8 cluster)*.

*IMP-1 cluster includes: IMP-1, IMP-4, IMP-3, IMP-6, IMP-10, IMP-25, IMP-30, IMP-34, IMP-40, IMP-42, IMP-52, IMP-55.

IMP-8 cluster includes: IMP-2, IMP-8, IMP-16, IMP-19, IMP-20, IMP-22, IMP-24.

3. Principle of the procedure

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit is designed for the qualitative detection and differentiation of carbapenemase-encoding genes (NDM, VIM, OXA, KPC and/or IMP) from bacterial isolates of clinical samples and rectal swabs. After DNA isolation, the identification of carbapenemases is performed by the amplification of a conserved region of the NDM gene, VIM gene, OXA-48 and OXA-48-like genes, KPC gene and IMP gene, using specific primers and a fluorescent-labelled probes.

VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal, which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. Each kit includes two kinds of multiplex reaction mix, and each one corresponds to one different assay. The first reaction mix detects *NDM* and/or *VIM* carbapenemase-encoding genes (*Carbapenemase-producing Enterobacteriaceae* 1) and contains the internal control (IC). The second reaction mix detects *OXA*, *KPC* and/or *IMP* carbapenemase-encoding genes (*Carbapenemase-producing Enterobacteriaceae* 2).

4. Reagents provided

VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products, and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).

VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit has been tested on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied BioSystem 7500 RT PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cyclers, DNA-Technology DTlite Real-Time PCR System.

Rotor-Gene® Q (Qiagen), Roche Molecular Diagnostics Cobas z480 Analyzer, Roche Molecular Diagnostics LightCycler 480 II and Mic Real Time PCR Cyclers.

When using the Applied Biosystems 7500 Fast with strips it is recommended to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set of exposure values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *Carbapenemase-producing Enterobacteriaceae* Reaction-Mix tube has been reconstituted, it may be kept at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended to store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-CPE136 and VS-CPE136E). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.

- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-CPE136 and VS-CPE136E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Make sure to use a well for determining *NDM* and *VIM* genes and another well for determining *OXA*, *KPC* and *IMP*. Be careful not to mix them throughout the process.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please see Annex 1 for Open and rotor-gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for Open and rotor-gene format with

extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit has been tested on bacterial isolates from different clinical samples and rectal swabs. Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, all specimens should be collected and labelled appropriately in clean containers. After collection, specimens should be placed in a biohazard bag and should be transported and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at Room Temperature (RT) for up to 2 hours, or at 2 to 8°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), we recommend shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 24 hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC guidelines (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>), the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94), and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from bacterial isolates from clinical samples and from rectal swabs, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- GXT NA Extraction kit in GenoXtract® system (Hain), using QIASymphony RGQ® (QIAGEN).
- MagDEA Dx SV Kit, using the magLEAD® 6gC or 12gC instrument (Precision System Science Co.).
- Invisorb® Spin Universal Kit (Invitex).

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for Carbapenemase-producing *Enterobacteriaceae* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	Carbapenemase-producing <i>Enterobacteriaceae</i> 1			Carbapenemase-producing <i>Enterobacteriaceae</i> 2			Interpretation of Controls
	NDM (FAM) ¹	VIM (ROX) ¹	Internal Control (HEX) ²	OXA (FAM) ¹	KPC (ROX) ¹	IMP (Cy5) ¹	
Positive Control (PC)	≤40	≤40	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≥40 or no signal	≤40	≥40 or no signal	≥40 or no signal	≥40 or no signal	Valid

Table 1. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid', and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40) in PC and NC control wells.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

Carbapenemase-producing <i>Enterobacteriaceae</i> 1			Carbapenemase-producing <i>Enterobacteriaceae</i> 2			Result
NDM (FAM)	VIM (ROX)	Internal Control (HEX)	OXA (FAM)	KPC (ROX)	IMP (Cy5)	
≤40	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≥40 or no signal	NDM-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≤40	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≥40 or no signal	VIM-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≤40	≥40 or no signal	≥40 or no signal	OXA-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≤40	≥40 or no signal	KPC-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≤40	IMP-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤35 ²	≥40 or no signal	≥40 or no signal	≥40 or no signal	Carbapenemase-producing <i>Enterobacteriaceae</i> not detected²
≥40 or no signal	≥40 or no signal	≥ 35 or no signal ²	≥40 or no signal	≥40 or no signal	≥40 or no signal	Invalid- Repeat Testing²

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

1 The Internal Control (IC) shows or not an amplification signal ($C_t \leq 40$ or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of Carbapenemase-producing *Enterobacteriaceae* target genes negative, IC must show an amplification signal with C_t less than 35. If there is an absence of signal or C_t value ≥ 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

Note: One individual patient sample simultaneously can contain several carbapenemase-encoding genes. Table 2 shows only the most representative results that can be expected with the VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 1).



Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System. (Carbapenemase-producing *Enterobacteriaceae* 2).

9.2. References with extraction control (references in Annex 3 and 4)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for Carbapenemase-producing *Enterobacteriaceae* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	Carbapenemase-producing <i>Enterobacteriaceae</i> 1			Carbapenemase-producing <i>Enterobacteriaceae</i> 2			Interpretation of Controls
	NDM (FAM) ¹	VIM (ROX) ¹	Extraction Control (HEX) ²	OXA (FAM) ¹	KPC (ROX) ¹	IMP (Cy5) ¹	
Positive Control (PC)	≤40	≤40	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≥40 or no signal	≤40	≥40 or no signal	≥40 or no signal	≥40 or no signal	Valid

Table 3. Expected Performance of Controls

¹ In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid', and retesting is required.

² The Extraction Control (EC) should show an amplification signal (Ct ≤40) in PC and NC control wells.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

Carbapenemase-producing <i>Enterobacteriaceae</i> 1			Carbapenemase-producing <i>Enterobacteriaceae</i> 2			Result
NDM (FAM)	VIM (ROX)	Extraction Control (HEX)	OXA (FAM)	KPC (ROX)	IMP (Cy5)	
≤40	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≥40 or no signal	NDM-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≤40	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≥40 or no signal	VIM-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≤40	≥40 or no signal	≥40 or no signal	OXA-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≤40	≥40 or no signal	KPC-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤40 or no signal ¹	≥40 or no signal	≥40 or no signal	≤40	IMP-expressing <i>Enterobacteriaceae</i>
≥40 or no signal	≥40 or no signal	≤35 ²	≥40 or no signal	≥40 or no signal	≥40 or no signal	Carbapenemase-producing <i>Enterobacteriaceae</i> not detected²
≥40 or no signal	≥40 or no signal	≥ 35 or no signal ²	≥40 or no signal	≥40 or no signal	≥40 or no signal	Invalid- Repeat Testing²

Table 2. Interpretation of individual patient sample results. Ct values, no signal = no amplification curve.


1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of Carbapenemase-producing *Enterobacteriaceae* target genes negative, EC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

Note: One individual patient sample simultaneously can contain several carbapenemase-encoding genes. Table 4 shows only the most representative results that can be expected with the VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 3. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 1).

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from bacterial isolates from clinical samples and rectal swabs.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from samples must be done.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Carbapenemase-producing Enterobacteriaceae* Positive Control, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Carbapenemase-producing Enterobacteriaceae* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the DNA during sample shipping/storage and/or processing.
 - A bacterial load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent drug-resistant infections or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable pathogens and does not imply that these pathogens are infectious or are the causative agents for clinical symptoms. However, a positive

result is indicative of the presence of targets Carbapenemase-producing *Enterobacteriaceae* sequences (NDM, VIM, KPC, OXA and IMP genes).

- Negative results do not preclude Carbapenemase-producing *Enterobacteriaceae* infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak pathogen levels during infections caused by Carbapenemase-producing *Enterobacteriaceae* have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogens.
- If diagnostic tests for other drug-resistant related illnesses are negative and the patient's clinical presentation and epidemiological information suggest that Carbapenemase-producing *Enterobacteriaceae* infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.
- A bioinformatics tool (NCBI BLAST) analysis of the primer alignments showed that VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit detects the following family members for each family of the carbapenemase genes:

Family	Family members
OXA	48, 162, 163, 181, 199, 204, 232, 244, 245, 247, 252, 370, 405, 416, 438, 439, 484, 505, 514, 515, 517, 519, 538, 546, 547, 566, 567, 788, 793, 833, 894
KPC	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 51, 52, 54, 56
VIM	1, 2, 3, 4, 6, 8, 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 23, 24, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 39, 40, 41, 42, 43, 44, 45, 46, 48, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 62, 63, 64, 65, 66, 67, 68, 70
NDM	1, 2, 3, 4, 5, 6, 7, 9, 11, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29
IMP	1, 4, 5, 6, 8, 9, 10, 11, 12, 13, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 28, 29, 30, 33, 34, 37, 38, 39, 40, 41, 42, 44, 45, 49, 52, 53, 55, 56, 58, 59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85

Table 5. Family members of each carbapenemase gene family detected by VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit.

11. Quality control

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit was tested using bacterial isolates from clinical samples and directly from clinical samples (rectal swabs). In

order to determine the clinical diagnostic accuracy, a multicenter evaluation has been conducted through collaboration with national entities. A summary of the sites, sample type and workflow is included in the following table:

	Site	Sample type	Workflow	Target
1.A	'Departamento de Microbiología Clínica' and the 'Àrea de Biología Molecular' of 'Laboratori de Referència de Catalunya'	Bacterial isolates from clinical samples	GXT NA Extraction kit in GenoXtract® system (Hain), following manufacturer's instructions + Applied BioSystem 7500 RT PCR System	Carbapenemase-producing <i>Enterobacteriaceae</i>
1.B	'Departamento de Microbiología Clínica' and the 'Àrea de Biología Molecular' of 'Laboratori de Referència de Catalunya'	Rectal swabs	GXT NA Extraction kit in GenoXtract® system (Hain), following manufacturer's instructions + Applied BioSystem 7500 RT PCR System	Carbapenemase-producing <i>Enterobacteriaceae</i>

Table 6. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following tables:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1.A	Phenotypical characterization + LightMix® modular carbapenemase assay	<u>Global:</u> Carbapenemase-producing <i>Enterobacteriaceae</i>	80	17	0	0	1 (0.95-1)	1 (0.80-1)	1 (0.95-1)	1 (0.80-1)
		OXA gene	33	64	0	0	1 (0.89-1)	1 (0.94-1)	1 (0.89-1)	1 (0.94-1)
		VIM gene	23	74	0	0	1 (0.85-1)	1 (0.95-1)	1 (0.85-1)	1 (0.95-1)
		NDM gene	19	78	0	0	1 (0.82-1)	1 (0.95-1)	1 (0.82-1)	1 (0.95-1)
		KPC gene	2	95	0	0	1 (0.16-1)	1 (0.96-1)	1 (0.16-1)	1 (0.96-1)
		IMP gene	7	90	0	0	1 (0.59-1)	1 (0.96-1)	1 (0.59-1)	1 (0.96-1)
1.B	Phenotypical characterization + LightMix® modular carbapenemase assay	<u>Global:</u> Carbapenemase-producing <i>Enterobacteriaceae</i>	34	16	0	0	1 (0.91-1)	1 (0.79-1)	1 (0.91-1)	1 (0.79-1)
		OXA gene	13	37	0	0	1 (0.75-1)	1 (0.90-1)	1 (0.75-1)	1 (0.90-1)
		VIM gene	11	39	0	0	1 (0.71-1)	1 (0.91-1)	1 (0.71-1)	1 (0.91-1)
		NDM gene	10	40	0	0	1 (0.69-1)	1 (0.91-1)	1 (0.69-1)	1 (0.91-1)
		KPC gene	5	45	0	0	1 (0.47-1)	1 (0.92-1)	1 (0.47-1)	1 (0.92-1)
		IMP gene	0	50	0	0	n.a*	1 (0.91-1)	n.a*	1 (0.91-1)

Table 7. True positive (TP) and true negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV) and Predictive Negative Values (NPV) for VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit.

* Due to all analyzed samples were negative for IMP gene, the analytical sensitivity of the test could not be performed.

Results show high agreement to detect Carbapenemase-producing *Enterobacteriaceae* using VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit has a detection limit of 10 DNA copies per reaction for carbapenemase-encoding genes *NDM*, *KPC* and *IMP*, 50 DNA copies per reaction for *OXA* and 100 DNA copies per reaction for *VIM* (Figure 4, 5, 6, 7 and 8).

Figure 4. Dilution series of *NDM* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 1, channel FAM).


Figure 5. Dilution series of *VIM* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 1, channel ROX).


Figure 6. Dilution series of OXA (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96TM Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 2, channel FAM).


Figure 7. Dilution series of KPC (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96TM Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 2, channel ROX).



Figure 8. Dilution series of IMP cluster (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96TM Real-Time PCR Detection System (Carbapenemase-producing *Enterobacteriaceae* 2, channel Cy5).

12.3. Analytical specificity

The specificity of the VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit assay was confirmed by testing a panel of different microorganisms which represents the most common drug resistant pathogens. No cross-reactivity of the VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit with genomic DNA of the selected pathogens was observed (Table 8), except for the targeted pathogens of each assay.

Cross-reactivity testing			
OXA-48 producing <i>Serratia marcescens</i> isolate	-/+	VanA-type <i>Enterococcus avium</i>	-
SHV-1 (non-ESBL), KPC-3, and OXA-48 producing <i>Klebsiella pneumonia</i> isolate	-/+	VanA- type <i>Enterococcus faecium</i> LMG16165 strain	-
TEM-1 (non-ESBL), SHV-1 (non-ESBL), CTX-M-2 (ESBL), and KPC-2 producing <i>Klebsiella pneumonia</i> isolate	-	VanA- type <i>Enterococcus faecium</i> IOWA 1	-
OXA-244 producing <i>Escherichia coli</i> isolate	-/+	VanB- type <i>Enterococcus faecium</i> IOWA 2	-
TEM-1 (non-ESBL) and IMP-1 producing <i>Escherichia coli</i> isolate	-/+	VanA-type <i>Enterococcus faecalis</i>	-
SHV-12 (ESBL), CTX-M-9 (ESBL) and OXA-48 producing <i>Enterobacter cloacae</i> isolate	-/+	VanB-type <i>E. faecalis</i> (Andrewes and Horder) Schleifer and Kilpper-Balz	-
TEM-1 (non-ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL) and NDM-1 producing <i>Enterobacter cloacae</i> isolate	+/-	VanC and VanB- types <i>Enterococcus gallinarum</i> ENT20120142	-
NDM-7 producing <i>Enterobacter cloacae</i> -complex isolate	+/-	VanC type- <i>Enterococcus gallinarum</i> (Bridge and Sneath 1982) Collins, Jones, Farrow, Kilpper-Balz and Schleifer 1984 VP	-
VIM-1 producing <i>Citrobacter braakii</i> isolate	+/-	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) N315	-
KPC-3 and VIM-4 producing <i>Citrobacter freundii</i> -complex isolate	+/+	Methicillin-resistant <i>Staphylococcus aureus</i> ST398	-
<i>H. pylori</i> Clarithromycin resistant (23S rRNA A2146G)	-	mecC Methicillin-resistant <i>Staphylococcus aureus</i>	-
<i>H. pylori</i> Clarithromycin resistant (23S rRNA A2147G)	-	cMRSA isolate (oxa ^R , PVL-positive, spa:t 310)	-

Table 8. Reference pathogenic microorganisms used in this study performed at CerTest facility. ESBL = Extended Spectrum β -lactamases.

12.4. Analytical reactivity

The reactivity of VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit was evaluated against TEM-1 (non-ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL) and NDM-1 producing *Enterobacter cloacae* isolate; NDM-7 producing *Enterobacter cloacae*-complex isolate; VIM-1 producing *Citrobacter braakii* isolate; KPC-3 and VIM-4 producing *Citrobacter freundii*-complex isolate; SHV-12 (ESBL), CTX-M-9 (ESBL) and OXA-48 producing *Enterobacter cloacae* isolate; OXA-244 producing *Escherichia coli* isolate; SHV-1 (non-ESBL), KPC-3, and OXA-48 producing *Klebsiella pneumonia* isolate; OXA-48 producing *Serratia marcescens* isolate; TEM-1 (non-ESBL), SHV-1 (non-ESBL), CTX-M-2 (ESBL), and KPC-2 producing *Klebsiella pneumonia* isolate and TEM-1 (non-ESBL) and IMP-1 producing *Escherichia coli* isolate, showing positive results.

ANNEX 1

OPEN AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCES
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-CPE106L
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-CPE106H
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-CPE112L
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-CPE112H
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-CPE136
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-CPE101L
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-CPE101H
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-CPE101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition.

Target	Channel	Gene
Carbapenemase-producing <i>Enterobacteriaceae</i> 1	NDM-producing <i>Enterobacteriaceae</i>	FAM
	VIM-producing <i>Enterobacteriaceae</i>	ROX
	Internal control (IC)	HEX, VIC or JOE *
Carbapenemase-producing <i>Enterobacteriaceae</i> 2	OXA-producing <i>Enterobacteriaceae</i>	FAM
	KPC-producing <i>Enterobacteriaceae</i>	ROX
	IMP-producing <i>Enterobacteriaceae</i>	Cy5

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.certest.es.

A1.2 Reagents provided

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3 and A1.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
Carbapenemase-producing <i>Enterobacteriaceae</i> 1 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1/3/6 x 8-well strip
Carbapenemase-producing <i>Enterobacteriaceae</i> 2 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Carbapenemase-producing <i>Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	2/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CPE101L, VS-CPE101H, VS-CPE106L, VS-CPE106H, VS-CPE112L and VS-CPE112H.

Reagent/Material	Description	Colour	Amount
Carbapenemase-producing <i>Enterobacteriaceae</i> 1 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	Transparent	2/9 x 4-well strip
Carbapenemase-producing <i>Enterobacteriaceae</i> 2 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Carbapenemase-producing <i>Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	4/18 x 4-cap strip

Table A1 4. Reagents and materials provided in VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CEL101 and VS-CPE136. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Carbapenemase-producing *Enterobacteriaceae* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized Carbapenemase-producing *Enterobacteriaceae* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A1 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*NDM* and *OXA* genes), ROX (*VIM* and *KPC* genes) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-CPE148T

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target		Channel	Gene
Carbapenemase-producing <i>Enterobacteriaceae</i> 1	NDM-producing <i>Enterobacteriaceae</i>	FAM	NDM
	VIM-producing <i>Enterobacteriaceae</i>	ROX	VIM
	Internal control (IC)	HEX, VIC or JOE *	-
Carbapenemase-producing <i>Enterobacteriaceae</i> 2	OXA-producing <i>Enterobacteriaceae</i>	FAM	OXA
	KPC-producing <i>Enterobacteriaceae</i>	ROX	KPC
	IMP-producing <i>Enterobacteriaceae</i>	Cy5	IMP

Table A2. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.cerest.es.

A2.2 Reagents provided

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
Carbapenemase-producing <i>Enterobacteriaceae</i> 1 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	2 vials
Carbapenemase-producing <i>Enterobacteriaceae</i> 2 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	2 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Carbapenemase-producing <i>Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CPE148T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

Carbapenemase-producing Enterobacteriaceae Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Carbapenemase-producing Enterobacteriaceae* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Carbapenemase-producing Enterobacteriaceae* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (NDM and OXA genes), ROX (VIM and KPC genes) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 3

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-CPE106LE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS- CPE106HE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS- CPE112LE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS- CPE112HE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-CPE136E
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS- CPE101LE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS- CPE101HE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-CPE101E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
Extraction control (EC)	HEX, VIC or JOE *	-
Carbapenemase-producing <i>Enterobacteriaceae</i> 1	NDM-producing <i>Enterobacteriaceae</i> FAM	NDM
	VIM-producing <i>Enterobacteriaceae</i> ROX	VIM
Carbapenemase-producing <i>Enterobacteriaceae</i> 2	OXA-producing <i>Enterobacteriaceae</i> FAM	OXA
	KPC-producing <i>Enterobacteriaceae</i> ROX	KPC
	IMP-producing <i>Enterobacteriaceae</i> Cy5	IMP

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.cerest.es.

A3.2 Reagents provided

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3 y A3.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips compatible devices.

Table A3.4 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es)

Reagent/Material	Description	Colour	Amount
<i>Carbapenemase-producing Enterobacteriaceae 1</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/3/6 x 8-well strip
<i>Carbapenemase-producing Enterobacteriaceae 2</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Carbapenemase-producing Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	2/6/12 x 8-cap strip

Table A3. 3. Reagents and materials provided in VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CPE101LE, VS-CPE101HE, VS-CPE106LE, VS-CPE106HE, VS-CPE112LE and VS-CPE112HE.

Reagent/Material	Description	Colour	Amount
<i>Carbapenemase-producing Enterobacteriaceae 1</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9 x 4-well strip
<i>Carbapenemase-producing Enterobacteriaceae 2</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Carbapenemase-producing Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	4/18 x 4-cap strip

Table A3. 4. Reagents and materials provided in VIASURE *Carbapenemase-producing Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CPE101E and VS-CPE136TE. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) by adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *Carbapenemase-producing Enterobacteriaceae* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

Carbapenemase-producing Enterobacteriaceae Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A3. 5. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (NDM and OXA genes), ROX (VIM and KPC genes) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Carbapenemase-producing <i>Enterobacteriaceae</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-CPE148TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target		Channel	Gene
Extraction control (EC)		HEX, VIC or JOE *	-
Carbapenemase-producing <i>Enterobacteriaceae</i> 1	NDM-producing <i>Enterobacteriaceae</i>	FAM	NDM
	VIM-producing <i>Enterobacteriaceae</i>	ROX	VIM
Carbapenemase-producing <i>Enterobacteriaceae</i> 2	OXA-producing <i>Enterobacteriaceae</i>	FAM	OXA
	KPC-producing <i>Enterobacteriaceae</i>	ROX	KPC
	IMP-producing <i>Enterobacteriaceae</i>	Cy5	IMP

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A4.2 Reagents provided

VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
Carbapenemase-producing <i>Enterobacteriaceae</i> 1 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	2 vials
Carbapenemase-producing <i>Enterobacteriaceae</i> 2 Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	2 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Carbapenemase-producing <i>Enterobacteriaceae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE Carbapenemase-producing *Enterobacteriaceae* Real Time PCR Detection Kit with Ref. VS-CPE148TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Carbapenemase-producing Enterobacteriaceae* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

Carbapenemase-producing Enterobacteriaceae Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Carbapenemase-producing Enterobacteriaceae* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Carbapenemase-producing Enterobacteriaceae* Reaction-Mix (white vial) into each tube.

2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Carbapenemase-producing Enterobacteriaceae* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tube in the thermocycler.

3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*NDM* and *OXA* genes), ROX (*VIM* and *KPC* genes) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

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- Miller, J. M., Binnicker, M. J., Campbell, S., Carroll, K. C., Chapin, K. C., Gilligan, P. H., Gonzalez, M. D., Jerris, R. C., Kehl, S. C., Patel, R., Pritt, B. S., Richter, S. S., Robinson-Dunn, B., Schwartzman, J. D., Snyder, J. W., Telford, S., Theel, E. S., Thomson, R. B., Weinstein, M. P., & Yao, J. D. (2018). A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2018 Update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1–e94. <https://doi.org/10.1093/cid/ciy381>
- Monteiro, J., Widen, R. H., Pignatari, A. C. C., Kubasek, C., & Silbert, S. (2012). Rapid detection of carbapenemase genes by multiplex real-time PCR. *Journal of Antimicrobial Chemotherapy*, 67(4), 906–909. <https://doi.org/10.1093/jac/dkr563>
- Pollett, S., Miller, S., Hindler, J., Uslan, D., Carvalho, M., & Humphries, R. M. (2014). Phenotypic and molecular characteristics of carbapenem-resistant Enterobacteriaceae in a health care system in Los Angeles, California, from 2011 to 2013. *Journal of Clinical Microbiology*, 52(11), 4003–4009. <https://doi.org/10.1128/JCM.01397-14>
- Singh-Moodley, A., & Perovic, O. (2016). Antimicrobial susceptibility testing in predicting the presence of carbapenemase genes in Enterobacteriaceae in South Africa. *BMC Infectious Diseases*, 16(1), 1–10. <https://doi.org/10.1186/s12879-016-1858-7>
- Teo, J. W. P., La, M. Van, & Lin, R. T. P. (2016). Development and evaluation of a multiplex real-time PCR for the detection of IMP, VIM, and OXA-23 carbapenemase gene families on the BD MAX open system. *Diagnostic Microbiology and Infectious Disease*, 86(4), 358–361. <https://doi.org/10.1016/j.diagmicrobio.2016.08.019>
- Van der Zee, A., Roorda, L., Bosman, G., Fluit, A. C., Hermans, M., Smits, P. H. M., van der Zanden, A. G. M., te Witt, R., Bruijnesteijn van Coppenraet, L. E. S., Cohen Stuart, J., & Ossewaarde, J. M. (2014). Multi-centre evaluation of real-time multiplex PCR for detection of carbapenemase genes OXA-48, VIM, IMP, NDM and KPC. *BMC Infectious Diseases*, 14(1), 1–5. <https://doi.org/10.1186/1471-2334-14-27>
- World Health Organization. (2017). Guidelines for the prevention and control of carbapenem-resistant Enterobacteriaceae, Acinetobacter baumannii and Pseudomonas aeruginosa in health care facilities. In WHO. <http://apps.who.int/bookorders>.

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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