VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

SARS-CoV-2

Handbook for the following references/

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VIASURE SARS-CoV-2 Real Time PCR Detection Kit 6 x 8-well strips, low profile

VS-NCO206LE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 6 x 8-well strips, high profile

VS-NCO206HE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 12 x 8-well strips, low profile

VS-NCO212LE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 12 x 8-well strips, high profile

VS-NCO212HE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 96-well plate, low profile

VS-NCO213LE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 96-well plate, high profile

VS-NCO213HE

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®

VS-NCO236E

VIASURE SARS-CoV-2 Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®

VS-NCO272E

ENGLISH

1. Intended use

VIASURE SARS-CoV-2 Real Time PCR Detection Kit is a real-time RT-PCR test designed for the qualitative detection of RNA from the SARS-CoV-2 in respiratory samples (nasopharyngeal swab, nasal swab and oropharyngeal swab) and saliva samples from individuals suspected of COVID-19 infection by their healthcare provider. This test is intended for use as an aid in the diagnosis of SARS-CoV-2 in combination with clinical and epidemiological risk factors. RNA is extracted from respiratory specimens and saliva, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for SARS-CoV-2.

2. Summary and Explanation

Coronavirus are enveloped non-segmented positive-sense RNA viruses and belong to *Coronaviridae* family. There are six coronavirus species known to cause human diseases. Four viruses (229E, OC43, NL63 and HKU1) cause common cold symptoms and the other two (severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV)) are zoonotic and producing more severe complications. SARS-CoV and MERS-CoV have caused more than 10,000 cumulative cases in the past two decades, with mortality rates of 34% MERS-CoV and 10% SARS-CoV.

In December 2019, some people that worked at or lived around the Huanan seafood market in Wuhan, Hubei Province, China, have presented pneumonia of unknown cause. Deep sequencing analysis of the respiratory samples indicated a novel coronavirus, which was named firstly 2019 novel coronavirus (2019-nCoV) and lately SARS-CoV-2.

Human-to-human transmission of the SARS-CoV-2 has been confirmed, even in the incubation period without symptoms, and the virus causes severe respiratory illness like those SARS-CoV produced. Although the pneumonia is the principal illness associated, a few patients have developed severe pneumonia, pulmonary edema, acute respiratory distress syndrome, or multiple organ failure and death. Centers of Disease Control and Prevention (CDC) believes that symptoms of SARS-CoV-2 may appear in as few as 2 days or as long as 14 days after exposure, being the most common fever, cough, myalgia and dyspnea. Less common symptoms are sore throat, headache, diarrhea and vomiting. It seems that people above 60 years old, males, and people with comorbidities most often have severe disease.

Diagnosis of SARS-CoV-2 is performed detecting conventional causes of pneumonia early and detected by next-generation sequencing or real-time RT-PCR methods. Several assays that detect the SARS-CoV-2 have been are currently available, such as China CDC (gene targets, *ORF1ab* and *N*), Charité – Germany (gene targets, *RdRP* and *E*) or US CDC (two targets in *N* gene).

CDC recommends upper respiratory tract specimens (nasopharyngeal (NP) swab, oropharyngeal (OP) swabs, nasal mid-turbinate swab, nasal swab, nasopharyngeal wash/aspirate or nasal wash/aspirate (NW) and a saliva specimens collected mainly by a healthcare provider) and/or lower respiratory specimens (sputum, endotracheal aspirate, or bronchoalveolar lavage in patients with more severe respiratory disease) for the identification of

SARS-CoV-2. In addition, other clinical specimens as blood, urine and stool may be collected to monitor the presence of the virus.

3. Principle of the procedure

VIASURE SARS-CoV-2 Real Time PCR Detection Kit is designed for the qualitative detection of RNA from the SARS-CoV-2 in respiratory and saliva samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of *ORF1ab* and *N* genes for SARS-CoV-2 using specific primers and a fluorescent-labeled probe.

VIASURE SARS-CoV-2 Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE SARS-CoV-2 Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in an stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity. *ORF1ab* gene is amplified and detected in FAM channel, *N* gene is amplified and detected in ROX channel and the Extraction control (EC) in HEX channel, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE SARS-CoV-2 Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 x 8-cap strip

Table 1. Reagents and materials provided in VIASURE SARS-CoV-2 Real Time PCR Detection Kit with Ref. VS-NCO206LE, VS-NCO206HE, VS-NCO212LE and VS-NCO212HE.

Reagent/Material	Description	Color	Amount
SARS-CoV-2 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table 2. Reagents and materials provided in VIASURE SARS-CoV-2 Real Time PCR Detection Kit with Ref VS-NCO213LE and VS-NCO213HE.

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 x 4-cap strip

Table 3. Reagents and materials provided in VIASURE SARS-CoV-2 Real Time PCR Detection Kit with Ref. VS-NCO236E and VS-NCO272E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE SARS-CoV-2 Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE SARS-CoV-2 Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Applied Biosystems QuantStudio™ 5 Real-Time PCR System, Applied Biosystems QuantStudio™ 12K Flex Real-Time PCR System, 96-well block, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), Roche Molecular Diagnostics Cobas z480 Analyzer, and geneLEAD VIII System ((Precision System Science Co., Ltd. (PSS)). When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.

- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-NCO213LE, VS-NCO213HE, VS-NCO236E and VS-NCO272E). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For VS- NCO236E and VS-NCO272E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local, state, and federal regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Specimen collection, transport and storage

The VIASURE SARS-CoV-2 Real Time PCR Detection kit has been validated on nasopharyngeal and oropharyngeal specimens collected with synthetic fiber swabs with plastic and placed immediately into a sterile transport tube containing Universal transport medium (UTM), Viral transport medium (VTM), or saline solution (0.9 % sodium chloride); nasal swabs collected with Aptima® Multitest Swab Specimen Collection Kit (HOLOGIC) and Aptima® Unisex Swab Specimen Collection Kit for Female Endocervical and Male Urethral Swab Specimens (HOLOGIC); and saliva samples collected in a sterile clean container (No preservative is required).

Patient samples must be collected, transport and storage according to appropriate laboratory guidelines. For details, refer to the CDC guidelines (Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19) (website <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>).

8.2. Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized SARS-CoV-2 Positive Control in an area away from the other components.

8.3. RNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

If the Extraction Control is used to monitor nucleic acid isolation and as PCR inhibition control, add 5 µL of the EC to the specimen and/or lysis buffer mixture. (clinical specimen, as well as, positive control and/or negative control). Close each tube and vortex for 10 seconds.

If the Extraction Control is used only as a PCR inhibition control, 1 µL of the EC should be added to the reconstituted Reaction-Mix.

For RNA extraction from respiratory samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).
- Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (ROCHE).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- MagCore® Viral Nucleic Acid Extraction kit, using the MagCore® HF16 automated Nucleic Acid Extractor System.
- EZ1 Virus Mini Kit, using EZ1 instrument (Qiagen).
- mSample Preparation Systems RNA, using the Abbott m2000 RealTime System (Abbott Molecular).
- QIAamp Viral RNA Mini Kit (Qiagen).
- QIAamp MinElute Virus Spin kit automated using the QIAcube instrument (Qiagen).
- Automatic extraction system using the chemagic™ 360-D instrument (PerkinElmer).
- RNEasy Mini Kit (Qiagen).
- epMotion® 5070/5075 (Eppendorf).

- MagNA Pure LC Total Nucleic Acid Isolation Kit (Roche Molecular Systems).
- MagNA Pure 96 External Lysis Buffer and MagNA Pure 96 DNA and Viral NA Small Volume Kit using the MagNA Pure 96 Instrument (using Viral NA Plasma ext lys SV protocol) (Roche Molecular Systems).
- Smart-32 Nucleic Acid Extraction Instrument (DaAn Gene).
- EMAG® Nucleic Acid Extraction System (bioMérieux, Inc).

For RNA extraction from saliva samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- QIAamp Viral RNA Mini Kit (Qiagen).
- Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).

8.4. Lyophilized positive control

SARS-CoV-2 Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized SARS-CoV-2 Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.5. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of RNA sample in different wells.

Add 5 µL of reconstituted SARS-CoV-2 Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*ORF1ab* gene), ROX (*N* gene) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

All the result of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions. Using the following tables 5 and 6 read and analyze the results.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for SARS-CoV-2 in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>ORF1ab</i> gene (FAM) ¹	<i>N</i> gene (ROX) ¹	Extraction Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≥40 or no signal	≤40	Valid

Table 5. Expected Performance of Controls

¹ In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

² The Extraction Control (EC) should show an amplification signal (Ct ≤40 in control wells (PC and NC)). Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted. For interpretation of patient sample results, use the following table:

ORF1ab gene (FAM)	N gene (ROX)	Extraction Control (HEX)	Interpretation for patients' samples	
≤40	≤40	≤40 or no signal ¹	Valid	SARS-CoV-2 RNA Detected
≤40	≥40 or no signal	≤40 or no signal ¹	Inconclusive	If only one SARS-CoV-2 target gene amplifies, repeat test depending on the available material: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat RT-qPCR with the same isolated RNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered SARS-CoV-2 positive.
≥40 or no signal	≤40	≤40 or no signal ¹	Inconclusive	
≥40 or no signal	≥40 or no signal	≤ 35 ²	Valid	SARS-CoV-2 RNA not Detected ²
≥40 or no signal	≥40 or no signal	≥ 35 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 6. Interpretation of patient sample results. Ct values. no signal = no amplification curve.


1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

2 In the case of SARS-CoV-2 target genes negative, EC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the RNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each RT-qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with RNA extracted from respiratory samples (nasopharyngeal swab, nasal swab, and oropharyngeal swab) and saliva samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, either the high number of cDNA template copies which contains each SARS-CoV-2 Positive Control vial, samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and SARS-CoV-2 Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.
- The specific primer and probe combinations for detection of the *ORF1ab* and *N* genes used in VIASURE SARS-CoV-2 Real Time PCR Detection Kit designed for the detection of SARS-CoV-2, do not show significant combined homologies with the human genome, human microflora, or other coronaviruses (with the exception of some *N* and/or *ORF1ab* sequences from SARS-CoV, and other coronaviruses identified in bats and pangolin), which might result in predictable false positive.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including RNA extraction).
 - Degradation of the viral RNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown SARS-CoV-2 variants.
 - A viral load in the specimen below the limit of detection for the assay.

- The presence of RT-qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent COVID-19 or used during the treatment of the infection have not been evaluated.
- Failure to follow instructions for use and the assay procedure.
- A single-gene amplification or even random (irregular) positive results is suggestive of slightly different amplification yield of the target regions. Samples with low viral load might usually result in *N* single-gene amplification. It has been described that *N* gene assay might be more sensitive than the *ORF1ab* gene assay detecting positive clinical specimens. In case of a doubt, it is recommended referring to a reference laboratory for further testing if clinically indicated.
- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets viral sequences (*ORF1ab* and/or *N* genes).
- Negative results or single target amplification do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by SARS-CoV-2 have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the virus.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that SARS-CoV-2 infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.

11. Quality control

VIASURE SARS-CoV-2 Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Extraction Control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE SARS-CoV-2 Real Time PCR Detection Kit was tested using 100 respiratory specimens (nasopharyngeal swab) from symptomatic patients with suspicion of COVID-19. These results were compared with those obtained with a molecular detection method by the Spanish National Reference Center (Institute of Health Carlos III (ISCIII)) (the protocol "2019-nCoV by real-time RT-PCR" suggested by Charité (Berlin), with modifications). The results were as follows:

VIASURE SARS-CoV-2 Real Time PCR Detection Kit	Reference Molecular detection method (Protocol suggested by Charité Universitätsmedizin Berlin)			
		+	-	Total
	+	2	0	2
	-	0	98	98
	Total	2	98	100

Table 7. Comparative results for SARS-CoV-2.

Sensitivity, specificity, Positive Predictive Value (PPV) and Negative Predictive Value (NPV) for VIASURE SARS-CoV-2 Real Time PCR detection kit compared with Charité protocol are >99.9%.

The clinical performance of VIASURE SARS-CoV-2 Real Time PCR Detection Kit was conducted using 80 anonymized respiratory specimens from the biobank of the National Center for Microbiology (NCM) (Spain). These results were compared with those obtained with a molecular detection method by the Spanish National Reference Center (ISCIII), (the protocol "2019-nCoV by real-time RT-PCR" suggested by Charité (Berlin), with modifications) (Corman et al., 2020). The results were as follows:

VIASURE SARS-CoV-2 Real Time PCR Detection Kit	Reference Molecular detection method (Protocol suggested by Charité Universitätsmedizin Berlin)			
		+	-	Total
	+	41	1*	42
	-	0	38	38
	Total	41	39	80

Table 8. Comparative results for SARS-CoV-2. *One sample resulted positive only for N gene.

Sensitivity and specificity for VIASURE SARS-CoV-2 Real Time PCR Detection Kit compared with Charité protocol are >99.9% and 97.5%, respectively.

The clinical performance of VIASURE SARS-CoV-2 Real Time PCR Detection Kit was retrospective and prospectively evaluated in comparison to two commercial molecular methods, Thermo-Fisher TaqPath™ COVID-19 Combo Kit or Aptima® SARS-CoV-2 Assay (Panther® System) (Hologic).

Firstly, a retrospective evaluation was carried out with leftover frozen nasopharyngeal samples collected in UTM or saline solution and nasal swabs collected with Aptima® Multitest Swab Specimen Collection Kit or Aptima® Unisex Swab Specimen Collection Kit for Female Endocervical and Male Urethral Swab Specimens (HOLOGIC), from individuals suspected of SARS-CoV-2 infection by the healthcare provider in the COVID-19 disease endemic

region(s). These samples have been previously characterized by Thermo-Fisher TaqPath™ COVID-19 Combo Kit as molecular reference methods. The results were as follows:

VIASURE SARS-CoV-2 Real Time PCR Detection Kit	Thermo-Fisher TaqPath™ COVID-19 Combo Kit			
		+	-	Total
	+	34	1	35
	-	0	37	37
	Total	34	38	72

Table 9. Comparative results for SARS-CoV-2.

The Positive Percent Agreement (PPA), Negative Percent Agreement (NPA) and Overall Percent Agreement (OPA) for VIASURE SARS-CoV-2 Real Time PCR detection kit were calculated in relation to the results from the TaqPath™ COVID-19 Combo Kit (Thermo Fisher Scientific) as shown in Table 10.

Microorganism	PPA (%)	NPA (%)	OPA (%)
SARS-CoV-2	100 (89.8 to 100)	98.9 (93.7 to >99.9)	99.3 (95.5 to >99.9)

Table 10. PPA, NPA and OPA values for VIASURE SARS-CoV-2 Real Time PCR detection kit (95% confidence interval).

A second similar prospective evaluation was also performed using leftover nasopharyngeal and nasal samples collected in UTM or saline solution from individuals suspected of SARS-CoV-2 infection. These samples were simultaneously characterized by Aptima® SARS-CoV-2 Assay (Panther® System) (Hologic). The results were as follows:

VIASURE SARS-CoV-2 Real Time PCR Detection Kit	Aptima® SARS-CoV-2 Assay (Panther® System) (Hologic)			
		+	-	Total
	+	7	0	7
	-	0	57	57
	Total	7	57	64

Table 11. Comparative results for SARS-CoV-2.

The Positive Percent Agreement (PPA), Negative Percent Agreement (NPA) and Overall Percent Agreement (OPA) for VIASURE SARS-CoV-2 Real Time PCR detection kit were calculated in relation to the results from the Aptima® SARS-CoV-2 Assay (Panther® System) (Hologic), as shown in Table 12.

Microorganism	PPA (%)	NPA (%)	OPA (%)
SARS-CoV-2	100 (56.1-100)	100 (92.1-100)	100 (93.2 – 100)

Table 12. PPA, NPA and OPA values for VIASURE SARS-CoV-2 Real Time PCR detection kit (95% confidence interval).

Besides, the clinical performance of VIASURE SARS-CoV-2 Real Time PCR Detection Kit was conducted testing in parallel 28 anonymized saliva specimens and oropharyngeal swabs from patients with signs and symptoms of COVID-19 infection. The results were as follows:

VIASURE SARS CoV-2 Real Time PCR Detection Kit (saliva specimens)	VIASURE SARS CoV-2 Real Time PCR Detection Kit (oropharyngeal swabs)			
		+	-	Total
	+	18	0	18
	-	1	9	10
	Total	19	9	28

Table 13. Comparative results for different types of clinical specimens.

The values obtained with the oropharyngeal samples were taken as the reference values and it was demonstrated that the results obtained when analyzing the saliva samples in parallel, were as expected.

The results show high agreement to detect SARS-CoV-2 using VIASURE SARS-CoV-2 Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE SARS-CoV-2 Real Time PCR Detection Kit has a detection limit of 20 genome copies per reaction for *ORF1ab* and *N* genes with a positive rate of $\geq 95\%$ (Figures 2 and 3).

Figure 2. Dilution series of *ORF1ab* gene (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (FAM channel).



Figure 3. Dilution series of *N* gene (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (ROX channel).

12.3. Analytical specificity

The specificity of the SARS-CoV-2 assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
Human Adenovirus types 1-5, 8, 14, 15, 31, 40 and 41	-	Influenza A/California/7/2009(H1N1)pdm09	-	<i>Legionella longbeachae</i>	-
Human Bocavirus	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella micdadei</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Virginia/ATCC2/2009 (H1N1)pdm09	-	Human metapneumovirus A, B and 8 Peru6-2003	-
<i>Bordetella parapertussis</i>	-	Influenza A/Victoria/210/2009 (H3N2)	-	<i>Moraxella catarrhalis</i>	-
<i>Bordetella pertussis</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus	-	<i>Mycoplasma pneumoniae</i>	-
<i>Candida albicans</i> 3147	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> not rifampin resistant	-
<i>Chlamydia caviae</i>	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Chlamydia psittaci</i> genotype A and C	-	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	<i>Pneumocystis jirovecii</i>	-
<i>Chlamydophila pneumoniae</i> CM-1 and TWAR strain 2023	-	Influenza A/Wisconsin/15/2009 (H3N2)	-	<i>Pseudomonas aeruginosa</i> PAO1-LAC	-
Human coronavirus 229E, OC43, NL63 and HKU1	-	Influenza A/DE-SH/Reihenente/AR8444/ 2016 (H5N8) virus	-	Human rhinovirus type C and 17	-
MERS Coronavirus	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
SARS Coronavirus Strain Frankfurt 1	-	Influenza B/Brisbane/60/2008	-	<i>Staphylococcus epidermidis</i> strain RP62A	-
Enterovirus 68, 71 and D68 US/MO/14-18947	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> 2022 and strain TIGR4	-
Enterovirus Echovirus 11 and 30	-	Influenza B/Phuket/3073/2013 virus	-	<i>Streptococcus pyogenes</i> typing strain T1	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Legionella bozemanii</i>	-	<i>Streptococcus salivarius</i> subsp. <i>salivarius</i> strain DSM 13084	-
<i>Haemophilus influenzae</i> MinnA and Rd	-	<i>Legionella dumoffii</i>	-	Respiratory syncytial virus (RSV) A and B	-
Influenza A/New Caledonia/20/99(H1N1) virus	-				

Table 14. Reference pathogenic microorganisms used in this study

12.4. Analytical reactivity

The reactivity of VIASURE SARS-CoV-2 Real Time PCR Detection Kit was evaluated against RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1, Human 2019-nCoV strain 2019-nCoV/Italy-INMI1, SARS-CoV-2 strain 2019nCoV/USA-WA1/2020, SARS-CoV-2 strain hCoV-19/Netherlands/Diemen_1363454/2020 and synthetic RNA controls for two variants of the SARS-CoV-2 virus: MT007544.1 (SARS-CoV2 isolate Australia/VIC01/2020) and MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1), showing positive result.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

VIASURE Real Time PCR Detection Kits are available in a ready-to-use lyophilized format placed inside wells with different dimensions, low or high profile. Depending on the thermal block of the equipment to be used, one measure or another will fit. Please, consult the table and check the specifications of your equipment. If the equipment does not appear in the list below, please contact with your supplier. This table is for guidance, it is recommended to check the equipment before running the (RT)-qPCR.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast / 7500 Fast Dx Real-Time PCR System ^{(1) (5)}
	QuantStudio™ 12K Flex 96-well Fast
	QuantStudio™ 6 Flex 96-well Fast
	QuantStudio™ 7 Flex 96-well Fast
	QuantStudio™ 3 Fast Real-Time PCR System ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
	StepOne Plus™ Real-Time PCR System ⁽²⁾
	StepOne™ Real-Time PCR System ⁽²⁾
	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96 Fast
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
	Mini Opticon™ Real-Time PCR Detection System ⁽³⁾
Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾
Cepheid	SmartCycler® ⁽⁴⁾
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
Qiagen	Rotor-Gene® Q ⁽⁴⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁵⁾
	LightCycler® 96 Real-Time PCR System ⁽⁵⁾
	Cobas z480 Analyzer ⁽⁵⁾

(1) Select Ramp Speed "Standard".

(2) No detection in Cy5 channel.

(3) Detection in FAM and HEX channels only

(4) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System tubes.

(5) Shell Frame grid plate which fits in these qPCR System is necessary.

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System ⁽⁵⁾
Applied Biosystems	7300 Real-Time PCR System ^{(2) (5)}
	7500 Real-Time PCR System ⁽⁵⁾
	7900 HT Real-Time PCR System ⁽²⁾
	ABI PRISM 7000 ⁽³⁾
	ABI PRISM 7700 ⁽²⁾
	QuantStudio™ 12K Flex 96-well
	QuantStudio™ 6 Flex 96-well
	QuantStudio™ 7 Flex 96-well
	QuantStudio™ 3 Real-Time PCR System ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
	iCycler iQ™ Real-Time PCR Detection System
	iCycler iQ™5 Real-Time PCR Detection System
	MyiQ™ Real-Time PCR Detection System ⁽³⁾
	MyiQ™2 Real-Time PCR Detection System ⁽³⁾
Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾
Cepheid	SmartCycler® ⁽⁴⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler
	DTlite Real-Time PCR System
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽⁴⁾
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be "none". Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Lightcycler®480II Roche	FAM	465/510	Colour Compensation is required for Roche Thermocyclers
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Cobas z 480 Roche	FAM	465/510	Colour Compensation is required for Roche Thermocyclers
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene/Agilent Technologies	FAM	FAM	Passive reference option for ROX must be "none"
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTLite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía

- Huang, C. *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30183-5.
- Zhu N. *et al.* A novel coronavirus from patients with pneumonia in China, 2019. *New England Journal of Medicine*, 2020. DOI : 10.1056/NEJMod2001017.

3. World Health Organization. MERS situation update. January 2020. Available from <https://applications.emro.who.int/docs/EMCSR254E.pdf?ua=1> Accessed November 2020.
4. Chen N. *et al.*. Epidemiological and Clinical Characteristics of 99 Cases of 2019-Novel Coronavirus (2019-nCoV) Pneumonia in Wuhan, China. *The Lancet*, 2020. DOI: 10.1016/S0140-6736(20)30211-7.
5. Lv D.F. *et al.* Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clinica Chimica Acta* 2020; 506: 172-175.
6. World Health Organization. Clinical management of COVID-19 disease" Interim guidance 27 May 2020. Available from <https://www.who.int/publications/i/item/clinical-management-of-covid-19> Accessed November 2020.
7. Lu R. *et al.* Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30251-8.
8. Rothe C. *et al.* Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*, 2020. DOI : 10.1056/NEJMc2001468.
9. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Symptoms of Coronavirus. Available from <https://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html> Accessed November 2020.
10. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Older Adults. Available from <https://www.cdc.gov/coronavirus/2019-ncov/need-extra-precautions/older-adults.html> Accessed November 2020.
11. World Health Organization. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases. Interim guidance. 19 March 2020. Available from <https://www.who.int/publications-detail/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed November 2020.
12. Yan Y *et al.* Laboratory testing of SARS-CoV, MERS-CoV, and SARS-CoV-2 (2019-nCoV): Current status, challenges, and countermeasures. *Reviews in Medical Virology* 2020; 30(3):e2106.
13. Centers of Disease Control and Prevention (CDC). 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes. <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed November 2020.
14. Chu D.K.W. *et al.* Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clinical Chemistry* 2020;66(4): 549-555.
15. Corman V.M. *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *European communicable disease bulletin* 2020;25(3).
16. European Centre for Disease Prevention and Control. Novel coronavirus disease 2019 (COVID-19) pandemic: increased transmission in the EU/EEA and the UK – sixth update – 12 March 2020. Stockholm: ECDC; 2020. Available from <https://www.ecdc.europa.eu/sites/default/files/documents/RRA-sixth-update-Outbreak-of-novel-coronavirus-disease-2019-COVID-19.pdf> Accessed November 2020.
17. Lim, Y. X., Ng, Y. L., Tam, J. P., & Liu, D. X. (2016). Human coronaviruses: a review of virus–host interactions. *Diseases*, 4(3), 26.
18. McBride R. *et al.* The coronavirus nucleocapsid is a multifunctional protein. *Viruses* 2014; 6(8):2991-3018.
19. Sheikh A. *et al.* Analysis of preferred codon usage in the coronavirus N genes and their implications for genome evolution and vaccine design. *Journal of Virological Methods* 2020; 277:113806.

20. World Health Organization. Public health surveillance for COVID-19. 7 August 2020. Available from [https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-\(2019-ncov\)](https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-(2019-ncov)). Accessed November 2020.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System ⁽⁵⁾
Applied Biosystems	7300 Real-Time PCR System ^{(2) (5)}
	7500 Real-Time PCR System ⁽⁵⁾
	7900 HT Real-Time PCR System ⁽²⁾
	ABI PRISM 7000 ⁽³⁾
	ABI PRISM 7700 ⁽²⁾
	QuantStudio™ 12K Flex 96-well
	QuantStudio™ 6 Flex 96-well
	QuantStudio™ 7 Flex 96-well
	QuantStudio™ 3 Real-Time PCR System ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
	iCycler iQ™ Real-Time PCR Detection System
	iCycler iQ™5 Real-Time PCR Detection System
	MyiQ™ Real-Time PCR Detection System ⁽³⁾
	MyiQ™2 Real-Time PCR Detection System ⁽³⁾
Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾
Cepheid	SmartCycler® ⁽⁴⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler
	DTlite Real-Time PCR System
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽⁴⁾
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.
- Magtration® and MagDEA® Dx are PSS patented technologies (Japanese Patent: No.3115501; US Patent: No.5,702,950; European Patent: No.687501)
- chemagic™ 360-D instrument is a registered trademark of PerkinElmer.
- epMotion® 5070/5075 is a registered trademark of Eppendorf.
- EMAG®Nucleic Acid Extraction System is a registered trademark of bioMérieux, Inc.

Revision: November 2020

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.cerTest.es

VIASURE online

F-362 rev01

VIASURE

CE

IVD

Real Time PCR Detection Kits

CerTest
BIOTEC