VIASURE

Real Time PCR Detection Kits


by CerTest
BIOTEC

Sexually Transmitted Diseases

Handbook for the following references/

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VIASURE <i>Sexually transmitted diseases</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-STD106L
VIASURE <i>Sexually transmitted diseases</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-STD106H
VIASURE <i>Sexually transmitted diseases</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-STD112L
VIASURE <i>Sexually transmitted diseases</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-STD112H
VIASURE <i>Sexually transmitted diseases</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-STD136

ENGLISH

1. Intended use

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and/or *Mycoplasma hominis*, in urogenital, endocervical and urine samples from patients with signs and symptoms of sexually transmitted diseases (STDs). This test is intended for use as an aid in the diagnosis of STDs in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*.

2. Summary and Explanation

Sexually transmitted infections (STIs) represent a group of diseases that affect the sexual and reproductive health of millions of people, being a public problem of interest. Etiological agents responsible for STIs include fungi, bacteria, parasites and viruses. Some of these microorganisms are eliminated after a period of time, while others are recurrent and some remain in the body asymptotically, allowing the progress of the disease and generating consequences such as inflammations of the genito-urinary tract, infertility and even the development of cancer.

Trichomonas vaginalis infection is one of the most common sexually transmitted diseases (STDs) in the world. *Trichomonas vaginalis* is a flagellated pathogen protozoan belonging to the order *Trichomonadida* located in both, male and female urogenital tract but it has also been isolated from the respiratory tract of infants and adults. In women it can be found in the vagina and in the urethra, while in men it can be found in the urethra, the prostate and the epididymis. *Trichomonas vaginalis* infection has been associated with vaginitis, cervicitis and urethritis, premature rupture of membranes and premature delivery in pregnant women. *Trichomonas vaginalis* infection has also been associated with an increased risk of HIV acquisition and transmission in women.

Among other sexually transmitted pathogens, four species of bacteria of the genus *Mycoplasma* can be found (*Mycoplasma hominis*, *Mycoplasma genitalium*, *Ureaplasma urealyticum* and *Ureaplasma parvum*) whose presence can also lead to various infections of the genitourinary tract in humans.

Mycoplasma hominis colonizes the lower urogenital tract and is associated with urogenital infections, particularly bacterial vaginosis and non-gonococcal urethritis. It is also involved in extra genital infections, such as postpartum or post-abortion fever, in post-cesarean wound infections or after a hysterectomy. In neonates, it can cause meningitis, brain abscesses and eye infections. In adults, bacteremia, septic arthritis, osteitis, endocarditis, mediastinitis, brain abscesses and respiratory infections have been described. Most patients have predisposing factors, including immunosuppression, trauma, respiratory problems or post-manipulation and/or surgery of the genitourinary tract. It is considered a bacterium present in the normal micro flora in 20% of the male population and in 40% of the female population.

Mycoplasma genitalium is a facultative anaerobic organism and a recognized cause of nongonococcal urethritis in men. In women, *M. genitalium* has been associated with cervicitis, endometritis, pelvic inflammatory disease (PID), infertility, susceptibility to human immunodeficiency virus (HIV), and adverse birth outcomes, indicating a consistent relationship with female genital tract pathology. The long-term reproductive consequences of *M. genitalium* infection in asymptomatic individuals need to be investigated further.

Ureaplasma spp. are often isolated from human genital mucosa of individuals with a lack of symptoms. In humans, two major species, namely *Ureaplasma parvum* and *Ureaplasma urealyticum* (family *Mycoplasmataceae*, genus *Ureaplasma*) are part of the genital flora of men and women and are present in almost 70% of sexually active population. These bacteria cause inflammation and lead to chorioamnionitis, preterm deliveries, and premature rupture of membranes. Recently, their pathogenic roles on spermatogenesis and subsequent male infertility have been also suggested. Although both male and female carriers are asymptomatic, *Ureaplasma* spp. is occasionally isolated from the neonatal lesions of bronchopulmonary dysplasia, intraventricular haemorrhage, and necrotizing enterocolitis.

Neisseria gonorrhoeae is an obligate human pathogen and is the etiological agent of gonorrhoea. Syndromes include cervicitis in women, and urethritis, pharyngitis and proctitis in both sexes. If untreated, women may experience severe sequelae of pelvic inflammatory disease, chronic pelvic pain, ectopic pregnancy and tubal infertility, while men may develop epididymitis, prostatitis and urethral stricture.

Chlamydia trachomatis is a small intracellular bacterium that requires living cells to multiply. There are 18 serotypes; D-K are the ones that cause sexually transmitted infections, as well as neonatal infections. The majority of genital chlamydial infections in both males and females are asymptomatic. When symptoms do occur, lower urogenital tract infection can manifest as cervicitis in females and urethritis in males and females. Whether symptomatic or asymptomatic, untreated chlamydia can ascend to the upper genital tract. In males, this can cause epididymitis, which is not thought to be an important cause of long-term sequelae. However, in females, upper tract infection can result in pelvic inflammatory disease (PID), a spectrum of clinical disorders involving infection and inflammation of the uterus, fallopian tubes, ovaries, or adjacent peritoneum.

Conventional diagnostic assays, lack sensitivity, require viable organisms and thus special shipment conditions and, sometimes, invasive sampling. As nucleic acid amplification tests, as real time PCR, allow us to overcome some of these limitations. Several molecular diagnostic assays have recently been commercialized to assist the syndromic diagnosis of STIs. In addition, their implicit multiplexing capacity allows for the detection of multiple pathogens in a single sample.

3. Principle of the procedure

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and/or *Mycoplasma hominis*, in urogenital, endocervical and urine samples from patients with signs and symptoms of sexually transmitted diseases (STDs). After DNA isolation, the identification of the STDs is performed by the amplification of a conserved region of the *T. vaginalis*-specific 2-kb repeated sequence, (*Trichomonas vaginalis*), *ureasa* gene (*Ureaplasma urealyticum* and *Ureaplasma parvum*),

yidC gene (*Mycoplasma hominis*), *porA* and *Opa* genes (*Neisseria gonorrhoeae*), a region within ORF2 of the chlamydial plasmid (*Chlamydia trachomatis*) and *MgPa adhesin* gene (*Mycoplasma genitalium*), using specific primers and a fluorescent-labelled probe.

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in an stabilized format, as well as an internal control to monitor PCR inhibition. Each kit includes two kinds of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of *Neisseria gonorrhoeae*, *Chlamydia trachomatis* and *Mycoplasma genitalium* (*Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium* 8/4-well strips). *Neisseria gonorrhoeae* DNA targets are amplified and detected in ROX channel, *Chlamydia trachomatis* DNA targets are amplified and detected in FAM channel and *Mycoplasma genitalium* DNA targets are amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2). The first strip (*Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium* 8/4-well strips) contains an internal control (IC) which is amplified and detected in Cy5 channel. The second strip contains the multiplex reaction mix for the detection of *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis* (*Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis* 8/4-well strips). *Trichomonas vaginalis* DNA targets are amplified and detected in FAM channel, *Ureaplasma urealyticum* DNA targets are amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2), *Ureaplasma parvum* DNA targets are amplified and detected in ROX channel and *Mycoplasma hominis* DNA targets are amplified and detected in Cy5 channel.

4. Reagents provided

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1 and 2. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>Neisseria gonorrhoeae</i> , <i>Chlamydia trachomatis</i> & <i>Mycoplasma genitalium</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	3/6 x 8-well strip
<i>Trichomonas vaginalis</i> , <i>Ureaplasma urealyticum</i> , <i>Ureaplasma parvum</i> & <i>Mycoplasma hominis</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	White	3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Sexual Transmitted Diseases Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit with Ref. VS-STD106L, VS-STD106H, VS-STD112L and VS-STD112H.

Reagent/Material	Description	Colour	Amount
<i>Neisseria gonorrhoeae</i> , <i>Chlamydia trachomatis</i> & <i>Mycoplasma genitalium</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9 x 4-well strip
<i>Trichomonas vaginalis</i> , <i>Ureaplasma urealyticum</i> , <i>Ureaplasma parvum</i> & <i>Mycoplasma hominis</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	Transparent	9 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Sexual Transmitted Diseases Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	18 X 4-cap strip

Table 2. Reagents and materials provided in VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit with Ref. VS-STD136. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.

- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid) and Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-STD136). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For VS-STD136 (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Make sure to use a well for determining *Neisseria gonorrhoeae*, *Chlamydia trachomatis* and *Mycoplasma genitalium* and another well for determining *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*. Be careful not to mix them throughout the process.

- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.
- The product VIASURE *Sexually Transmitted diseases* Real Time PCR Detection kit has only been validated with the equipment mentioned in Section 5 of these Instructions for Use. The rest of the compatible Real Time PCR instrument indicated in Annex 1 is based on bibliographic data.

8. Test procedure

8.1. Sample preparation

Clinical specimens (urogenital specimens, urine samples, rectal and endocervical samples) should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples should be frozen at -20°C . In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenize sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

8.1.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from clinical samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions.

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- US121 NX-48 Urine/Swab DNA Kit using the Nextactor® NX-48 (Genolution Inc.).

- MagDEA® Dx SV, using magLEAD® 6gC (Precision System Science)).
- NucleoSpin® RNA Virus (MACHEREY-NAGEL).
- NucleoMag® Pathogen (MACHEREY-NAGEL).
- STARMag 96x4 Universal Cartridge Kit (Seegene, Korea)
- Microlab STAR Let automatic extraction system (Hamilton, Switzerland)
- Invisorb® Spin Universal Kit (Stratec).

8.2. Lyophilized positive control

Sexually Transmitted Diseases Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Sexually Transmitted Diseases* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *Sexually Transmitted Diseases* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 1) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Trichomonas vaginalis* and *Chlamydia trachomatis*), HEX, JOE or VIC (*Ureaplasma urealyticum* and *Mycoplasma genitalium*), ROX (*Ureaplasma parvum* and *Neisseria gonorrhoeae*), and Cy5 channels (*Mycoplasma hominis* and Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in negative control well and the presence of signal for *Sexually Transmitted Diseases* Positive Control well. The analysis of the samples is done by the software itself of the used real time PCR equipment according to manufacturer's instructions.

Interpretation of results for *Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium* 4/8-well strips:

- A sample is considered positive for ***Chlamydia trachomatis*** if there is an amplification signal in FAM channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for ***Neisseria gonorrhoeae*** if there is an amplification signal in ROX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for ***Mycoplasma genitalium*** if there is an amplification signal in HEX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.
- In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

Interpretation of results for *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis* 4/8-well strips:

- A sample is considered positive for *Trichomonas vaginalis* if there is an amplification signal in FAM channel and the Ct value obtained is less than 40.
- A sample is considered positive for *Ureaplasma urealyticum* if there is an amplification signal in HEX/VIC/JOE channel and the Ct value obtained is less than 40.
- A sample is considered positive for *Ureaplasma parvum* if there is an amplification signal in ROX channel and the Ct value obtained is less than 40.
- A sample is considered positive for *Mycoplasma hominis* if there is an amplification signal in Cy5 channel and the Ct value obtained is less than 40.
- A sample is considered negative, if the sample shows no amplification signal in the detection system.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium*).



Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*).

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from urogenital specimens and endocervical specimens for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*, rectal specimens (rectal-vaginal specimens and rectal samples) for *Neisseria gonorrhoeae*, *Ureaplasma urealyticum* and *Mycoplasma hominis*; and urine samples for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*.
- ON 0318 only intervenes in the conformity assessment of the test for *Chlamydia trachomatis*. The scope of the CE certification covers the detection of *Chlamydia trachomatis* from urine samples and urogenital and endocervical swabs. The detection of *Chlamydia trachomatis* in pharyngeal, semen, serum and rectal, urethral and vaginal swabs are out of the scope of certification by ON 0318. The rest of the pathogens have self-certification for CE marking.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by contamination by Sexually Transmitted Diseases, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Sexually Transmitted Diseases* Real Time PCR Detection Kit was tested using 948 different specimens first-void urines, rectal swabs, urethral swabs, endocervical swabs/exudate, vaginal swabs and pharyngeal swabs from symptomatic patients. Of all the analyzed samples, 46 discordant samples were obtained. These discrepancies were mainly observed in samples close to the limit of detection.

These results of both Multiplex reactions were compared with those obtained with a molecular detection method (Allplex™ STI Essential Assay (Seegene)). Of the 46 discordant samples, it was possible to solve the discrepancy using FTD Urethritis plus (Fast Track Diagnostics) in 35 samples, the rest could not be re-analysed for lack of sample.

The results were as follows:

VIASURE <i>Sexually Transmitted Diseases</i> Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	19	1	20
	-	1	927	928
	Total	20	928	948

Table 4. Comparative results for *Neisseria gonorrhoeae*.

VIASURE <i>Sexually Transmitted Diseases</i> Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	68	3	71
	-	2	875	877
	Total	70	878	948

Table 5. Comparative results for *Chlamydia trachomatis*.

VIASURE <i>Sexually Transmitted Diseases</i> Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	28	1	29
	-	2	917	918
	Total	30	918	948

Table 6. Comparative results for *Mycoplasma genitalium*.

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	20	3	23
	-	0	925	925
	Total	20	928	948

Table 7. Comparative results for *Trichomonas vaginalis*.

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	166	6	172
	-	7	769	776
	Total	173	775	948

Table 8. Comparative results for *Ureaplasma urealyticum*.

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	402	13	415
	-	5	528	533
	Total	407	541	948

Table 9. Comparative results for *Ureaplasma parvum*.

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit	Allplex™ STI Essential Assay (Seegene)			
		+	-	Total
	+	183	2	185
	-	0	763	763
	Total	183	765	948

Table 10. Comparative results for *Mycoplasma hominis*.

Sensitivity, specificity, PPV and NPV values for VIASURE Sexually Transmitted Diseases Real Time PCR detection kit (CerTest Biotec) compared with Allplex™ STI Essential Assay (Seegene) are shown in next table:

Microorganism	SE (%)	SP (%)	VPP (%)	VPN (%)
<i>C. trachomatis</i>	97.14	99.6	95.77	99.77
<i>M. genitalium</i>	93.33	99.89	96.55	99.78
<i>N. gonorrhoeae</i>	95	99.89	95	99.79
<i>T. vaginalis</i>	100	99.68	86.96	100
<i>U. urealyticum</i>	95.95	99.23	96.51	99.1
<i>U. parvum</i>	98.77	97.6	96.87	99.06
<i>M. hominis</i>	100	99.74	98.92	100

Table 11. Sensitivity, specificity, PPV and NPV values for VIASURE Sexually Transmitted Diseases Real Time PCR detection kit.

In conclusion, the results show a high sensitivity and specificity to detect sexually transmitted diseases using VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction (Figure 3, 4, 5, 6, 7, 8 and 9).

Figure 3. Dilution series of *Trichomonas vaginalis* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*, channel FAM)


Figure 4 Dilution series of *Ureaplasma urealyticum* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*, channel HEX)


Figure 5. Dilution series of *Ureaplasma parvum* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*, channel ROX).


Figure 6. Dilution series of *Mycoplasma hominis* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*, channel Cy5).


Figure 7. Dilution series of *Chlamydia trachomatis* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium*, channel FAM).


Figure 8. Dilution series of *Mycoplasma genitalium* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium*, channel HEX).



Figure 9. Dilution series of *Neisseria gonorrhoeae* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium*, channel ROX).

12.3. Analytical specificity

The specificity of the Sexually Transmitted Diseases assay was confirmed by testing a panel consisting of different microorganisms representing the most common Sexually Transmitted Diseases pathogens. No cross-reactivity was detected against any of the following microorganisms tested:

Cross-reactivity testing					
<i>Candida parapsilosis</i>	-	<i>Enterococcus faecium</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Candida tropicalis</i>	-	<i>Escherichia coli</i> 0.1285; O18:H7:K1	-	<i>Serratia marcescens</i>	-
<i>Candida glabrata</i>	-	<i>Gardnerella vaginalis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
<i>Candida krusei</i>	-	<i>Haemophilus influenza</i> MinnA	-	<i>Stenotrophomonas maltophilia</i>	-
<i>Candida dubliniensis</i>	-	<i>Haemophilus ducreyi</i> class 1	-	<i>Streptococcus agalactiae</i>	-
<i>Candida albicans</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Streptococcus pneumoniae</i>	-
<i>Aspergillus fumigatus</i>	-	<i>Klebsiella pneumoniae</i>	-	<i>Treponema pallidum</i>	-
<i>Bacteroides fragilis</i>	-	<i>Listeria ivanovii</i>	-	<i>Trichomonas vaginalis</i>	-/+
<i>Acinetobacter baumannii</i>	-	<i>Listeria monocytogenes</i>	-	<i>Ureaplasma parvum</i>	-/+
<i>Chlamydia trachomatis</i> (LGV)	-/+	<i>Listeria innocua</i>	-	<i>Ureaplasma urealyticum</i>	-/+
<i>Chlamydia trachomatis</i> (SW)	-/+	<i>Mycoplasma genitalium</i>	-/+	Cytomegalovirus	-
<i>Chlamydia trachomatis</i> (Genovar F)	-/+	<i>Mycoplasma hominis</i>	-/+	Herpes simplex virus 1	-
<i>Enterobacter cloacae</i>	-	<i>Neisseria gonorrhoeae</i>	-/+	Herpes simplex virus 2	-
<i>Enterobacter aerogenes</i>	-	<i>Neisseria meningitidis</i>	-	Human papillomavirus 16	-
<i>Enterococcus faecalis</i>	-	<i>Proteus mirabilis</i>	-	Human papillomavirus 18	-

Table 12. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Neisseria gonorrhoeae* was evaluated against *Neisseria gonorrhoeae* strains St 49226 and Lvl Ng PorA, showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Chlamydia trachomatis* was evaluated against *Chlamydia trachomatis* (SW and LGV), Swedish strain, and Genovar F, and Serovars D, E, I, K and J; showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Mycoplasma genitalium* was evaluated against *Mycoplasma genitalium* showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Trichomonas vaginalis* was evaluated against *Trichomonas vaginalis* showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Ureaplasma urealyticum* was evaluated against *Ureaplasma urealyticum* showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Ureaplasma parvum* was evaluated against *Ureaplasma parvum* showing positive result.

The reactivity of VIASURE Sexually Transmitted Diseases Real Time PCR Detection Kit for *Mycoplasma hominis* was evaluated against *Mycoplasma hominis* showing positive result.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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

13. Bibliography/Bibliografía

1. D. León *et al.* Molecular detection of sexually transmitted agents in a symptomatic group of men and its relationship with sexual behaviour. *Revista Chilena de Infectología* 2016; 33(5):505-512.
2. A. K. Wendel *et al.* Trichomonas vaginalis Polymerase Chain Reaction Compared with Standard Diagnostic. *Clinical Infectious Diseases* 2020; 35:576-80.
3. G. Kusdian *et al.* The biology of Trichomonas vaginalis in the light of urogenital tract infection. *Molecular and Biochemical Parasitology Journal* (2015).
4. T. Posse *et al.* Bacteriemia por Mycoplasma hominis: un agente etiológico subestimado. *Revista Argentina Microbiología*. 2017;02.009.
5. F. Kazumasa *et al.* Rapid and simple detection of Ureaplasma species from vaginal swab samples using a loop-mediated isothermal amplification method. *American Journal of Reproductive Immunology* 2017;e12771.
6. S. Ona *et al.* Mycoplasma genitalium: An Overlooked Sexually Transmitted Pathogen in Women? *Infectious Diseases in Obstetrics and Gynecology* 2016; 2016: 4513089.
7. N. Lai-King *et al.* The laboratory diagnosis of Neisseria gonorrhoeae. *Canadian Journal of Infectious Diseases and Medical Microbiology*. 2005 Jan-Feb; 16(1): 15–25.
8. CDC Grand Rounds: Chlamydia prevention: challenges and strategies for reducing disease burden and sequelae. *Centers for Disease Control and Prevention (CDC)*.
9. G. Fernández *et al.* Usefulness of a novel multiplex real-time PCR assay for the diagnosis of sexually-transmitted infections. *Enfermedades Infecciosas y Microbiología Clínica* 2016.
10. S. Alexander *et al.* Evaluation of strategies for confirming Neisseria gonorrhoeae nucleic acid amplification tests. *Journal of Medical Microbiology* 2011; 60: 909-912.
11. J. Schirm *et al.* Trichomonas vaginalis detection using real-time TaqMan PCR. *Journal of Microbiological Methods* 2007; 68(2): 243-247.
12. J. Yi *et al.* Detection and biovar discrimination of Ureaplasma urealyticum by real-time PCR. *Molecular and Cellular Probes* 2005; 19(4): 255-260.
13. J.S. Jensen *et al.* Use of TaqMan 5' nuclease real-time PCR for quantitative detection of Mycoplasma genitalium DNA in males with and without urethritis who were attendees at a sexually transmitted disease clinic. *Journal of Clinical Microbiology* 2004; 42(2): 683-692.
14. T.A. Halse *et al.* A multiplexed real-time PCR assay for rapid detection of Chlamydia trachomatis and identification of serovar L-2, the major cause of Lymphogranuloma venereum in New York. *Molecular and Cellular Probes* 2006 Oct;20(5):290-7. Epub 2006 Mar 6.
15. D. Férandon *et al.* Development of a real-time PCR targeting the yidC gene for the detection of Mycoplasma hominis and comparison with quantitative culture. *Clinical Microbiology and Infection* 2011; 17(2): 155-159.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: September 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC

*El Organismo Notificado 0318 solo interviene en la identificación de *Chlamydia trachomatis* incluida en NMT/
The Notified Body 0318 is only involved in the identification of *Chlamydia trachomatis* included in NMT.