VIASURE

Real Time PCR Detection Kit

Aspergillus differentiation

CE

IVD

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PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101

Table A 1. References for Open and Rotor-Gene format with internal control products. / Referencias para productos Open and Rotor-Gene Format con control interno.

TUBE FORMAT WITH INTERNAL CONTROL (SEE ANNEX 2) / FORMATO TUBO CON CONTROL INTERNO (VER ANEXO 2)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-ASP196T

Table A 2. References for Tube format with internal control products. / Referencias para productos formato Tubo con control interno.

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL (SEE ANNEX 3) / OPEN AND ROTOR-GENE FORMAT CON CONTROL DE EXTRACCIÓN (VER ANEXO 3)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101E

Table A 3. References for Open and Rotor-Gene format with extraction control products. / Referencias para productos Open and Rotor-Gene Format con control de extracción.

TUBE FORMAT WITH EXTRACTION CONTROL (SEE ANNEX 4) / FORMATO TUBO CON CONTROL DE EXTRACCIÓN (VER ANEXO 4)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-ASP196TE

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EN For download IFUS from other languages, please enter in **certest.es/viasure/labeling**. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use	7
2.	Summary and Explanation	7
3.	Principle of the procedure	8
4.	Reagents provided	8
5.	Reagents and equipment to be supplied by the user	8
6.	Transport and storage conditions	9
7.	Precautions for users	9
8.	Test procedure	11
8.1.	Specimen collection, transport and storage	11
8.2.	DNA extraction	11
9.	Result interpretation	12
9.1.	References with internal control (references in Annex 1 and 2)	12
9.2.	References with extraction control (references in Annex 3 and 4)	14
10.	Limitations of the test	16
11.	Quality control	17
12.	Performance characteristics	17
12.1.	Clinical sensitivity and specificity	17
12.2.	Analytical sensitivity	18
12.3.	Analytical specificity	19
12.4.	Analytical reactivity	20
	ANNEX 1	21
A1.1	Principle of the procedure	21
A1.2	Reagents provided	21
A1.3	Test procedure	22
A1.3.1	Lyophilized positive control	22
A1.3.2	PCR protocol	23
	ANNEX 2	24
A2.1	Principle of the procedure	24
A2.2	Reagents provided	24
A2.3	Test procedure	24
A2.3.1	Lyophilized positive control	24
A2.3.2	Lyophilized reaction mix tube	25
A2.3.3	PCR protocol	25
	ANNEX 3	27
A3.1	Principle of the procedure	27
A3.2	Reagents provided	27
A3.3	Test procedure	29
A3.3.1	Lyophilized extraction control	29
A3.3.2	Lyophilized positive control	29
A3.3.3	PCR protocol	29
	ANNEX 4	31
A4.1	Principle of the procedure	31
A4.2	Reagents provided	31

A4.3 Test procedure	32
A4.3.1 Lyophilized extraction control	32
A4.3.2 Lyophilized positive control.....	32
A4.3.3 Lyophilized reaction mix tube.....	32
A4.3.4 PCR protocol	33

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ENGLISH

1. Intended use

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit is a real time PCR assay designed for the qualitative detection and differentiation of DNA from *Aspergillus fumigatus*, *Aspergillus flavus* and/or *Aspergillus terreus*, from bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS), and sputum clinical samples from individuals suspected of aspergillosis by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of aspergillosis caused by *Aspergillus fumigatus*, *Aspergillus flavus* and/or *Aspergillus terreus*, in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using real time PCR, and detected using fluorescent reporter dye probes specific for *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus*.

2. Summary and Explanation

The genus *Aspergillus* encompasses more than 250 species based on morphology and is one of the largest genera of filamentous fungi causing human diseases. *Aspergillus* species cause a wide spectrum of diseases in humans. Depending on the underlying immune status of the host, *Aspergillus* diseases can be roughly classified into three groups with distinct pathogenic mechanisms, clinical manifestations, and overlapping features, based on the immune response: 1) hypersensitivity-allergy, such as Severe Asthma with Fungal Sensitization (SAFS) or Allergic Bronchial Pulmonary Aspergillosis (ABPA) (in atopic patients); 2) structural diseases with more/less inflammation, such as Aspergilloma or Chronic Pulmonary Aspergillosis (CPA) (in immunocompetent patient); 3) severe immunodeficiency, such as Invasive Pulmonary Aspergillosis (IPA) or Invasive Bronchial Aspergillosis (IBA) (in immunocompromised patients). *Aspergillus* has emerged as one of the most common causes of infectious death in severely immunocompromised patients, with mortality rates up to 40%-50% in patients with acute leukemia and recipients of hematopoietic stem cells transplantation (HSCTs).

There are clinical classifications of aspergillosis: allergic (extrinsic asthma, extrinsic allergic alveolitis, and allergic bronchopulmonary aspergillosis), saprophytic colonizing (pulmonary, extrapulmonary), and invasive (pulmonary and extrapulmonary). Most human infections arise through inhalation of fungal spores, which are present in indoor and outdoor environments. Accordingly, isolation of *Aspergillus* spp. from cultures of the respiratory tract of asymptomatic patients with no evidence of invasive or allergic disease is common. Even though it is not always associated with infection, *Aspergillus* colonization is associated with a substantial increase in risk for development of invasive infection across a wide range of immunocompromised individuals.

Aspergillus fumigatus is a saprotrophic fungus which has been considered the most ubiquitous fungal species in the environment. It has been described as responsible for the one of the most devastating lung infections in terms of morbidity and mortality. *A. flavus* is the second-leading pathogen causing invasive and noninvasive aspergillosis, next to *Aspergillus fumigatus*. Infection due to *A. flavus* is predominant in Asia, the Middle East and Africa, possible due to its better ability to survive in hot and arid climatic conditions compared to other *Aspergillus* spp. In contrast, *A. terreus* is an emerging opportunistic fungus associated with disseminated infections and high mortality rates when treated with polyenes.

The most common diagnostic methods are its isolation from culture combined with histopathologic characteristics or radiography finding compatible with disease. Biomarkers such as 1,3- β -D-glucan and

galactomannan are also used. However, these have some limitations, showing low sensitivity in certain population groups, such as patients treated with antifungals. Molecular diagnostic based on qPCR is proposed as a future alternative for the detection of this pathogen in a more sensitive and specific way.

3. Principle of the procedure

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit is designed for the qualitative detection and differentiation of DNA from *Aspergillus fumigatus*, *Aspergillus flavus* and/or *Aspergillus terreus*, from bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS) and sputum clinical samples. After DNA isolation, the identification of *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* is performed by the amplification of a conserved region of the *ITS1* region of *A. flavus*, the *ITS1* region of *A. terreus*, and the *23S rRNA* gene of *A. fumigatus*, using specific primers and a fluorescent-labelled probes.

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal, which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to verify the correct functioning of the amplification mix.

4. Reagents provided

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products, and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).
- Lysis buffer and Proteinase K.

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit has been tested on the following equipments: 7500 Fast Real-Time PCR Systems (Applied Biosystems), CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad), AriaMx Real-Time PCR System (Agilent Technologies), DTprime Real-time Detection Thermal Cycler (DNA-Technology), DTLite Real-Time PCR System (DNA-Technology), Rotor-Gene® Q (Qiagen), Cobas z480 Analyzer (Roche Molecular Diagnostics), LightCycler® 480 Instrument II (Roche Molecular Diagnostics), QuantGene 9600 (Bioer Technology) and VIASURE V-Lab96 Cycler (CerTest Biotec S.L.). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTLite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *Aspergillus differentiation* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.

-
- Do not remove desiccant from reagent pouches once is open.
 - Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-ASP113L, VS-ASP113H, VS-ASP136, VS-ASP172, VS-ASP113LE, VS-ASP113HE, VS-ASP136E and VS-ASP172E). Remove any air excess in the pouches prior to closing.
 - Do not use reagents if the foil has been broken or damaged.
 - Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
 - Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
 - For references VS-ASP101, VS-ASP136, VS-ASP172, VS-ASP101E, VS-ASP136E and VS-ASP172E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
 - An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
 - Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
 - In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
 - Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
 - Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
 - Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
 - Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
 - In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
 - Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please see Annex 1 for open and rotor-gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for open and rotor-gene format with extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit has been tested in respiratory samples (bronchoalveolar lavages -BAL-, bronchoalveolar aspirates -BAS-, and sputum). Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, all specimens should be collected and labelled appropriately in clean containers. After collection, specimens should be placed in a biohazard bag and should be transported and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at Room Temperature (RT) for up to 2 hours, or at 2°-8°C for up to 5 days, following the local and national regulations for the transport of pathogen material. For long term transport (more than 5 days), it is recommend shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 2°-8°C for up to 5 days, or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94).

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

Due to VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit is also available with extraction control (references include in Annexes 3 and 4), if the Extraction Control (EC) is used to monitor nucleic acid isolation and as PCR inhibition control, add 5µl of the reconstituted EC to the specimen and/or lysis buffer mixture (clinical specimen, as well as positive control and/or negative control). Close each tube and vortex for 10 seconds. If the Extraction Control is used only as a PCR inhibition control, 1µl of the reconstituted EC should be added to the reconstituted Reaction-Mix.

For DNA extraction from clinical samples, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

- Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).

NOTE: One of the major hurdles of using PCR to diagnose aspergillosis is extraction of the DNA. The presence of a complex, sturdy cell wall that is resistant to lysis impairs extraction of the DNA by conventional methods employed for bacteria. Several fungal DNA extraction protocols have been described in the literature (Griffiths *et al.*). To improve *Aspergillus* detection, a pre-treatment step consisting of a previous 10 minutes' incubation at 56 °C in Lysis Buffer and proteinase K mix showed stress importance.

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>A. flavus</i> (FAM) ¹	<i>A. fumigatus</i> (HEX) ¹	<i>A. terreus</i> (ROX) ¹	Internal Control (Cy5) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	> 40 or no signal	> 40 or no signal	> 40 or no signal	≤40	Valid

Table 1. Expected Performance of Controls. Ct values no signal = no amplification curve.

1. In cases where either or both control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid', and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40) in PC and NC control wells.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

<i>A. flavus</i> (FAM)	<i>A. fumigatus</i> (HEX)	<i>A. terreus</i> (ROX)	Internal Control (Cy5)	Interpretation for patients' individual samples	
≤40	> 40 or no signal	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. flavus</i> DNA Detected
> 40 or no signal	≤40	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. fumigatus</i> DNA Detected
> 40 or no signal	> 40 or no signal	≤40	≤40 or no signal ¹	Valid	<i>A. terreus</i> DNA Detected
≤40	≤40	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. flavus</i> and <i>A. fumigatus</i> DNA Detected
≤40	> 40 or no signal	≤40	≤40 or no signal ¹	Valid	<i>A. flavus</i> and <i>A. terreus</i> DNA Detected
> 40 or no signal	≤40	≤40	≤40 or no signal ¹	Valid	<i>A. fumigatus</i> and <i>A. terreus</i> DNA Detected
≤40	≤40	≤40	≤40 or no signal ¹	Valid	All targets' DNA Detected
> 40 or no signal	> 40 or no signal	> 40 or no signal	≤35 ²	Valid	Targets not Detected
> 40 or no signal	> 40 or no signal	> 40 or no signal	> 35 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 2. Interpretation of individual patient sample results. Ct values no signal = no amplification curve.


1 The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* target genes negative, IC must show an amplification signal with Ct ≤ 35. If there is an absence of signal or Ct value > 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- repeat qPCR with the same isolated DNA sample, or
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

9.2. References with extraction control (references in Annex 3 and 4)

All the results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>A. flavus</i> (FAM) ¹	<i>A. fumigatus</i> (HEX) ¹	<i>A. terreus</i> (ROX) ¹	Extraction Control (Cy5) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	> 40 or no signal	> 40 or no signal	> 40 or no signal	≤40	Valid

Table 3. Expected Performance of Controls. Ct values no signal = no amplification curve.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid', and retesting is required.

2. The Extraction Control (EC) should show an amplification signal (Ct ≤40) in PC and NC control wells.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

<i>A. flavus</i> (FAM)	<i>A. fumigatus</i> (HEX)	<i>A. terreus</i> (ROX)	Extraction Control (Cy5)	Interpretation for patients' individual samples	
≤40	> 40 or no signal	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. flavus</i> DNA Detected
> 40 or no signal	≤40	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. fumigatus</i> DNA Detected
> 40 or no signal	> 40 or no signal	≤40	≤40 or no signal ¹	Valid	<i>A. terreus</i> DNA Detected
≤40	≤40	> 40 or no signal	≤40 or no signal ¹	Valid	<i>A. flavus</i> and <i>A. fumigatus</i> DNA Detected
≤40	> 40 or no signal	≤40	≤40 or no signal ¹	Valid	<i>A. flavus</i> and <i>A. terreus</i> DNA Detected
> 40 or no signal	≤40	≤40	≤40 or no signal ¹	Valid	<i>A. fumigatus</i> and <i>A. terreus</i> DNA Detected
≤40	≤40	≤40	≤40 or no signal ¹	Valid	All targets' DNA Detected
> 40 or no signal	> 40 or no signal	> 40 or no signal	≤35 ²	Valid	Targets not Detected
> 40 or no signal	> 40 or no signal	> 40 or no signal	> 35 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 4. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.


1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction control between the controls and the clinical samples, due to the extraction process.

2 In the case of *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* target genes negative, EC must show an amplification signal with Ct ≤ 35. If there is an absence of signal or Ct value > 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- repeat qPCR with the same isolated DNA sample, or
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

Figure 2. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from respiratory samples (bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS), and sputum).
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus*, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Aspergillus differentiation* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction). Please note that a pre-treatment step consisting of a previous 10 minutes' incubation at 56 °C in Lysis Buffer and proteinase K mix showed stress importance for a correct *Aspergillus* DNA extraction.
 - Degradation of the DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus*.
 - A fungal load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics, or immunosuppressant drugs, used to prevent *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* infections or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable fungi and does not imply that these fungi are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* target sequences.
- Negative results do not preclude *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak fungal levels during infections caused by *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the fungi.

- If diagnostic tests for other respiratory diseases are negative and the patient's clinical presentation and epidemiological information suggest that *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit was tested using respiratory samples (bronchoalveolar lavages -BAL-, bronchoalveolar aspirates -BAS-, and sputum), from donors with signs and symptoms of lower respiratory infection compatible with *Aspergillus* spp.. In order to determine the clinical diagnostic accuracy, an evaluation has been conducted through collaboration with a national entity. A summary of the sites, sample type and workflow is included in the following table:

Site	Sample type	Workflow	Target
1 CerTest Biotec, in collaboration with the Biobanco del Sistema de Salud de Aragón, (Zaragoza, Spain)	BAL	MagDEA® Dx SV Kit using the magLEAD® 12gC instrument (Precision System Science Co.) + Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit using the Maxwell® RSC 16 Instrument (Promega) + LightCycler® 480 Instrument II (Roche Molecular Diagnostics)	<i>A. fumigatus</i>
			<i>A. terreus</i>
	BAS	MagDEA® Dx SV Kit using the magLEAD® 12gC instrument (Precision System Science Co.) + Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit using the Maxwell® RSC 16 Instrument (Promega) + LightCycler® 480 Instrument II (Roche Molecular Diagnostics)	<i>A. flavus</i>
			<i>A. fumigatus</i>
			<i>A. terreus</i>
	Sputum	MagDEA® Dx SV Kit using the magLEAD® 12gC instrument (Precision System Science Co.) + Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit using the Maxwell® RSC 16 Instrument (Promega) + LightCycler® 480 Instrument II (Roche Molecular Diagnostics)	<i>A. fumigatus</i>
			<i>A. terreus</i>

Table 5. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	Sample type	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Culture + MALDI-TOF	<i>A. fumigatus</i>	BAL	14	123	0	0	1(0.76–1)	1 (0.97 – 1)	1(0.76–1)	1 (0.97 – 1)
		<i>A. terreus</i>		2	135	0	0	1 (0.15-1)	1 (0.97 – 1)	1 (0.15-1)	1 (0.97 – 1)
		<i>A. flavus</i>	BAS	3	680	0	0	1(0.54-1)	1 (0.99 – 1)	1(0.54-1)	1 (0.99 – 1)
		<i>A. fumigatus</i>		23	297	0	0	1 (0.85 – 1)	1 (0.99 – 1)	1 (0.85 – 1)	1 (0.99 – 1)
		<i>A. terreus</i>		3	316	0	1	0.75 (0.2 – 1)	1 (0.98 – 1)	1 (0.3 – 1)	0.99 (0.97 – 0.99)
		<i>A. fumigatus</i>	Sputum	4	219	0	0	1(0.39–1)	1 (0.98 – 1)	1(0.39–1)	1 (0.98 – 1)
		<i>A. terreus</i>		4	219	0	0	1 (0.4– 1)	1 (0.98 – 1)	1 (0.4– 1)	1 (0.98 – 1)
		<i>A. fumigatus</i>	Global	41	639	0	0	1(0.91 – 1)	1 (0.99 – 1)	1(0.91 – 1)	1 (0.99 – 1)
		<i>A. terreus</i>		9	670	0	1	0.9 (0.55– 0.99)	1 (0.99 – 1)	1(0.62 – 0.98)	0.99 (0.99 – 0.99)

Table 6. True positive (TP) and true negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV) and Predictive Negative Values (NPV) for VIASURE *Aspergillus* differentiation Real Time PCR Detection Kit.

Results show high agreement to detect *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* using VIASURE *Aspergillus* differentiation Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Aspergillus* differentiation Real Time PCR Detection Kit has a detection limit of 4×10^{-3} CFU per reaction for *Aspergillus flavus*, 4×10^{-2} CFU per reaction for *Aspergillus fumigatus*, and 8×10^{-2} CFU per reaction for *Aspergillus terreus*, with a positive rate of $\geq 95\%$, on serum samples, and 8×10^{-2} CFU per reaction for *Aspergillus flavus*, 0.1 CFU per reaction for *Aspergillus fumigatus* and 4×10^{-2} CFU per reaction for *Aspergillus terreus*, with a positive rate of $\geq 95\%$, on BAL samples.

Figure 3. Dilution series of *Aspergillus flavus* (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).


Figure 4. Dilution series of *Aspergillus fumigatus* (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel HEX).



Figure 5. Dilution series of *Aspergillus terreus* (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel ROX).

12.3. Analytical specificity

The specificity of the VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit assay was confirmed by testing a panel of different microorganisms associated to respiratory diseases. No cross-reactivity of the VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with any of the following microorganisms tested was observed.

Cross-reactivity testing (wet testing)					
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	<i>Haemophilus influenzae</i> MinnA	-	<i>Legionella bozemanii</i> serovar 1	-
<i>Aspergillus niger</i>	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Legionella dumoffii</i>	-
Bocavirus	-	Influenza A/California/7/2009(H1N1)pdm09 virus	-	<i>Legionella longbeachae</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella micdadei</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i> serogroup 1	-
<i>Bordetella parapertussis</i>	-	Influenza A/Victoria/210/2009 (H3N2) virus	-	Human metapneumovirus A1, A2 and B	-
<i>Bordetella pertussis</i>	-	Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	<i>Moraxella catarrhalis</i>	-
<i>Chlamydia caviae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycoplasma pneumoniae</i>	-

Cross-reactivity testing (wet testing)					
<i>Chlamydia psittaci</i> genotypes A and C	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> not rifampin resistant	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
Human coronavirus 229E, OC43, NL63	-	Influenza A/DE-SH/Reihenente/AR8444/ 2016 (H5N8) virus	-	<i>Pneumocystis jirovecii</i> Type A1 and g885652	-
MERS-CoV	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	Human rhinovirus type C	-
Enterovirus 68 and 71	-	Influenza B/Brisbane/60/2008 virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Enterovirus Echovirus 30	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> Z022	-
Enterovirus Coxsackievirus A24, A9 and B3	-	Influenza B/Phuket/3073/2013 virus	-	Respiratory syncytial virus (RSV) A and B	-

Table 7. Reference pathogenic microorganisms used in this study performed at CerTest facility.

12.4. Analytical reactivity

The reactivity of VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit was evaluated against DNA extracted from cultures of *Aspergillus flavus* Link 1809 strain MUCL 19006, *Aspergillus fumigatus* Fresenius 1863 strain DSM 819 and *Aspergillus terreus* var. *terreus* Thom 1918 strain DSM 1958 (as templates), showing positive results.

ANNEX 1

OPEN AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCES
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Aspergillus flavus</i>	FAM	<i>ITS1</i> region
<i>Aspergillus fumigatus</i>	HEX, VIC or JOE *	<i>23S rRNA</i> gene
<i>Aspergillus terreus</i>	ROX	<i>ITS1</i> region
Internal Control (IC)	Cy5	-

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.ceratest.es.

A1.2 Reagents provided

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.ceratest.es).

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP101L, VS-ASP101H, VS-ASP106L, VS-ASP106H, VS-ASP112L and VS-ASP112H.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref VS-ASP113L and VS-ASP113H.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP101, VS-ASP136 and VS-ASP172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Aspergillus differentiation Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized

Aspergillus differentiation Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Aspergillus differentiation* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the strips or plate in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*A. flavus*), HEX, JOE or VIC (*A. fumigatus*), ROX (*A. terreus*) and Cy5 (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-ASP196T

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Aspergillus flavus</i>	FAM	<i>ITS1</i> region
<i>Aspergillus fumigatus</i>	HEX, VIC or JOE *	23S rRNA gene
<i>Aspergillus terreus</i>	ROX	<i>ITS1</i> region
Internal Control (IC)	Cy5	-

Table A2. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.certest.es.

A2.2 Reagents provided

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

Aspergillus differentiation Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized

Aspergillus differentiation Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Aspergillus differentiation* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Aspergillus differentiation* Reaction-Mix (white vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Aspergillus differentiation* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.cerTEST.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*A. flavus*), HEX, JOE or VIC (*A. fumigatus*), ROX (*A. terreus*) and Cy5 (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 3

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Aspergillus flavus</i>	FAM	<i>ITS1</i> region
<i>Aspergillus fumigatus</i>	HEX, VIC or JOE *	23S rRNA gene
<i>Aspergillus terreus</i>	ROX	<i>ITS1</i> region
Extraction Control (EC)	Cy5	-

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.ceratest.es.

A3.2 Reagents provided

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3, A3.4 and A.3.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A3.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A3.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.ceratest.es).

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A3. 3. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP101LE, VS-ASP101HE, VS-ASP106LE, VS-ASP106HE, VS-ASP112LE and VS-ASP112HE.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A3. 4. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref VS-ASP113LE and VS-ASP113HE.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A3. 5. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP101E, VS-ASP136E and VS-ASP172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) by adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *Aspergillus differentiation* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

Aspergillus differentiation Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Aspergillus differentiation* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Aspergillus differentiation* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the strips or plate in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A3. 6. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*A. flavus*), HEX, JOE or VIC (*A. fumigatus*), ROX (*A. terreus*) and Cy5 (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-ASP196TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Aspergillus flavus</i>	FAM	<i>ITS1</i> region
<i>Aspergillus fumigatus</i>	HEX, VIC or JOE *	<i>23S rRNA</i> gene
<i>Aspergillus terreus</i>	ROX	<i>ITS1</i> region
Extraction Control (EC)	Cy5	-

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A4.2 Reagents provided

VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
<i>Aspergillus differentiation</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Aspergillus differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE *Aspergillus differentiation* Real Time PCR Detection Kit with Ref. VS-ASP196TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 μ L of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *Aspergillus differentiation* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

Aspergillus differentiation Positive Control contains high copies of the template; the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Aspergillus differentiation* Positive Control (red vial) by adding 100 μ L of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Aspergillus differentiation* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 μ L of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C $\pm$ 5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Aspergillus differentiation* Reaction-Mix (white vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Aspergillus differentiation* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*A. flavus*), HEX, JOE or VIC (*A. fumigatus*), ROX (*A. terreus*) and Cy5 (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

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Lugar	Kit comparador	Diana	Tipo muestra	TP	TN	FP	FN	Sensibilidad	Especificidad	PPV	NPV
1	Cultivo + MALDI-TOF	<i>A. fumigatus</i>	BAL	14	123	0	0	1 (0.76–1)	1 (0.97–1)	1 (0.76–1)	1 (0.97–1)
		<i>A. terreus</i>		2	135	0	0	1 (0.15–1)	1 (0.97–1)	1 (0.15–1)	1 (0.97–1)
		<i>A. flavus</i>		3	680	0	0	1 (0.54–1)	1 (0.99–1)	1 (0.54–1)	1 (0.99–1)
		<i>A. fumigatus</i>	BAS	23	297	0	0	1 (0.85–1)	1 (0.99–1)	1 (0.85–1)	1 (0.99–1)
		<i>A. terreus</i>		3	316	0	1	0.75 (0.2–1)	1 (0.98–1)	1 (0.3–1)	0.99 (0.97–0.99)
		<i>A. fumigatus</i>	Espudo	4	219	0	0	1 (0.39–1)	1 (0.98–1)	1 (0.39–1)	1 (0.98–1)
		<i>A. terreus</i>		4	219	0	0	1 (0.4–1)	1 (0.98–1)	1 (0.4–1)	1 (0.98–1)
		<i>A. fumigatus</i>	Global	41	639	0	0	1 (0.91–1)	1 (0.99–1)	1 (0.91–1)	1 (0.99–1)
		<i>A. terreus</i>		9	670	0	1	0.9 (0.55–0.99)	1 (0.99–1)	1 (0.62–0.98)	0.99 (0.99–0.99)

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PRODUCTO	REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101L
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101H
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101

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PRODUCTO	REFERENCIAS
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-ASP106LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-ASP106HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-ASP112LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-ASP112HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-ASP113LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-ASP113HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-ASP136E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-ASP172E
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-ASP101LE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-ASP101HE
VIASURE <i>Aspergillus differentiation</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-ASP101E

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.