CHORUS VITAMIN EXTRACTION

REF 86500

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INSTRUCTIONS FOR USE

CHORUS VITAMIN EXTRACTION

Kit for extraction of samples to be tested with the CHORUS VIT B12 and CHORUS FOLATE kits

For *in vitro* diagnostic use only

1. INTENDED USE

Kit for the extraction of samples to be tested with the CHORUS VIT B12 (REF. 86400) and CHORUS FOLATE (REF. 86402) tests.

It must only be used by professional laboratory personnel.

2. INTRODUCTION

The kit allows the correct treatment for the extraction of the samples needed to perform the CHORUS VIT B12 and CHORUS FOLATE tests.

3. PRINCIPLE OF THE METHOD

The kit consists of three reagents to be used for the extraction of the samples on which to dose vitamin B12 and/or folate.

The extraction process allows these vitamins to free themselves from the carrier proteins to which they are normally bound for their absorption and functionality so that they can be correctly detected.

After performing the sample extraction procedure, it is ready to be dosed.

For analysis of the extracted sample, refer to the instructions for use of the CHORUS VIT B12 (REF. 86400) and/or CHORUS FOLATE (REF. 86402) test.

4. PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY.

Wait for the serum samples to reach room temperature (18-30°C) before performing the extraction procedure.

The used, deteriorated or expired product must be treated as potentially infected residue, therefore, disposed of in accordance with the provisions of current laws.

Personnel safety warnings

- Any material of human origin must be considered potentially infected. All samples must be handled in accordance with the safety rules normally adopted in the laboratory.

- Use single-use gloves and eye protection when handling samples.
- Wash your hands thoroughly after use.
- Refer to the Safety Data Sheet (www.diesse.it) for the safety characteristics of the reagents contained in the kit.
- Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried. All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste. Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

- Collect the sample taking care not to contaminate the product.
- Do not use the reagents after the expiry date.
- Take all the necessary precautions to avoid contamination of the product during its use.

5. KIT COMPOSITION

The package is sufficient for 72 extractions.

RELEASING AGENT

RELEASING BUFFER

1 x 6,0 ml

Content: Solution containing sodium hydroxide (<1.5%) and potassium cyanide (<0.03%).

Liquid, to be mixed with the stabilizing buffer.

STABILIZING AGENT

STABILIZING BUFFER

1 x 0,220 ml

Content: Solution containing tris-2 carboxyethyl phosphine (TCEP) (<15%).

Liquid, to be mixed with the releasing buffer.

NEUTRALIZING AGENT

NEUTRALIZING BUFFER

1 x 5,0 ml

Content: Monobasic sodium phosphate salt buffer (<8.5%) able to reduce the pH of the extracted sample, with red dye (0.1%).

Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- Vortex
- Test tubes in which to do sample extraction
- Test tube shaker
- Centrifuge
- Micropipettes capable of accurately withdrawing volumes of 1-1000 µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C.

The expiry date is indicated on the outer label of the package.

The reagents have limited stability after opening:

RELEASING AGENT	8 weeks
STABILIZING AGENT	8 weeks
NEUTRALIZING AGENT	8 weeks

7. SAMPLE TYPE AND STORAGE

The sample is serum obtained from blood collected by normal venous blood sampling and handled with the precautions required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C; for longer storage periods, freeze at temperatures < -20°C.

The sample can be thawed for a maximum of 3 times. Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

The extracted sample can be used within 8 hours from extraction if kept at 4°C.

It is strongly recommended to avoid any intensive sample processing.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

1. Dispense 100 µl of the sample to be extracted into a test tube (not provided).
2. Mix the stabilizing buffer (Stabilizing Agent) and the releasing buffer (Releasing Agent) to prepare a 1/40 extracting solution (Stabilizing/Releasing).
For example, for the extraction of 6 samples, mix 10 µl of stabilizing buffer and 390 µl of releasing buffer.
The 1/40 solution is stable at room temperature for 8 hours after preparation.
3. Dispense the previously prepared extraction reagent into the 50 µl test tube.
4. Using the vortex, mix 3 times for 1-2 seconds **immediately after adding the extraction reagent.**
5. Incubate for 15 minutes at room temperature.
6. Dispense the neutralizing buffer into the 50 µl test tube.
7. Using the vortex, mix 3 times for 1-2 seconds **immediately after adding the neutralizing buffer.**
8. Use the extracted sample as indicated in the IFU of the kit to be used.

Should you have any problems, contact Scientific Support:

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

9. INCIDENT REPORTING

If a serious incident has occurred in relation to this device in the European Union market territory, please report it to the manufacturer and the competent authority of your Member State without delay.

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diesse.it

FOLATE & VITAMIN B12 SAMPLE PREPARATION

PURPOSE:

The following illustrated steps were created to help train technicians and simplify the process of sample preparation

SAMPLES:

Prepare all patients, calibrators and controls before use in the assay

SAMPLE STABILITY:

Fresh serum can be kept up to 5 days at 2-8°C; for longer periods, please freeze at -20°C or lower. Extracted sample can be used within 8 hours since the extraction, when maintained at 4°C.

EXTRACTION STEPS

STEP 1

- Add 100µl of sample (patient, calibrator or control) into a tube
- Hint: touch side of tube with pipette while dispensing at an angle near the bottom of the tube
- Hint: prepare all samples

STEP 2

- Preparation of Extraction Reagent: mix the Stabilizing Agent and the Releasing Agent to prepare a 1/40 solution (Stabilizing/Releasing).
Ex. For 6 samples, mix 10µl of Stabilizing Agent and 390µl of Releasing Agent
- Add 50µl of Extraction Reagent into sample tubes
- Hint: touch side of tube with pipette while dispensing at an angle slightly above the sample (do not touch the sample to avoid contamination)
- Mix on vortex 1-2 seconds using multiple touch technique (3 x 1-2 seconds), immediately after the addition of the Extraction Reagent
- Hint: mix before moving onto the next tube

STEP 3

- Incubate 15 minutes at room temperature

STEP 4

- Add 50µl of Neutralizing Agent into sample tubes
- Hint: touch side of tube with pipette while dispensing at an angle slightly above the sample (do not touch the sample to avoid contamination)
- Mix on vortex 1-2 seconds using multiple touch technique (3 x 1-2 seconds), immediately after adding the Neutralizing Agent
- Hint: mix before moving onto the next tube

STEP 5

- Use prepared samples accordingly to IFU procedure

diesse.it

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CS	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricației
	EL	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusque
	CS	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi până la
	EL	Ημερομηνία λήξης		
	IT	Nu reutilizați	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CS	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	EL	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CS	Pozor, čtěte přiložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atenție, consultați documentele însoțitoare
	EL	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CS	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Producător
	EL	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CS	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Conținut suficient pt <n> teste
	EL	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CS	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limită de temperatură
	EL	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CS	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru utilizare consultați instrucțiunile
	EL	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CS	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risc biologic
	EL	Βιολογικοί κίνδυνοι		
	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CS	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Număr de catalog
	EL	Αριθμός καταλόγου		
	IT	Dispositiv medical pentru diagnosticare in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CS	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dispositiv medical pentru diagnosticare in vitro
	EL	In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CS	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	EL	Αριθμός Παρτίδας		
	IT	Marcajul de conformitate CE	ES	Marcado CE de conformidad
	EN	CE marking of conformity	FR	Marquage de conformité CE
	CS	Oznakowanie zgodności CE	PT	Marcação CE de conformidade
	DE	CE-Konformität Skennzeichnung	RO	Marcajul de conformitate CE
	EL	Σημάνση συμμορφωσής CE		