Coagulation Control Plasmas

- REF 5186
- REF 5187
- REF 5183
- REF 5482
- Routine Control N
- Routine Control A
- Routine Control SA
- Routine Coagulation Control Set

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0482P 2014/02 (13)

Coagulation Control Plasmas

Instructions for use

en

INTENDED USE

Routine Control N, Routine Control A and Routine Control SA are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assayed for Fibrinogen, TCT and ATIII, and are prepared from normal human plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

REF	Component	Content	Description
5186	Routine Control - N	10 x 1 mL	Prepared from pooled normal plasma.
5187	Routine Control - A	10 x 1 mL	Prepared from adsorbed human plasma.
5183	Routine Control - SA	10 x 1 mL	Prepared from adsorbed human plasma.
5482	Routine Coagulation Control Set:		
	Routine Control - N	4 x 1 mL	
	Routine Control - A	3 x 1 mL	
	Routine Control - SA	3 x 1 mL	

Each kit contains instructions for use.
Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Coagulation Control Plasmas may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with all suitable, commercial reagents.

STORAGE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 8 hours when kept at 2 – 8°C or 4 weeks at -20°C when flash frozen. Keep covered.

SAMPLE COLLECTION AND PREPARATION

Not applicable.

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS

Routine Control N should give values within the laboratory normal range for PT, aPTT and fibrinogen assays. Routine Control A and

Routine Control SA have been standardised to give prolonged and markedly prolonged PT and aPTT times respectively. Lot and instrument specific expected values are provided with each pack of controls.

LIMITATIONS

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{1,2,3}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility			
Intra-assay precision			
Sample	n	aPTT CV (%)	PT CV (%)
Routine Control N	5	2.83	1.01
Routine Control A	5	2.76	1.71
Routine Control SA	5	1.72	1.03

REFERENCES

- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.

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