# VIASURE

## Real Time PCR Detection Kits

by CerTest  
BIOTEC

### Yellow Fever

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## ENGLISH

### 1. Intended use

VIASURE Yellow Fever Real Time PCR Detection Kit is designed for specific identification of Yellow fever in clinical samples from patients with signs and symptoms of Yellow virus infection. This test is intended for use as an aid in the diagnosis of the Yellow fever in combination with clinical and epidemiological risk factors. RNA is extracted from specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for Yellow fever.

### 2. Summary and Explanation

Yellow fever virus is an RNA virus that belongs to the *Flaviviridae* family. The "yellow" in the name refers to the jaundice caused by liver involvement that affects some patients.

Yellow fever is an acute viral haemorrhagic zoonotic disease transmitted between humans and from monkeys to humans by mosquitoes, primarily by *Aedes (aegypti)* and *Haemogogus*, in 1 of 3 cycles: urban yellow fever, jungle (sylvatic) yellow fever, and intermediate (savannah) yellow fever. More than 900 million inhabitants and travellers are at risk in the 44 endemic countries of Latin America and Africa.

The incubation period is generally 3–6 days from infection until illness. Many people do not experience symptoms, but when these do occur, the most common are an acute febrile phase occurs with myalgia, headache, back pain, anorexia, nausea, and sometimes vomiting; symptoms typically resolve within 1 week. Small proportions of patients enter a second, more toxic phase within 24 hours of recovering from initial symptoms, and develop severe symptoms (with high fever, jaundice, bleeding, and kidney damage). Half of the patients who enter this toxic phase die; due to there is currently no specific anti-viral drug. Although, specific care to treat dehydration, liver and kidney failure, and fever improves outcomes. Therefore, vaccination and mosquito control are the most important means of preventing yellow fever.

Yellow fever is difficult to diagnose, especially during the early stages. It can be confused with severe malaria, leptospirosis, viral hepatitis (especially the fulminating forms of hepatitis B and D), other hemorrhagic fevers, infection with other flaviviruses (Dengue, Zika and West Nile viruses), and poisoning. Usually, the detection of yellow fever virus is carried out by virus cultivation on susceptible cells followed by immunofluorescence, ELISA, immunohistochemical examination or PCR. For quantifying the load of infectious particles in cell culture or serum samples, plaque assay is still the commonly used method. However, this is cumbersome and time-consuming technique. With the recent development of Real-Time PCR, a rapid and accurate alternative to quantify viral load in blood is available.

### 3. Principle of the procedure

VIASURE Yellow Fever Real Time PCR Detection Kit is designed for the diagnosis of the Yellow fever in clinical samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed

generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of the 5' UTR gene using specific primers and a fluorescent-labelled probe.

VIASURE Yellow Fever Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE Yellow Fever Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase in an stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity. Yellow fever RNA targets are amplified and detected in FAM channel and the Extraction control (EC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

#### 4. Reagents provided

VIASURE Yellow Fever Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

| Reagent/Material              | Description                                                                          | Colour      | Amount              |
|-------------------------------|--------------------------------------------------------------------------------------|-------------|---------------------|
| Yellow Fever 8-well strips    | A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format | White       | 6/12 x 8-well strip |
| Rehydration Buffer            | Solution to reconstitute the stabilized product                                      | Blue        | 1 vial x 1.8 mL     |
| Yellow Fever Positive Control | Non-infectious synthetic lyophilized cDNA                                            | Red         | 1 vial              |
| Extraction Control            | Non-infectious nucleic acid lyophilized                                              | Green       | 1 vial              |
| Negative control              | Non template control                                                                 | Violet      | 1 vial x 1 mL       |
| Water RNase/DNase free        | RNase/DNase free water                                                               | White       | 1 vial x 1 mL       |
| Tear-off 8-cap strips         | Optical caps for sealing wells during thermal cycling                                | Transparent | 6/12 X 8-cap strip  |

Table 1. Reagents and materials provided in VIASURE Yellow Fever Real Time PCR Detection Kit with Ref. VS-FEV106LE, VS-FEV106HE, VS-FEV112LE and VS-FEV112HE.

| Reagent/Material              | Description                                                                          | Color       | Amount           |
|-------------------------------|--------------------------------------------------------------------------------------|-------------|------------------|
| Yellow Fever 96-well plate    | A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format | White       | 1 plate          |
| Rehydration Buffer            | Solution to reconstitute the stabilized product                                      | Blue        | 1 vial x 1.8 mL  |
| Yellow Fever Positive Control | Non-infectious synthetic lyophilized cDNA                                            | Red         | 1 vial           |
| Extraction Control            | Non-infectious nucleic acid lyophilized                                              | Green       | 1 vial           |
| Negative control              | Non template control                                                                 | Violet      | 1 vial x 1 mL    |
| Water RNase/DNase free        | RNase/DNase free water                                                               | White       | 1 vial x 1 mL    |
| Tear-off 8-cap strips         | Optical caps for sealing plate during thermal cycling                                | Transparent | 12 X 8-cap strip |

Table 2. Reagents and materials provided in VIASURE Yellow Fever Real Time PCR Detection Kit with Ref VS-FEV113LE and VS-FEV113HE.

| Reagent/Material              | Description                                                                          | Colour      | Amount              |
|-------------------------------|--------------------------------------------------------------------------------------|-------------|---------------------|
| Yellow Fever 4-well strips    | A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format | Transparent | 9/18 x 4-well strip |
| Rehydration Buffer            | Solution to reconstitute the stabilized product                                      | Blue        | 1 vial x 1.8 mL     |
| Yellow Fever Positive Control | Non-infectious synthetic lyophilized cDNA                                            | Red         | 1 vial              |
| Extraction Control            | Non-infectious nucleic acid lyophilized                                              | Green       | 1 vial              |
| Negative control              | Non template control                                                                 | Violet      | 1 vial x 1 mL       |
| Water RNase/DNase free        | RNase/DNase free water                                                               | White       | 1 vial x 1 mL       |
| 4-cap strips                  | Optical caps for sealing wells during thermal cycling                                | Transparent | 9/18 X 4-cap strip  |

Table 3. Reagents and materials provided in VIASURE Yellow Fever Real Time PCR Detection Kit with Ref. VS-FEV136E and VS-FEV172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Yellow Fever Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE Yellow Fever Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-FEV113LE, VS-FEV113HE, VS-FEV136E and VS-FEV172E). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-FEV136E and VS-FEV172E (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.

- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

## 8. Test procedure

### 8.1. Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Yellow Fever* Positive Control in an area away from the other components.

### 8.2. RNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

If the Extraction Control is used to monitor nucleic acid isolation and as PCR inhibition control, add 5 µL of the EC to the specimen and/or lysis buffer mixture. (clinical specimen, as well as, positive control and/or negative control). Close each tube and vortex for 10 seconds.

If the Extraction Control is used only as a PCR inhibition control, 1 µL of the EC should be added to the reconstituted Reaction-Mix.

For RNA extraction from blood samples you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- QIAamp Viral RNA Mini Kit (Qiagen).
- ZP02003 MagPurix Viral Nucleic Acid Extraction Kit and ZP02013 MagPurix Viral RNA Extraction Kit.

### 8.3. Lyophilized positive control

Yellow Fever Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized Yellow Fever Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.4. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of RNA sample in different wells.

Add 5 µL of reconstituted Yellow Fever Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

| Cycles | Step                                   | Time   | Temperature |
|--------|----------------------------------------|--------|-------------|
| 1      | Reverse transcription                  | 15 min | 45°C        |
| 1      | Initial denaturation                   | 2 min  | 95°C        |
| 45     | Denaturation                           | 10 seg | 95°C        |
|        | Annealing/Extension (Data collection*) | 50 seg | 60°C        |

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (Yellow fever) and HEX, JOE or VIC channels (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne Plus™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for Yellow fever in the positive control well. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

| Yellow fever (FAM) | Extraction control (HEX) | Negative Control | Positive Control | Interpretation        |
|--------------------|--------------------------|------------------|------------------|-----------------------|
| +                  | +/-                      | -                | +                | Yellow fever Positive |
| -                  | +                        | -                | +                | Yellow fever Negative |
| -                  | -                        | -                | +                | Experiment fail       |
| +                  | +                        | +                | -                | Experiment fail       |

Table 5. Sample interpretation

+: Amplification curve

-: No amplification curve


A sample is considered positive if the Ct value obtained is less than 40 and the Extraction control shows or not an amplification signal. Sometimes, the detection of Extraction control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the Extraction control is positive. A failure in the extraction procedure and/or an inhibition of the PCR reaction can be excluded by the amplification of Extraction Control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of signal in both, extraction control and target in sample wells, we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of nucleic acids purification and/or inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

Differences can be observed in the values of Ct in the extraction controls between the control samples (negative control and positive control) and the clinical samples, due to the extraction process.

For samples at the limit of detection (LOD) or negative with a value of Ct for the extraction control  $\geq 35$ , it is recommended to filter the sample, or dilute the extracted (1:10 and 1:100) to avoid possible interference and inhibitors in the amplification reaction. It is recommended to review the instructions for using the extraction process used by the user.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with RNA extracted from blood samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.

- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Yellow fever either by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and Yellow fever Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.

## 11. Quality control

VIASURE Yellow fever Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Extraction control (EC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE Yellow fever Real Time PCR Detection Kit was evaluated using Tropical Virus panels (Dengue and West Nile Viruses) from QCMD for several years (2014, 2016, 2017 and 2018). These QCMD panels included 5 positive samples for yellow fever virus. The results were compared with the programs final reports. Other Tropical Virus including Dengue Virus Type 1, Dengue Virus Type 2, Dengue Virus Type 3, Dengue Virus Type 4, Zika Virus, Chikungunya Virus, Japanese Encephalitis virus, West Nile Virus (strains NY99 and Heja), and Tick borne Encephalitis virus, were not detected.


In conclusion, the results show a high sensitivity and specificity to detect Yellow fever using VIASURE Yellow fever Real Time PCR Detection Kit.

### 12.2. Analytical sensitivity

VIASURE Yellow fever Real Time PCR Detection Kit has a detection limit of  $\geq 10$  RNA copies per reaction. (Figure 2).

Figure 2. Dilution series of Yellow fever ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

### 12.3. Analytical specificity

The specificity of the yellow fever assay was confirmed by testing a panel consisting of different microorganisms representing the most common flavivirus, bloodstream pathogens or that cause acute febrile illness. No cross-reactivity was detected between almost any of the following microorganisms tested.

| Cross-reactivity testing                   |   |                                              |   |
|--------------------------------------------|---|----------------------------------------------|---|
| <i>Plasmodium falciparum</i>               | - | <i>Coxiella burnetii</i>                     | - |
| <i>Enterotoxigenic Escherichia coli</i>    | - | <i>Rickettsia conorii</i>                    | - |
| <i>Enteropathogenic Escherichia coli</i>   | - | Zika Virus Asian strain                      | - |
| <i>Salmonella enterica subsp. enterica</i> | - | Zika Virus French Polynesian strain 11474/16 | - |
| <i>Salmonella typhi</i>                    | - | Zika Virus African strain                    | - |
| <i>Pseudomonas aeruginosa</i>              | - | Chikungunya virus strain S27 Petersfield     | - |
| <i>Klebsiella oxytoca</i>                  | - | Dengue 1 virus strain Hawaii                 | - |
| <i>Listeria monocytogenes</i>              | - | Dengue 2 virus strain New Guinea C           | - |
| <i>Staphylococcus aureus subsp. aureus</i> | - | Dengue 3 virus strain H87                    | - |
| <i>Streptococcus pneumoniae</i>            | - | Dengue 4 virus strain H241                   | - |
| <i>Streptococcus pyogenes</i>              | - | St Louis Encephalitis virus strain 17D       | - |
| <i>Streptococcus agalactiae</i>            | - | Japanese encephalitis virus                  | - |
| <i>Haemophilus influenzae</i>              | - | West Nile virus strain NY99                  | - |
| <i>Bartonella henselae</i>                 | - | West Nile virus Heja                         | - |
| <i>Trypanosoma cruzi</i>                   | - | West Nile virus Ug37                         | - |

Table 6. Reference pathogenic microorganisms used in this study.

### 12.4. Analytical reactivity

The reactivity of VIASURE Yellow Fever Real Time PCR Detection Kit was evaluated against Yellow Fever virus strain 17D showing positive results.

## ANNEX 1

## COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

| Table A.1 LOW PROFILE BLOCK THERMOCYCLERS |                                                             |
|-------------------------------------------|-------------------------------------------------------------|
| Manufacturer                              | Model                                                       |
| Agilent Technologies                      | AriaMx/AriaDx Real-Time PCR System                          |
| Applied Biosystems                        | 7500 Fast Real-Time PCR System <sup>(1)</sup>               |
| Applied Biosystems                        | 7500 Fast Dx Real-Time PCR System <sup>(1)</sup>            |
| Applied Biosystems                        | QuantStudio™ 12K Flex 96-well Fast                          |
| Applied Biosystems                        | QuantStudio™ 6 Flex 96-well Fast                            |
| Applied Biosystems                        | QuantStudio™ 7 Flex 96-well Fast                            |
| Applied Biosystems                        | QuantStudio™ 3 Fast Real-Time PCR System <sup>(5)</sup>     |
| Applied Biosystems                        | QuantStudio™ 5 Fast Real-Time PCR System                    |
| Applied Biosystems                        | StepOne Plus™ Real-Time PCR System <sup>(5)</sup>           |
| Applied Biosystems                        | StepOne™ Real-Time PCR System <sup>(5)</sup>                |
| Applied Biosystems                        | ViiA™ 7 Fast Real-Time PCR System                           |
| BIONEER                                   | Exicycler™ 96                                               |
| Bio-Rad                                   | CFX96™ / CFX96™ IVD Real-Time PCR Detection System          |
| Bio-Rad                                   | Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup> |
| Cepheid                                   | SmartCycler® <sup>(3)</sup>                                 |
| Qiagen                                    | Rotor-Gene® Q <sup>(3)</sup>                                |
| Roche                                     | LightCycler® 480 Real-Time PCR System <sup>(4)</sup>        |
| Roche                                     | LightCycler® 96 Real-Time PCR System <sup>(4)</sup>         |
| Roche                                     | Cobas z480 Analyzer <sup>(4)</sup>                          |

| Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS |                                                                        |
|--------------------------------------------|------------------------------------------------------------------------|
| Manufacturer                               | Model                                                                  |
| Abbott                                     | Abbott m2000 RealTime System                                           |
| Applied Biosystems                         | 7300 Real-Time PCR System <sup>(5)</sup>                               |
| Applied Biosystems                         | 7500 Real-Time PCR System                                              |
| Applied Biosystems                         | 7900 HT Real-Time PCR System <sup>(5)</sup>                            |
| Applied Biosystems                         | ABI PRISM 7000 <sup>(6)</sup>                                          |
| Applied Biosystems                         | ABI PRISM 7700 <sup>(5)</sup>                                          |
| Applied Biosystems                         | QuantStudio™ 12K Flex 96-well                                          |
| Applied Biosystems                         | QuantStudio™ 6 Flex 96-well                                            |
| Applied Biosystems                         | QuantStudio™ 7 Flex 96-well                                            |
| Applied Biosystems                         | QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>                     |
| Applied Biosystems                         | QuantStudio™ 5 Real-Time PCR System                                    |
| Applied Biosystems                         | ViiA™ 7 Real-Time PCR System                                           |
| Analytik Jena Biometra                     | TOptical                                                               |
| Analytik Jena Biometra                     | qTOWER 2.0                                                             |
| BIONEER                                    | Exicycler™ 96                                                          |
| Bio-Rad                                    | CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System |
| Bio-Rad                                    | iCycler iQ™ Real-Time PCR Detection System                             |
| Bio-Rad                                    | iCycler iQ™5 Real-Time PCR Detection System                            |
| Bio-Rad                                    | MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>                    |
| Bio-Rad                                    | MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>                   |
| Cepheid                                    | SmartCycler® <sup>(3)</sup>                                            |
| DNA-Technology                             | DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>              |
| DNA-Technology                             | DTlite Real-Time PCR System <sup>(2)</sup>                             |
| Eppendorf                                  | Mastercycler™ep realplex                                               |
| Qiagen                                     | Rotor-Gene® Q <sup>(3)</sup>                                           |
| Stratagene / Agilent Technologies          | Mx3000P™ Real Time PCR System                                          |
| Stratagene / Agilent Technologies          | Mx3005P™ Real Time PCR System                                          |
| VIASURE                                    | VIASURE 48 Real Time PCR System <sup>(2)</sup>                         |
| VIASURE                                    | VIASURE 96 Real Time PCR System <sup>(2)</sup>                         |

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

## DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

| REAL-TIME PCR THERMOCYCLER    | VIASURE CHANNEL | DETECTION CHANNEL | OBSERVATIONS                                                                                                                                                                                                                                                                                                                                                    |
|-------------------------------|-----------------|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Bio-Rad CFX96™                | FAM             | FAM               | Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.                                                                                      |
|                               | HEX             | HEX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| ABI 7500 Applied Biosystems   | FAM             | FAM               | Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.         |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Roche Lightcycler®480II       | FAM             | 465/510           | Colour Compensation is required                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | 533/580           |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | 533/610           |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | 618/660           |                                                                                                                                                                                                                                                                                                                                                                 |
| Smartcycler® Cepheid          | FAM             | Channel 1         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | Channel 2         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Channel 3         |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Channel 4         |                                                                                                                                                                                                                                                                                                                                                                 |
| Abbott m2000rt                | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Mx3000P™ Mx 3005P™ Stratagene | FAM             | FAM               | Passive reference option for ROX must be none                                                                                                                                                                                                                                                                                                                   |
|                               | HEX             | VIC               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| AriaMx Agilent                | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | HEX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |
| Rotor-Gene®Q Qiagen           | FAM             | Green             | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".                                                                                                   |
|                               | HEX             | Yellow            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Orange            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Red               |                                                                                                                                                                                                                                                                                                                                                                 |
| Mic Real Time PCR Cycler bms  | FAM             | Green             | In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10) |
|                               | HEX             | Yellow            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | Orange            |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Red               |                                                                                                                                                                                                                                                                                                                                                                 |
| Exicycler™ 96 BIONEER         | FAM             | FAM               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | HEX             | JOE               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | ROX             | ROX               |                                                                                                                                                                                                                                                                                                                                                                 |
|                               | Cy5             | Cy5               |                                                                                                                                                                                                                                                                                                                                                                 |

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

**OPTICAL MEASUREMENT EXPOSURE SETTING**

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cyclers (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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### 13. Bibliography/Bibliografía

1. Drosten C. *et al.* Rapid detection and quantification of RNA of Ebola and Marburg viruses, Lassa virus, Crimean-Congo hemorrhagic fever virus, Rift Valley fever virus, dengue virus, and yellow fever virus by real-time reverse transcription-PCR. *Journal of Clinical Microbiology* 2002; 40(7):2323-30.
2. A.G. Fernandes-Monteiro *et al.* New approaches for the standardization and validation of a real-time qPCR assay using TaqMan probes for quantification of yellow fever virus on clinical samples with high quality parameters. *Human Vaccines & Immunotherapeutics* 2015; 11(7): 1865-1871.
3. World Health Organization. Yellow fever. <http://www.who.int/mediacentre/factsheets/fs100/en/>
4. World Health Organization. Emergencies preparedness, response. Yellow fever. <http://www.who.int/csr/disease/yellowfev/en>
5. Centers for Disease Control and Prevention. Yellow fever. <https://www.cdc.gov/yellowfever/>.

### 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: December 2019

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

**VIASURE**

CE

IVD

Real Time PCR Detection Kits

**CerTest**  
BIOTEC