Universal Turbilatex Sample Collection Vial

IU_MST-0019U

Cat.: IU_MST-0019U: Universal Turbilatex Sample Collection Vial (100x2.4mL)

INTENDED USE

Universal Turbilatex Sample Collection Vial is a single use vial containing diluent for faecal samples dilution. This product is intended to be used with FOB Turbilatex, Transferrin Turbilatex, Calprotectin Turbilatex and Pancreatic Elastase Turbilatex Combo reagents for the quantitative detection of human haemoglobin, human transferrin, human calprotectin and human pancreatic elastase in human stool samples, respectively.

For professional *in vitro* diagnostic use only.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- A trained person in Turbidimetric technique and autoanalyzer use is required.
- Do not use after expiration date. Sample diluent in unopened vials are stable up to the expiry date indicated on the package.
- Do not use the sample diluent if pack or vial is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.
- Sample diluent contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

STORAGE AND STABILITY

Kit components must be stored at temperature (2-30°C). Do not freeze.

Sample diluents are stable up to the expiration date printed on the label, always considering that vials must be properly closed to avoid any contamination, must be kept away from the sunlight and conserved at 2-30°C.

REAGENTS

Materials provided

- MST-0019U Universal Turbilatex Sample Collection Vial contains buffer and sodium azide <0.1% as preservative, (100x2.4mL).

Materials and instruments required but not provided

- Specimen collection container
- Shaker
- Microcentrifuge (10000g)
- Freezer (-20°C)
- Disposable gloves and laboratory equipment

SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of faeces (1-2 g or mL for liquid sample). Stool samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 3 days prior to testing. If not immediately tested, freeze the stored samples at -20°C. If not immediately tested, freeze the stored samples at -20°C maximum 6 months. In this case, the sample will be totally thawed, and brought to room temperature (15-30°C) before testing. Homogenise stool samples as thoroughly as possible prior to preparation.


Sample preparation

1. Collect a stool sample in a sterile container.
2. Pull out the shaped stick by turning the screw-cap (white) counter-clockwise. Do not remove the blue plastic.
3. Dip the shaped stick up to 4 different zones in the previously collected stool sample. and repeatedly rotate the stick in order to make the sample as homogenous as possible and make sure all grooves are filled with faecal material.
4. Before placing the stick into the tube, remove the stool excess by rotating the stick on the internal wall of the stool container.
5. Fully insert the stick with faecal material into the tube containing the extraction solution through the separator (blue). While turning the screw cap clockwise until completely closed, keep the separator (blue) by the fingers in order to avoid its opening. The sample dilution vial with diluted sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

Extraction Procedure

1. Bring the samples stored in a refrigerator at room temperature.
2. Shake each device by Vortex for 30-60 seconds in order to properly homogenate the content. If the grooves of the devices are visually empty of any stool material, proceed to step 3, otherwise repeat the vortexing for other 30-60 seconds. After a maximum vortexing time of 120 seconds proceed to step 3 regardless if there is some residual stool sample in the grooves.
3. Centrifuge the tube for 10 minutes at 1000-3000 x g or allow it to stand upright for at least 15 minutes in order to remove the possible residuals of faecal material, which might damage the aspirating and dispensing systems of the automated clinical chemistry instruments. After centrifugation, remove the device from the centrifuge.

4. After the completion of the previous steps, remove the screw cap (white) and the separator (blue) from devices by rotating the separator clockwise.
5. The faecal extract is now ready to be tested. Insert the opened device directly into the clinical chemistry analyser or transfer the supernatant to an adaptor cup (not provided).
6. Close the collection device by the provided additional screw cap and store or discard.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

IVD	For <i>in vitro</i> diagnostic use only	Keep dry
REF	Catalogue number	Temperature limitation
Σ	Use by	Lot number
n	Contains sufficient for <n> test	Manufacturer
CE	Do not reuse	DIL
	CE Marking	Sample diluent
		Keep away from sunlight

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