



CDV Ab Instructions For Use


Fast result in 10 minutes

For use with Canine Distemper Virus Antibody Test


Veterinary in-vitro diagnostic use only

Store between 4~30℃ (40~86 °F) **DO NOT FREEZE**


STEP 1 CHECK THE KIT CONTENTS BEFORE USE




STEP 2 TEST PROCEDURE



Use a sampling pipe to draw some whole blood/serum/plasma specimens (approximately 10 µl) through the capillary effect.




Peel off aluminum foil seal from the top of the extraction tube containing the extraction buffer, insert sampling end into pipe and cover tightly, mixed completely.



Holding the sampling pipe upright, carefully take off the cap of sampling pipe, transfer 3 drops (approximately 90µl) to the specimen well(S) of the test device.

STEP 3 INTERPRETATION OF TEST RESULT



Positive:

Two colored lines appear. One line should always appear in the control line region(C), and another one apparent colored line should appear in the test line region(T).



Negative:

Only one colored line appears in the control line region (C), and no colored line appears in the test line region (T).



Invalid:

No colored line appears in the control line region (C), indicating that the test result is ineffective. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. In this case, read the package insert carefully and test again with a new test device.




LIMITATIONS

Although the Canine Distemper Virus Antibody Test is very accurate in detecting canine distemper virus antibodies, a low incidence of false results may be occurred. Other clinically or laboratory tests might be required if questionable results are obtained. As other diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be diagnosed by the veterinarian after all clinical and laboratory findings have been evaluated.

INTERFERING SUBSTANCE

Analytes	Conc.	Analytes	Conc.
Whole Blood	20ul/ml	Oxymetazoline	0.6mg/ml
Mucin	50ug/ml	Phenylephrine	12mg/ml
Budesonide	200ul/ml	Rebetol	4.5ug/ml
Dexamethasone	0.8mg/ml	Relenza	282ng/ml
Flunisolide	6.8ng/ml	Tamiflu	1.1ug/ml
Mupirocin	12mg/ml		

SYMBOL



CLINICAL PERFORMANCE

The clinical performance was evaluated using different samples.

Clinical Result For CDV				
	Antigen	ELISA	Sensitivity	Specificity
Positive	115	120	95.83%	/
Negative	145	150	/	96.67%
95% Confidence interval			90.37%-98.46%	92.23%-98.78%

Zhejiang PushKang Biotechnology Co., Ltd.
Add: C408, Science and Technology Innovation Park NO.398, Mahuan Road, Binhai new Area, 312366 Shaoxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA.
WEB: www.pushkangbio.com

INTENDED USE

Canine Distemper Virus Antibody Test is a lateral flow immunochromatographic assay for the qualitative detection of canine distemper virus antibodies in dog's whole blood/serum/plasma.

PRECAUTIONS AND WARNINGS

- **Important:** Do not remove device from the pouch until ready for use.
- Do not use the components after expiration date.
- Humidity and temperature can influence the results.
- Do not mix components from kits with different lot numbers.
- All wastes should be properly decontaminated prior to disposal.
- Use a separate transfer pipette for each patient.
- Decontaminate and dispose of all samples, used kits and potentially contamination materials safely in accordance with national and local regulations.

KIT COMPONENTS

1. Test devices
2. Buffer
3. Sampling pipe
4. Package Insert

STORAGE AND STABILITY

- The test device is sealed and stored away from light at room temperature or refrigerated 4~30 °C (40~86°F). Do not freeze.
- The test kit should be used before the expiration date marked on the package insert.

DIRECTION FOR USE