VIASURE

Real Time PCR Detection Kit

Leishmania

CE IVD

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Content

1.	Intended use.....	6
2.	Summary and Explanation	6
3.	Principle of the procedure	6
4.	Reagents provided	7
5.	Reagents and equipment to be supplied by the user	7
6.	Transport and storage conditions.....	8
7.	Precautions for users	8
8.	Test procedure	10
9.	Result interpretation	11
10.	Limitations of the test	14
11.	Quality control.....	15
12.	Performance characteristics.....	16
	ANNEX 1	20
A1.1	Principle of the procedure.....	20
A1.2	Reagents provided.....	20
A1.3	Test procedure	21
	ANNEX 2	23
A2.1	Principle of the procedure.....	23
A2.2	Reagents provided.....	23
A2.3	Test procedure	23
	ANNEX 3	26
A3.1	Principle of the procedure.....	26
A3.2	Reagents provided.....	26
A3.3	Test procedure	28
	ANNEX 4	30
A4.1	Principle of the procedure.....	30
A4.2	Reagents provided.....	30
A4.3	Test procedure	30

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ENGLISH

1. Intended use

VIASURE *Leishmania* Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative detection of DNA from species of *Leishmania* in skin biopsy, blood and bone-marrow aspirate samples from individuals suspected of *Leishmania* spp. infection (or leishmaniasis) by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of visceral leishmaniasis (VL) and cutaneous leishmaniasis (CL) in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using qPCR and detected using fluorescent reporter dye probes specific for *Leishmania* species.

2. Summary and Explanation

Leishmaniasis is a vector-borne disease caused by a protozoan parasite from the family *Trypanosomatidae*, order *Kinetoplastida*, genus *Leishmania*. This parasite can infect both humans and other mammals after being bitten by a female phlebotomine sandfly, a tiny 2–3 mm long insect vector. Although asymptomatic infection is the most common outcome after parasite inoculation, more than 20 species of *Leishmania* cause clinical manifestations (disease development), which are classified in three main forms: cutaneous leishmaniasis (CL), visceral leishmaniasis (VL), and mucocutaneous leishmaniasis (MCL). CL is the most common variety of leishmaniasis worldwide, consisting of ulcers ranging from small-localized lesions to large ulcers all over the body. There may be many lesions, leading to disfiguring and/or incapacitating scars. *Leishmania* parasites causing CL are commonly divided into Old World species (like *L. major*, *L. tropica*, and *L. aethiopica*, that are prevalent around the Mediterranean Basin, the Middle East, the Horn of Africa, or the Indian subcontinent) and New World species (like *L. amazonensis*, *L. mexicana*, *L. braziliensis*, and *L. guyanensis*, endemic in Central and South America). These parasites are classified in one of two subgenera, *Leishmania* or *Viannia*, based on the site of development and attachment of the protozoa in the sandfly gut. MCL, often considering as a form of CL, is a potentially life-threatening and highly disfiguring condition due to the late-stage destruction of the oronasopharyngeal mucosa and cartilage, occasionally affecting the larynx, leading to aspiration pneumonia. 1-10% of MCL cases are caused after infection by a strain from the *Vianna* subgenera (*L. braziliensis* and *L. panamensis*, and, sometimes, *L. guyanensis* and *L. amazonensis*). VL comprises irregular bouts of fever, significant weight loss, swelling of the spleen and liver and serious anemia. This leishmaniasis form is caused by *L. donovani* in Asia and Africa, and *L. infantum* in the Mediterranean Basin, the Middle East, central Asia, South America and Central America, and its fatality rate can be as high as 100% within 2 years if it is not treated.

Leishmaniasis has been historically widespread in tropical climates across multiple territories. It is endemic in Asia, the Middle East, Northern Africa, the Mediterranean and South and Central America, being found in 89 countries. 1.5 to 2 new million new cases occur worldwide annually, and 70,000 deaths per year are attributed to this disease. According to the WHO, the majority of new VL cases reported worldwide in 2018 were in Brazil, Sudan, South Sudan, Ethiopia, India, and Kenya. Likewise, the corresponding to CL reported worldwide in 2018 were in Brazil, Bolivia, Peru, Colombia, Algeria, Tunisia, Iraq, Syria, Iran, Turkmenistan, and Afghanistan.

The diagnosis of leishmaniasis has been carried out traditionally by combination of direct and indirect diagnostic methods. Histopathologic examination under the microscope of fixed tissues, biopsies, or *in vitro* culture to detect amastigote parasite stage are considered the gold standard of diagnosis. However, the sensitivity of direct observation depends on the tissue, making necessary the utilization of alternative techniques. Several serological

assays are available, including direct agglutination test, ELISA, immunofluorescence and western blot, but should be interpreted in the context of clinical history. Molecular methods based on amplification of nuclear or kinetoplast DNA are very sensitive and allow for the identification of the *Leishmania* species. Molecular tests are especially important where simpler techniques fail (eg., in mucosal lesions where parasites are sporadic, and in chronic lesions). PCR platforms show sensitivity values between 92 and 100% and specificity of 100%. Other targets have been selected for PCR assays, such as Heat Shock Protein 70 kDa (HSP70), the internal transcribed spacer 1 (ITS-1), or the small 18S ribosomal subunit. The 18S rRNA region, part of the eukaryotic ribosome small subunit (SSU), is commonly utilized for the diagnosis of different *Leishmania* species, reflecting its highly conserved nature.

3. Principle of the procedure

VIASURE *Leishmania* Real Time PCR Detection Kit is designed for the qualitative detection of DNA from *Leishmania* spp. in skin biopsy, blood and bone-marrow aspirate samples. After DNA isolation, the detection of *Leishmania* spp. is performed by the amplification of a conserved region of the 18S rRNA gene of *Leishmania* spp. using specific primers and a fluorescent-labeled probe.

VIASURE *Leishmania* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Leishmania* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an **internal control** to monitor PCR inhibition.

4. Reagents provided

VIASURE *Leishmania* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Leishmania* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tube format (Annex 2 and 4).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.

- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).

VIASURE *Leishmania* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cyclers, DNA-Technology DTlite Real-Time PCR System, Roche Molecular Diagnostics LightCycler® 480, Roche Molecular Diagnostics Cobas z480 Analyzer, Applied Biosystems StepOne™ Real-Time PCR System, Linear NEOS-96 Real Time PCR System, Bioer Technology QuantGene 9600, Applied Biosystems QuantStudio®6 instrument, and Qiagen Rotor-Gene® Q. When using the Applied Biosystems 7500 Fast with strips it is recommended to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set of exposure values:

- DTprime Real-time Detection Thermal Cyclers (DNA-Technology): FAM channel -500*, HEX channel - 1000, ROX channel - 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel - 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *Leishmania* Reaction-Mix tube has been reconstituted, it may be kept at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended to store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.

- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-LEI113L, VS-LEI113H, VS-LEI136, VS-LEI172, VS-LEI113LE, VS-LEI113HE, VS-LEI136E and VS-LEI172E). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-LEI136, VS-LEI172, VS-LEI136E and VS-LEI172E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Leishmania* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification

criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.

- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for Open and Rotor-Gene format with extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *Leishmania* Real Time PCR Detection Kit has been tested in skin biopsy, blood and bone-marrow aspirate samples. Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at the following temperature conditions recommended by CDC following the local and national regulations for the transport of pathogen material.

- Skin biopsies: in sterile buffer medium at RT < 48h or at 2-8°C >48h and in sterile container at 2-8°C < 24 h.
- Blood samples and bone-marrow aspirates: <2h at RT or >2h at -80°C.

For long term transport (more than 2 hours), it is recommended shipping at -20°C or lower.

It is also recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 2 hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

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8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from clinical samples, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- EZ1 DNA tissue kit using the EZ1 Advanced automated instrument (QIAGEN).
- EZ1 DNA blood kit using the EZ1 Advanced automated instrument (QIAGEN).
- High Pure PCR purification kit for tissue samples (Roche Diagnostics).
- Speedtools DNA Extraction kit (Biotools).

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Leishmania spp.* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>Leishmania spp.</i> (FAM) ¹	Internal control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≤40	Valid

Table 1. Expected Performance of Controls. Ct no signal = no amplification curve.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

<i>Leishmania</i> spp. (FAM)	Internal control (HEX)	Interpretation for patients' individual samples
≤40	≤40 or no signal ¹	<i>Leishmania</i> spp. DNA Detected
≥40 or no signal	≤ 35 ²	<i>Leishmania</i> spp. DNA not Detected ²
≥40 or no signal	≥ 35 or no signal ²	Test Failure – Repeat Testing ²


Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

1 The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of *Leishmania* spp. target gene negative, IC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

9.2. References with extraction control (references in Annex 3 and 4)

All the results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Leishmania* spp. in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>Leishmania</i> spp. (FAM) ¹	Extraction Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	≥40 or no signal	≤40	Valid

Table 3. Expected Performance of Controls. Ct no signal = no amplification curve.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The Extraction Control (EC) should show an amplification signal (Ct ≤40) in control wells (PC and NC). Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following criteria for read and analyze the results:

<i>Leishmania</i> spp. (FAM)	Extraction control (HEX)	Interpretation for patients' individual samples
≤40	≤40 or no signal ¹	<i>Leishmania</i> spp. DNA Detected
≥40 or no signal	≤ 35 ²	<i>Leishmania</i> spp. DNA not Detected ²
≥40 or no signal	≥ 35 or no signal ²	Test Failure – Repeat Testing ²


Table 4. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction control between the controls and the clinical samples, due to the extraction process.

2 In the case of *Leishmania* spp. target gene negative, EC must show an amplification signal with Ct less than 35. If there is an absence of signal or Ct value ≥ 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from skin biopsy, blood and bone-marrow aspirate samples.
- The quality of the test depends on the quality of the sample; properly extracted nucleic acid from clinical samples must be done.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Leishmania spp.*, either the high number of DNA template copies which contains each *Leishmania* Positive Control vial, samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Leishmania* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.
- The specific primer and probe combinations for detection of the *18S rRNA* gene used in VIASURE *Leishmania* Real Time PCR Detection Kit designed for the detection of *Leishmania spp.*, do not show significant combined homologies with the human genome or human microflora, which might result in predictable false positive.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the protozoan DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown *Leishmania* variants.

- A protozoan load in the specimen below the limit of detection for the assay.
- The presence of qPCR inhibitors or other types of interfering substances. The impacts of any kind of drugs used to prevent leishmaniasis or used during the treatment of the infection have not been evaluated.
- Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable protozoa and does not imply that these protozoa are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of protozoan target sequences (18S rRNA gene).
- Negative results do not preclude leishmaniasis and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak protozoan levels during infections caused by *Leishmania* spp. have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other vector-borne illnesses are negative and the patient's clinical presentation and epidemiological information suggest that leishmaniasis is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.
- Analysis *in silico* of the primers and probe alignments showed that VIASURE *Leishmania* Real Time PCR Detection Kit detects the following species of *Leishmaniinae* subfamily, with no clinical relevance for leishmaniasis diagnosis:

Number of analyzed sequences	983
	<i>L. aethiopica</i> <i>L. infantum</i> <i>L. amazonensis</i> <i>L. major</i> <i>L. braziliensis</i> <i>L. mexicana</i> <i>L. chagasi</i> <i>L. panamensis</i> <i>L. donovani</i> <i>L. tropica</i> <i>L. gerbilli</i> <i>L. turanica</i> <i>L. guyanensis</i>
Max Fw mismatches	0
Max Rv mismatches	0
Max probe mismatches	0
Species detected without mismatches	<i>Crithidia expoeki</i> <i>Crithidia bombi</i> <i>Crithidia brevicula</i> <i>Crithidia mellificae</i>
Species detected with mismatches	<i>Lotmaria passim</i>

Table 5. Analysis *in silico* of primers and probe designed for *Leishmania* spp. detection.

11. Quality control

VIASURE *Leishmania* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Leishmania* Real Time PCR Detection Kit was tested using clinical samples (skin biopsy samples and visceral samples) from patients with suspicion of leishmaniasis. In order to determine the clinical diagnostic accuracy, multicenter evaluations have been conducted through collaboration with national entities. A summary of the sites, sample type and workflow is included in the following table:

	Site	Sample type	Workflow	Target
1	Certest Biotec and Instituto de Salud Carlos III (ISCIII) (Madrid, Spain)	All type of clinical samples	Speedtools DNA Extraction kit (Biotools) + Applied Biosystems 2720 thermal cycler + DNA Technologies DTlite Real-Time PCR System + Bio-Rad CFX96 Touch™ Real-Time PCR Detection System	<i>Leishmania spp.</i>
1.A		Skin biopsy samples		
1.B		Peripheral blood samples		
1.C		Bone-marrow aspirate samples		
2	'Hospital de la Sta. Creu i Sant Pau' and University of Barcelona, (Barcelona, Spain)	All type of clinical samples	EZ1 DNA tissue kit and EZ1 DNA blood kit (QIAGEN) using the EZ1 Advanced automated instrument + High Pure PCR purification kit for tissue samples (Roche Diagnosis) + Applied Biosystems QuantStudio®6 instrument	<i>Leishmania spp.</i>
2.A		Skin biopsy samples		
2.B		Visceral clinical samples		

Table 6. Summary of the sites, sample type and workflow carried out during the clinical evaluation.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *Leishmania* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Nested-PCR: Ln-PCR as described by Cruz et al.2002 ¹ , 2006 ²	<i>Leishmania</i> spp.	88	79	0	0	1 (0.95-1)	1 (0.95-1)	1 (0.95-1)	1 (0.95-1)
1.A			35	26	0	0	1 (0.90-1)	1 (0.86-1)	1 (0.90-1)	1 (0.86-1)
1.B			26	26	0	0	1 (0.86-1)	0.96 (0.86-1)	1 (0.86-1)	0.96 (0.86-1)
1.C			27	27	0	0	1 (0.87-1)	1 (0.87-1)	1 (0.87-1)	1 (0.87-1)
2	Reference algorithm *	<i>Leishmania</i> spp.	27	21	1	0	1 (0.87 – 1)	0.95 (0.77 – 0.99)	0.96 (0.79 – 0.99)	1 (0.80 – 1)
	<i>Leishmania</i> spp. Real-TM PCR kit (Sacace biotechnologies)		28	16	0	5	0.84 (0.63 – 0.94)	1 (0.79 – 1)	1 (0.84 – 1)	0.76 (0.52 – 0.90)
2.A	Reference algorithm *		23	20	1	0	1 (0.85-1)	0.95 (0.76-0.99)	0.95 (0.76-0.99)	1 (0.79-1)
	<i>Leishmania</i> spp. Real-TM PCR kit (Sacace biotechnologies)		24	16	0	4	0.85 (0.67 – 0.96)	1 (0.79 – 1)	1 (0.82 – 1)	0.80 (0.55-0.93)
2.B	Reference algorithm *		4	1	0	0	1 (0.39 – 1)	1 (0.02 – 1)	1 (0.39 – 1)	1 (0.02 – 1)
	<i>Leishmania</i> spp. Real-TM PCR kit (Sacace biotechnologies)		4	0	0	1	0.8 (0.28-0.99)	n.a.#	1 (0.47-1)	n.a.#

Table 7. True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *Leishmania* Real Time PCR Detection Kit.

Due to the lack of true negative samples, the calculation of the clinical specificity of the test could not be performed.

* Reference algorithm: Light microscopy + sample culture using the NNN (Novy-MacNeal-Nicolle) and the Schneider *Drosophyla* media + In-house real-time PCR as described by Mary et al. 2004³ with variations according to Martín-Ezquerro et al., 2009⁴ and Velez et al., 2020⁵.

1 Cruz, I., Canavate, C., Rubio, J.M., Morales, M.A., Chicharro, C., Laguna, F., Jimenez-Mejias, M., Sirera, G., Videla, S., Alvar, J., Spanish HIV-Leishmania Study Group. (2002). A nested polymerase chain reaction (Ln-PCR) for diagnosing and monitoring *Leishmania* infantum infection in patients co-infected with human immunodeficiency virus. *Trans R Soc Trop Med Hyg*, 96, S1/185-S1/189.

2 Cruz, I., Chicharro, C., Nieto, J., Bailo, B., Cañavate, C., Figueras, M.C., Alvar, J. (2006) Comparison of new diagnostic tools for management of pediatric Mediterranean visceral leishmaniasis. *J Clin Microbiol*, 44(7), 2343-7.

3 Mary, C., Faraut, F., Lascombe, L., et al. (2004) Quantification of *Leishmania* infantum DNA by a rea-time PCR assay with high sensitivity. *J. Clin Microbiol*, 42, 5249-55. doi: 10.1128/JCM.42.11.5249-5255.2004.

4 Martín-Ezquerro, G., Fisa, R., Riera, C., Rocamora, V., Fernández-Casado, A., Barranco, C., Serra, T., Baró, T., Pujol, R.M. (2009) Role of *Leishmania* spp. infestation in nondiagnostic cutaneous granulomatous lesions: Report of a series of patients from a Western Mediterranean area. *British Journal of Dermatology* 161, 320–325. doi: 0.1111/j.1365-2133.2009.09282.x.


5 Velez, R., Domenech, E., Rodríguez-Cortés, A., Barrios, D., Tebar, S., Fernández-Arévalo, A., Aguilar, R., Dobaño, C., Alberola, J., Cairó, J., Gállego, M. (2020). Evaluation of canine leishmaniosis vaccine CaniLeish® under field conditions in native dog populations from an endemic Mediterranean area – a randomized controlled trial. *Acta Tro*, doi: 10.1016/j.actatropica.2020.105387.

The results show high agreement to detect *Leishmania* spp. using VIASURE *Leishmania* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Leishmania* Real Time PCR Detection Kit has a detection limit of 10 genome copies per reaction for *Leishmania* spp. (18S rRNA gene) with a positive rate of $\geq 95\%$ on blood samples.

Figure 3. Dilution series of *Leishmania* spp. (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

12.3. Analytical specificity

The specificity of *Leishmania* spp. assay was confirmed by testing a panel consisting of different microorganisms associated to blood-borne diseases. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Anaplasma marginale</i>	-	Chikungunya virus strain S27 Petersfield	-	St Louis Encephalitis virus	-
<i>Bartonella henselae</i> strain Houston-1	-	Chikungunya virus F24	-	<i>Theileria annulata</i>	-
<i>Borrelia hermsii</i>	-	<i>Coxiella burnetii</i> strain Nine Mile Q	-	Tick-Borne Encephalitis virus strain Neudorfl	-
<i>Borrelia lusitanae</i>	-	Dengue virus serotype 1 strain Hawaii A	-	<i>Treponema phagedenis</i>	-
<i>Borrelia valaisiana</i>	-	Dengue virus serotype 2 strain New Guinea C	-	<i>Trypanosoma cruzi</i>	-
<i>Borrelia azfeli</i> strain P-Ko/1984	-	Dengue virus serotype 3 strain H87	-	West Nile virus strain NY99	-
<i>Borrelia bavariensis</i>	-	Dengue virus serotype 4 strain H241	-	West Nile virus strain Heja	-
<i>Borrelia bisetti</i>	-	Japanese Encephalitis virus	-	West Nile virus strain Ug37	-
<i>Borrelia burgdorferi</i> strain IRS	-	Japanese Encephalitis virus strain Nakayama	-	Yellow Fever virus strain 17D	-
<i>Borrelia burgdorferi sensu stricto</i> strain B31	-	<i>Leptospira</i>	-	Yellow Fever virus strain French Neurotropic	-
<i>Borrelia garinii</i> Type strain	-	<i>Plasmodium falciparum</i> strain 3D7	-	Zika FB-GWUH-2016	-
<i>Borrelia japonica</i>	-	<i>Rickettsia conorii</i> strain Moroccan	-	Zika virus African strain	-
<i>Borrelia miyamotoi</i>	-	Rift Valley Fever Virus strain AR21229	-	Zika virus Asian strain PF13/251013-18	-
<i>Borrelia spielmanii</i>	-	Rift Valley Fever Virus strain MP12	-	Zika Virus French Polynesian strain 11468/16	-

Table 8. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Leishmania* Real Time PCR Detection Kit was evaluated against synthetic DNA, and against DNA extracted from clinical samples infected by different species of *Leishmania*: *L. infantum*, *L. major*, *L. mexicana*, *L. tropica*, and subgenus *L. Viannia* (*L. braziliensis* and *L. panamensis*), showing positive results.

ANNEX 1

OPEN AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-LEI106L
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-LEI106H
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-LEI112L
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-LEI112H
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-LEI113L
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-LEI113H
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-LEI136
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-LEI172
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-LEI101L
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-LEI101H
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS-LEI101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *Leishmania* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition.

Target	Channel	Gene
<i>Leishmania</i> spp.	FAM	18S rRNA gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.certest.es.

A1.2 Reagents provided

VIASURE *Leishmania* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI101L, VS-LEI101H, VS-LEI106L, VS-LEI106H, VS-LEI112L and VS-LEI112H.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref VS-LEI113L and VS-LEI113H.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI101, VS-LEI136 and VS-LEI172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Leishmania Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Leishmania* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Leishmania* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (18S *rRNA* gene) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-LEI196T

Table A2. 1.References.

A2.1 Principle of the procedure

VIASURE *Leishmania* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to monitor the PCR inhibition.

Target	Channel	Gene
<i>Leishmania</i> spp.	FAM	18S rRNA gene
Internal control (IC)	HEX, VIC or JOE *	-

Table A2. 2.Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.certest.es.

A2.2 Reagents provided

VIASURE *Leishmania* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

Leishmania Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Leishmania* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Leishmania* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Leishmania* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Leishmania* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (18S rRNA gene) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR

System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 3

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-LEI106LE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-LEI106HE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-LEI112LE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-LEI112HE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-LEI113LE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-LEI113HE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-LEI136E
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-LEI172E
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-LEI101LE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-LEI101HE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit 1 x 4-well strips, Rotor-Gene®	VS-LEI101E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE *Leishmania* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Leishmania</i> spp.	FAM	18S rRNA gene
Extraction control (EC)	HEX, VIC or JOE *	-

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A3.2 Reagents provided

VIASURE *Leishmania* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3, A3.4 and A3.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A3.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A3.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A3. 3. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI101LE, VS-LEI101HE, VS-LEI106LE, VS-LEI106HE, VS-LEI112LE and VS-LEI112HE.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A3. 4. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref VS-LEI113LE and VS-LEI113HE.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A3. 5. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI101E, VS-LEI136E and VS-LEI172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) by adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *Leishmania* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

Leishmania Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Leishmania* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Leishmania* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A3. 6. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (18S *rRNA* gene) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Leishmania</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-LEI196TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE *Leishmania* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Leishmania</i> spp.	FAM	18S rRNA gene
Extraction control (EC)	HEX, VIC or JOE *	-

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.ceratest.es.

A4.2 Reagents provided

VIASURE *Leishmania* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
<i>Leishmania</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Leishmania</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE *Leishmania* Real Time PCR Detection Kit with Ref. VS-LEI196TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *Leishmania* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

Leishmania Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Leishmania* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Leishmania* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Leishmania* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *Leishmania* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µl of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (18S rRNA gene) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

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Bibliography/Bibliografía

1. Bruschi, F., & Gradoni, L. (2018). The leishmaniasis: Old neglected tropical diseases. In *The Leishmaniasis: Old Neglected Tropical Diseases*. <https://doi.org/10.1007/978-3-319-72386-0>
2. Burza, S., Croft, S. L., & Boelaert, M. (2018). Leishmaniasis. *The Lancet*, 392(10151), 951–970. [https://doi.org/10.1016/S0140-6736\(18\)31204-2](https://doi.org/10.1016/S0140-6736(18)31204-2)
3. Conter, C. C., Mota, C. A., dos Santos, B. A., de Souza Braga, L., de Souza Terron, M., Navasconi, T. R., Fernandes, A. C. B. S., Demarchi, I. G., de Castro, K. R. R., Aristides, S. M. A., Lonardoni, M. V. C., Teixeira, J. J. V., & Silveira, T. G. V. (2019). PCR primers designed for new world Leishmania: A systematic review. *Experimental Parasitology*, 207, 107773. <https://doi.org/10.1016/j.exppara.2019.107773>
4. Gállego, M., & Riera, C. (n.d.). LAS LEISHMANIOSIS HUMANAS: LEISHMANIOSIS AUTÓCTONA POR *Leishmania infantum*.
5. León, C. M., Muñoz, M., Hernández, C., Ayala, M. S., Flórez, C., Teherán, A., Cubides, J. R., & Ramírez, J. D. (2017). Analytical performance of Four Polymerase Chain Reaction (PCR) and real time PCR (qPCR) assays for the detection of six *Leishmania* species DNA in Colombia. *Frontiers in Microbiology*, 8, 1907. <https://doi.org/10.3389/fmicb.2017.01907>
6. Maxfield, L., & Crane, J. S. (2020). Leishmaniasis. In *StatPearls [Internet]*. Treasure Island (FL): StatPearls Publishing; 2021 Jan-. <https://www.ncbi.nlm.nih.gov/books/NBK531456/>
7. Moreira, O. C., Yadon, Z. E., & Cupolillo, E. (2017). The applicability of real-time PCR in the diagnostic of cutaneous leishmaniasis and parasite quantification for clinical management: Current status and perspectives. *Acta Tropica*. <https://doi.org/10.1016/j.actatropica.2017.09.020>
8. Steverding, D. (2017). The history of leishmaniasis. *Parasites and Vectors*, 10(82), 1–10. <https://doi.org/10.1186/s13071-017-2028-5>
9. World Health Organization. (2021). *Leishmaniasis*. https://www.who.int/health-topics/leishmaniasis#tab=tab_1

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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