VIASURE

Real Time PCR Detection Kit

M. tuberculosis complex

CE IVD

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PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101L
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101H
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106L
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106H
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112L
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112H
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113L
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113H
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172

Table A 1. References for Open and Rotor-Gene format with internal control products. / Referencias para productos Open y Rotor-Gene Format con control interno.

TUBE FORMAT WITH INTERNAL CONTROL (SEE ANNEX 2) / FORMATO TUBO CON CONTROL INTERNO (VER ANEXO 2)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE M. tuberculosis complex Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTC196T

Table A 2. References for Tube format with internal control products. / Referencias para productos formato Tubo con control interno.

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL (SEE ANNEX 3) / OPEN Y ROTOR-GENE FORMAT CON CONTROL DE EXTRACCIÓN (VER ANEXO 3)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101LE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101HE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106LE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106HE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112LE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112HE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113LE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113HE
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101E
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136E
VIASURE M. tuberculosis complex Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172E

Table A 3. References for Open and Rotor-Gene format with extraction control products. / Referencias para productos Open y Rotor-Gene Format con control de extracción.

TUBE FORMAT WITH EXTRACTION CONTROL (SEE ANNEX 4) / FORMATO TUBO CON CONTROL DE EXTRACCIÓN (VER ANEXO 4)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
VIASURE M. tuberculosis complex Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTC196TE

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EN For download IFUS from other languages, please enter in **certest.es/viasure/labeling**. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use.....	7
2.	Summary and Explanation	7
3.	Principle of the procedure	8
4.	Reagents provided	8
5.	Reagents and equipment to be supplied by the user	8
6.	Transport and storage conditions.....	9
7.	Precautions for users	9
8.	Test procedure	10
8.1.	Specimen collection, transport and storage	11
8.2.	DNA extraction.....	11
9.	Result interpretation	12
9.1.	References with internal control (references in Annex 1 and 2).....	12
9.2.	References with extraction control (references in Annex 3 and 4).....	13
10.	Limitations of the test	15
11.	Quality control.....	16
12.	Performance characteristics.....	16
12.1.	Clinical sensitivity and specificity	16
12.2.	Analytical sensitivity	18
12.3.	Analytical specificity.....	18
12.4.	Analytical reactivity	19
	ANNEX 1	20
A1.1	Principle of the procedure.....	20
A1.2	Reagents provided.....	20
A1.3	Test procedure	21
A1.3.1	Lyophilized positive control.....	21
A1.3.2	PCR protocol.....	22
	ANNEX 2	23
A2.1	Principle of the procedure.....	23
A2.2	Reagents provided.....	23
A2.3	Test procedure	23
A2.3.1	Lyophilized positive control.....	23
A2.3.2	Lyophilized reaction mix tube.....	24
A2.3.3	PCR protocol.....	24
	ANNEX 3	25
A3.1	Principle of the procedure.....	25
A3.2	Reagents provided.....	25
A3.3	Test procedure	27
A3.3.1	Lyophilized extraction control	27
A3.3.2	Lyophilized positive control.....	27
A3.3.3	PCR protocol.....	27
	ANNEX 4	29
A4.1	Principle of the procedure.....	29

A4.2 Reagents provided.....	29
A4.3 Test procedure	29
A4.3.1 Lyophilized extraction control	29
A4.3.2 Lyophilized positive control.....	30
A4.3.3 Lyophilized reaction mix tube.....	30
A4.3.4 PCR protocol.....	30

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ENGLISH

1. Intended use

VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative detection of DNA from *Mycobacterium tuberculosis* complex (MTBC) species in respiratory samples (sputum) from individuals suspected of Tuberculosis by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of the infection by species of *Mycobacterium tuberculosis* complex, in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using real-time PCR, and detected using fluorescent reporter dye probes specific for species of *Mycobacterium tuberculosis* complex.

2. Summary and Explanation

Tuberculosis is a contagious, chronic, and granulomatous disease caused by *Mycobacterium tuberculosis*. This disease was declared in 1993 as a "global health emergency" due to its magnitude as a public health problem. Infection occurs through the inhalation of aerosols that contain the pathogen and are transmitted by people with active pulmonary tuberculosis. After inhalation, the bacteria deposit in the alveoli and spread through the lymphatic circulation. Further dissemination to other parts of the lung and occasionally to other organs is achieved through hematogenous circulation. The most common form of the disease is pulmonary tuberculosis, although tuberculous meningitis, miliary tuberculosis (disseminated), intestinal tuberculosis, lymphadenitis, osteomyelitis and Pott's disease (affected bones) also occur. Primary infection leads to active disease in approximately 10% of infected people, and in 80% of these cases, in the period of two years. In the remaining 90%, the immune system controls the infection, and the individual is not infectious or asymptomatic. It is estimated that one third of the world population has latent tuberculosis; that is, these people are infected with the bacillus, but (still) have not become ill or can transmit the infection. In this clinical state, TB bacilli can remain inactive for years (latent TB). However, when the immune system weakens, the latent infection can be reactivated. In a person infected with HIV, the risk of reactivation of latent TB is more than 10% per year, compared to a lifetime risk of 10-20% for HIV negative people.

Mycobacteria are a group of microorganisms that constitute one of the most serious health problems worldwide. Three groups can be defined within the genus *Mycobacterium*: 1) *Mycobacterium tuberculosis* complex that produces tuberculosis and is formed by the species *M. tuberculosis*, *M. bovis* (including *M. bovis* BCG), *M. africanum*, *M. microti*, *M. caprae*, *M. canettii*, *M. pinnipedii*, *M. mungi* and *M. suricattae*; 2) *M. leprae* that produces leprosy; 3) Other non-tuberculous mycobacteria (NTM) that are opportunistic and produce non-tuberculous symptoms with less pathogenic power.

When the active form of the disease occurs, the symptoms (cough, fever, night sweats, weight loss, etc.) can be mild for many months. As a result, patients sometimes take time to seek medical attention and transmit the bacteria to others. Over a year, a tubercular patient can infect 10 to 15 people through close contact. If they do not receive the appropriate treatment, up to two thirds of the tubercular patients die.

The application of molecular techniques for diagnosis and typing, more sensitive, specific and quick than traditional tests, allow to improve the knowledge of the epidemiology of the infection and facilitate decisions for its control.

3. Principle of the procedure

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is designed for the qualitative detection of DNA from strains of *M. tuberculosis complex* in respiratory samples (sputum). After DNA isolation, the detection of strains of *M. tuberculosis complex* is performed by the amplification of a conserved region of the insertion sequences *IS6110* and *IS1081*, using specific primers and fluorescent-labelled probes.

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer, polymerase) in a stabilized format, as well as an internal control to verify the correct functioning of the amplification mix.

4. Reagents provided

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit has been validated on the following equipments: 7500 Fast Real-Time PCR System (Applied Biosystems), StepOne™ Real-Time PCR System (Applied Biosystems), CFX96™ Real-Time PCR Detection System (Bio-Rad), AriaMx Real-Time PCR System (Agilent Technologies), DTprime

Real-time Detection Thermal Cycler (DNA-Technology), DTLite Real-Time PCR System (DNA-Technology), Rotor-Gene® Q Cycler (Qiagen), Cobas z480 Analyzer (Roche Molecular Diagnostics), ViiA™ 7 Real-Time PCR System (Applied Biosystems) and VIASURE V-Lab96 Cycler (CerTest Biotec S.L). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTLite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *M. tuberculosis complex* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at – 20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-MTC113L, VS-MTC113H, VS-MTC136, VS-MTC172, VS-MTC113LE, VS-MTC113HE, VS-MTC136E and VS-MTC172E). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.

- For references VS-MTC101, VS-MTC136, VS-MTC172, VS-MTC101E, VS-MTC136E and VS-MTC172E (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

Please, see Annex 1 for Open and Rotor-Gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for Open and Rotor-Gene format with extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit has been tested on respiratory samples (sputum). Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at room temperature for up to 2 hours, or at 4°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), we recommend shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 4°C for up to 24 hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>), the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94) and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cencenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

For sputum samples, the recommended amount of sputum is approximately 200µL that can be pipetted. * In order to improve the yield and quality of bacterial DNA from sputum specimens, it is recommended a sample pre-treatment with N-acetyl-L-cysteine–sodium hydroxide (NALC-NaOH). Consult DNA Extraction section for more information.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

Due to VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is also available with extraction control (references include in Annexes III and IV), if the Extraction Control (EC) is used to monitor nucleic acid isolation and as PCR inhibition control, add 5µl of the reconstituted EC to the specimen and/or lysis buffer mixture (clinical specimen, as well as positive control and/or negative control). Close each tube and vortex for 10 seconds. If the Extraction Control is used only as a PCR inhibition control, 1µl of the reconstituted EC should be added to the reconstituted Reaction-Mix.

For DNA extraction from respiratory samples (sputum*), you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- Invisorb® Spin Universal kit (Invitek).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

* In order to improve the yield and quality of bacterial DNA from sputum specimens, it is recommended a sample pre-treatment with N-acetyl-L-cysteine–sodium hydroxide (NALC-NaOH). After decontamination, the cell pellets should be resuspended in a maximum final volume of 1 mL of phosphate buffer. In addition, usage of small elution volumes (50-100 µL) may raise the DNA concentration.

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Mycobacterium tuberculosis* complex in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>M. tuberculosis</i> complex (FAM) ¹	Internal Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	≤40	Valid

Table 1. Expected Performance of Controls.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The Internal Control (IC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read, and analyse the results:

<i>M. tuberculosis</i> complex (FAM)	Internal Control (HEX)	Interpretation for patient's individual samples	
≤40	≤40 or no signal ¹	Valid	DNA of <i>M. tuberculosis</i> complex strains Detected
>40 or no signal	≤35 ²	Valid	Targets not Detected ²
>40 or no signal	>35 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curves.

1 The internal control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.


2 In case of *M. tuberculosis* complex target genes negative, IC must show an amplification signal with Ct ≤35. If there is an absence of signal in Ct value >35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate.

Depending on the available material:

- Repeat qPCR with the same isolated DNA sample, or
- Re-extract and retest another aliquot of the same specimen, or
- Obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

9.2. References with extraction control (references in Annex 3 and 4)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for *Mycobacterium tuberculosis* complex in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>M. tuberculosis</i> complex (FAM) ¹	Extraction Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	≤40	Valid

Table 3. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The Extraction Control (EC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following criteria for read and analyse the results:

<i>M. tuberculosis</i> complex (FAM)	Extraction Control (HEX)	Interpretation for patient's individual samples	
≤40	≤40 or no signal ¹	Valid	DNA of <i>M. tuberculosis</i> complex strains Detected
>40 or no signal	≤35 ²	Valid	Targets not Detected²
>40 or no signal	>35 or no signal ²	Invalid	Test Failure – Repeat Testing²

Table 4. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.


2 In the case of *M. tuberculosis* complex target gene negative, EC must show an amplification signal with Ct ≤ 35. If there is an absence of signal or Ct value > 35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- Repeat qPCR with the same isolated DNA sample, or
- Re-extract and retest another aliquot of the same specimen, or

c) Obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with DNA extracted from respiratory samples (sputum).
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Mycobacterium tuberculosis* complex, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *M. tuberculosis* complex Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the bacterial DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown *Mycobacterium tuberculosis* complex strains.
 - A bacterial load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent infection of *Mycobacterium tuberculosis* complex strains or used during the treatment of the infection have not been evaluated.

- Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable bacteria and does not imply that these bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targeted bacterial sequences (insertion sequences *IS6110* and *IS1081*).
- Negative results do not preclude *Mycobacterium tuberculosis* complex strains infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak bacterial levels during infections caused by *Mycobacterium tuberculosis* complex strains have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that *Mycobacterium tuberculosis* complex strain/s infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit was tested using respiratory samples (sputum) from patients with suspicion of tuberculosis. In order to determine the clinical diagnostic accuracy, different multicenter evaluations have been conducted through collaboration with national entities. A summary of the sites, sample type and workflow are included in the following table. The results were as follows:

	Site	Sample type	Workflow	Target
1	"Instituto Aragonés de Ciencias de la Salud" (Hospital Universitario Miguel Servet - HUMS-, Zaragoza, Spain)	DNA extracted from culture (Sputum specimens)	No information is available on the DNA extraction process of the clinical sample panel + StepOne™ Real-Time PCR System (Applied Biosystems)	<i>Mycobacterium tuberculosis</i> complex
	"Instituto Aragonés de Ciencias de la Salud" (Hospital Universitario Miguel Servet - HUMS-, Zaragoza, Spain)	DNA extracted from culture (Sputum specimens)	No information is available on the DNA extraction process of the clinical sample panel + ViiA™ 7 Real-Time PCR System (Applied Biosystems)	<i>Mycobacterium tuberculosis</i> complex
2	Hospital Ernest Lluch (HEL, Calatayud, Zaragoza, Spain)	bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS), sputum	MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.) + DTlite Real-Time PCR System (DNA-Technology)	<i>Mycobacterium tuberculosis</i> complex
3	Certest Biotec S.L (San Mateo de Gállego, Spain)	bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS), sputum	MagDEA® Dx Nucleic Acid Extraction kit using the magLEAD® 12gC (Precision System Science Co) + LightCycler® 480 Real-Time PCR System (Roche).	<i>Mycobacterium tuberculosis</i> complex

Table 5. Site, sample type, workflow, and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit in sputum samples were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Culture, RFLP-IS6110 and spoligotyping	<i>Mycobacterium tuberculosis</i> complex	61	4	1	1	0.98 (0.91-1)	0.8 (0.28-0.99)	0.98 (0.9-0.99)	0.8 (0.28-0.99)
	Culture, RFLP-IS6110 and spoligotyping	<i>Mycobacterium tuberculosis</i> complex	60	5	0	0	1 (0.94-1)	1 (0.47-1)	1 (0.94-1)	1 (0.47-1)
2	Culture and bacilloscopic	<i>Mycobacterium tuberculosis</i> complex (sputum)	10	35	0	0	1 (0.65-1)	1 (0.87-1)	1 (0.65-1)	1 (0.87-1)
		<i>Mycobacterium tuberculosis</i> complex (BAL)	1	6	0	0	1 (0.054-1)	1 (0.51-1)	1 (0.054-1)	1 (0.51-1)
		<i>Mycobacterium tuberculosis</i> complex (BAS)	2	6	0	0	1 (0.19-1)	1 (0.51-1)	1 (0.19-1)	1 (0.51-1)
3	Amplified <i>Mycobacterium Tuberculosis</i> Direct Test (Hologic) + smear and culture	<i>Mycobacterium tuberculosis</i> complex	2	325	0	0	1 (0.15-1)	1 (0.98 – 1)	1 (0.15-1)	1 (0.98 – 1)

Table 6. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit.


In addition to this, 10 pleural fluid samples, as well as 6 urines and 1 faecal specimen were also analysed with VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit, obtaining concordant results comparing to reference methods used.

Results show high agreement to detect *Mycobacterium tuberculosis* complex strains using VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit has a detection limit of 10 DNA copies per reaction for the insertion sequences IS6110 and IS108, with a positive rate of $\geq 95\%$ (Figure 3).

Figure 3. Dilution series of insertion sequences IS6110 and IS1081 (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (FAM channel).

12.3. Analytical specificity

The specificity of the VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit was confirmed by testing a panel consisting of different microorganisms associated to respiratory diseases. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	Influenza A/South Australia/55/2014, IVR-175)	-	<i>Mycobacterium gordonae</i>	-
<i>Bordetella parapertussis</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycobacterium kansasii</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus (clade 3C2a.1)	-	<i>Mycobacterium marinum</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus	-	<i>Mycobacterium peregrinum</i>	-
<i>Chlamydia caviae</i>	-	Influenza B/Brisbane/60/2008 virus	-	<i>Mycobacterium smegmatis</i>	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza B/Florida/04/06 virus	-	<i>Mycoplasma pneumoniae</i>	-
<i>Chlamydia psittaci</i> genotypes A and C	-	Influenza B/Phuket/3073/2013 virus	-	Human Adenovirus 1-5, 8, 15, 31, 40 and 41	-
<i>Haemophilus influenzae</i> MinnA	-	<i>Legionella bozemanii</i>	-	Human Bocavirus	-
Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Legionella micdadei</i>	-	Human coronavirus 229E, OC43 and NL63	-
Influenza A/California/7/2009(H1N1)pdm09 virus	-	<i>Legionella dumoffii</i>	-	Human metapneumovirus A and B	-
Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Legionella pneumophila</i>	-	Human parainfluenza 1, 2, 3 and 4 viruses	-

Influenza A/Hong Kong/4801/2014, NYMC X-263B	-	<i>Legionella longbeache</i>	-	Human rhinovirus	-
Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Moraxella catarrhalis</i>	-	MERS Coronavirus	-
Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Mycobacterium abscessus</i>	-	<i>Streptococcus pneumoniae</i> Z022	-
Influenza A/Victoria/210/2009(H3N2) virus	-	<i>Mycobacterium avium</i> subsp. <i>avium</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Mycobacterium fortuitum</i>	-	Respiratory Syncytial virus (RSV) type A and B	-

Table 7. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit for *M. tuberculosis* complex strains was evaluated using DNA extracted from reference /clinical strains (as templates): *Mycobacterium bovis* AF 2122/97, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *M. tuberculosis* L2 sublineage Beijing, *M. tuberculosis* L4 Euro-American sublineages, *M. tuberculosis* T, *M. tuberculosis* Haarlem, *M. tuberculosis* LAM, *M. tuberculosis* X, *Mycobacterium africanum* (L5 and L6) and *Mycobacterium microti*, showing positive results .

ANNEX 1

OPEN AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Mycobacterium tuberculosis complex</i>	FAM	Insertion sequences <i>IS6110</i> and <i>IS1081</i>
Internal control (IC)	HEX, VIC or JOE *	-

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.cerTEST.es.

A1.2 Reagents provided

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.cerTEST.es).

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref. VS-MTC101L, VS-MTC101H, VS-MTC106L, VS-MTC106H, VS-MTC112L and VS-MTC112H.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref VS-MTC113L and VS-MTC113H.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref. VS-MTC101, VS-MTC136 and VS-MTC172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 mL (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

M. tuberculosis complex Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized

M. tuberculosis complex Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis complex*) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis</i> complex Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTC196T

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Mycobacterium tuberculosis</i> complex	FAM	Insertion sequences <i>IS6110</i> and <i>IS1081</i>
Internal control (IC)	HEX, VIC or JOE *	-

Table A2. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.certes.es.

A2.2 Reagents provided

VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref. VS-MTC196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

M. tuberculosis complex Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis* complex Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *M. tuberculosis complex* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *M. tuberculosis complex* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.cerTEST.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis complex*) and HEX, JOE or VIC (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.cerTEST.es).

ANNEX 3

OPEN AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101E
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136E
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Mycobacterium tuberculosis complex</i>	FAM	Insertion sequences <i>IS6110</i> and <i>IS1081</i>
Extraction Control (EC)	HEX, VIC or JOE *	-

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A3.2 Reagents provided

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3, A3.4 and A3.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A3.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A3.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es)

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A3. 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref. VS-MTC101LE, VS-MTC101HE, VS-MTC106LE, VS-MTC106HE, VS-MTC112LE and VS-MTC112HE.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A3. 4. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref VS-MTC113LE and VS-MTC113HE.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A3 5. Reagents and materials provided in VIASURE *M. tuberculosis* complex Real Time PCR Detection Kit with Ref. VS-MTC0101E, VS-MTC136E and VS-MTC172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 mL (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *M. tuberculosis complex* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

M. tuberculosis complex Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *M. tuberculosis complex* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A3.6. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis* complex) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTC196TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>Mycobacterium tuberculosis complex</i>	FAM	Insertion sequences <i>IS6110</i> and <i>IS1081</i>
Extraction Control (EC)	HEX, VIC or JOE *	-

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A4.2 Reagents provided

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis complex</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis complex</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit with Ref. VS-MTC196TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *M. tuberculosis complex* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

M. tuberculosis complex Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *M. tuberculosis complex* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *M. tuberculosis complex* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *M. tuberculosis* complex Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tubes in the thermocycler.

3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis* complex) and HEX, JOE or VIC (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

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PRODUCTO	REFERENCIA
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172

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PRODUCTO	REFERENCIA
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTC101LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTC101HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113LE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113HE
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTC101E
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136E
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172E

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.