VIASURE

Real Time PCR Detection Kit

*M. tuberculosis complex +
non-tuberculous mycobacteria*

CE IVD

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EN For download IFUS from other languages, please enter in certest.es/viasure/labeling. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use.....	8
2.	Summary and Explanation	8
3.	Principle of the procedure	9
4.	Reagents provided	9
5.	Reagents and equipment to be supplied by the user	9
6.	Transport and storage conditions.....	10
7.	Precautions for users	10
8.	Test procedure	12
8.1.	Specimen collection, transport and storage	12
8.2.	DNA extraction	12
9.	Result interpretation	13
9.1.	References with internal control (references in Annex 1 and 2).....	13
9.2.	References with extraction control (references in Annex 3 and 4)	15
10.	Limitations of the test	18
11.	Quality control	19
12.	Performance characteristics.....	19
12.1.	Clinical sensitivity and specificity	19
12.2.	Analytical sensitivity	21
12.3.	Analytical specificity	22
12.4.	Analytical reactivity	23
	ANNEX 1	24
A1.1	Principle of the procedure.....	24
A1.2	Reagents provided.....	24
A1.3	Test procedure	26
A1.3.1	Lyophilized positive control.....	26
A1.3.2	PCR protocol.....	26
	ANNEX 2	27
A2.1	Principle of the procedure.....	27
A2.2	Reagents provided.....	27
A2.3	Test procedure	28
A2.3.1	Lyophilized positive control.....	28
A2.3.2	Lyophilized reaction mix tube.....	28
A2.3.3	PCR protocol.....	28
	ANNEX 3	30
A3.1	Principle of the procedure.....	30
A3.2	Reagents provided.....	30
A3.3	Test procedure	32
A3.3.1	Lyophilized extraction control	32
A3.3.2	Lyophilized positive control.....	32
A3.3.3	PCR protocol.....	32
	ANNEX 4	34
A4.1	Principle of the procedure.....	34
A4.2	Reagents provided.....	34

A4.3 Test procedure	35
A4.3.1 Lyophilized extraction control	35
A4.3.2 Lyophilized positive control.....	35
A4.3.3 Lyophilized reaction mix tube.....	35
A4.3.4 PCR protocol.....	35

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ENGLISH

1. Intended use

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative detection of DNA from the genus *mycobacterium*, the *M. tuberculosis* complex (MTBC) and the specie *M. tuberculosis* in clinical isolates of mycobacteria and decontaminated sputum (with positive or negative acid-fast bacilli (AFB) smear) from individuals suspected of Tuberculosis infection by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of the genus *mycobacterium*, the *M. tuberculosis* complex (MTBC) and/or the specie *M. tuberculosis* infection, in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, amplified using real time PCR, and detected using specific primers and a fluorescent reporter dye probe specific for the genus *mycobacterium*, the *M. tuberculosis* complex (MTBC) and/or the specie *M. tuberculosis*.

2. Summary and Explanation

Tuberculosis (TB) is a contagious, chronic and granulomatous disease caused mainly by the specie *M. tuberculosis*. This disease was declared in 1993 as a "global health emergency" due to its magnitude as a public health problem. Tuberculosis usually affects the lungs, but it can also affect other parts of the body such as the brain, kidneys or spine. Only people who suffer from pulmonary tuberculosis transmit the infection. Primary infection leads to active disease in approximately 10% of infected people, and the 80% of these cases, in the period of two years. In the remaining 90%, the immune system controls the infection, and the individual is not infectious or asymptomatic. It is estimated that one third of the world population has latent tuberculosis; that is, these people are infected with the bacillus, but (still) have not become ill or can transmit the infection. In this clinical state, TB bacilli can remain inactive for years (latent TB). However, when the immune system weakens, the latent infection can be reactivated. In a person infected with HIV, the risk of reactivation of latent TB is more than 10% per year, compared to a lifetime risk of 10-20% for HIV negative people.

Mycobacteria are a group of microorganisms that constitute one of the most serious health problems worldwide. Three groups can be defined within the genus *Mycobacterium*: 1) *Mycobacterium tuberculosis* complex that produces tuberculosis and is formed by the species *M. tuberculosis*, *M. bovis* (including *M. bovis* BCG), *M. africanum*, *M. microti*, *M. caprae*, *M. canettii*, *M. pinnipedii*, *M. mungi* and *M. suricattae*; 2) *M. leprae* that produces leprosy; 3) Other non-tuberculous mycobacteria (NTM) that are opportunistic and produce non-tuberculous symptoms with less pathogenic power. The isolation of NTM is increasingly frequent and its differentiation from the *M. tuberculosis* complex is of great clinical and public health importance, since it defines the isolation of patients in special rooms of health centres and the study of patient contacts.

The symptoms of active tuberculosis will depend on where mycobacteria develop within the body. The general symptoms of tuberculosis include fatigue, weight loss, fever and night sweats. Active pulmonary tuberculosis can cause breathing difficulties, chest pain and bloody expectorations. The symptoms of tuberculosis in other parts of the body depend on the area or organ affected. Over a year, a tubercular patient can infect 10 to 15 people through close contact. If they do not receive the appropriate treatment, up to two thirds of the tubercular patients die.

The early diagnosis of tuberculosis is essential for the control of the disease in order to interrupt the chain of transmission. The application of molecular techniques for diagnosis and typing, more sensitive, specific and quick than traditional tests, allow to improve the knowledge of the epidemiology of the infection and facilitate decisions for its control.

3. Principle of the procedure

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is designed for the qualitative detection of DNA from the genus *mycobacterium*, the *M. tuberculosis* complex (MTBC) and/or the specie *M. tuberculosis* in clinical strains and respiratory samples (decontaminated sputum). After DNA isolation, the detection of the genus *mycobacterium*, *M. tuberculosis* complex and/or *M. tuberculosis* specie is performed by the amplification of a conserved region of the *16S rRNA* gene (genus *mycobacterium*), the insertion sequences IS6110 and IS1081 (*M. tuberculosis* complex), and/or a fragment of the *TbD1* deletion region (*M. tuberculosis* specie), using specific primers and fluorescent-labeled probes.

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal, which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer, polymerase) in a stabilized format, as well as an internal control to verify the correct functioning of the amplification mix.

4. Reagents provided

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products, Annex 2 for tube format with internal control products, Annex 3 for open and rotor-gene format with extraction control products and Annex 4 for tube format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1 and 3).

- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or $\leq -70^{\circ}\text{C}$.
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).

VIASURE *M. tuberculosis complex + non-tuberculous mycobacteria* Real Time PCR Detection Kit has been validated on the following equipments: 7500 Fast Real-Time PCR System (Applied Biosystems), CFX96™ Real-Time PCR Detection System (Bio-Rad), AriaMx Real-Time PCR System (Agilent Technologies), DTprime Real-time Detection Thermal Cycler (DNA-Technology), DTlite Real-Time PCR System (DNA-Technology), Rotor-Gene® Q (Qiagen), Cobas z480 Analyzer (Roche Molecular Diagnostics) and VIASURE V-Lab96 Cycler (CerTest Biotec S.L). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *M. tuberculosis complex + non-tuberculous mycobacteria* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival

- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-MTD113L, VS-MTD113H, VS-MTD136, VS-MTD172, VS-MTD113LE, VS-MTD113HE, VS-MTD136E and VS-MTD172E). Remove any excess air in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-MTD101, VS-MTD136, VS-MTD172, VS-MTD101E, VS-MTD136E and VS-MTD172E (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *M. tuberculosis complex + non-tuberculous mycobacteria* Real Time PCR Detection Kit and any additional reagents or equipments required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions, and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format with internal control products Test Procedure, Annex 2 for Tube format with internal control products Test Procedure, Annex 3 for Open and Rotor-Gene format with extraction control products Test Procedure and Annex 4 for Tube format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *M. tuberculosis complex + non-tuberculous mycobacteria* Real Time PCR Detection Kit has been tested in mycobacterial culture, clinical isolates of mycobacteria and decontaminated sputum (with positive or negative acid-fast bacilli (AFB) smear). Other types of samples must be validated by the user.

Collection, storage and transport of specimens should be maintained per the conditions validated by the user. Overall, clinical samples should be collected and labelled appropriately in clean containers and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at Room Temperature (RT) for up to 2 hours, or at 4°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), It is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 4°C for up to 24 hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>), the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94) and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cencenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

For sputum samples, the recommended amount of sputum is approximately 200µL that can be pipetted. * In order to improve the yield and quality of bacterial DNA from sputum specimens, it is recommended a sample pre-treatment with N-acetyl-L-cysteine–sodium hydroxide (NALC-NaOH). Consult DNA Extraction section for more information.

8.2. DNA extraction

Prepare the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

Due to VIASURE *M. tuberculosis complex + non-tuberculous mycobacteria* Real Time PCR Detection Kit is also available with extraction control (references include in Annexes III and IV), if the Extraction Control (EC) is used to monitor nucleic acid isolation and as PCR inhibition control, add 5µl of the reconstituted EC to the specimen

and/or lysis buffer mixture (clinical specimen, as well as positive control and/or negative control). Close each tube and vortex for 10 seconds. If the Extraction Control is used only as a PCR inhibition control, 1 µl of the reconstituted EC should be added to the reconstituted Reaction-Mix.

For DNA extraction from mycobacterial culture, clinical isolates of mycobacteria and decontaminated sputum*, you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

* In order to improve the yield and quality of bacterial DNA from sputum specimens, it is recommended a sample pre-treatment with N-acetyl-L-cysteine–sodium hydroxide (NALC-NaOH). After decontamination, the cell pellets should be resuspended in a maximum final volume of 1 mL of phosphate buffer. In addition, usage of small elution volumes (50-100 µL) may raise the DNA concentration.

9. Result interpretation

9.1. References with internal control (references in Annex 1 and 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control (IC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for the genus *mycobacterium*, the *M. tuberculosis* complex (MTBC) and the specie *M. tuberculosis* in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>M. tuberculosis</i> complex (FAM) ¹	Genus <i>mycobacterium</i> (ROX) ¹	<i>M. tuberculosis</i> (Cy5) ¹	Internal Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	>40 or no signal	>40 or no signal	≤40	Valid

Table 1. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

Note that to consider a positive result, the fluorescence level must be at least 20% of the total fluorescence value.

For interpretation of individual patient sample results, use the following table, read and analyze the results:

<i>M. tuberculosis</i> complex (FAM)	Genus <i>mycobacterium</i> (ROX)	<i>M. tuberculosis</i> (Cy5)	Internal control (HEX)	Interpretation for patient's individual samples	
≤40	≤40	≤40	≤40 or no signal ¹	Valid	<i>M. tuberculosis</i> specie Positive (If Ct value FAM < Ct value ROX) ³ or NTM Positive (If Ct value ROX channel < Ct value FAM and Cy5 channels) ³
≤40	>40 or no signal [#]	≤40	≤40 or no signal ¹	Valid	<i>M. tuberculosis</i> specie Positive (If Ct value FAM > 35) ³ * If Ct value FAM < 35, please repeat the assay: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample. After retesting one time, if the result holds, the sample should be considered <i>M. tuberculosis</i> specie positive.
≤40	>40 or no signal [#]	>40 or no signal [#]	≤40 or no signal ¹	Valid	MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct value FAM > 35) ³ * If Ct value FAM < 35, please repeat the assay: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered MTBC strains different from <i>M. tuberculosis</i> positive.
≤40	≤40	>40 or no signal	≤40 or no signal ¹	Valid	MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct value FAM = Ct value ROX and/or Ct value FAM < Ct value ROX) ³ or NTM Positive (If Ct value FAM > Ct value ROX) ³
>40 or no signal	≤40 [†]	>40 or no signal	≤40 or no signal ¹	Valid	NTM positive (In case Cy5 amplifications appear with a Ct > 37, the sample should be considered as NTM Positive).
>40 or no signal	>40 or no signal	>40 or no signal	≤ 35 ²	Valid	Targets not Detected ²
>40 or no signal	>40 or no signal	>40 or no signal	> 35 or no signal ²	Invalid	Test Failure – Repeat Testing²

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve. MTBC= *M. tuberculosis* complex, NTM = non-tuberculous mycobacteria.

[#]: Because of the amplified targets, FAM signal is stronger than ROX or ROX and Cy5 signals. It could be observed preferential amplification curve in this channel and absence in ROX or ROX and Cy5 channels. This result also indicates the presence of MTBC strain different from *M. tuberculosis* or *M. tuberculosis* specie.

[†]: Please be aware that amplification curves with later Ct may appear in ROX channel due to environmental NTM.

1 The internal control (IC) shows or not an amplification signal ($C_t \leq 40$ or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of *genus mycobacterium*, the *M. tuberculosis* complex (MTBC) and the specie *M. tuberculosis* target genes negative, IC must show an amplification signal with $C_t \leq 35$. If there is an absence of signal or C_t value > 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.


3 Only the integer value of C_t will be considered. If the C_t value given by the PCR equipment has decimals, these will not be considered for the interpretation.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate.

Depending on the available material:

- repeat qPCR with the same isolated DNA sample, or,
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

9.2. References with extraction control (references in Annex 3 and 4)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run validates the reaction by checking the absence of signal in the negative control well and the presence of signal for the *genus mycobacterium*, the *M. tuberculosis* complex

(MTBC) and the specie *M. tuberculosis* in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	<i>M. tuberculosis</i> complex (FAM) ¹	Genus <i>mycobacterium</i> (ROX) ¹	<i>M. tuberculosis</i> (Cy5) ¹	Extraction Control (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	>40 or no signal	>40 or no signal	≤40	Valid

Table 3. Expected Performance of Controls

1. In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.
2. The Extraction Control (EC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

Note that to consider a positive result, the fluorescence level must be at least 20% of the total fluorescence value. For interpretation of individual patient sample results, use the following criteria for read and analyze the results:

<i>M. tuberculosis</i> complex (FAM)	Genus <i>mycobacterium</i> (ROX)	<i>M. tuberculosis</i> (Cy5)	Extraction Control (HEX)	Interpretation for patient's individual samples	
≤40	≤40	≤40	≤40 or no signal ¹	Valid	<i>M. tuberculosis</i> specie Positive (If Ct value FAM < Ct value ROX) ³ or NTM Positive (If Ct value ROX channel < Ct value FAM and Cy5 channels) ³
≤40	>40 or no signal [#]	≤40	≤40 or no signal ¹	Valid	<i>M. tuberculosis</i> specie Positive (If Ct value FAM > 35) ³ * If Ct value FAM < 35, please repeat the assay: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered <i>M. tuberculosis</i> specie positive.
≤40	>40 or no signal [#]	>40 or no signal [#]	≤40 or no signal ¹	Valid	MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct value FAM > 35) ³ * If Ct value FAM < 35, please repeat the assay: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another

					aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered <i>MTBC</i> strains different from <i>M. tuberculosis</i> positive.
≤40	≤40	>40 or no signal	≤40 or no signal ¹	Valid	MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct value FAM = Ct value ROX and/or Ct value FAM < Ct value ROX)) ³ or NTM Positive (If Ct value FAM > Ct value ROX) ³
>40 or no signal	≤40 [†]	>40 or no signal	≤40 or no signal ¹	Valid	NTM positive (In case Cy5 amplifications appear with a Ct > 37, the sample should be considered as NTM Positive).
>40 or no signal	>40 or no signal	>40 or no signal	≤ 35 ²	Valid	Targets not Detected ²
>40 or no signal	>40 or no signal	>40 or no signal	> 35 or no signal ²	Invalid	Test Failure – Repeat Testing²

Table 4. Interpretation of individual patient sample results. Ct values, no signal = no amplification curve. MTBC= *M. tuberculosis* complex, NTM = non-tuberculous mycobacteria.

#: Because of the amplified targets, FAM signal is stronger than ROX or ROX and Cy5 signals. It could be observed preferential amplification curve in this channel and absence in ROX or ROX and Cy5 channels. This result also indicates the presence of MTBC strain different from *M. tuberculosis* or *M. tuberculosis* specie.

†: Please be aware that amplification curves with later Ct may appear in ROX channel due to environmental NTM.

1 The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

2 In the case of *genus mycobacterium*, the *M. tuberculosis* complex (MTBC) and the specie *M. tuberculosis* target gene negative, EC must show an amplification signal with Ct ≤35. If there is an absence of signal or Ct value >35 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.


3 Only the integer value of Ct will be considered. If the Ct value given by the PCR equipment has decimals, these will not be considered for the interpretation.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate.

Depending on the available material:

- a) repeat qPCR with the same isolated DNA sample, or
- b) re-extract and retest another aliquot of the same specimen or,
- c) obtain a new specimen and retest.

Figure 2. Correct run of negative and positive control run on the CFX96™ Real-Time PCR Detection System (Bio-Rad).

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from mycobacterial culture, clinical isolates of mycobacteria and decontaminated sputum (with positive or negative acid-fast bacilli (AFB) smear).
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or the specie *M. tuberculosis*, either by samples containing high concentrations of target DNA, or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *M. tuberculosis* complex + non-tuberculous mycobacteria Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in separate laboratory areas.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the bacterial DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown species of the genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or *M. tuberculosis*.
 - A bacterial load in the specimen below the limit of detection for the assay.

- The presence of qPCR inhibitors or other types of interfering substances. The impact of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection by the genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or *M. tuberculosis* or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable bacteria and does not imply that these bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets sequences (16S rRNA gene, the insertion sequences IS6110 and IS1081, and/or a fragment of the *TbD1* deletion region).
- Negative results do not preclude genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or *M. tuberculosis* infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak bacterial levels during infections caused by genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or *M. tuberculosis* have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that genus *mycobacterium*, *M. tuberculosis* complex (MTBC) and/or *M. tuberculosis* infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or the Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit was tested using clinical strains and the remnant of respiratory specimens (decontaminated sputum) from patients with symptoms of respiratory infection. In order to determine the clinical diagnostic accuracy, an evaluation has been conducted through collaboration with a national entity. A summary of the site, sample type and workflow is included in the following table.

	Site	Sample type	Workflow	Target
1	Institut d'Investigació en Ciències de la Salut Germans Trias i Pujol (IGTP), Badalona, Barcelona, Spain	clinical strains	Seegene DNA extraction protocol + CFX96™ Touch Real-Time PCR Detection System (Bio-Rad).	<i>Mycobacterium</i> spp. (NTM)
				<i>M. tuberculosis</i> complex (MTBC)
				<i>M. tuberculosis</i> (MTB)
	Institut d'Investigació en Ciències de la Salut Germans Trias i Pujol (IGTP), Badalona, Barcelona, Spain	sputums	Seegene DNA extraction protocol + CFX96™ Touch Real-Time PCR Detection System (Bio-Rad).	<i>Mycobacterium</i> spp. (NTM)
				<i>M. tuberculosis</i> complex (MTBC)
				<i>M. tuberculosis</i> (MTB)
2	Certest Biotec S.L (San Mateo de Gállego, Spain)	bronchoalveolar lavages (BAL)	MagDEA® Dx Nucleic Acid Extraction kit using the magLEAD® 12gC (Precision System Science Co) + LightCycler® 480 Real-Time PCR System (Roche).	<i>M. tuberculosis</i> complex (MTBC)
		bronchoalveolar aspirates (BAS)		<i>Mycobacterium</i> spp. (NTM)
		sputum		<i>M. tuberculosis</i> (MTB)

Table 5. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following tables:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Culture	NTM	14	38	0	0	1 (0.76-1)	1 (0.9-1)	1 (0.76-1)	1(0.9-1)
		MTBC	15	37	0	0	1(0.78-1)	1(0.9-1)	1(0.78-1)	1(0.9-1)
		MTB	24	27	0	1	0.96 (0.79-0.99)	1 (0.87-1)	1 (0.87-1)	0.96 (0.79-0.99)
	Reference method*	NTM	22	91	9	8	0.73 (0.54-0.87)	0.91(0.83-0.95)	0.7(0.51-0.85)	0.91(0.84-0.96)
		MTBC	52	77	0	1	0.98(0.89-1)	1(0.95-1)	1(0.91-1)	0.98 (0.92-0.99)
		MTB	47	78	1	4	0.92(0.81-0.97)	0.98 (0.93-1)	0.97 (0.87-0.99)	0.95(0.87-0.98)
2	Amplified Mycobacterium Tuberculosis Direct Test (Hologic) + smear and culture	MTBC	2	325	0	0	1 (0.15-1)	1 (0.98 – 1)	1 (0.34-1)	1 (0.99 – 1)
		NTM	5	322	0	0	1 (0.47 – 1)	1 (0.98 – 1)	1 (0.65 – 1)	(0.99 – 1)
		MTB	0	327	0	0	NA	1 (0.98 – 1)	NA	1 (0.99 – 1)

Table 6. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit. MTBC = *M. tuberculosis* complex, MTB = *M. tuberculosis*, NTM = non-tuberculosis mycobacteria.

* Characterization performed with culture, and in addition, in patients with respiratory NTM infection, NTM species were identified with INNO-LiPA MYCOBACTERIA v2 (Fujirebio), and patients with TB diagnosis, MTBC identification was performed with the SD BIOLINE TB Ag MPT64 Rapid Test (Abbott).

Results show high agreement to detect targets using VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit has a detection limit of 10 DNA copies per reaction for the insertion sequences IS6110 and IS1081, the 16S rRNA gene and the *TbD1* deletion region, with a positive rate of $\geq 95\%$.

Figure 3. Dilution series of insertion sequences IS6110 and IS1081 (10^7 - 10^1 copies/rxn) templates run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel FAM).


Figure 4. Dilution series of 16S rRNA gene (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel ROX).



Figure 5. Dilution series of *TbD1* deletion region (10^7 - 10^1 copies/rxn) template run on the CFX96™ Real-Time PCR Detection System (Bio-Rad) (channel Cy5)

12.3. Analytical specificity

The specificity of the VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit was confirmed by testing a panel consisting of different microorganisms associated to respiratory diseases. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	Influenza A/Victoria/210/2009(H3N2) virus	-	<i>Legionella longbeache</i>	-
<i>Bordetella parapertussis</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Moraxella catarrhalis</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/South Australia/55/2014, IVR-175)	-	<i>Mycoplasma pneumoniae</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Human Adenovirus 1-5, 8, 15, 31, 40 and 41	-
<i>Chlamydia caviae</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus (clade 3C2a.1)	-	Human Bocavirus	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus	-	Human coronavirus 229E, OC43 and NL63	-
<i>Chlamydia psittaci</i> genotypes A and C	-	Influenza B/Brisbane/60/2008 virus	-	Human metapneumovirus A and B	-
<i>Haemophilus influenzae</i> MinnA	-	Influenza B/Florida/04/06 virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
Influenza A/Anhui/1/2013 (H7N9) virus	-	Influenza B/Phuket/3073/2013 virus	-	Human rhinovirus	-
Influenza A/California/7/2009(H1N1)pdm09 virus	-	<i>Legionella bozemanii</i>	-	MERS Coronavirus	-
Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Legionella micdadei</i>	-	<i>Streptococcus pneumoniae</i> Z022	-
Influenza A/Hong Kong/4801/2014, NYMC X-263B	-	<i>Legionella dumoffii</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i>	-	Respiratory Syncytial virus (RSV) type A and B	-
Influenza A/New Caledonia/20/99(H1N1) virus	-				

Table 7. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit for the specie *M. tuberculosis* was evaluated using DNA extracted from reference strains (*Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331 and *Mycobacterium tuberculosis* X004439) and from clinical strains containing *M. tuberculosis* lineages 1, 2, 3, and lineage 4 (T sublineage, LAM sublineage, and X sublineage) as templates, showing positive results.

The reactivity of VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit for the *Mycobacterium* genus was evaluated using DNA extracted from reference strains (*Mycobacterium bovis* AF 2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium africanum*, *Mycobacterium microti*, *Mycobacterium abscessus*, *Mycobacterium avium* subsp. *avium*, *Mycobacterium chelonae*, *Mycobacterium fortuitum*, *Mycobacterium gordonae*, *Mycobacterium kansasii*, *Mycobacterium marinum*, *Mycobacterium peregrinum*, *Mycobacterium smegmatis*, *Mycobacterium scrofulaceum*, *Mycobacterium xenopi*), from clinical strains (*Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), *Mycobacterium africanum* lineage 5 and 6, *Mycobacterium avium* complex (MAC), *Mycobacterium avium* subsp. *paratuberculosis*, *Mycobacterium avium* subsp. *silvaticum*, *Mycobacterium chelonae* complex (group III-*M. abscessus*), *Mycobacterium intracellulare*) and from clinical samples (*Mycobacterium celatum*,, *Mycobacterium chimaera-intracellulare* as templates, showing positive results.

The reactivity of VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit for *M. tuberculosis* complex strains was evaluated using DNA extracted from reference strains (*Mycobacterium bovis* AF 2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium africanum* and *Mycobacterium microti*, as well as from clinical strains (*Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), and *Mycobacterium africanum* lineage 5 and 6) as templates, showing positive results.

ANNEX 1

OPEN FORMAT AND ROTOR-GENE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTD106L
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTD106H
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTD112L
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTD112H
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 96-well plate, low profile	VS-MTD113L
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 96-well plate, high profile	VS-MTD113H
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTD136
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTD172
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTD101L
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTD101H
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTD101

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>M. tuberculosis</i> complex	FAM	Insertion sequences IS6110 and IS1081
<i>Mycobacterium</i> spp.	ROX	16S rRNA gene
<i>Mycobacterium tuberculosis</i>	Cy5	TbD1 deletion region
Internal control (IC)	HEX, VIC or JOE *	-

Table A1 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.cerfetest.es.

A1.2 Reagents provided

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible

devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

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Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref VS-MTD113L and VS-MTD113H.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	2/9/18 X 4-cap strip

Table A1 5. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD101, VS-MTD136 and VS-MTD172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

M. tuberculosis complex + non-tuberculous mycobacteria Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.cerTEST.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 6. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis complex*), ROX (*Mycobacterium spp.*), Cy5 (*M. tuberculosis*) and HEX, JOE or VIC channels (Internal Control (IC)). (Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.cerTEST.es)).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTD196T

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>M. tuberculosis</i> complex	FAM	Insertion sequences IS6110 and IS1081
<i>Mycobacterium</i> spp.	ROX	16S rRNA gene
<i>Mycobacterium tuberculosis</i>	Cy5	<i>TbD1</i> deletion region
Internal control (IC)	HEX, VIC or JOE *	-

Table A2. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.cerest.es.

A2.2 Reagents provided

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD196T.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

M. tuberculosis complex + non-tuberculous mycobacteria Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *M. tuberculosis complex + non-tuberculous mycobacteria* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *M. tuberculosis complex + non-tuberculous mycobacteria* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate, the strips, or tubes in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis* complex), ROX (*Mycobacterium* spp.), Cy5 (*M. tuberculosis*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es).

ANNEX 3

OPEN FORMAT AND ROTOR-GENE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTD106LE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTD106HE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTD112LE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTD112HE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 96-well plate, low profile	VS-MTD113LE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 96-well plate, high profile	VS-MTD113HE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTD136E
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTD172E
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-MTD101LE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-MTD101HE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit 2 x 4-well strips, Rotor-Gene®	VS-MTD101E

Table A3. 1. References.

A3.1 Principle of the procedure

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>M. tuberculosis</i> complex	FAM	Insertion sequences IS6110 and IS1081
<i>Mycobacterium</i> spp.	ROX	16S rRNA gene
<i>Mycobacterium tuberculosis</i>	Cy5	<i>TbD1</i> deletion region
Extraction Control (EC)	HEX, VIC or JOE *	-

Table A3. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.certest.es.

A3.2 Reagents provided

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A3.3, A3.4 and A3.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A3.3 includes materials and reagents to be used with 8-well strips

compatible devices. Table A3.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A3.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.cerTEST.es).

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

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Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A3. 4. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref VS-MTD113LE and VS-MTD113HE.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis complex + non-tuberculous mycobacteria</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis complex + non-tuberculous mycobacteria</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A3. 5. Reagents and materials provided in VIASURE *M. tuberculosis complex + non-tuberculous mycobacteria* Real Time PCR Detection Kit with Ref. VS-MTD101E, VS-MTD136E and VS-MTD172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A3.3 Test procedure

A3.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control in an area away from the other components.

A3.3.2 Lyophilized positive control

M. tuberculosis complex + non-tuberculous mycobacteria Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A3.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *M. tuberculosis* complex + *non-tuberculous mycobacteria* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A3. 6. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis* complex), ROX (*Mycobacterium* spp.), Cy5 (*M. tuberculosis*) and HEX, JOE or VIC channels (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es).

ANNEX 4

TUBE FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-MTD196TE

Table A4. 1. References.

A4.1 Principle of the procedure

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Target	Channel	Gene
<i>M. tuberculosis</i> complex	FAM	Insertion sequences IS6110 and IS1081
<i>Mycobacterium</i> spp	ROX	16S rRNA gene
<i>Mycobacterium tuberculosis</i>	Cy5	<i>TbD1</i> deletion region
Extraction Control (EC)	HEX, VIC or JOE *	-

Table A4. 2. Target, channel and genes.

*Depending on the equipment used select the proper detection channel, to check most common detection channels consult website www.ceratest.es.

A4.2 Reagents provided

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A4.3.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A4. 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD196TE.

A4.3 Test procedure

A4.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used to reconstitute the lyophilized *M. tuberculosis complex* + *non-tuberculous mycobacteria* Positive Control in an area away from the other components.

A4.3.2 Lyophilized positive control

M. tuberculosis complex + *non-tuberculous mycobacteria* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex* + *non-tuberculous mycobacteria* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A4.3.3 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *M. tuberculosis complex* + *non-tuberculous mycobacteria* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A4.3.4 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *M. tuberculosis* complex + *non-tuberculous mycobacteria* Reaction-Mix (white vial) into each tube.

2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of DNA sample in different wells.

Add 5 µL of reconstituted *M. tuberculosis* complex + *non-tuberculous mycobacteria* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tubes in the thermocycler.

3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A4. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*M. tuberculosis* complex), ROX (*Mycobacterium spp.*), Cy5 (*M. tuberculosis*) and HEX, JOE or VIC channels (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (to check the most common detection channels consult website www.certest.es).

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.