

γ -GT (GGT)

L- γ -Glutamyltransferase Manual

INTENDED USE

For the quantitative *in vitro* determination of γ -glutamyl transferase (γ -GT) in serum or plasma. This product is suitable for Manual use.

Cat. No.

GT 2750	R1a. Buffer/Glycylglycine	1 x 70 ml
20 x 3 ml	R1b. Substrate	20 x 3 ml

GTIN: 05055273203431

CLINICAL SIGNIFICANCE⁽¹⁾


Gamma-glutamyltransferase (γ -GT) in serum originates primarily from the hepatobiliary system. Therefore γ -GT is elevated in all forms of liver disease and has been shown to be more sensitive than alkaline phosphatase in detecting obstructive jaundice, cholangitis and cholecystitis. High levels of γ -GT are also seen in patients with primary or secondary liver cancer. Increased levels are also observed in cases of alcohol abuse and in alcoholic liver cirrhosis. In patients receiving anticonvulsant drugs such as phenytoin and phenobarbital, increased levels of the enzyme in serum may reflect induction of new enzyme activity and the toxic effects of alcohol and other drugs on the microsomal structures in liver cells. Although γ -GT is the most sensitive enzymatic indicator of hepatobiliary disease it can not be used to differentiate between different types of hepatobiliary disease. However, γ -GT can be used in combination with other biochemical markers to discriminate between different types of hepatobiliary disease.

PRINCIPLE⁽²⁾

COLORIMETRIC METHOD

The substrate L- γ -glutamyl-3-carboxy-4-nitroanilide, in the presence of glycylglycine is converted by γ -GT in the sample to 5-amino-2-nitrobenzoate which can be measured at 405nm.

L- γ -glutamyl-3-carboxy-4-nitroanilide + glycylglycine



L- γ -glutamylglycylglycine + 5-amino-2-nitrobenzoate

SAMPLE COLLECTION AND PREPARATION

Serum: use only non-haemolysed serum.

Plasma: use only EDTA plasma that is free from haemolysis.

Other anticoagulants interfere with this test.

γ -GT in serum and plasma is stable for 7 days at +2 to +8°C or for 3 months at -20°C.

REAGENT COMPOSITION

Contents	Concentration in the Test
R1a. Buffer/Glycylglycine	
Tris buffer	100 mmol/l, pH 8.25
Glycylglycine	100 mmol/l
R1b. Substrate	
L- γ -glutamyl-3-carboxy-4-nitroanilide	2.9 mmol/l

SAFETY PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Solution R1a and R1b contain Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

Health and Safety Data Sheets are available on request.

Please dispose of all biological and chemical materials according to local guidelines.

The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.

STABILITY AND PREPARATION OF REAGENTS

R1a. Buffer/Glycylglycine

Contents stable as supplied up to the expiry date when stored at +2 to +8°C.

R1b. Substrate

Reconstitute one vial of Substrate R1b with 3.0 ml of Buffer/Glycylglycine R1a.

Stable for 21 days at +2 to +8°C or 5 days at +15 to +25°C.

PROCEDURE

MATERIALS PROVIDED

Buffer/Glycylglycine
Substrate

MATERIALS REQUIRED BUT NOT PROVIDED

Randox Assayed Multi-sera Level 2 (Cat. No. HN 1530) and Level 3 (Cat. No. HE 1532).
Pipetting device to deliver 100 μ l and 1 ml.
Timing device and water bath or heating block to maintain temperature at 25, 30 or 37°C.
Spectrophotometer with wavelength capability of 400 - 420 nm.

TEST METHOD

Wavelength:	Hg 405 nm (400 - 420 nm)
Cuvette:	1 cm light path
Temperature:	25°C, 30°C, 37°C
Measurement:	against air

Pipette into cuvette:

Sample	0.10 ml
Reagent (25°C, 30°C, 37°C)	1.00 ml

Mix, read initial absorbance and start timer simultaneously.
Read again after 1, 2 and 3 min.

CALCULATION

To calculate the GGT activity use the following formula.

$$U/L = 1158 \times \Delta A \text{ 405 nm/min}$$

STANDARDISATION

Randox γ -GT Standard included in the kit is traceable to γ -GT reference materials AD452 (IFCC) and JSCC TS01.

QUALITY CONTROL

Randox Assayed Multi-sera, Level 2 and Level 3 are recommended for daily quality control. Two levels of controls should be assayed at least once a day. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error, the following steps should be taken:

1. Check instrument settings and light source.
2. Check cleanliness of all equipment in use.
3. Check water, contaminants i.e. bacterial growth may contribute to inaccurate results.
4. Check reaction temperature.
5. Check expiry date of kit and contents.
6. Contact Randox Laboratories Customer Technical Support, Northern Ireland (028) 94422413.

INTERFERENCE

Avoid haemolysis as it interferes with the assay.

NORMAL VALUES IN SERUM⁽³⁾

	25°C	30°C	37°C
Men	6-28 U/l	8-38 U/l	11-50 U/l
Women	4-18 U/l	5-25 U/l	7-32 U/l

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

SPECIFIC PERFORMANCE CHARACTERISTICS

The following performance characteristics were obtained using a Cobas Mira analyser at 37°C. Users are advised to establish their own performance data for other automated analysers.

LINEARITY

The method is linear to 400 U/l. If the absorbance change per minute exceeds 0.2 dilute 0.1 ml of sample with 0.5 ml of 0.9% (w/v) NaCl solution and reassay. Multiply the result by 6.

SENSITIVITY

The minimum detectable level has been determined, by replicate dilution of a suitable control, as 7U/l.

PRECISION

Within Run

	Level 2	Level 3
Mean (U/l)	52.3	185.1
SD	0.66	1.39
CV (%)	1.26	0.95
n	20	20

Between Run

	Level 2	Level 3
Mean (U/l)	54.7	189.5
SD	1.26	2.74
CV (%)	2.3	1.5
n	20	20

METHOD COMPARISON

The Randox method (Y) was compared to another commercially available method (X). Forty patient samples, with values spanning the range 8 to 154 U/l were tested. Linear regression analysis of the data resulted in the following equation:

$$Y = 0.979x + 0.765$$

with a correlation coefficient of $r = 0.998$.

REFERENCES

1. Teitz, N. N., Fundamentals of Clinical Chemistry ed 3 Philadelphia, W B Saunders Co. 1987, pg 391.
2. Szasz, G., Clin Chem. 1969; **22**: 124-136.
3. Szasz, G., Bergmeyer H.U., ed. Method of Enzymatic analysis. Weinheim Verlag Chemie, 1974.

Revised 05 May 16 bi
Rev. 002