VIASURE

Real Time PCR Detection Kit

Gastrointestinal Panel I

CE

IVD

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EN For download IFUS from other languages, please enter in certest.es/viasure/labeling. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use.....	5
2.	Summary and Explanation	5
3.	Principle of the procedure	7
4.	Reagents provided	8
5.	Reagents and equipment to be supplied by the user	8
6.	Transport and storage conditions.....	8
7.	Precautions for users	9
8.	Test procedure	10
9.	Result interpretation	11
10.	Limitations of the test	18
11.	Quality control.....	19
12.	Performance characteristics.....	19
	ANNEX 1	23
A1.1	Principle of the procedure.....	23
A1.2	Reagents provided.....	24
A1.3	Test procedure.....	24
	ANNEX 2	26
A2.1	Principle of the procedure.....	26
A2.2	Reagents provided.....	27
A2.3	Test procedure	27

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ENGLISH

1. Intended use

VIASURE Gastrointestinal Panel I Real Time PCR Detection Kit is a real-time RT-PCR test designed for the specific and qualitative detection of DNA/RNA from *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica*, Norovirus GI and/or Norovirus GII (NoV), Rotavirus (RV), Adenovirus (AdV), Astrovirus (AstV) and/or Sapovirus (SaV) in human faeces samples from patients with signs and symptoms of gastrointestinal infection by their healthcare professional (HCP). The test is intended for use as an aid in the diagnosis of *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV in combination with clinical and epidemiological risk factors.

DNA/RNA is extracted from clinical samples, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV.

2. Summary and Explanation

Viruses, bacteria and/or parasites can cause gastroenteritis, an inflammation of the gastrointestinal tract that affects both the stomach and the small intestine. Gastrointestinal infections are self-limited and resolved within a few days. However, in a healthcare setting and in specific populations (newborns/infants, immunocompromised patients, or older populations), are potentially serious. Therefore, rapid diagnosis, appropriate treatment, and infection control measures are particularly important in these settings.

The genus *Salmonella* is taxonomically divided into 6 subspecies of *Salmonella enterica* and *Salmonella bongori*. *Salmonella* causes two types of illnesses: enteric fever (typhoid) and acute gastroenteritis, commonly called salmonellosis. *Salmonella* transmission occurs through contaminated food (poultry, poultry products, beef, pork, eggs, milk, and shellfish), water, or through contact with infected animals. *Salmonella*-infected patients frequently experience nausea, vomiting, abdominal pain, muscle aches, fever, severe diarrhoea, and/or headaches.

The genus *Campylobacter* consists of 26 species, 2 provisional species, and 9 subspecies. The *Campylobacter* species most commonly associated with human diseases are *C. jejuni* and *C. coli*, but there are other species than can also cause human infections. Poultry is the largest source of transmission of *Campylobacter* to humans. Other risk factors include consumption of animal products and water, contact with animals, and even person-to-person transmission (faecal-oral route through vomiting). *Campylobacter* infection causes gastroenteritis characterised by fever, vomiting, headaches, head and abdominal pain with watery or bloody diarrhoea, averaging 6 days.

The genus *Yersinia* includes three pathogenic species in humans and animals: *Yersinia enterocolitica*, *pestis* and *pseudotuberculosis*. There are six biotypes of *Y. enterocolitica*; five of which are considered human pathogens (biotypes 1B, 2, 3, 4 and 5). In addition, there are 60 to 70 serotypes, among which O:3, O:9, O:8, O:5. 27 are associated with disease in humans. *Yersinia enterocolitica* is a foodborne pathogen and its clinical manifestation often include nausea, vomiting, abdominal pain, diarrhoea, and fever. There are clear indications that animal foods, especially pork and dairy products, are responsible for human infections.

Shigellosis is a type of bacillary dysentery characterised by severe diarrhoea that contains blood and mucus. The disease is caused by any of the four *Shigella* species (*S. dysenteriae*, *S. flexneri*, *S. boydii* and *S. sonnei*), by enteroinvasive *Escherichia coli* (EIEC). In fact, the differentiation of *Shigella* and *E. coli* enteroinvasive is complicated due to the ability to cause dysentery using the same invasive method. Shigellosis is a frequently foodborne illness, especially in foods that require processing and / or are prepared from raw or previously cooked products without reheat. The low infective dose (10 cells) allows the disease to spread effectively through infected food or water, and through human-to-person contact.

Cryptosporidium can be transmitted by several routes: person-to-person contact and with companion and / or farm animals, as well as ingestion of contaminated food, drinking water and water for recreational use. *Cryptosporidium* can be present asymptotically or cause diarrhea and abdominal discomfort with weight loss and malabsorption. *Giardia lamblia* infection occurs by fecal-oral transmission, through direct contact or by ingestion of contaminated food or water. Diarrhea is the main symptom of this disease, but the clinical manifestations vary, from asymptomatic infections to acute diarrhea, abdominal cramps, bloating and flatulence often accompanied by nausea and weight loss. *Entamoeba histolytica* disease is spread mainly by water or food contaminated with cysts but can also be transmitted from person to person through the fecal-oral route. The clinical features of amoebiasis range from asymptomatic colonization (90% of cases) to amoebic dysentery (stomach pain, nausea or vomiting, and severe hemorrhagic diarrhea and fever) or illness invasive extraintestinal, most commonly in the form of liver abscesses.

Norovirus (NoV) belongs to the *Caliciviridae* family and is considered the leading cause of acute non-bacterial gastroenteritis in all age groups worldwide. NoVs are classified into six distinct genogroups (GI-GVI) based on the phylogenetic tree inferred by comparison of the complete peptide sequence of VP1. In particular, the GI, GII, and GIV genogroups infect humans. After an incubation period of 12 to 48 hours, Norovirus disease is characterised by projectile vomiting, non-bloody diarrhoea, nausea, abdominal cramps, and low-grade fever. Noroviruses can infect humans through multiple routes, including the oral route, transmitted through contact with fecal matter or by aerosolization of the vomit of infected people, as well as contaminated surfaces, food, or water.

Rotavirus (RV) infection is the world's leading cause of severe gastroenteritis in infants and children under 5 years of age. The genus *Rotavirus* belongs to the *Reoviridae* family. Rotaviruses can be divided into 7 distinct groups (A-G). Rotavirus A-C serogroups are human pathogens, although Rotavirus A produces more than 90% of infections. The main route of transmission is fecal-oral due to contamination of food, water, or objects. Characteristic symptoms are vomiting, watery diarrhoea, and abdominal pain.

Adenoviruses belong to the *Adenoviridae* family of unencapsulated double-stranded DNA (dsDNA) viruses. There are more than 50 immunologically differentiated serotypes of human Adenoviruses grouped into 6 species (HADV-A to HADV-F) that can cause infections in humans ranging from respiratory diseases (mainly HADV-B and C), and conjunctivitis (HADV-B. and D), up to gastroenteritis (AdVH-F serotypes 40 and 41). In addition to the latter, AdVH serotypes 31, 12, 18, 1, 2, 5 and 6 have also been associated with acute diarrhea. Although the epidemiological characteristics of Adenovirus vary by type, all are transmitted by direct contact or by the fecal-oral route and are considered pathogens for all age groups of the disease. population regardless of age, especially for children. In fact, these viruses are highly infectious, and a few viral particles are enough to cause infection.

Human Astroviruses (HAstVs) belong to the *Astroviridae* family. The 8 classic HAstVs serotypes represent 2-9% of all acute non-bacterial gastroenteritis in children worldwide, although infections have also been reported in healthy adults, the elderly, and immunocompromised patients. Typically, HAstV infection induces mild, watery diarrhea that lasts 2 to 3 days, along with vomiting, fever, anorexia, and abdominal pain. The infection is transmitted primarily through the fecal-oral route, even though food and water can act as vehicles for the transmission of human enteric viruses.

Sapoviruses (SaV), formerly known as "Sapporo-like viruses," belong to the *Caliciviridae* family and cause acute gastroenteritis in humans and pigs. Fifteen genotypes grouped into four genogroups (GI.1 to GI.8, GII.1 to GII.5, GIV, and GV) have been described as causing infection in humans. SaV is considered an important cause of gastroenteritis in children under 5 years of age, while in adults it is less important. The clinical symptoms of sapovirus infection are believed to be milder than symptoms associated with norovirus, including mild or acute watery diarrhea, stomach cramps, nausea, vomiting, and occasionally fever. SaVs can be transmitted through the fecal-oral route from contaminated food and water, as well as through person-to-person contact.

Culture, antigen detection, and microscopy have been considered the gold standard methods for the diagnosis of these pathogens; however, it is time consuming and not very sensitive. The real-time PCR assay is less labor-intensive and has higher sensitivity and specificity, making it an attractive alternative.

3. Principle of the procedure

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit is designed for the diagnosis of RNA/DNA *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV in faeces samples. After DNA isolation, the identification of AdV, *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica* is performed by the amplification of a conserved region of the specific genes using specific primers and fluorescent-labelled probes. After RNA isolation, the detection of RNA viruses NoV GI, NoV GII, RV, AtV and SaV is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of specific genes using specific primers and a fluorescent-labelled probes.

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in a stabilized format, as well as an internal control to monitor PCR inhibition (if present).

4. Reagents provided

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open format with internal control products and Annex 2 for open format with extraction control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates).
- DNA/RNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or $\leq -70^{\circ}\text{C}$.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 μL , 20-200 μL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit has been validated on the following equipments: 7500 Fast Real-Time PCR System (Applied Biosystems), CFX96™ Real-Time PCR Detection System (Bio-Rad), AriaMx Real-Time PCR System (Agilent Technologies), Cobas z480 Analyzer (Roche Molecular Diagnostics), Rotor-Gene® Q (Qiagen), DTprime Real-time Detection Thermal Cycler (DNA-Technology), DTlite Real-time Detection Thermal Cycler (DNA-Technology). To check thermocycler compatibility and most common detection channels consult website www.certest.es.

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with VIASURE Real Time PCR Detection Kits. This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.

- Keep components away from light.

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- **Make sure to put the strip in the correct direction. The first well is marked with a hole on the top right corner. Be careful not to invert the strip throughout PCR process.**
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit and any additional reagents or equipment required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

- In accordance with Regulation (EC) No 1907/2006 (REACH), "VIASURE Real Time PCR Detection Kits" does not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because it does not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions, and procedures.

8. Test procedure

Please see Annex 1 for Open format with internal control products Test Procedure and Annex 2 for Open format with extraction control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit has been validated in faecal samples. Other types of samples must be validated by the user.

Collection, storage and transport specimens should be maintained per the conditions validated by the user. Overall, faecal samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at 2 to 8°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), It is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 24 hours or frozen at -20°C or ideally at -80°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The faecal specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf> and the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94).

8.2. RNA/DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

Due to VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit is also available with extraction control (References include in Annex 2), if the Extraction Control (EC) is used to monitor nucleic acid isolation and as PCR inhibition control, **add 5µL** of the reconstituted EC to the specimen and/or lysis buffer mixture (clinical specimen, as well as positive control and/or negative control). Close each tube and vortex for 10 seconds.

If the Extraction Control is used only as a PCR inhibition control, 1µL of the reconstituted EC should be added to the reconstituted Reaction-Mix.

For RNA/DNA extraction from faecal samples, you can use your manual or automatic routine optimized system, or any commercially available RNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- Invisorb® Spin Universal Kit (Invitek).
- QIAamp DNA Stool Mini Kit (QIAGEN).
- RIDA® Xtract (r-Biopharm).
- Maxwell® RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).
- MagDEA Dx SV kit, using the magLEAD® 6gC or 12gC instrument (Precision System Science Co.)

9. Result interpretation

9.1. References with internal control (references in Annex 1)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Gastrointestinal Panel I* positive control of each well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	Target detected ¹	Internal Control ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal ³	≤40	Valid

Table 1. Expected Performance of Controls

1 In cases where one or more controls fail (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2. The Internal Control (IC) should show an amplification signal (Ct ≤40) in control wells (PC and NC).

3. Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter* (SCY well). Therefore, a cut-off Ct of 37 is set for this target.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

In case that all target genes of the different pathogens results are negative, IC must show an amplification signal with Ct less than 40. If there is an absence of signal or Ct value > 40 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA/RNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.


For interpretation of patient sample results, use the following table, read and analyse the results:

	Reaction Mix	Pathogens	Channels			
			FAM	HEX	ROX	Cy5
1	SCY	<i>Salmonella</i>	Positive			
		Internal Control		Positive/Negative		
		<i>Yersinia enterocolitica</i>			Positive	
		<i>Campylobacter</i>				Positive*
2	SHY	<i>Shigella/EIEc</i>	Positive			
		Internal Control		Positive/Negative		
3	KGE	<i>Giardia lamblia</i>	Positive			
		Internal Control		Positive/Negative		
		<i>Entamoeba histolytica</i>			Positive	
		<i>Cryptosporidium</i>				Positive
4	NOR	Norovirus GII	Positive			
		Internal Control		Positive/Negative		
		Norovirus GI				Positive
5	RTV	Rotavirus	Positive			
		Internal Control		Positive/Negative		
6	ADV	Adenovirus	Positive			
		Internal Control		Positive/Negative		
7	ATV	Astrovirus	Positive			
		Internal Control		Positive/Negative		
8	SAV	Sapovirus	Positive			
		Internal Control		Positive/Negative		

Table 2. Interpretation of patient sample results. Positive: Amplification curve. Empty: No amplification curve.

* Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter* (SCY well). Therefore, a cut-off Ct of 37 is set for this target.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.


Table 3. Correct run of negative and positive controls run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

9.2. References with extraction control (references in Annex 2)

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the extraction procedure and correct functioning of the amplification mix (1 KGE reaction mix). The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Gastrointestinal Panel I* positive control of each well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	Targets detected ¹	Extraction/Internal Control ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal ³	≤40	Valid

Table 4. Expected Performance of Controls.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The Extraction/Internal Control (EC) should show an amplification signal (Ct ≤40) in control wells (PC and NC). Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

3. Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter* (SCY well). Therefore, a cut-off Ct of 37 is set for this target.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

In case that all target genes of the different pathogens results are negative, EC must show an amplification signal with Ct less than 40. If there is an absence of signal or Ct value >40 of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA/RNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

The Extraction Control (EC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction control between the controls and the clinical samples, due to the extraction process.

For interpretation of patient sample results, use the following table, read, and analyse the results:


	Reaction Mix	Pathogens	Channels			
			FAM	HEX	ROX	Cy5
1	SCY	<i>Salmonella</i>	Positive			
		Internal Control		Positive/Negative		
		<i>Yersinia enterocolitica</i>			Positive	
		<i>Campylobacter</i>				Positive*
2	SHY	<i>Shigella/EIEc</i>	Positive			
		Internal Control		Positive/Negative		
3	1KGE	<i>Giardia lamblia</i>	Positive			
		Extraction Control		Positive/Negative		
		<i>Entamoeba histolytica</i>			Positive	

		<i>Cryptosporidium</i>				Positive
4	NOR	Norovirus GII	Positive			
		Internal Control		Positive/Negative		
		Norovirus GI				Positive
5	RTV	Rotavirus	Positive			
		Internal Control		Positive/Negative		
6	ADV	Adenovirus	Positive			
		Internal Control		Positive/Negative		
7	ATV	Astrovirus	Positive			
		Internal Control		Positive/Negative		
8	SAV	Sapovirus	Positive			
		Internal Control		Positive/Negative		

Table 5. Interpretation of patient sample results. Positive: Amplification curve. Empty: No amplification curve.

* Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter* (SCY well). Therefore, a cut-off Ct of 37 is set for this target.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user, to verify the correct performance of each qPCR steps and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.


Table 6. Correct run of negative and positive controls run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with RNA/DNA extracted from faecal samples. See section 12 to check the matrix specimens validated.
- The quality of the test depends on the quality of the sample; properly extracted nucleic acid from clinical samples must be extracted.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV, either the high number of DNA or cDNA template copies which contains each *Gastrointestinal Panel I* Positive Control vial, samples containing high concentrations of target RNA/DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Gastrointestinal Panel I* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory area.
- The specific primer and probe combinations for detection of the genes used in VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit designed for the detection of *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV, do not show significant combined homologies with the human genome or human microflora, which might result in predictable false positive.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including RNA/DNA extraction).
 - Degradation of the viral RNA/DNA and/or viral/bacterial DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV variants.
 - A viral/bacterial load in the specimen below the limit of detection for the assay.
 - Failure to follow instructions for use and the assay procedure.

If in doubt, refer to section 9 to check the correct interpretation of the results.

- A positive test result does not necessarily indicate the presence of viable viruses/bacteria and does not imply that these viruses/bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets viral/bacterial sequences.

- Negative results or single target amplification do not preclude gastrointestinal infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral/bacterial levels during infections caused by *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogens.
- If diagnostic tests for other gastrointestinal illnesses are negative and the patient's clinical presentation and epidemiological information suggest that *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*, *Shigella/EIEC*, *Cryptosporidium*, *G. lamblia* and/or *E. histolytica*, NoV GI and/or NoV GII, RV, AdV, AstV and/or SaV is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.

11. Quality control

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) or Extraction control (EC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit was tested using stool clinical samples from patients with clinical signs and symptoms of gastrointestinal infection. To determine the clinical diagnostic accuracy, a multicentre evaluation has been conducted through collaboration with national entities. A summary of the sites, sample type and workflow are included in the following table:

Site	Sample type	Workflow	Master Mix Code	Target
1 CerTest Biotec S.L., Parasitology Reference and Research Laboratory, Spanish National Centre for Microbiology (SNCM), Health Institute Carlos III (ISCIII), Hospital Clínico Universitario Lozano Blesa (HCULB)	Stool samples	Extracted and purified stool samples + DT Prime Real-Time PCR System (DNA Technologies)	GP01	<i>Salmonella</i>
				<i>Y. enterocolitica</i>
				<i>Campylobacter</i>
				<i>Shigella/EIEC</i>
				<i>G. lamblia</i>
				<i>E. histolytica</i>
				<i>Cryptosporidium</i> spp.
				Norovirus GI
				Norovirus GII
				Rotavirus
				Adenovirus

				Astrovirus
				Sapovirus

Table 7. Site, sample type, workflow and target

Besides, true positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values of VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Master Mix Code	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Routine analysis *	GP01	<i>Salmonella</i>	12	164	0	0	1 (0.73-1)	1 (0.97-1)	1 (0.69-1)	1 (0.97-1)
			<i>Y. enterocolitica</i>	12	164	0	0	1 (0.73-1)	1 (0.97-1)	1 (0.69-1)	1 (0.97-1)
			<i>Campylobacter</i>	14	161	1	0	1 (0.76-1)	0.99 (0.96-1)	0.93 (0.66-0.99)	1 (0.97-1)
			<i>Shigella</i> /EIEC	10	166	0	0	1 (0.69-1)	1 (0.97-1)	1 (0.65-1)	1 (0.97-1)
			<i>G. lamblia</i>	23	153	0	0	1 (0.85-1)	1 (0.97-1)	1 (0.82-1)	1 (0.96-1)
			<i>E. histolytica</i>	20	156	0	0	1 (0.83-1)	1 (0.97-1)	1 (0.79-1)	1 (0.97-1)
			<i>Cryptosporidium</i> spp.	24	152	0	0	1 (0.85-1)	1 (0.97-1)	1 (0.82-1)	1 (0.96-1)
			Norovirus GI	3	172	1	0	1 (0.29-1)	0.99 (0.96-1)	0.75 (0.21-0.98)	1 (0.97-1)
			Norovirus GII	17	159	0	1	0.94 (0.72-0.99)	1 (0.97-1)	1 (0.77-1)	1 (0.97-1)
			Rotavirus	18	155	3	0	1 (0.81-1)	0.98 (0.94-0.99)	0.85 (0.62-0.96)	1 (0.96-1)
			Adenovirus	18	155	2	1	0.94 (0.74-0.99)	0.98 (0.95-0.99)	0.94 (0.74-0.99)	0.98 (0.95-0.99)
			Astrovirus	9	167	0	0	1 (0.66-1)	1 (0.97-1)	1 (0.62-1)	1 (0.97-1)
			Sapovirus	9	167	0	0	1 (0.66-1)	1 (0.97-1)	1 (0.62-1)	1 (0.97-1)

Table 8. True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit.

* Routine analysis: Culture method, MALDI-TOF mass spectrometry (if colonies were not easily identifiable) and immunochromatography for diagnosis of gastrointestinal viruses.

12.2. Analytical sensitivity

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA/RNA copies per reaction for *Y. enterocolitica* (*ystB* gene); *Shigella*/EIEC; *Giardia lamblia* and *Entamoeba histolytica*; Norovirus GI and Norovirus GII; Rotavirus; Adenovirus; Astrovirus; and Sapovirus; a detection limit of ≥ 50 DNA copies per reaction for *Cryptosporidium*, *Salmonella*, *Y. enterocolitica* (*ail* gene) and *Campylobacter*, with a positive rate of $\geq 95\%$.

12.3. Analytical specificity

The specificity of the *Gastrointestinal Panel I* assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter coli</i>	-/+	<i>Enterococcus faecalis</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-/+	<i>Yersinia enterocolitica</i> O:3	-/+
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter upsaliensis</i>	-/+	<i>Yersinia enterocolitica</i> O:9	-/+
<i>Helicobacter heilmannii</i>	-	<i>Proteus vulgaris</i>	-	<i>Bacteroides fragilis</i>	-
<i>Shigella flexneri</i>	-/+	<i>Citrobacter freundii</i>	-	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	-
<i>Shigella dysenteriae</i>	-/+	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Aeromonas caviae</i>	-
<i>Salmonella typhi</i>	-/+	<i>Serratia liquefaciens</i>	-	Adenovirus serotypes 40/41	-/+
<i>Salmonella paratyphi</i> A	-/+	<i>Vibrio parahaemolyticus</i>	-	Rotavirus A	-/+
<i>Salmonella paratyphi</i> B	-/+	<i>Clostridium difficile</i>	-	Astrovirus Genotype I-VIII	-/+
<i>Salmonella typhimurium</i>	-/+	<i>Clostridium perfringens</i>	-	Sapovirus	-/+
<i>Salmonella bongori</i>	-/+	Enteropathogenic <i>Escherichia coli</i>	-	<i>Entamoeba histolytica</i>	-/+
<i>Salmonella enteritidis</i>	-/+	Enterotoxigenic <i>Escherichia coli</i>	-	<i>Entamoeba dispar</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i>	-/+	<i>Klebsiella oxytoca</i>	-	<i>Dientamoeba fragilis</i>	-
<i>Salmonella pullorum</i>	-/+	<i>Listeria monocytogenes</i>	-	<i>Cryptosporidium parvum/hominis</i>	-/+
<i>Salmonella gallinarum</i>	-/+	<i>Candida albicans</i>	-	<i>Giardia intestinalis</i>	-/+
<i>Campylobacter lari</i>	-/+	<i>Arcobacter butzleri</i>	-	<i>Blastocystis hominis</i>	-
<i>Campylobacter fetus</i>	-/+	<i>Pseudomonas aeruginosa</i>	-	Norovirus Genotypes I and II	-/+

Table 9. Reference pathogenic microorganisms used in this study

12.4. Analytical reactivity

The reactivity of VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit for ***Salmonella*, *Campylobacter* & *Y. enterocolitica* Reaction mix**, was evaluated against *Salmonella* (*Salmonella bongori* strain 1224.72, serotype Brookfield (66:z41:-); *Salmonella enterica* subsp. *enterica* serovar *enteritidis*; *Salmonella enterica* subsp. *enterica* strain LT2, serovar *typhimurium*, serotype 4,5,12:i:1,2; *Salmonella enterica* subsp. *enterica* serovar *gallinarum*; *Salmonella enterica* subsp. *enterica* serovar *pullorum*; *Salmonella enterica* subsp. *enterica* serovar *paratyphi* A; *Salmonella enterica* subsp. *enterica* serovar *paratyphi* B and *Salmonella enterica* subsp. *enterica* serovar *typhi*); *Campylobacter* (*Campylobacter coli*; *Campylobacter hyointestinalis* strain 80-4577-4; *Campylobacter jejuni* subsp. *jejuni* strain LMG 18455; *Campylobacter lari* and *Campylobacter upsaliensis* strain CCUG 24803) and *Yersinia enterocolitica* (*Yersinia enterocolitica* O:3 and *Yersinia enterocolitica* O:9), showing positive result.

The reactivity of VIASURE ***Shigella*/EIEC Reaction mix**, was evaluated against *Shigella dysenteriae*, *Shigella flexneri*, *Shigella boydii* and *Shigella sonnei*, showing positive result.

The reactivity of VIASURE ***Cryptosporidium*, *Giardia* & *E. histolytica* Reaction mix**, was evaluated against *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica* strain DS4-868, showing positive result.

The reactivity of VIASURE ***Norovirus GI + GII Reaction mix***, was evaluated against *Norovirus GI* (genotypes GI.1, GI.3, GI.7, GI.8 and GI.Pb) and *Norovirus GII* (genotypes GII.1, GII.2, GII.3, GII.4, GII.4 New Orleans 2009, GII.4 Sydney 2012, GII.5, GII.6, GII.7, GII.10, GII.12, GII.13, GII.16, GII.17, GII.21, GII.b, GII.P4, GII.P17, GII.P21, GII.P16_GII.4 2012 recombinant, GII.P16_GII.2 recombinant, GII.Pe, G II.Pe_G II.4 and GII.Pg_GII.1), showing positive result.

The reactivity of VIASURE **Rotavirus Reaction mix**, was evaluated against human Rotavirus A strain ATCC® VR-2018TM (VP7 serotype 1 (G1)) and Rotavirus genotypes G1P[8], G2P[4], G2P[8], G3P[8], G4P[8], G1+G3P[8], G6P[8], G9P[8], G10+G12P[8], G12P[8] and G-UDP[8] strains, showing positive result.

The reactivity of VIASURE **Adenovirus Reaction mix**, was evaluated against Human Adenovirus 40 strain Dugan, Human Adenovirus 41 strain Tak, Human Adenovirus Type 31 strain 1315/63, Human Adenovirus type 1, Human Adenovirus 2 strain Adenoid 6 and Human Adenovirus 5, Human Adenovirus type 6, showing positive result.

The reactivity of VIASURE **Astrovirus Reaction mix**, was evaluated against human Astrovirus serotypes I-VIII, showing positive result.

The reactivity of VIASURE **Sapovirus Reaction mix**, was evaluated against Sapovirus genotypes GI.1, GI.2, GI.3, GII.1, GII.2, GII.3 and OH8021/2008/JP-like (putative GI.5), showing positive result.

ANNEX 1

OPEN FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Gastrointestinal Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-GP0112L
VIASURE <i>Gastrointestinal Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-GP0112H

Table A1 1. References

A1.1 Principle of the procedure

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in a stabilized format, as well as an internal control to monitor PCR inhibition.

The reaction mixes places in each well of ***Gastrointestinal Panel I*** 8-well strip (GP01) allow the detection of the specific target pathogens through specific genes in the following channels (Table A1.2):


		Code	Target	Channel	Gene
	1	SCY	<i>Salmonella</i>	FAM	<i>invA</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
			<i>Y. enterocolitica</i>	ROX	<i>ail</i> and <i>ystB</i>
			<i>Campylobacter</i>	Cy5	<i>16S rRNA</i>
	2	SHY	<i>Shigella/EIEC</i>	FAM	<i>ipaH</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
	3	KGE	<i>Giardia lamblia</i>	FAM	<i>18S rRNA</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
			<i>E. histolytica</i>	ROX	<i>18S rRNA</i>
			<i>Cryptosporidium</i>	Cy5	<i>18S rRNA</i>
	4	NOR	<i>Norovirus GII</i>	FAM	<i>ORF1-ORF2</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
			<i>Norovirus GI</i>	Cy5	<i>ORF1-ORF2</i>
	5	RTV	Rotavirus	FAM	<i>NSP3</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
	6	ADV	Adenovirus	FAM	<i>Hexon</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
	7	ATV	Astrovirus	FAM	<i>ORF1b</i>
			Internal Control (IC)	HEX, VIC or JOE *	-
	8	SAV	Sapovirus	FAM	<i>ORF1</i>
			Internal Control (IC)	HEX, VIC or JOE *	-

Table A1 2. Target, channel and genes. *Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.cerfess.es. Note that the first well is marked with a hole on the top right corner.

A1.2 Reagents provided

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A1.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.4 includes materials and reagents to be used with 8-well strips compatible devices. (Consult the thermocycler compatibility on CerTest's website www.cer-test.es).

Reagent/Material	Description	Colour	Amount
<i>Gastrointestinal Panel I</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Gastrointestinal Panel I</i> Positive Control	Non-infectious synthetic lyophilized cDNA/DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit with Ref. VS-GP0112L and VS-GP0112H.

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Gastrointestinal Panel I Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Gastrointestinal Panel I* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips. **Make sure to put the strip in the correct direction (Table A1 2). The first well is marked with a hole on the top right corner.**

- 1) Reconstitute the number of strips you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA/DNA sample, reconstituted *Gastrointestinal Panel I* Positive Control (red vial) or Negative Control (violet vial) in the 8 wells of each strip and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A1 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM, ROX, Cy5 and/or HEX (JOE or VIC) channels in the following Table A1.2. Depending on the equipment used select the proper detection channel consult website www.certest.es.

ANNEX 2

OPEN FORMAT WITH EXTRACTION CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE <i>Gastrointestinal Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-GP0112LE
VIASURE <i>Gastrointestinal Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-GP0112HE

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity.

The reaction mixes places in each well of ***Gastrointestinal Panel I*** 8-well strip (GP01) allow the detection of the specific target pathogens through specific genes in the following channels (Table A1.2):

	Code	Target	Channel	Gene
1	SCY	<i>Salmonella</i>	FAM	<i>invA</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
		<i>Y. enterocolitica</i>	ROX	<i>ail</i> and <i>ystB</i>
		<i>Campylobacter</i>	Cy5	<i>16S rRNA</i>
2	SHY	<i>Shigella/EIEC</i>	FAM	<i>ipaH</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
3	1KGE	<i>Giardia lamblia</i>	FAM	<i>18S rRNA</i>
		Extraction Control (EC)	HEX, VIC or JOE *	-
		<i>E. histolytica</i>	ROX	<i>18S rRNA</i>
		<i>Cryptosporidium</i>	Cy5	<i>18S rRNA</i>
4	NOR	Norovirus <i>GII</i>	FAM	<i>ORF1-ORF2</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
		Norovirus <i>GI</i>	Cy5	<i>ORF1-ORF2</i>
5	RTV	Rotavirus	FAM	<i>NSP3</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
6	ADV	Adenovirus	FAM	<i>Hexon</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
7	ATV	Astrovirus	FAM	<i>ORF1b</i>
		Internal Control (IC)	HEX, VIC or JOE *	-
8	SAV	Sapovirus	FAM	<i>ORF1</i>
		Internal Control (IC)	HEX, VIC or JOE *	-

Table A2. 2. Target, channel and genes. *Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website www.certest.es. Note that the first well is marked with a hole on the top right corner.

A2.2 Reagents provided

VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A2.4 includes materials and reagents to be used with 8-well strips compatible devices. (Consult the thermocycler compatibility on CerTest's website www.cerTEST.es)

Reagent/Material	Description	Colour	Amount
<i>Gastrointestinal Panel I</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	12x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Gastrointestinal Panel I</i> Positive Control	Non-infectious synthetic lyophilized cDNA/DNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	12 x 8-cap strip

Table A2. 3. Reagents and materials provided in VIASURE *Gastrointestinal Panel I* Real Time PCR Detection Kit with Ref. VS-GP0112LE and VS-GP0112HE.

A2.3 Test procedure

A2.3.1 Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Gastrointestinal Panel I* Positive Control in an area away from the other components.

A2.3.2 Lyophilized positive control

Gastrointestinal Panel I Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Gastrointestinal Panel I* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.3 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips. **Make sure to put the strip in the correct direction (Table A2 2). The first well is marked with a hole on the top right corner.**

- 1) Reconstitute the number of strips you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the 8 wells of each strip for negative control.

Add 5 µL of RNA/DNA sample in different wells.

Add 5 µL of reconstituted *Gastrointestinal Panel I* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A2. 4 PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM, ROX, Cy5 and/or HEX (JOE or VIC) channels in the following Table A1.2. Depending on the equipment used select the proper detection channel consult website www.certest.es.

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Bibliography/Bibliografía

1. J. Hoorfar, et al. Automated 59 Nuclease PCR Assay for Identification of *Salmonella enterica*. *Journal of Clinical Microbiology* 2000; 38: 3429–3435.
2. D. Rodríguez-Lázaro et al. A rapid and direct real time PCR-based method for identification of *Salmonella* spp. *Journal of Microbiological Methods* 2003; 54: 381– 390.
3. F. Barletta et al. Multiplex real-time PCR for detection of *Campylobacter*, *Salmonella*, and *Shigella*. *Journal of Clinical Microbiology* 2013; 51: 2822-2829.
4. N.O. Kaakoush et al. Global Epidemiology of *Campylobacter* Infection. *Clinical Microbiology Reviews* 2015; 28(3): 687-720.
5. S.M. Man. The clinical importance of emerging *Campylobacter* species. *Nature Reviews Gastroenterology & Hepatology* 2011; 8(12): 669-685.
6. R.F. de Boer et al. Detection of *Campylobacter* species and *Arcobacter butzleri* in stool samples by use of real-time multiplex PCR. *Journal of Clinical Microbiology* 2013; 51(1): 253-259.
7. N. Botteldoorn et al. Quantification of *Campylobacter* spp. in chicken carcass rinse by real-time PCR. *Journal of Applied Microbiology* 2008; 105(6):1909-1918.
8. M. Arrausi-Subiza, et al. Evaluation of different enrichment methods for pathogenic *Yersinia* species detection by real time PCR. *BMC Veterinary Research* 2014; 10(1): 192.
9. A.D. Jourdan et al. Development of a fluorogenic 5' nuclease PCR assay for detection of the *ail* gene of pathogenic *Yersinia enterocolitica*. *Applied and Environmental Microbiology* 2000; 66(9):3750-3755.
10. S.T. Lambertz et al. Real-time PCR method for detection of pathogenic *Yersinia enterocolitica* in food. *Applied and Environmental Microbiology* 2008;74(19):6060-6067.
11. J. Liang et al. Two novel *ail*-positive biotype 1A strains of *Yersinia enterocolitica* isolated in China with unequal adhesion and invasion properties. *Infection, Genetics and Evolution* 2014;27:83-88.
12. A. Ub-Din et al. Relationship among *Shigella* spp. and enteroinvasive *Escherichia coli* (EIEC) and their differentiation. *Brazilian Journal of Microbiology* 2014; 45(4): 1131-1138.
13. V. D. Thiem et al. Detection of *Shigella* by a PCR assay targeting the *ipaH* gene suggests increased prevalence of shigellosis in Nha Trang, Vietnam. *Journal of Clinical Microbiology* 2004; 42(5): 2031-2035.
14. W. S. Lin et al. A quantitative PCR assay for rapid detection of *Shigella* species in fresh produce. *Journal of Food Protection* 2010; 73(2): 221-233.
15. W. Mokhtari et al. Real-Time PCR using SYBR Green for the detection of *Shigella* spp. in food and stool samples. *Molecular and Cellular Probes* 2013, 27; 53-59.
16. S. Ghosh et al. Genetic characterization of *Shigella* spp. isolated from diarrhoeal and asymptomatic children. *Journal of Medical Microbiology* 2014, 63; 903-910.
17. J.J.R. Limor. Detection and differentiation of *Cryptosporidium* parasites that are pathogenic for humans by real-time PCR. *Journal of Clinical Microbiology* 2002; 40(7): 2335-2338.
18. S.E. Stroup et al. Real-time PCR detection and speciation of *Cryptosporidium* infection using Scorpion probes. *Journal of Medical Microbiology* 2006; 55: 1217-1222.
19. S.J. Hadfield et al. Detection and differentiation of *Cryptosporidium* spp. in human clinical samples by use of real-time PCR. *Journal of Clinical Microbiology* 2011; 49(3): 918-924.
20. D.V. den Bossche et al. Comparison of four rapid diagnostic tests, ELISA, microscopy and PCR for the detection of *Giardia lamblia*, *Cryptosporidium* spp. and *Entamoeba histolytica* in feces. *Journal of microbiological methods* 2015; 110: 78-84.

21. R.D. Adam. Biology of *Giardia lamblia*. *Clinical Microbiology Reviews* 2001; 14(3): 447-475.
22. J.J. Verweij et al. Real-time PCR for the detection of *Giardia lamblia*. *Molecular and Cellular Probes* 2003; 17(5): 223-225.
23. R. Fotedar et al. PCR Detection of *Entamoeba histolytica*, *Entamoeba dispar*, and *Entamoeba moshkovskii* in Stool Samples from Sydney, Australia. *Journal of Clinical Microbiology* 2007; 45 (3): 1035-1037.
24. R. Haque et al. Multiplex real-time PCR assay for detection of *Entamoeba histolytica*, *Giardia intestinalis*, and *Cryptosporidium* spp. *The American Journal of Tropical Medicine and Hygiene* 2007; 76(4): 713-717.
25. Y. Lau et al. Real-time PCR assay in differentiating *Entamoeba histolytica*, *Entamoeba dispar*, and *Entamoeba moshkovskii* infections in Orang Asli settlements in Malaysia. *Parasites & Vectors* 2013; 6(1): 250.
26. C. Fuentes et al. Standardized multiplex one-step qRT-PCR for hepatitis A virus, norovirus GI and GII quantification in bivalve mollusks and water. *Food Microbiology* 2014; 40: 55-63.
27. T. Kageyama et al. Broadly reactive and highly sensitive assay for Norwalk-like viruses based on real-time quantitative reverse transcription-PCR. *Journal of Clinical Microbiology* 2003; 41(4): 1548-1557.
28. N. Jothikumar et al. Rapid and sensitive detection of noroviruses by using TaqMan-based one-step reverse transcription-PCR assays and application to naturally contaminated shellfish samples. *Applied and Environmental Microbiology* 2005; 71(4): 1870-1875.
29. J. Vinjé. Advances in laboratory methods for detection and typing of norovirus. *Journal of Clinical Microbiology* 2015; 53(2): 373-381.
30. S. Mijatovic et al. Sensitive and Specific Quantitative Detection of Rotavirus A by One-Step Real-Time Reverse Transcription-PCR Assay without Antecedent Double-Stranded-RNA Denaturation. *Journal of Clinical Microbiology* 2013; 51(9):3047-3054.
31. XL. Pang et al. Increased detection of rotavirus using a real time reverse transcription-polymerase chain reaction (RT-PCR) assay in stool specimens from children with diarrhea. *Journal of Medical Virology* 2004; 72 (3): 496-501.
32. S Q Zeng. et al. One-step quantitative RT-PCR for the detection of rotavirus in acute gastroenteritis. *Journal of virological methods* 2008; 153(2): 238-240.
33. I Wilhelmi et al. Viruses causing gastroenteritis. *Clinical Microbiology and infection* 2003; 9: 247-262.
34. A Heim. et al. Rapid and Quantitative Detection of Human Adenovirus DNA by Real-Time PCR. *Journal of Medical Virology* 2003; 70: 228–239.
35. N Jothikumar. et al. Quantitative Real-Time PCR Assays for Detection of Human Adenoviruses and Identification of Serotypes 40 and 41. *Applied and environmental microbiology* 2005; 71: 3131–3136.
36. A Allard. et al. Rapid typing of human adenoviruses by a general PCR combined with restriction endonuclease analysis. *Journal of Clinical Microbiology* 2001; 39(2): 498-505.
37. S. Guix et al. Human Astrovirus diagnosis and typing: current and future prospects. *Letters in Applied Microbiology* 2005; 41 (2): 103-105.
38. 2. D. Sano et al. Detection of oxidative damages on viral capsid protein for evaluating structural integrity and infectivity of human norovirus. *Environmental Science & Technology* 2010; 44(2): 808-812.
39. 3. A. Bosch et al. Human Astroviruses. *Clinical Microbiology Reviews* 2014; 27(4): 1048-1074.
40. M. Okada et al. The detection of human sapoviruses with universal and genogroup-specific primers. *Archives of Virology* 2006; 151(12): 2503-2509.

41. 2. T. Oka et al. Detection of human sapovirus by real-time reverse transcription-polymerase chain reaction. *Journal of Medical Virology* 2006; 78(10): 1347-1353.
42. 3. M. Kitajima et al. Detection and genetic analysis of human sapoviruses in river water in Japan. *Applied and Environmental Microbiology* 2010; 76(8): 2461-2467.
43. 4. C. Logan et al. Real-time reverse transcription PCR detection of norovirus, sapovirus and astrovirus as causative agents of acute viral gastroenteritis. *Journal of Virological Methods* 2007; 146(1-2): 36-44.
44. 5. F. Rovida et al. Molecular detection of gastrointestinal viral infections in hospitalized patients. *Diagnostic Microbiology and Infectious Disease* 2013; 77(3): 231-235.
45. 6. D. Sano et al. Quantification and genotyping of human sapoviruses in the Llobregat river catchment, Spain. *Applied and Environmental Microbiology* 2011; 77(3): 1111-1114.
46. 7. P.J. Johansson. A nosocomial sapovirus-associated outbreak of gastroenteritis in adults. *Scandinavian Journal of Infectious Diseases* 2005; 37(3): 200-204.

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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