

# Anti-CCP

Anti-Cyclic Peptide Containing Citrulline Antibody Test Kit  
(Rate Scattering Turbidimetric Method)


## Instructions for Use

Version: A/8

REF HP-Anti-CCP-25

### Manufacturer

 Shijiazhuang Hipro Biotechnology Co., Ltd.  
No. 3 Building, Block C, Fangyi Science Park, No. 365 Huai'an  
East Road, Hi-tech Zone, Shijiazhuang, 050000 Hebei P.R. China  
After sale service: 400-0191-606  
www.hipro.us

  Lotus NL B.V.

Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague,  
Netherlands.  
Tel: +31644168999

### Product Name

General Name:Anti-Cyclic Peptide Containing Citrulline  
Antibody Test Kit (Rate Scattering Turbidimetric Method)

### Specification

Package Specification  
25 Tests/ Kit.

### Intended Use

This product is used to determine the content of anti-cyclic  
peptide containing citrulline antibody (Anti-CCP) in human  
serum.

Anti-cyclic peptide antibody is a polypeptide fragment of annu-  
lar polyfilament protein. It is a kind of autoantibody based on  
IgG, and can be used in the auxiliary diagnosis of rheumatoid  
arthritis (RA).

### Test Principle

The cyclic citrulline peptide antigen is coated on the latex  
surface. The anti-cyclic peptide containing citrulline antibody  
in the sample and the antibody become to immune complexes  
by Latex agglutination reaction. The immune complexes will  
produce the phenomenon of light scattering which is propor-  
tional to the intensity of scattered light and samples of  
Anti-CCP levels. Using specific protein analyzer to measure the  
intensity of scattered light, the concentration of Anti-CCP is  
determined by comparing the turbidity of samples to the stan-  
dard concentration.

### Component

The Anti-CCP test kit consists of two reagents R1 and R2, as  
shown on Figure 1.




Figure 1

Name	Content	Quantity
Reagent 1 (R1)	Phosphate buffer(pH 6.5)	0.1mol/L
	Sodium azide	0.1%
	Polyethylene glycol 6000	3.5%
Reagent 2 (R2)	Phosphate buffer(pH 8.0)	0.1mol/L
	Cyclic citrullin peptide antigen with latex	4.3mL/L
IC card (optional)	/	1

Do not mix different batches of reagents.

### Storage&stability

Store the test kit at 2°C-8°C until the expiration date indicated  
on the label. The test kit is stable for one year when unopened.  
Use up the test kit within one month after opening the pack-  
age.

Do not freeze the test kit.

Do not mix different lots of the test kit.

### Special Instrument Requirements

HP-083/4-I POCT Immunoassay System,  
HP-083/4-II POCT Immunoassay System,  
HP-AFS/1 Automatic Immunoassay System,  
HP-AFS/3 Automatic Immunoassay System.

### Specimen type

Serum, avoid hemolysis. Fasting blood collection and separa-  
tion of serum as soon as possible. The sampl e store at 2-8°C  
for 3 days, -20°C for 1 month. Avoid repeated freezing. Before  
test, ensure fully mixed.


### Procedures

#### HP-083/4-I&HP-083/4-II POCT Immunoassay System


#### Note:

- Please read user manual of HP-083/4-I and HP-083/4-II  
before use.
- The analyzer will finish self check after start-up.
- Insert the IC card of Anti-CCP test kit to let analyzer read the  
parameter.
- The analyzer calibration can be done with app. It is recom-  
mended that analyzer calibration should be done for each  
new lot of test kit.


#### Step 1 Sample Preparation



- 1a** Allow the test kit back to  
room temperature for 30  
minutes before use.



- 1b** Use the R2 with quantitative  
capillary to collect sample.



- 1c** Insert the R2 into R1.

#### Note:

- The parameter is built in the IC card.
- Please insert the corresponding IC card into analyzer to let  
the analyzer read the parameter before each assay test.
- The capillary of the R2 should be fully filled.

#### Step 2 Testing




- 2a** Hold the narrow side of R1  
and shake from left to right  
for 12 times to let the sample  
completely mix with R1.


**Note:**Do not hold the wide sides of R1.



- 2b** Press the piston on R2.



- 2c** Hold the narrow sides of R1  
and shake for 3-5 seconds to  
let R1 and R2 mix well.




- 2d** Insert the R1 into test channel  
of HP-083/4-II and the analyz-  
er will finish the test and print  
out results automatically.

#### HP-AFS/1&HP-AFS/3 Automatic Immunoassay System


#### Note:

- Please read user manual of HP-AFS/1 and HP-AFS/3 before  
use.
- The analyzer will finish the self check after start-up.
- It is recommend to do analyzer calibration monthly and for  
each new lot of test kit.

#### Step 1 Sample Preparation



- 3a** Allow the test kit back to  
room temperature for 30  
minutes before use.



- Note:**
- Please update the standard curve with the barcode on the R1 cuvette if a new lot test kit is to be used.
  - The capillary of the R2 should be fully filled.

**Step 2 Testing**

**4a** Insert the R1 into test channel of analyzer.

**4b** The analyzer will mix the sample from R2 capillary with R1 automatically.

**4c** The analyzer will mix the R2 and R1 automatically.

**4d** The analyzer will test and print the results automatically.

**Calibration**

The calibration values for the different lots of the kits are stored on the calibration IC card or the two-dimensional code on the cuvette. Before test the new lot of kits, read the calibration card parameters first. Or the instrument automatically scan the two-dimensional code on the cup to obtain the corresponding calibration curve during testing.

**Quality control**

3- level calibration system guarantee the results' reliability for each lot of test kits, including the instrument calibration, remote reagent calibration and the third party calibration.

The third party calibration applicable for:

1. The daily indoor quality control test.
2. New lots of reagent.
3. New operator training.
4. The results can not match the clinical symptoms.
5. The first use of the reagent.

If still can not be calibrated, contact the manufacture for further technical support.

**Reference Value**

<45U/mL. Recommended that each laboratory establish its own reference range.

**Interpretation**

The test results  $\geq 45$ U/mL indicate that patients have rheumatoid arthritis, it is recommended to conduct further examinations.

The result only for clinical reference, comprehensive consideration should be combined with the clinical management of patients with symptoms / signs, medical history, other laboratory tests and treatment response.

All laboratory tests depend on random errors. If the test results are in doubt, or if they do not match the clinical symptoms, re-test the sample or confirm the results with other methods.

**Limitations**

Hemoglobin>5g/L, triglyceride>5mmol/L, bilirubin>97 $\mu$ mol/L will affect the test result.

**Performance Characteristics**

1. Linearity range: 5U/mL ~ 400U/mL.
2. Precision  
Test the control material by Anti-Cyclic Peptide Containing Citrulline Antibody Test Kit 2 times per day for 20 days (n=80)  
The data as below:

**a.**

HP-083/4-II POCT Immunoassay System					
Sample	Mean U/mL	Within-Run		Between-Run	
		S.D.	%C.V.	S.D.	%C.V.
Control 1	22.67	1.38	6.1	1.41	6.2
Control 3	72.63	2.91	4.0	2.98	4.1

**b.**

HP-AFS/3 Automatic Immunoassay System					
Sample	Mean U/mL	Within-Run		Between-Run	
		S.D.	%C.V.	S.D.	%C.V.
Control 1	22.68	1.41	6.2	1.43	6.3
Control 3	72.71	3.05	4.2	3.20	4.4

**c.**

HP-AFS/1 Automatic Immunoassay System					
Sample	Mean U/mL	Within-Run		Between-Run	
		S.D.	%C.V.	S.D.	%C.V.
Control 1	22.58	1.35	6.0	1.42	6.3
Control 3	71.49	4.00	5.6	4.29	6.0

**3.Methodology comparison**

Compared to Anti-CCP LIA(x) by test the same serum sample, the relative data as below:

HP-AFS/3 Automatic Immunoassay System				
Site No.	Sample Type	No.of Assays	Regression Line	Coefficient correlation
1	Serum	50	Y= 1.04X+0.10	0.97

The concentration of sample is about 5U/mL ~ 400U/mL.

**Precautions**

**⚠ Attention:**

Only for in vitro diagnostic.

Only for professional use.

All samples and reactive wastes are treated as sources of infection.

Do not use the kits beyond shelf life.

Do not mix different batches of reagents.

**⚠ Warning:**

To avoid error, do not forced to take out the cuvette from the device. Follow the device operation manual strictly, If the problem cannot be solved, contact the manufacturer for further technical support.

**SYMBOLS USED ON LABELS**

Symbol	Usage	Symbol	Usage
	Use-By date		Do not freeze
	Batch code		Biological risks
	Manufacturer		Do Not Reuse
	Temperature Limit		
	Contains sufficient for <n> tests		
	Do not use if package is damaged		
	Consult Instructions for use		
	Keep Away from Sunlight		
	In Vitro Diagnostic Medical device		
	Authorized Representative in the European Community		

**References**

- 1、Sehellekens GA, Visser H, De Jong BA, et al. The diagnostic properties of rheumatoid arthritis antibodies recognizing a cyclic citrullinated peptide. Arthritis Rheum, 2000, 43: 155-163.
- 2、Liao KP, Batra KL, Chibnik L, et al. Anti-CCP revised criteria for the classification of rheumatoid arthritis. Ann Rheum Dis, 2008, 30.
- 3、Yamane T, Hashiramoto A, Tanaka Y, et al. Easy and Accurate Diagnosis of Rheumatoid Arthritis Using Anti-Cyclic Citrullinated Peptide 2 Antibody, Swollen Joint Count, and C-Reactive Protein/Rheumatoid Factor. J Rheumatol, 2008, 35(3): 414-420.

**Approval Date&Revision Date**

Approval Date: Sep 9,2015  
Revision Date: May 6,2016  
Revision Date: May 1, 2017  
Revision Date: Feb 1,2019  
Revision Date: Apr 1,2019  
Revision Date: Jan 1, 2021  
Revision Date: Apr 1, 2021  
Revision Date: Jan 1, 2023  
Revision Date: Dec 22,2023

尺寸:24\*25cm展开尺寸,横向三折页再垂直方向两次对折