

TT3

Total Triiodothyronine Test Kit (Fluorescence Immunoassay Method) Instructions for Use

Version:A/3

REF HP-Palm-TT3-25

Manufacturer

Shijiazhuang Hipro Biotechnology Co.,Ltd.
No. 3 Building, Block C, Fangyi Science Park, No. 365 Huai'an East Road, Hi-tech Zone, Shijiazhuang, 050000 Hebei P.R. China
After sale service: 400-0191-606
www.hipro.us

EC REP Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Tel: +31644168999

Product Name

General Name: Total Triiodothyronine Test Kit
(Fluorescence Immunoassay Method)

Package Specification

25 Tests/Kit

Note: other specifications are available, 10Tests/Kit, 50 Tests/Kit, 100 Tests/Kit.

Intended Use

This Total Triiodothyronine (TT3) Test Kit is intended for in vitro diagnostic use in the quantitative determination of total triiodothyronine (TT3) in human serum, and plasma (EDTA and lithium heparin) or whole blood, using the fluorescence immunoassay analyzer. Measurements of the total triiodothyronine are used in diagnosis of thyroid or pituitary disorders.

For in vitro diagnostic use only .

For prescription use only.

Test Principle

The TT3 test kit is a rapid quantitative assay based on the principle of double antigens competitive fluorescence immunoassay.


Pad contains mouse monoclonal anti-human total T3 antibody.

Membrane coats with total T3 recombinant antigen on T line and goat polyclonal anti-mouse IgG antibodies on C line.

The dissociated total T3 antigen in the sample competes with the total T3 recombinant antigen on T line ,then the total T3 antigen reacts with mouse monoclonal anti-human total T3 antibody in the pad to form antigen-antibody complexes. The total T3 recombinant antigen on T line also reacts with mouse monoclonal anti-human total T3 antibody in the pad to form the fluorescent antigen-antibody complexes. The concentration between total T3 in the sample and fluorescent complexes is in reverse correlation .The rest mouse monoclonal anti-human total T3 antibody reacts with goat polyclonal anti-mouse IgG antibodies to form fluorescent antibody complexes.The fluorescence immunoassay analyzers will calculate the total T3 concentration by detecting the fluorescence intensity of T and C line.

Material Provided

Component	25Kits
Sample Diluent A	25
Sample Diluent B	25
Calibration Cassette	2
Test Cassette	25
Instructions for Use	1



Material Needed but Not Provided

1. Timing device
2. Pipette


Storage&Stability

Store test kit at: 2°C-30°C until the expiration date indicated on the label.

Avoid exposure to direct sunlight. The test kit is stable for 18months when unopened.

Perform the test within 1 hour when opened.

Applicable Instrument



**Palm F Fluorescence
Immunoassay Analyzer**

Specimen type

The specimen type of Hipro TT3 test kit is human serum, plasma (EDTA or lithium heparin) or whole blood

-Take blood samples in observance of the standard precautions for the withdrawal of biological fluids.

-Do not use samples that have remained at room temperature for more than 8 hours.

-Freeze the sample at -20°C or lower temperature if dosage is not performed within 72 hours.

-Due to possible evaporation effects, samples on the analyzer should be measured within 2 hours.

Storage

Fresh collected specimens are stable if stored at 2°C-8°C for up to 24 hours. For longer storage, aliquot, cap tightly, and freeze at -20°C for up to 3months.

Avoid repeated freezing and thawing.

Avoid hemolysate.

Required volume

30μL is required for serum, plasma or whole blood, used for dilution. 80μL is required of diluted sample for each determination. This volume does not include the dead volume(unusable volume in the sample container) or the additional volume required to make replicates or other tests to be performed on the same sample.


TT3 Test Kit Instructions Guide

Before You Begin

- 1a** Please read the TT3 IFU and applicable instrument User Manual carefully before operation. Please follow the instructions strictly.




- 1b** A timing device (phone, clock or timer) and a pipette are required but not provided.



Before You Begin


- 1c** Allow the sample back to room temperature




- 1d** Ensure the test kit is at room temperature for at least 30 minutes prior to use.



- 1e** Check on the expiration date on labels of different component.




- 1f** Make sure the lot number on calibration cassette and test cassette are same.
Do not misuse the cassettes with different lot.




Test Procedures


2a Turn on the applicable instrument according to the User Manual.




2b Insert the calibration cassette into the analyzer as the direction shows, to let the built-in parameter be read.




2c Use the pipette to add 30μL specimen(serum/plasma/whole blood) into the sample diluent A.



2d Use the pipette to add 45μL sample diluent B into sample diluent A tube and wait for 5-10seconds.




2e Wait for 15minutes to let sample and diluent completely mix.



Test Procedures


3a Use the pipette to add 80μL mixed liquid from sample diluent A tube onto sample well of test cassette.



3b Wait for 10 minutes.



3c Insert the test cassette into the analyzer and the analyzer will display the results automatically.



! Note: Avoid bubbles when collecting samples.

Reference Value

The normal reference range is 0.92-2.79nmol/L.

159 healthy patients were measured and the results (95%) meet the linearity range 0.92-2.79nmol/L.

The value is indicative only and may differ from other published values as a result of differences in methods and in the population being studied. It is recommended that each laboratory establish its own reference range.

Performance Characteristics

- Strip width**
The width of strip should be $\geq 2.5\text{mm}$.
- Liquid flow speed**
The speed should be $\geq 10\text{mm/min}$.
- Linearity**
Within range 0.61-9.22 nmol/L, the correlation coefficient (r) should be ≥ 0.990 .
- Within-lot precision**
The coefficient of variation (CV) should be $\leq 15\%$.

- Between-lot precision**
The coefficient of variation (CV) should be $\leq 20\%$
- Accuracy**
The relative deviation should be $\leq 15\%$.
- Limit of blank(LoB)**
The LoB should be $\leq 0.61\text{nmol/L}$.

Precautions and limitations

For in vitro diagnostic use

Only experienced laboratory personnel should use this test and handling all laboratory reagents should be based on normal precautions required.

Safety precautions

- Do not smoke, eat or apply cosmetics in areas in which patients' samples or kit reagents are handled.
- Laboratory gloves should be worn while handling patients' samples or disposing of solid or liquid wastes.
- Disposal of all waste material should be in accordance with local guidelines.
- Avoid electromagnetic or vibration environment for analyzer.
- Do not insert the polluted cassette into analyzer.

Potential biohazard warning

Some reagents of Hipro test kits contain material of animal origin, even if they are certified as deriving from healthy animals, it is recommended to handle them with the same precaution used for potentially infectious samples.

Limitations-interference

As with all immunoassays, the results of this test can be influenced by factors present in some patients' specimens. The factors include hemoglobin $>6\text{mg/mL}$, triglyceride $>10\text{mg/mL}$, bilirubin $>0.6\text{mg/mL}$, For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Safety Phases

- S 37 Wear suitable gloves
- S 60 This material and /or its container must be disposed of as hazardous waste.

Symbols used on labels

Symbol	Usage
	Use-By date
	Batch code
	Manufacturer
	Keep Away from Sunlight
	Temperature Limit
	In Vitro Diagnostic Medical device
	Authorized Representative in the European Community
	CE Mark
	Consult Instructions for use
	Do not freeze
	Contains sufficient for <n> tests
	Biological risks
	Do not use if package is damaged
	Date of manufacture
	Do Not Reuse

References

- Zhang Haiyan, Dong Lei, Wu Caijun, et al. [J]. Chin J emergency med, 2013,22 (10) : 1,132-1,135.
- Wu Guanhui, Kong Fanzhen, Cheng Qingzhang, et al. [J]. Chin J neuromed, 2014,13 (11) : 1 139-1 142. (in Chinese)

Approval Date&Revision Date

Approval Date: Sept 9,2019

Revision Date: Jul 15, 2021

Revision Date: Jul 25,2023

Revision Date: Dec 22, 2023

工艺：70g双胶纸，双面印刷，横向4折页，竖向对折 展开尺寸：32*25cm