CHORUS

ANTI-CCP

REF 86094

REF 86094/12

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17. BIBLIOGRAFIA

1. The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Pruijnet et al., Arthritis Research & Therapy 2010, 12:203
2. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273–281
3. 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
4. Anti-CCP antibodies: the past, the present and the future Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

INSTRUCTIONS FOR USE

CHORUS ANTI-CCP

For the semiquantitative determination of IgG antibodies anti-CCP

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of IgG class antibodies against cyclic citrullinated peptide (anti-CCP) in human serum, using a disposable device applied on the chorus and Chorus TRIO instruments.

2. INTRODUCTION

Rheumatoid Arthritis (RA) is the most common inflammatory joints disease with a prevalence ranging between 0.5% and 1% worldwide.

Antibodies against cyclic citrullinated peptide (anti-CCP) are highly specific for Rheumatoid Arthritis (RA) and their detection is one of the ACR-criteria for classification of RA (ACR : American College of Rheumatology) .

3. PRINCIPLE OF THE METHOD

The anti-CCP device is ready to use for the detection of antibodies against CCP in the CHORUS instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen is bound to the solid phase.

The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human immunoglobulins conjugated to horse radish peroxidase.

The unbound conjugate is eliminated and the TMB substrate is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus instruments.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-

approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.

7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
11. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
12. Pipette the test serum (50 µl) in well 1 of the device (see figure).
13. Do not use the device after the expiry date.
14. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86094).


The kit is sufficient for 12 tests. (REF 86094/12).

DD DEVICES

6 packages each containing 6 devices (REF 86094).

2 packages each containing 6 devices (REF 86094/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL coated with highly purified synthetic peptides containing citrulline residues

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: saline proteic solution with Proclin (0.1%)

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 ml

Contents: Diluted human serum, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: Diluted human serum, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Heat-inactivation can rise to erroneous results. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results. Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used.

The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense (add) 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result in AU/ml calculated on the basis of a lot-dependent graph stored in the instrument. The results can be interpreted as follows:

POSITIVE: when the concentration in the sample is > 18 AU/ml

NEGATIVE: when the concentration in the sample is < 12 AU/ml.

DOUBTFUL/EQUIVOCAL: when the concentration in the sample is between 12 and 18 AU/ml

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

11. LIMITATIONS

All the positive results require careful interpretation. The test can not be used alone for a clinical diagnosis. A negative result does not preclude the possibility of disease. The test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range of 3.0-300 AU/ml.

For samples > 300AU/ml retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

13. ANALYTICAL SPECIFICITY

5 samples were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)
Bilirubin (4.5-45 mg/dl)
Triglycerides (10-250 mg/dl)
Hemoglobin (5-30 mg/ml)

The presence in the serum sample of interfering substances described above do not affect the test result.

14. CROSS REACTION

232 samples from the sera collection were tested, their clinical characteristics are listed in the following table:

Samples Group	n° TOT	Neg < 12AU/ml	Dub 12-18 AU/ml	Pos > 18AU/ml
Healthy Donors	90	88	0	2
Routine samples	58	58	0	0
ANA	14	13	0	1
Thyroiditis	20	20	0	0
Rubella IgM Positive	5	5	0	0
Legionella IgM Positive	4	4	0	0
Chlamydia Positive	5	5	0	0
Salmonella Positive	3	3	0	0
Toxoplasma g. IgM Positive	6	6	0	0
Mycoplasma p. IgM Positive	9	9	0	0
Anti-streptolysin Positive	3	3	0	0
EBV IgM Positive	4	4	0	0
Parvovirus IgM Positive	5	5	0	0
LUE IgM Reactive	6	5	1	0

The clinical specificity is:

Healthy Donors 88/90 = 97.8% (CI_{95%}: 92.3-99.4%)

Routine Samples 58/58 = 100% (CI_{95%}: 93.8-100%)

Autoimmune Disease 33/34 = 97.1% (CI_{95%}: 85.1-99.5%)

Infectious Disease 50/50 = 100% (CI_{95%}: 92.7-100%)

15. METHOD COMPARISON

174 samples have been tested with Diesse kit and with a competitor kit. Data are summarized in the following table :

		Reference		
		+	-	Tot.
Diesse	+	72	2	74
	-	9	91	100
	Tot.	81	93	174

Percent Positive Agreement (~Diagnostic Sensitivity): 88.9%
CI_{95%}: 80.2 – 94.0

Percent Negative Agreement: (~Diagnostic Specificity): 97.8%
CI_{95%}: 92.5 – 99.4

The agreement between the two methods is high with a Cohen's Kappa of 0.88.

16. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Average (AU/ml)	CV%	Average (AU/ml)	CV%
1	64.9	2.8	65.1	4.1
2	26.3	2.5	26.5	13.4
3	27.0	3.6	26.7	12.4
4	15.9	7.5	15.6	13.3
5	16.8	5.8	15.8	11.5
6	8.1	5.6	7.5	12.8

Sample	Between lots		Between Instruments	
	Average (AU/ml)	CV%	Average (AU/ml)	CV%
1	63.3	3.3	63.3	3.5
2	25.2	7.6	25.2	7.1
3	25.5	7.7	25.5	7.4
4	15.7	7.1	15.7	8.2
5	15.8	6.3	15.9	6.4
6	7.2	6.3	7.2	7.1

17. REFERENCES

1. The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Puijnet et al., Arthritis Research & Therapy 2010, 12:203
2. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273–281
3. 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
4. Anti-CCP antibodies: the past, the present and the future
Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

Manufactured by
DIESSE Diagnostica Senese S.p.A.
Strada dei laghi 39
53035 Monteriggioni (Siena)
Italy

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- The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Pruijnet et al., Arthritis Research & Therapy 2010, 12:203
- Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273-281
- 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
- Anti-CCP antibodies: the past, the present and the future Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

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Percent Positive Agreement (~Sensibilidad de Diagnóstico):
88.9% CI_{95%}: 80.2 – 94.0

Percent Negative Agreement: (~Especificidad de Diagnóstico):
97.8% CI_{95%}: 92.5 – 99.4

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- 1.The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Pruijnet et al. Arthritis Research & Therapy 2010, 12:203
2. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies.G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273–281
- 3.2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
4. Anti-CCP antibodies: the past, the present and the future

Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

Fabricado por
DIESSE Diagnostica Senese S.p.A.
Strada dei laghi 39
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAPHIE

- The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Pruijnet et al., Arthritis Research & Therapy 2010, 12:203
- Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies.G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273–281
- 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
- Anti-CCP antibodies: the past, the present and the future Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

Fabriqué par
DIESSE Diagnostica Senese S.p.A.
Strada dei laghi 39
53035 Monteriggioni (Siena)
Italy

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Percent Positive Agreement (~Sensibilitatea Diagnosticului):
88.9% CI_{95%}: 80.2 – 94.0

Percent Negative Agreement: (~Specificitatea Diagnosticului):
97.8% CI_{95%}: 92.5 – 99.4

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17. BIBLIOGRAFIE

1. The use of citrullinated peptides and proteins for the diagnosis of rheumatoid arthritis Ger JM Pruijnet et al., Arthritis Research & Therapy 2010, 12:203
2. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies.G A Schellekens et al. J Clin Invest. 1998 January 1; 101(1): 273–281
3. 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative D. Aletaha et al. Ann Rheum Dis 2010;69:1580-1588 doi:10.1136/ard.2010.138461
4. Anti-CCP antibodies: the past, the present and the future
Walther J. van Venrooij et al. Nature Reviews Rheumatology 7, 391-398 (July 2011)

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