

BLOOD GAS CONTROL - LEVEL 2 (BG CONTROL 2)

CAT. NO. BG5002 **LOT NO.** 291BG
SIZE: 30 x 1.8 ml **EXPIRY:** 2025-06-28
GTIN: 05055273227116

INTENDED USE

This product is intended for in vitro diagnostic use, in the quality control of Blood Gas analysis.

DEVICE DESCRIPTION

The Blood Gas Controls are supplied at 3 levels, 1, 2 and 3. Target values and ranges are supplied for the following analytes: Calcium, Chloride, Glucose, Lactate, pCO₂, pH, pO₂, Potassium, Sodium and Total CO₂.

SAFETY PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. Health and Safety Data Sheets are available on request.

Not suitable for instruments that do not recommend products with artificial dye.

STORAGE AND STABILITY

UNOPENED: The product is stable to expiration date when stored at +2°C to +8°C. Avoid exposure to freezing and temperatures greater than +30°C.

OPENED: For pH/blood gas values, the control should be analysed within 1 minute of opening. For electrolyte measurements, the control should be analysed within 1 hour after opening.

PREPARATION FOR USE

The Blood Gas Control should be brought to +20°C to +23°C before use. Allow at least 4 hours for ampoules to equilibrate to this temperature, prior to testing. Before use, hold the ampoule at the top and bottom (with forefinger and thumb) and shake 15 - 20 times to mix the solution. Tap the ampoule to restore the liquid to the bottom of the ampoule. Open the ampoule by snapping off the tip at the score. Use gauze, tissue, gloves or an appropriate ampoule opener to protect fingers from cuts. Immediately introduce the liquid from the ampoule to the analyser.

MATERIALS PROVIDED

Blood Gas Control - Level 2 30 x 1.8 ml

ASSIGNED VALUES

Due to the variation caused by test equipment, test reagents and laboratory technique, the quoted ranges are provided for guidance. It is recommended that these ranges are used until each laboratory has established its own ranges, based on individual laboratory requirements.

Each batch of Blood Gas Control is submitted to a number of external laboratories and values are assigned from a consensus of results obtained by these laboratories.

EC	REP
----	-----

Randox Teoranta, Meenmore,
Dungloe, Donegal,
F94 TV06, Ireland

17 Oct '24 me

Table of Content

Method

2

Method		Blood Gas Control - Level 2 (BG CONTROL 2)						
Lot. No: 291BG Cat. No: BG5002 Expiry: 2025-06-28								
Size: 30 x 1.8 ml		Range						
Analyte	Unit	Target	Low	High	1SD	2SD	Method	
Calcium	mg/dl	4.09	3.68	4.50	0.205	0.410	Colorimetric	
	mmol/l	1.02	0.918	1.12	0.050	0.100		
	mg/dl	4.17	3.75	4.59	0.210	0.420	Ion Selective Electrode	
	mmol/l	1.04	0.936	1.14	0.050	0.100		
Chloride	mmol/l	91.9	84.5	99.3	3.70	7.40	ISE Indirect	
Glucose	mg/dl	90.5	76.9	104	6.75	13.5	Colorimetric	
	mmol/l	5.02	4.27	5.77	0.375	0.750		
	mg/dl	90.5	76.9	104	6.75	13.5	Enzymatic Electrode	
	mmol/l	5.02	4.27	5.77	0.375	0.750		
	mg/dl	90.3	76.8	104	6.85	13.7	Glucose Oxidase	
	mmol/l	5.01	4.26	5.76	0.375	0.750		
	mg/dl	90.6	77.0	104	6.70	13.4	Hexokinase	
	mmol/l	5.03	4.28	5.78	0.375	0.750		
Lactate	mg/dl	30.5	25.0	36.0	2.75	5.50	Colorimetric	
	mmol/l	3.38	2.77	3.99	0.305	0.610		
	mg/dl	33.4	27.4	39.4	3.00	6.00	Enzymatic Electrode	
	mmol/l	3.71	3.04	4.38	0.335	0.670		
	mg/dl	30.5	25.0	36.0	2.75	5.50	Optical Fluorescence	
	mmol/l	3.39	2.78	4.00	0.305	0.610		
	pCO2	kPa	5.36	4.29	6.43	0.535	1.07	Ion Selective Electrode
		kPa	5.15	4.12	6.18	0.515	1.03	Optical Fluorescence
pO2	kPa	17.3	13.8	20.8	1.75	3.50	Ion Selective Electrode	
	kPa	16.9	13.5	20.3	1.70	3.40	Optical Fluorescence	
Potassium	mmol/l	3.89	3.58	4.20	0.155	0.310	ISE method - direct	
Sodium	mmol/l	140	133	147	3.50	7.00	ISE method - direct	