

TOTAL ANTIOXIDANT STATUS (TAS CONTROL) CONTROL SERUM

CAT. NO. NX 2331
LOT NO. 394NX
SIZE: 10 x 5 ml
EXPIRY: 2025-03-28
GTIN: 05055273204728

This Control Serum is human in origin.

INTENDED USE

This product is intended for use in the control of the Randox Total Antioxidant Status assay, Catalogue Number: NX2332.

PREPARATION

Open the vial very carefully, avoiding any loss of material and reconstitute in exactly 5 ml of freshly distilled water at +20°C to +25°C. Replace the rubber stopper, close vial and leave to stand for 30 minutes before use.

STABILITY

UNOPENED: Stable to expiry, when stored at +2°C to +8°C.

OPENED: Once reconstituted, the control is stable for 2 days at +2°C to +8°C or 12 hours at +15°C to +25°C.

SAFETY

Human source material, from which this product has been derived, has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests.

However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

ASSAYED VALUES

Target value (mmol/l)	Range (mmol/l)
1.81	1.45 – 2.17

The presence of a vertical bar in the margin indicates a technical update from the previous revision

EC REP Randox Teoranta, Meenmore,
Dungloe, Donegal,
F94 TV06, Ireland

Rev. 04 Nov '24 Id

THIS PAGE IS INTENTIONALLY BLANK