VIASURE

Real Time PCR Detection Kit

Monkeypox virus

RUO

For Research Use Only (RUO)

This product has no declared clinical intended purpose and is not for clinical diagnostic use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

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EN For download IFUS from other languages, please enter in **certest.es/viasure/labeling**. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended use.....	6
2.	Summary and Explanation	6
3.	Principle of the procedure	7
4.	Reagents provided	7
5.	Reagents and equipment to be supplied by the user	7
6.	Transport and storage conditions.....	8
7.	Precautions for users	8
8.	Test procedure	9
8.1.	Specimen collection, transport and storage	9
8.2.	DNA extraction	9
9.	Result interpretation	10
10.	Limitations of the test	11
11.	Quality control.....	12
	ANNEX 1	13
A1.1	Principle of the procedure.....	13
A1.2	Reagents provided.....	13
A1.3	Test procedure	15
A1.3.1	Lyophilized positive control.....	15
A1.3.2	PCR protocol.....	15
	ANNEX 2	16
A2.1	Principle of the procedure.....	16
A2.2	Reagents provided.....	16
A2.3	Test procedure	17
A2.3.1	Lyophilized positive control.....	17
A2.3.2	Lyophilized reaction mix tube.....	17
A2.3.3	PCR protocol.....	17

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ENGLISH

1. Intended use

VIASURE Monkeypox virus Real Time PCR Detection Kit is a real-time PCR test designed for the qualitative identification of DNA from Monkeypox virus in clinical samples from individuals suspected of Monkeypox virus infection by their healthcare professional (HCP). This test is intended to be used for research purposes, without any medical objective is not regarded as devices for performance evaluation. DNA is extracted from clinical specimens, amplified using real-time PCR and detected using fluorescent reporter dye probes specific for Monkeypox virus.

2. Summary and Explanation

Monkeypox virus (MPXV) is an enveloped double-stranded DNA virus. It belongs to the *Orthopoxvirus* genus of the *Poxviridae* family. The *Orthopoxvirus* genus also includes Vaccinia virus, Cowpox virus, Variola virus and several other, animal-related, poxviruses. Two phylogenetically distinct clades of MPXV have been identified through genomic sequencing: the Central African (Congo Basin) clade and the West African clade. Typically, the Central African MPXV is associated with more severe disease, higher mortality, and more frequent human-to-human transmission.

MPXV is transmitted to humans through contact with an infected animal or human, or with material contaminated with the virus. The virus enters the body through broken skin, the respiratory tract or the mucous membranes. The incubation period of monkeypox is usually from 6 to 13 days but can range from 5 to 21 days. Virus transmission through direct or indirect contact with live or dead animals is assumed to be the main factor for human MPX infections. This may occur by bite or scratch, bush meat preparation, direct contact with body fluids or lesions from an infected animal or contaminated material (indirect contact). Eating inadequately cooked meat of an infected animal is an additional possible risk factor. Similar to smallpox, human-to-human transmission of MPXV occurs mostly through large respiratory droplets during direct and prolonged face-to-face contact. In addition, MPXV can be transmitted by direct contact with body fluids of an infected person or with contaminated objects, such as bedding or clothing.

Human monkeypox often begins with a combination of the following symptoms: fever, headache, chills, exhaustion, asthenia, lymph node swelling, back pain and muscle aches. Commonly, within one to three days after onset of fever, the patient develops a rash, which tends to first appear on the face and then spreads to other parts of the body, including hands and feet. The cutaneous lesions often first present as macules, evolving successively to papules, vesicles, pustules, crusts and scabs. The main difference between smallpox and MPX is that MPXV causes lymphadenopathy (e.g. in the cervical or inguinal region) while smallpox virus and chickenpox virus usually do not. The study of Immune markers has provided evidence of asymptomatic MPX infections in individuals vaccinated against smallpox and others that had not been vaccinated. For most people, MPX is a self-limited disease, typically lasting two to four weeks and resulting in complete recovery.

The key objectives of surveillance and case investigation for monkeypox in the current context are to rapidly identify cases, clusters, and the sources of infection as soon as possible in order to provide optimal clinical care, isolate cases to prevent further transmission, identify and manage contacts and tailor effective control and prevention

methods based on most commonly identified routes of transmission; in this sense, Polymerase chain reaction (PCR) is the preferred laboratory test given its accuracy and sensitivity.

3. Principle of the procedure

VIASURE *Monkeypox virus* Real Time PCR Detection Kit is designed for the identification of Monkeypox virus in clinical samples. After DNA isolation, the detection of Monkeypox virus is performed by the amplification of a conserved region of *G2R* and *F3L* genes for Monkeypox virus using specific primers and fluorescent-labelled probes.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Monkeypox virus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer and polymerase) in a stabilized format, as well as an **endogenous internal control** to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an **endogenous internal control (EIC)** (hemoglobin subunit beta gene, *HBB* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

4. Reagents provided

VIASURE *Monkeypox virus* Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open and rotor-gene format with internal control products and Annex 2 for tube format with internal control products.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Monkeypox virus* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2).
- DNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments). Only for Rotor Gene format (Annex 1).

To check thermocycler compatibility and most common detection channels consult website www.certest.es.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the *Monkeypox virus* Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

7. Precautions for users

- This VIASURE *Monkeypox virus* Real Time PCR Detection Kit is for Research Use Only
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-MPX113LRUO, VS-MPX113HRUO, VS-MPX136RUO and VS-MPX172RUO). Remove any excess air in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-MPX101RUO, VS-MPX136RUO and VS-MPX172RUO (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Monkeypox virus* Real Time PCR Detection Kit and any additional reagents or equipment required for testing are not contaminated. Always avoid microbial and ribonuclease

(RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.

- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink, or smoke in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions, and procedures.

8. Test procedure

Please see Annex 1 for Open and Rotor-Gene format with endogenous internal control products Test Procedure and Annex 2 for Tube format with endogenous internal control products Test Procedure.

8.1. Specimen collection, transport and storage

The VIASURE *Monkeypox virus* Real Time PCR Detection Kit has been designed to test clinical samples.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For details, refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>), the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94), and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

8.2. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from respiratory samples, you can use your manual or automatic routine optimized system, or any commercially available DNA extraction kit and follow the manufacturer's instructions

9. Result interpretation

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (EIC) signal to verify the correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

It is recommended to set the threshold values for each channel (target) independently by the end-user. Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	Monkeypox virus (FAM) ¹	Endogenous Internal Control (EIC) (HEX) ²	Interpretation of Controls
Positive Control (PC)	≤40	≤40	Valid
Negative Control (NC)	>40 or no signal	>40 or no signal*	Valid

Table 1. Expected Performance of Controls

* It is possible that non-specific amplifications in late Ct values (Ct>35) may appear due to human manipulation of the process, since human HBB gene is a human housekeeping gene and is expected to be present in all nucleated human cells.

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The positive template control includes human housekeeping HBB gene target; therefore, amplification signals are observed in all target channels, including the Endogenous Internal Control (EIC).

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following table, read, and analyze the results:

Monkeypox virus (FAM)	Endogenous Internal Control (EIC) (HEX)	Interpretation for patient's individual samples	
≤40	≤40 or no signal ¹	Valid	Monkeypox virus DNA Detected
>40 or no signal	≤35 ²	Valid	Targets not Detected ²
>40 or no signal	>35 or no signal ²	Invalid	Test Failure – Repeat Testing ²

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curves.

1 The Endogenous Internal Control (EIC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 In the case of Monkeypox virus target sites negative, EIC must show an amplification signal with Ct ≤35. The Ct value could be very variable due to the endogenous internal control is a human housekeeping gene that should be present in all human nucleated cells in the original samples. If there is an absence of signal or Ct value > 35 of the endogenous internal control, the result is considered as 'invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use; the extraction process used by the user; to verify the correct performance of each qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- repeat qPCR with the same isolated DNA sample, or
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

10. Limitations of the test

- The test must be used for RUO purposes (Exclusive Use in Research), without any medical objective.
- The quality of the test depends on the quality of the sample: nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Monkeypox virus, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the viral DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown Monkeypox virus variants.
 - A viral load in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics, or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.

If in doubt, refer to section 9 to check the correct interpretation of the results.

- A positive test result does not necessarily indicate the presence of viable virus and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets virus sequences.
- Negative results do not preclude Monkeypox virus infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by Monkeypox virus have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other diseases are negative and the patient's clinical presentation and epidemiological information suggest that Monkeypox virus infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.
- Some samples may fail to exhibit *HBB* amplification curves due to low human cell numbers in the original clinical sample. A negative IC signal does not preclude the presence of Monkeypox virus in a clinical specimen.

11. Quality control

VIASURE Monkeypox virus Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the endogenous internal control (EIC) in each well confirms the correct performance of the technique.

ANNEX 1

OPEN FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 1 x 8-well strips, low profile	VS-MPX101LRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 1 x 8-well strips, high profile	VS-MPX101HRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 6 x 8-well strips, low profile	VS-MPX106LRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 6 x 8-well strips, high profile	VS-MPX106HRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 12 x 8-well strips, low profile	VS-MPX112LRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 12 x 8-well strips, high profile	VS-MPX112HRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 96-well plate, low profile	VS-MPX113LRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 96-well plate, high profile	VS-MPX113HRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 2 x 4-well strips, Rotor-Gene®	VS-MPX101RUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 9 x 4-well strips, Rotor-Gene®	VS-MPX136RUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO 18 x 4-well strips, Rotor-Gene®	VS-MPX172RUO

Table A1 1. References

A1.1 Principle of the procedure

VIASURE Monkeypox virus Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an Endogenous Internal Control (EIC) (hemoglobin subunit beta gene, *HBB* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

Target	Channel	Gene
Monkeypox virus	FAM	G2R and F3L genes
Endogenous Internal control (EIC)	HEX	HBB gene (β-globin)

Table A1 2. Target, channel and genes.

A1.2 Reagents provided

VIASURE Monkeypox virus Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3, A1.4 and A1.5. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. Table A1.5 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips. (Consult the thermocycler compatibility on CerTest's website www.certest.es).

Reagent/Material	Description	Colour	Amount
<i>Monkeypox virus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	1/6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	1/6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX01LRUO, VS-MPX101HRUO, VS-MPX106LRUO, VS-MPX106HRUO, VS-MPX112LRUO and VS-MPX112HRUO.

Reagent/Material	Description	Colour	Amount
<i>Monkeypox virus</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref VS-MPX113LRUO and VS-MPX113HRUO.

Reagent/Material	Description	Colour	Amount
<i>Monkeypox virus</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	Transparent	2/9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Caps for sealing wells during thermal cycling	Transparent	2/9/18 x 4-cap strip

Table A1 4. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX101RUO, VS-MPX136RUO and VS-MPX172RUO. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

A1.3 Test procedure

A1.3.1 Lyophilized positive control

Monkeypox virus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Monkeypox virus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Monkeypox virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A1 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Monkeypox virus*) and HEX (Endogenous Internal Control (EIC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

ANNEX 2

TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO, 4 tubes x 24 reactions	VS-MPX196TRUO
VIASURE Monkeypox virus Real Time PCR Detection Reagents RUO, 2 tubes x 24 reactions	VS-MPX148TRUO

Table A2. 1. References.

A2.1 Principle of the procedure

VIASURE Monkeypox virus Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an Endogenous Internal Control (EIC) (hemoglobin subunit beta gene, *HBB* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

Target	Channel	Gene
Monkeypox virus	FAM	G2R and F3L genes
Endogenous Internal control (EIC)	HEX	<i>HBB</i> gene (β-globin)

Table A2. 2. Target, channel and genes. *Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website www.ceratest.es.

A2.2 Reagents provided

VIASURE Monkeypox virus Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3 and A2.4.

Reagent/Material	Description	Colour	Amount
Monkeypox virus Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Monkeypox virus Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE Monkeypox virus Real Time PCR Detection Kit with Ref. VS-MPX196TRUO.

Reagent/Material	Description	Colour	Amount
<i>Monkeypox virus</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format	White	2 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Monkeypox virus</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table A2. 4. Reagents and materials provided in VIASURE *Monkeypox virus* Real Time PCR Detection Kit with Ref. VS-MPX148TRUO.

A2.3 Test procedure

A2.3.1 Lyophilized positive control

Monkeypox virus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Monkeypox virus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Monkeypox virus* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *Monkeypox virus* Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Monkeypox virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website www.certest.es).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 sec	95°C
	Annealing/Extension (Data collection*)	50 sec	60°C

Table A2. 5. PCR protocol.

Fluorogenic data should be collected during the extension step (*) through the FAM (*Monkeypox virus*) and HEX (Endogenous Internal Control (EIC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website www.certest.es).

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Bibliography/Bibliografía

1. Yu Li *et al.* Real-time PCR assays for the specific detection of monkeypox virus West African and Congo Basin strain DNA. *Journal of Virological Methods*, 2010.169 (2010) 223–227.
2. Rinat A. Maksyutov *et al.* Species-specific differentiation of variola, monkeypox, and varicella-zoster viruses by multiplex real-time PCR assay. *Journal of Virological Methods*, 2106.236 (2016) 215–220.
3. Li D *et al.* Evaluation of the GeneXpert for Human Monkeypox Diagnosis. *American Journal of Tropical Medicine and Hygiene*. 2017 Feb 8;96(2):405-410.
4. Shchelkunov SN *et al.* Species-specific identification of variola, monkeypox, cowpox, and vaccinia viruses by multiplex real-time PCR assay. *Journal of Virological Methods*. 2011 Aug;175(2):163-9.
5. European Centre for Disease Prevention and Control. Available from <https://www.ecdc.europa.eu/en/all-topics-z/monkeypox/factsheet-health-professionals> Accessed May 2022.
6. World Health Organization. Monkeypox. Available from https://www.who.int/health-topics/monkeypox#tab=tab_1. Accessed May 2022.

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Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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