Factor Deficient Plasmas

(PT Based - Factors II, V, VII, X)

- REF 5191
- REF 5192
- REF 5195
- REF 5790
- REF 5791
- REF 5792
- REF 5795
- Factor V Deficient Plasma (Congenital)
- Factor VII Deficient Plasma (Congenital)
- Factor X Deficient Plasma (Congenital)
- Factor II Deficient Plasma (Immunodepleted)
- Factor V Deficient Plasma (Immunodepleted)
- Factor VII Deficient Plasma (Immunodepleted)
- Factor X Deficient Plasma (Immunodepleted)

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0599P 2012/03 (12)

Factor Deficient Plasmas

(PT Based - Factors II, V, VII, X)

Instructions for use

en

INTENDED USE

The Factor Deficient Plasmas are intended for the quantitative determination of the respective factor in patients suspected of having a congenital or acquired deficiency of this coagulation protein. Quantitative measurement of individual coagulation factors by the one-stage method requires substrate plasma lacking the factor to be measured. A dilution of the test plasma is mixed with the factor deficient plasma and the clot time of the mixture determined. The degree of clot time correction with the patient plasma is compared to the correction with a reference material, allowing the % activity of the patient plasma to be determined¹. The Helena Biosciences Europe Factor Deficient Plasmas can be used on any instrument capable of performing PT-based factor assay testing. Refer to the instrument Operators' Manual for appropriate instructions.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample. All immunodepleted plasmas are HCV negative.

COMPOSITION			
Component	Content	Description	Preparation
Factor Deficient Plasmas	10 x 1 mL	All the Factor Deficient Plasmas listed above are derived from human plasma and contain less than 1% residual factor activity.	Reconstitute each vial with 1 mL of purified water. Swirl gently and allow to stand for 15 minutes. Mix well before use (do not shake).

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

- REF 5185
- REF 5375
- REF 5265H / 5265 / 5267 / 5269
- Calibration Plasma
- Owren's Buffer
- Thromboplastin LI

STORAGE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Once reconstituted, the reagent is stable for 8 hours when kept at *2 –*8°C. The lyophilised product should appear as a dry, straw coloured plug or pieces. Any deviation from this appearance may indicate signs of product deterioration.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes*.

PROCEDURE

Manual Method

Prepare all reagents as instructed with each pack. Pre-warm the recalcified thromboplastin mixture to *37°C before use.

1. Standard Curve Preparation:
- a. Prepare the following dilutions in Owren's Buffer:

Tube	Calibration Plasma (mL)	Owren's Buffer (mL)	Activity (%)
1	0.1	0.4	100
2	0.1	0.9	50
3	0.1	1.9	25
4	0.1	3.9	12.5

- b. Mix without shaking.

2. Patient Sample Preparation:

- a. Prepare a 1+ 4 dilution of the patient plasma or control plasma in Owren's Buffer.
- b. Mix without shaking.

3. Testing:

- a. Pipette, in duplicate, 0.1 mL of factor deficient plasma into a reaction tube.
- b. Add 0.1 mL of standard, patient or control plasma dilution and incubate at *37°C for 2 minutes.
- c. Add 0.2 mL of recalcified thromboplastin reagent while simultaneously starting a stopwatch.
- d. Determine the clot time for each of the standard, control or patient dilutions.
- e. Plot % Activity (X-axis) versus Mean Clot Time (Y-axis) for the standards on 2 cycle log-log graph paper.
- f. A straight line should be obtained.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Interpolate patient values from the curve. Calculate exact patient or control values by correcting for differences in calibration plasma values as follows:

Correct Activity (%) =

(Calibration plasma Reference Value / 100)

x

Interpolated Patient or Control Value

LIMITATIONS

The results obtained with Factor Deficient Plasmas depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{3,4,5}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

- REF 5301
- REF 5302
- Speciality Assayed Control N (SAC-N)
- Speciality Assayed Control A (SAC-A)

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Expected values for factor activity are 50–150%⁶.

PERFORMANCE CHARACTERISTICS

Each laboratory should establish its own performance data. Helena Biosciences Europe intrinsic factor assays are designed to give a linear standard curve from 10–150%. Within run and between run precisions are expected to be <5%, using a range of automated instruments.

REFERENCES

1. Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
2. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
3. Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
4. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
5. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.
6. Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path, Chicago, 36

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RIFERIMENTI

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, 37:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, 55:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, 59:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path., Chicago, 36

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