CHORUS SPECIFIC IgE CAPTURE

DIESSE
DIESSE

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Disponibili su richiesta.

14. BIBLIOGRAFIA

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111

8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M at al. Nature 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

INSTRUCTIONS FOR USE

CHORUS SPECIFIC IGE CAPTURE

For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS SPECIFIC IgE CAPTURE is an immunoassay kit for the quantitative automated determination of allergen-specific IgE class antibodies.

The test is performed on human serum using a disposable device applied on instruments CHORUS TRIO/CHORUS EVO.

The kit is intended to detect the presence of specific IgE antibodies as an aid in screening, confirmation/determination or monitoring of atopic allergic diseases. The specified target population to be tested with the CHORUS SPECIFIC IgE CAPTURE product is represented by patients with suspected presence of IgE class antibodies in serum.

It is to be used only by professional laboratory personnel.

2. INTRODUCTION

Immunoglobulins of class E (IgE) contain 2 types of polypeptide chains, their molecular weight is about 200.000 Da and they bind to the surfaces of mast cells and basophil granulocytes. The binding of allergens to cell-bound IgE causes the release of histamine and other vasoactive substances with a consequent series of events known as allergic reaction. It has been shown that most patients with atopic allergic diseases, such as atopic asthma, atopic dermatitis, and hay fever, show increased levels of IgE in the blood.

Allergen-specific IgE antibodies circulate in the blood of allergic patients, like subjects suffering from hay fever, atopic eczema, or asthma.

The research of allergen-specific IgE antibodies, in combination with other clinical data, is useful to ensure a proper identification of the substances responsible for allergies.

3. PRINCIPLE OF THE METHOD

The device is ready to use for the detection of allergen-specific IgE class antibodies, in the CHORUS TRIO/CHORUS EVO instruments.

The test is based on the ELISA principle (Enzyme Linked ImmunoSorbent Assay).

Anti-human IgE monoclonal antibody is bound to the solid phase. Specific immunoglobulins bind to the antibody following incubation with diluted human serum. After washing to remove unreacted proteins, incubation is performed with the specific biotinylated antigen. After washing to remove unreacted proteins, the Streptavidin-Peroxidase conjugate is added. The unbound conjugate is eliminated and the peroxidase substrate is added. The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow. The

colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test.

The results are expressed in Units (KU/L).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the instruments.
4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use, bring the package containing the devices to room temperature (18-30°C) for at least 30 minutes.

Use the device within 60 minutes.

1. Discard devices which show the substrate (well 4) blue colored.
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO/CHORUS EVO instrument; the Instructions for Use must be carefully followed, and the Instrument Operating Manual must be consulted.
5. Check that the instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.


7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the guidance in the chapter "Analytical Specificity").
11. Do not use the device after the expiry date.
12. **Make sure that the CHORUS TRIO instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: ALLERGEN-BIOTIN

Contents: allergen extract conjugated with biotin in containing Proclin 0.029% and gentamicin sulfate 0.0005%.

Position 7: MICROPLATE WELL

Coated with anti-human IgE monoclonal antibodies (Maximum concentration: 10 µg/ml).

Position 6: Empty

Position 5: SAMPLE DILUENT

Contents: Proteic solution containing Proclin 0.029% and gentamicin sulfate 0.0005%.

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: STOP SOLUTION

Contents: Sulfuric acid solution 0.3 M

Position 2: CONJUGATE

Contents: Streptavidin-Peroxidase in a buffer solution containing Proclin 0.029%.

Position 1: EMPTY WELL (CHORUS TRIO)

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.250 ml

Contents: Horse serum containing human IgE antibodies and preservative (Proclin 0.029%, gentamicin sulfate 0.0005%). Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.450 ml

Contents: Horse serum containing human IgE antibodies and preservative (Proclin 0.029%, gentamicin sulfate 0.0005%). Liquid, ready for use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced from a sample with known antigen concentration and diluted to reach a specific concentration, which is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different antigens levels.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- CHORUS EVO SANITIZING SOLUTION **REF** 82052
- CHORUS DECONTAMINATION BUFFER **REF** 82050 (CHORUS EVO)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at temperature ≤ -20°C for a least 2 months, can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

To assess the possibility of performing the test at the same time as other tests in the CHORUS line, please contact Customer Care.

CHORUS TRIO

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense in well no. 1 of each device:

SAMPLE	75 µl/device
CALIBRATOR	75 µl/device
POSITIVE CONTROL	75 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the Instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

CHORUS EVO

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Place the devices in the instrument. At each change of batch, use a device for the calibrator.
4. Place the primary test tube containing the sample, calibrator (if necessary), and the positive control in the designated rack.
5. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact Customer Care.

Tel: 0039 0577 319554
email: scientificsupport@diesse.it;
customer care@diesse.it

10. INTERPRETATION OF THE RESULTS

The instrument expresses the result in Units (KU/L) calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

KU/L	Interpretation
< 0.10	Absent or not detectable
0.10 – 0.28	Doubtful/Equivocal or Very Low
> 0.28	Sensitization/Positive

11. LIMITATIONS

The product should be used only by professional laboratory personnel.

The test is not suitable for samples different from human serum. All values obtained need careful interpretation that does not prescind from other indicators related to the same patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

Because there is no international standardization, both for the specific IgE units and for the reference allergens, the specific IgE concentrations obtained with various methods can be quite different. Therefore, it is advisable to take care when comparing or seeking to correlate different methods.

The specific IgE affinity constants may differ from allergen to allergen, as a consequence be aware that identical results for different allergens may not necessarily involve an equivalent clinical situation.

12. CALIBRATION RANGE

Calibration range 0.10 – 50.0 KU/L.

For samples > 50.0 KU/L retest the diluted (ex. 1:3) sample in the Negative Control/Sample Diluent (REF 83607-not supplied with the kit).

13. PERFORMANCES CHARACTERISTICS

Available on request.

14. REFERENCES

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Buckley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capron A, Dessaint JP, Capron M et al. Nature 1975; 253:474
14. Johansson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

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14. BIBLIOGRAFIA

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M at al. Nature 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

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