

PDM签审页

PDM版本：

PDM编码：

产品名称	出口_甲功试剂_甲状腺球蛋白抗体试剂盒 A-TG_说明书_英文		
库存编码	1041028	版本号	20200330
成品尺寸	210×297mm	单位	mm
印刷色	单色	允差	±2mm
材质	80g胶版纸，双面印刷		
备注			
设计			
审核			
批准			

# Antibodies to Thyroglobulin Detection Kit (Chemiluminescence Immunoassay)

## 【Product Name】

Antibodies to Thyroglobulin Detection Kit (Chemiluminescence Immunoassay)

## 【Package Specification】

Package Specification	Reagent Kit Composition
1×50 Tests/kit	1×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
1×50 Tests/kit (without Calibrator and Control)	1×Reagent
2×50 Tests/kit	2×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
2×50 Tests/kit (without Calibrator and Control)	2×Reagent
1×100 Tests/kit	1×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
1×100 Tests/kit (without Calibrator and Control)	1×Reagent
2×100 Tests/kit	2×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
2×100 Tests/kit (without Calibrator and Control)	2×Reagent
4×100 Tests/kit	4×Reagent, 2×Calibrator (High), 2×Calibrator (Low), 2×Control (Level 1), 2×Control (Level 2)
4×100 Tests/kit (without Calibrator and Control)	4×Reagent
1×200 Tests/kit	1×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
1×200 Tests/kit (without Calibrator and Control)	1×Reagent
2×200 Tests/kit	2×Reagent, 1×Calibrator (High), 1×Calibrator (Low), 1×Control (Level 1), 1×Control (Level 2)
2×200 Tests/kit (without Calibrator and Control)	2×Reagent

## 【Intended Use】

For quantitative determination of antibodies to thyroglobulin (A-TG) in human serum or plasma in vitro.

Thyroglobulin antibodies were discovered in 1956 by Roitt et al. in a serological study of autoimmune thyroiditis. A-TG target antigen Thyroglobulin (TG) is a soluble iodized glycoprotein synthesized and secreted by thyroid epithelial cells with a molecular weight of 660 kD consisting of 2748 amino acids. It is a biosynthesis precursor of triiodothyronine (T3) and thyroxine (T4), which is mainly stored in the thyroid follicular cavity in the form of colloid. Its content is very low in normal human serum. Thyroglobulin antibody is a disease indicator that is not pathogenic and is used to help diagnose whether a patient has an autoimmune thyroid disease. Thyroglobulin antibodies are generally considered to have no effect on the thyroid. Thyroglobulin antibodies bind to thyroglobulin and activate NK cells by interacting with the bound antibody via Fc receptors to attack target cells, resulting in destruction of thyroid cells. Thyroglobulin antibodies also affect the uptake, processing, and TG hydrolysis of TG antigens, thus affecting the autoimmune response of non-significant T cell epitopes, leading to the deterioration of autoimmune thyroid disease. 80% to 90% of patients with chronic thyroiditis and 60% of patients with hyperthyroidism can be strongly positive. Due to the heterogeneity of thyroglobulin, thyroglobulin antibodies are also detectable in elderly patients with other disorders, in patients with type I diabetes, and in patients with normal clinical and thyroid functions.

## 【Test Principle】

A-TG kit is based on the use of chemiluminescence immunoassay of indirect method for detection. The reagent kit comprises three parts of R1, R2 and R3. R1 is TG-coated magnetic particles, R2 is acridinium

ester-labeled antibodies, and R3 is PBS Buffer. Antibodies to thyroglobulin in samples to be tested combine with TG-coated particles. After rinse, combine with acridinium ester-labeled antibodies to human IgG, and form the antigen-antibody compound. The A-TG content in a sample is proportional to the relative light unit (RLU) detected by the system.

The system automatically performs the following steps:

- ◆ Put sample and reagent into the cuvette, incubate at 37°C;
- ◆ Separate magnetic particles, and rinse with washing buffer;
- ◆ Add acid trigger reagent and alkaline trigger reagent to stimulate chemiluminescence reaction.

## 【Main Components】

Composition	Main Components	Content
Reagent	R1	TG-coated magnetic particles
	R2	Acridinium ester-labeled antibodies to human IgG
	R3	PBS buffer
Calibrator (High, Low)	PBS buffer added A-TG	See the label
Control (Level 1, Level 2)	PBS buffer added A-TG	See the label

Note 1: The components in different batches of kits are not interchangeable.

Note 2: Different batches of Calibrators and Controls have different contents and batch specificity. The Calibrator fixed values and the target values of Controls are detailed in the bottle label.

Note 3: The necessary materials not provided are Acid Trigger Reagent, Alkaline Trigger Reagent, and Washing Buffer. The tests are carried out according to the instrument user manual and the instructions of the above reagents.

Note 4: Calibrators can be traced to national standard material.

## 【Storage Conditions & Shelf Life】

- The reagent kit shall be stored at 2°C ~ 8°C, away from sunlight, kept airtight and upright. For the shelf life refer to the label.
- After being used for the first time, the reagent can be stable for 28 days if stored on the instrument or at 2°C ~ 8°C. The calibrator and control after being opened for the first time can be stable for 28 days at 2°C ~ 8°C.

【Date of Manufacture & Expiry Date】 See the label.

## 【Applicable Instrument】

CM Series Chemiluminescence Immunoassay Analyzer and CSM Series Integrated System

## 【Sample Requirements】

- Specimens for tests are serum or plasma.
- Adopt correct medical technology to collect samples.
- Serious hemolysis, lipemia and turbid samples cannot be used for tests.
- Samples can be stable for 48h at 2°C ~ 8°C. If a test is not finished within 48h, freeze samples at -20°C or lower temperature.
- Samples can only be frozen once. Mix well after thawing.
- Before putting a sample in the system, ensure that the sample is without fibrous protein or other particles, and bubbles.

## 【Test Method】

- Reagent preparing

R1, R2 and R3 are all ready-to-use reagents, which can be used directly. Mix the reagents before loading them into the system. Visual inspection of the bottom of the kit ensures that all magnetic particles have been dispersed and re-suspended to avoid bubbles. Calibrator and control are ready-to-use and can be used directly. Before use, mix calibrator and control, balance them to room temperature and use them.

- Test procedure

Before loading reagents on the system, mix all reagents by hand. Visually inspect the reagent bottle bottom, to guarantee magnetic particles divided or resuspended. For detailed operation steps refer to the instrument user manual.

### 3. Calibration

When using new batches of reagents, the A-TG determination item needs to be re-calibrated and the calibration information registration card scanned (support manual input registration). By measuring low and high Calibrators, each calibration point on the pre-input main calibration curve is adjusted to a new calibration curve.

In the following cases, it should be calibrated again:

- ◆ Use the reagent kit of a new batch number.
- ◆ Replace trigger reagent with that of a new batch number.
- ◆ When the QC repeated results are not within the prescribed range.

### 4. QC

- 1) Two levels of Controls are determined on the day of testing a sample each time.
- 2) Controls must be tested when performing calibration. All calibrators and controls should be disposed of regarding as samples.
- 3) Take the following measures when control results are not within the lab regulated acceptable range.

- ◆ Ensure the reagent used has not expired.
- ◆ Ensure required maintenance is executed.
- ◆ Ensure test procedures are performed strictly following the instructions.

- ◆ Use new control to re-test.
- ◆ Use new calibrator to re-calibrate.
- ◆ Ask local technicians or distributor for help if necessary.

### 5. Calculation on test results

The instrument can calculate each sample concentration automatically, the unit is IU/mL.

### 【Reference Range】

Reference range is < 6.96IU/mL.

The lab should study the above reference range. The lab is suggested to set its own reference range due to geographical, patient dietary habit and environmental factors.

### 【Interpretation of Test Results】

1. Test results are not the only one as diagnosis index of clinical indications. Clinical significance is analyzed specifically combined with other test indices and clinical manifestation.
2. There is no direct comparability between sample A-TG concentration test by other ways and product test results.
3. Test results outside the reagent kit linear range should be tested again after diluting samples to the linear range.
4. When A-TG concentration is up to 100000IU/mL, there is no high dose hook effect.
5. A-TG measurement is different in test method, site identification, specificity and interfering factors, thus, A-TG test results are different for a specified sample; Inspectors should indicate the test method when supplying a laboratory test report to doctors. No direct comparability between test results obtained from different test methods. Direct cross use may lead to misinterpretation of its clinical significance; in the continuous monitoring of the efficacy of patients, before the method can be changed halfway, it is necessary to go through a full parallel experiment between the old and new methods and confirm its feasibility.

### 【Limitations of Test Method】

1. Patients of frequent exposure to animals and animal serum products and those who have used antibodies for in vivo diagnosis and

treatment may contain heterophilic antibodies, which may lead to false positive or false negative.

2. In the sample, when hemoglobin > 1000mg/dL, triglyceride > 1000mg/dL, bilirubin > 20mg/dL, total protein > 10g/dL test results may be affected.

3. RF in the sample may cause test result false positive or false negative.

### 【Product Performance Indices】

1. Precision: test national standard material, the relative deviation between concentration test value and nominal value should be within ± 10%.
2. Minimum detection limit: < 1.0IU/mL.
3. Linearity: linear range is 1.0IU/mL ~ 1000IU/mL, linear correlation coefficient  $r \geq 0.9900$ .
4. Repeatability:  $CV \leq 8.0\%$ .
5. Between-batch difference:  $CV \leq 15.0\%$ .














### 【Matters Needing Attention】

1. This product is only used for in vitro diagnosis.
2. Considering the possible evaporation effect, the samples, calibrators and controls on board should be analyzed/measured within 2 hours.
3. Please treat the samples as dangerous substances that may be infected with HIV, HBV, HCV, etc. To avoid or reduce the risk of infection, disposable gloves and eye/face protective items should be worn.
4. If the reagent enters the eye or mouth by mistake, or touches the skin, please rinse it with water quickly and receive medical treatment if necessary.
5. Samples and waste liquids are potentially biologically contagious. Operators should abide by laboratory safety regulations and treat waste liquids in accordance with local medical wastes, infectious wastes and industrial wastes.
6. Clinical samples should be treated as infectious samples, and operate according to the relevant laboratory specifications and requirements promulgated by the Health and Planning Commission, the Ministry of Science and Technology, and National Medical Products Administration and other relevant departments.
7. Avoid freezing the reagents.

### 【References】

1. National Committee for Clinical Laboratory Standards. How to define and determine reference intervals in the clinical laboratory. Approved Guideline. NCCLS Document C28-A2.Wayne(PA):NCCLS; 2000 June.
2. Nordyke RA, Gilbert FI Jr, Miyamoto LA, et al. The superiority of antimicrobial over antithyroglobulin antibodies for detecting Hashimoto's thyroiditis. Arch Intern Med 1993; 153:862-5.
3. Ruf J, Feldt-Rasmussen U, Hegedüs L, et al. Bispecific thyroglobulin and thyroperoxidase autoantibodies in patients with various thyroid and autoimmune diseases. J ClinEndocrinolMetab1994; 79(5): 1404-9.
4. Kricka, L. Interferences in Immunoassays still a threat. ClinChem2000; 46: 1037.
5. Clinical and Laboratory Standards Institute (formerly NCCLS). Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition, Wayne, PA: Clinical and Laboratory Standards Institute; 2005.CLSIEP7-A2.

【Instruction Approved & Modified Date】 03/2020

Symbol Explanation									
	Temperature limit		Batch code		Keep away from sunlight		Consult instructions for use		R22 S28 Danger level
	Authorised Representative		Catalogue Number		In vitro diagnostic medical device		Use-by date		This Way Up
	Biological Risk		Manufacturer		Comply with In Vitro Diagnostic Devices Directive (98/79/EC)				

**DIRUI®**

**DIRUI INDUSTRIAL CO., LTD.**  
95 Yunhe Street, New & High Tech.  
Development Zone Changchun, Jilin 130012 P.R.China

Tel: +86 431 85100409  
Fax: +86 431 85172581  
E-mail: dirui@dirui.com.cn  
http://www.dirui.com.cn

**CE**  
**IVD**

**EC REP**  
EMERGO EUROPE Prinsessegracht 20  
2514 AP The Hague The Netherlands