Thromboplastin LI

REF 5265H
REF 5265
REF 5267
REF 5269

REF 5265H
REF 5265
REF 5267
REF 5269

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0439P 2012/03 (13)

Thromboplastin LI

Instructions for use

INTENDED USE

The first standardised one-stage prothrombin time test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X)^{1,2}. It is also used for the induction and monitoring of oral anticoagulant therapy³⁻⁴ and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders. Thromboplastin LI is of rabbit brain origin but resembles human BCT in its low International Sensitivity Index (ISI). The ISI of Thromboplastin LI is approximately 1.1 (using manual calibration techniques) and is calibrated against the WHO international reference preparation⁵. Thromboplastin LI is particularly suited to the monitoring of oral anticoagulant therapy and, in conjunction with the appropriate factor deficient plasma, the measurement of factor activity in the extrinsic pathway. Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to normal citrated plasma, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

Composition	Content	Description	Preparation
Thromboplastin LI	4 x 2,5 mL (REF 5265H) 10 x 2,5 mL (REF 5265) 10 x 5 mL (REF 5267) 10 x 10 mL (REF 5269)	Each vial contains rabbit brain thromboplastin suspension. The reagent is a fine suspension of rabbit brain tissue.	The contents of the vial should be mixed well before and during use.
Calcium Chloride: 0.025M	4 x 2,5 mL (REF 5265H) 10 x 2,5 mL (REF 5265) 10 x 5 mL (REF 5267) 10 x 10 mL (REF 5269)	Each vial contains 0.025 M calcium chloride solution. The reagent is a clear, colourless solution.	The reagent is ready for use as packaged.

Each kit contains Instructions For Use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5519	ISI Calibrant Plasma Set
REF 5490	INR Reference Set
REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA

STORAGE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thromboplastin LI	Opened vials are stable for 1 month at *2 –8°C and 5 days at *15 –30°C. The working reagent (mixed with calcium chloride) is stable for 10 days at *2 –8°C or 5 days at *15 –30°C. DO NOT FREEZE. Large clumps of particles or changes in expected values may indicate product deterioration.
Calcium Chloride: 0.025M	Opened vials are stable for 1 month at *2 –8°C or 5 days at *15 –30°C. Particulate contamination, cloudiness or a change in expected performance may indicate product deterioration.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –8°C or *18 –24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes*.

PROCEDURE

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Helena Biosciences Europe ISI Calibrant Plasma Set (REF 5519) is recommended for this purpose⁵. This should be performed for each new reagent batch. The Helena Biosciences Europe INR Reference Set (REF 5490) should be used to check for shifts in the local system ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variances.

Manual Method

- Mix sufficient Thromboplastin LI and calcium chloride solution 1 + 1, to complete the anticipated testing for the day and incubate at *37°C for no more than 8 hours.
- Pre-warm 0.1 mL of the test plasma at *37°C for 2 minutes.
- Add 0.2 mL of freshly mixed Thromboplastin reagent to the plasma while simultaneously starting a stopwatch.
- Note the time for clot formation to the nearest 0.1 seconds.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other. %PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time), which should be a straight line when plotted on log-log graph paper.

INR values can be calculated using the following formula: INR = (PT Time Patient / Mean Normal PT Time)^{ISI}

For clear guidance on the indications for and management of patients on warfarin, please refer to The British Society for Haematology, for their most current edition of 'Guidelines on oral anticoagulation with warfarin'. At time of printing this is the 2011 fourth edition⁹.

LIMITATIONS

The use of serial dilutions of a reference plasma for the %PT curve is not recommended as this can lead to discrepancies caused by the low fibrinogen in the reference plasma dilutions which are not reflected in patient samples having predominantly normal fibrinogen levels.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA
REF 5490	INR Reference Set

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration. Using the manual (tilt-tube) technique, normal values ranging from 12.5-16.5 seconds are typical.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility

	Intra-assay precision				Inter-assay precision				Inter-batch precision			
Sample	n	Clot formation (seconds)	CV (%)		n	Clot formation (seconds)	CV (%)		n	Clot formation (seconds)	CV (%)	
Normal	10	17.0	0.92	100	13.7	2.09	3	13.1	2.0			
Abnormal	10	26.9	1.09	100	23.0	2.38	3	27.3	3.3			
Abnormal				100	25.6	2.15	3	44.0	2.2			

Helena Thromboplastin LI is insensitive to Heparin levels of up to 2 U/mL.

Method Comparison:

Comparison of INR values determined using Thromboplastin LI and two WHO Reference Thromboplastins yielded the following correlation equations:

Thromboplastin LI = 0.9936 (RTF95) + 0.0129	r ² = 0.9834	n = 79
Thromboplastin LI = 0.9589 (BCT441) + 0.1543	r ² = 0.9500	n = 79

REFERENCES

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci*, **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol*, **103**: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol*, **103**: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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REFERENZEN

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci*, **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*,

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RIFERIMENTI

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci.* **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol*, **103**: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol*, **103**: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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