CERTEST

SARS-CoV-2 + Flu A + Flu B

ONE STEP SARS-CoV-2, Influenza
type A and Influenza type B Ag
COMBO CARD TEST

CERTEST BIOTEC S.L.

ENGLISH

INTENDED USE

CerTest SARS-CoV-2 + Flu A + Flu B one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of nucleoprotein antigen of SARS-CoV-2, *Influenza* type A and *Influenza* type B from nasopharyngeal swab samples from patients suspected of COVID-19 infection and/or *Influenza A* and/or *Influenza B* infection.

CerTest SARS-CoV-2 + Flu A + Flu B combo card test offers a simple and highly sensitive screening assay as aid in the diagnosis of SARS-CoV-2 and/or *Influenza* type A and/or *Influenza* type B infection.

INTRODUCTION

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus.

The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Coronaviruses (CoVs) are RNA, enveloped viruses, with a large genome (29.9 kb), amongst which β -CoVs and α -CoVs can infect mammals. Viruses belonging to the Coronaviridae family have the most abundant membrane protein (M) among other proteins, i.e. spike glycoprotein (S), nucleocapsid protein (N) and an envelope protein (E). Spike glycoprotein (S) is one of the targets of T cell response in the immune system. The S protein also promotes the binding of the virus envelope to the ACE2 receptor and the entry of the virus into the cell target. ACE2 receptors are present in cells of the arteries, veins, smooth muscles, small intestine, alveoli of the lungs, hair follicles, cardiac myofibroblasts, skin, brain, and kidney, thus SARS-CoV-2 could potentially infect these tissues.

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of the viral infection include a sore throat, aches, pains, and shortness of breath. In some cases, runny nose, nausea, and diarrhea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks. Due to the high infectious rate of SARS-CoV-2, detection of asymptomatic positive patients is probably one of the key points to controlling the outbreak.

Influenza is caused by a virus that attacks mainly the upper respiratory tract - the nose, throat and bronchi and rarely also the lungs. The infection usually lasts for about a week. It is characterized by sudden onset of high fever, myalgia, headache and severe malaise, non-productive cough, sore throat, and rhinitis. Most people recover within one to two weeks without requiring any medical treatment. In the very young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, *Influenza* poses a serious risk. In these people, the infection may lead to severe complications of underlying diseases, pneumonia and death.

The currently circulating influenza viruses that cause human disease are divided into two groups: A and B. *Influenza A* has 3 subtypes which are important for humans: A (H3N2), A (H1N1) and A (H5N1), of which the former is currently associated with most deaths. *Influenza* viruses are defined by 2 different protein components, known as antigens, on the surface of the virus. They are spike-like features called haemagglutinin (H) and neuraminidase (N) components.

TEST PRINCIPLE

CerTest SARS-CoV-2 + Flu A + Flu B one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of nucleoprotein antigen of SARS-CoV-2, *Influenza* type A and *Influenza* type B from nasopharyngeal swab samples from patients suspected of COVID-19 infection and/or *Influenza A* and/or *Influenza B* infection.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza* type A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Influenza* type A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip C consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza* type B and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Influenza* type B) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

CERTEST SARS-CoV-2 + Flu A + Flu B

*One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette*

If the sample is SARS-CoV-2 positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres) in the strip A and if the sample is *Influenza A* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*Influenza* type A monoclonal antibodies-red polystyrene microspheres) in the strip B, which were previously pre-dried on the absorbent pad and if the sample is *Influenza* type B positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*Influenza* type B monoclonal antibodies-red polystyrene microspheres) in the strip C, which were previously pre-dried on the absorbent pad. The mixtures then move upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane of strip A (test line), the anti-*Influenza* type B antibodies present on the membrane of strip B (test line) and the anti-*Influenza* type B antibodies present on the membrane of strip C (test line) capture the coloured conjugate and a red line will be visible in the strips. These bands are used to interpret the result.

If the sample is negative, there is no SARS-CoV-2, *Influenza* type A, *Influenza* type B antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-SARS-CoV-2, the anti-*Influenza* type A and the anti-*Influenza* type B antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the three strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on the three membranes will capture control green-conjugate complex and the three control lines will always appear. The presence of these green lines serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors. Please remember that COVID-19 infection is very dangerous.
- The tests should be discarded in a proper biohazard container after testing.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest SARS-CoV-2 + Flu A + Flu B test. Do not use any other commercial kit component.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All positive results should be processed following local laws and regulations.

SPECIMEN COLLECTION AND PREPARATION


Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the wall.
(to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.

Follow the test procedure.

Specimen collection

Samples previously extracted by transport media:

Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTM, UTM or Saline Buffer.

1. Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0 mL. Using higher volumes will affect sensitivity of the system.
2. Dilute the extracted sample 1:1 in Reagent SARS-CoV-2 + Flu A + Flu B (sample diluent) provided.

Follow immediately with the point 3 of test procedure.

MATERIALS

MATERIALS PROVIDED

- CerTest SARS-CoV-2 + Flu A + Flu B combo card tests
- Reagent SARS-CoV-2 + Flu A + Flu B (sample diluent)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- Laboratory rack
- CerTest SARS-CoV-2 Positive Control swab + Instructions for use
- CerTest Influenza A Positive Control swab + Instructions for use
- CerTest Influenza B Positive Control swab + Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED



- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Add 15 drops (approx. 500µL) (1) Reagent SARS-CoV-2 + Flu A + Flu B and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest SARS-CoV-2 + Flu A + Flu B combo card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense 3 drops (approx. 120µL) from the testing tube, into the circular window marked with the letter A (3) and 3 drops (approx. 120µL), using the same tube, into the circular window marked with the letter B (4), and 3 drops (approx. 120µL) using same tube, into the circular window marked with the letter C (5).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent SARS-CoV-2 + Flu A + Flu B until seeing the liquid running through the reaction zone.


CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

SARS-CoV-2 Strip A-procedure


Add 3 drops in A-circular window

(3)

Influenza type A Strip B-procedure


Add 3 drops in B-circular window

(4)

Influenza type B Strip C-procedure

Add 3 drops in C-circular window

(5)

INTERPRETATION OF THE RESULTS (please refer to the illustration below)


- A: Green→Negative SARS-CoV-2
B: Green→Negative Influenza A
C: Green→Negative Influenza B
- A: Green/Red→Positive SARS-CoV-2
B: Green/Red→Positive Influenza A
C: Green/Red→Positive Influenza B
- A: Green/Red→Positive SARS-CoV-2
B: Green/Red→Positive Influenza A
C: Green→Negative Influenza B
- A: Green/Red→Positive SARS-CoV-2
B: Green→Negative Influenza A
C: Green/Red→Positive Influenza B
- A: Green/Red→Positive SARS-CoV-2
B: Green→Negative Influenza A
C: Green→Negative Influenza B
- A: Green→Negative SARS-CoV-2
B: Green/Red→Positive Influenza A
C: Green/Red→Positive Influenza B
- A: Green→Negative SARS-CoV-2
B: Green/Red→Positive Influenza A
C: Green→Negative Influenza B
- A: Green→Negative SARS-CoV-2
B: Green→Negative Influenza A
C: Green/Red→Positive Influenza B

INVALID RESULTS

-

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

	A (SARS-CoV-2)	B (Influenza A)	C (Influenza B)	Interpretation of the results
1.	-	-	-	There is no SARS-CoV-2, <i>Influenza</i> type A, <i>Influenza</i> type B antigens presence in sample.
	GREEN	GREEN	GREEN	
2.	+	+	+	There is SARS-CoV-2, <i>Influenza</i> type A, <i>Influenza</i> type B antigens presence in sample.
	GREEN-RED	GREEN-RED	GREEN-RED	
3.	+	+	-	There is SARS-CoV-2 and <i>Influenza</i> type A antigens presence in the sample. There is no <i>Influenza</i> type B antigens presence in the sample.
	GREEN-RED	GREEN-RED	GREEN-RED	
4.	+	-	+	There is SARS-CoV-2 and <i>Influenza</i> type B antigens presence in the sample. There is no <i>Influenza</i> type A antigens presence in the sample.
	GREEN-RED	GREEN-RED	GREEN	
5.	+	-	-	There is SARS-CoV-2 antigens presence in the sample. There is no <i>Influenza</i> type A and type B antigens presence in the sample.
	GREEN-RED	GREEN	GREEN-RED	
6.	-	+	+	There is <i>Influenza</i> type A and type B antigens presence in the sample. There is no SARS-CoV-2 antigens presence in the sample.
	GREEN	GREEN-RED	GREEN-RED	
7.	-	+	-	There is <i>Influenza</i> type A antigens presence in the sample. There is no SARS-CoV-2 and <i>Influenza</i> type B antigens presence in the sample.
	GREEN-RED	GREEN-RED	GREEN	
8.	-	-	+	There is <i>Influenza</i> type B antigens presence in the sample. There is no SARS-CoV-2 and <i>Influenza</i> type A antigens presence in the sample.
	GREEN	GREEN	GREEN-RED	
9.	Invalid result: either A, B or C, we recommend repeating the assay using the same sample with another test.			

Whatever result detailed in the above-mentioned table (2 to 8) should be followed up with additional confirmatory diagnostic procedures.

Whatever results detailed between 2 and 6 are not very common. Dual simultaneous respiratory infections are more frequent than triple respiratory virus infection.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

RECOMMENDATIONS

Recommendations World Health Organization: *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 September 2020.*

1. To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.
2. Where possible, all positive samples giving positive results should be transported to laboratories with NAAT (Nucleic acid amplification tests) capability for confirmatory testing.
3. This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT- confirmed COVID-19 outbreaks.
4. To monitor trends in disease incidence in communities.
5. For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
6. A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
7. Even CerTest SARS-CoV-2 + Flu A + Flu B test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so CerTest SARS-CoV-2 + Flu A + Flu B could detect as positive.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest SARS-CoV-2 + Flu A + Flu B should be used only with nasopharyngeal swab samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of SARS-CoV-2 and/or Influenza type A and/or Influenza type B antigens. A diagnosis of SARS-CoV-2, Influenza type A and Influenza type B infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. Positive results do not rule out co-infections with other pathogens.
6. Only VTM, UTM and Saline buffer transport media has been validated with CerTest SARS-CoV-2 + Flu A + Flu B test. When using transport media the test sensitivity can be reduced due to excessive dilution of sample.
7. A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique.

EXPECTED VALUES

On January 30, 2020, the World Health Organization declared COVID-19 a public health emergency of international concern. Since the onset of the pandemic, SARS-CoV-2 has affected millions of people, causing thousands of victims worldwide.

Recent studies have shown that the primary viral infection found in symptomatic patients was identical to that seen in asymptomatic patients, suggesting that COVID-19 patients with asymptomatic infections are capable of transmitting SARS-CoV-2. During the incubation period, asymptomatic patients are also likely to transmit the virus to healthy people.

Although SARS-CoV-2 infection can affect people of any age, children and the elderly with asthma are in a higher risk category for the disease. Those with underlying health problems such as cancer, chronic respiratory diseases, cardiovascular diseases, and diabetes are also reported to have a higher chance of developing serious illnesses.

Yearly Influenza epidemics can seriously affect all age groups, but the highest risk of complications occur among children younger than age two, adults age 65 or older and people of any age with certain medical conditions, such as chronic heart, lung, kidney, liver, blood or metabolic diseases (such as diabetes) or weakened immune systems.

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

Influenza affects between 20% and 30% of the United States population annually. Three fourths of the infected individuals develop an acute respiratory illness and one third of these seek medical attention. In a typical year, more than 20,000 to 40,000 Americans die and more than 100,000 are hospitalized because of complications related to influenza.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value (typical) of CerTest SARS-CoV-2 + Flu A + Flu B (strip A: SARS-CoV-2) is 1.00 ng/mL of SARS-CoV-2 recombinant protein or 1×10^3 TCID₅₀/mL of 2019nCoV/USA-WA1/2020.

Detection limit value (typical) of CerTest SARS-CoV-2 + Flu A + Flu B (strip B: Influenza A) is 6.25 ng/mL of Influenza A recombinant nucleoprotein.

Detection limit value (typical) of CerTest SARS-CoV-2 + Flu A + Flu B (strip C: Influenza B) is 25.00 ng/mL of Influenza B recombinant nucleoprotein.

Clinical sensitivity and specificity

Evaluations for SARS-CoV-2 + Flu A + Flu B (strip A: SARS-CoV-2)

An evaluation, with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2 + Flu A + Flu B, CerTest (SARS-CoV-2 strip) vs qPCR technique (VIASURE). The results were as follows:

IC test: CerTest SARS-CoV-2 + Flu A + Flu B (SARS-CoV-2)		qPCR technique		
		+	-	Total
+		26	1	27
-		2	233	235
Total		28	234	262

CerTest SARS-CoV-2 + Flu A + Flu B (SARS-CoV-2) vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	92.9%	76.5–99.1%
Specificity	99.6%	97.6–100.0%
PPV	96.3%	81.0–99.9%
NPV	99.1%	97.0–99.9%

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays* (11 September 2020) from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

Evaluations for SARS-CoV-2 + Flu A + Flu B (strip B: Influenza A)

Results showed in tables below have been obtained from several evaluations developed using an immunochromatographic (CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A), CerTest) and others commercial test kit.

1574 respiratory samples have been used for Influenza A evaluation (130 positive samples, 11 false positives, 1402 negatives samples and 31 false negatives). The results were as follows:

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

	Evaluation number	True Positives	False Negatives	True Negatives	False Positives
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs XPECT FLU A&B (Remel)	I	7	0	85	1
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs Instand EQAS program (March 2016)	II	101	26	955	3
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs Binax NOW® Influenza A&B Card (Alere)	III	11	2	87	0
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs SD Bioline, Influenza Ag A/B/A (H1N1) (SD Bioline)	IV	1	0	97	0
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs MariPOC® respi (ArcDia International)	V	10	3	178	7
	Total	130	31	1402	11

CerTest SARS-CoV-2 + Flu A + Flu B (Influenza A) vs other test used in evaluation (I, II, III, IV, V) for Influenza A		
	Mean Value	95% confidence interval
Sensitivity	80.7%	73.8 – 86.5%
Specificity	99.2%	98.6 – 99.6%
PPV	92.2%	86.5 – 96.0%
NPV	97.8%	96.9 – 98.5%

Evaluations for SARS-CoV-2 + Flu A + Flu B (strip C: Influenza B)

Results showed in tables below have been obtained from several evaluations developed using an immunochromatographic (CerTest SARS-CoV-2 + Flu A + Flu B (Influenza B), CerTest) and others commercial test kit.

1276 respiratory samples have been used for Influenza B evaluation (75 positive samples, 6 false positives, 1181 negatives samples and 14 false negatives). The results were as follows:

	Evaluation number	True Positives	False Negatives	True Negatives	False Positives
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza B) vs XPECT FLU A&B (Remel)	I	7	2	83	1
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza B) vs Instand EQAS program (March 2016)	II	63	12	1005	5
CerTest SARS-CoV-2 + Flu A + Flu B (Influenza B) vs SD Bioline, Influenza Ag A/B/A (H1N1) (SD Bioline)	III	5	0	93	0
	Total	75	14	1181	6

CerTest SARS-CoV-2 + Flu A + Flu B (Influenza B) vs other test used in evaluation (I, II, III) for Influenza B		
	Mean Value	95% confidence interval
Sensitivity	84.3%	75.0 – 91.1%
Specificity	99.5%	98.9 – 99.8%
PPV	92.6%	84.6 – 97.2%
NPV	98.8%	98.0 – 99.4%

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

The results showed a high sensitivity and specificity to detect SARS-CoV-2, Influenza type A and type B using CerTest SARS-CoV-2 + Flu A + Flu B.

Hook Effect

CerTest SARS-CoV-2 + Flu A + Flu B test does not show inhibitory hook effect at concentration higher than 100,000xLoD for Strip A: SARS-CoV-2, at concentrations higher than 30,000xLoD for Strip B: Influenza A and at concentrations higher than 8,000xLoD for Strip C: Influenza B.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest SARS-CoV-2+Flu A + Flu B; no cross reactivity against organisms, pathogens, substances that cause infections:

For Strip A. SARS-CoV-2

Adenovirus	Entamoeba histolytica	Legionella pneumophila	Salmonella typhimurium
Astrovirus	Escherichia coli O157	Listeria monocytogenes	Shigella boydii
Bovine Haemoglobin	Giardia	Norovirus GI	Shigella dysenteriae
Campylobacter jejuni	Helicobacter pylori	Norovirus GII	Shigella flexneri
Clostridium difficile antigen GDH	Human Calprotectin	Pig Haemoglobin	Shigella sonnei
Clostridium difficile Toxin A	Human Haemoglobin	Rotavirus	Streptococcus pneumococcal
Clostridium difficile Toxin B	Human Lactoferrin	Respiratory Syncytial Virus	Streptococcus pyogenes
Coronavirus strain 229E	Human Transferrin	Salmonella enteritidis	Yersinia Enterocolitica O:3
Coronavirus strain NL63	Influenza type A	Salmonella paratyphi A	Yersinia Enterocolitica O:9
Coronavirus strain OC43	Influenza type B	Salmonella typhi	

CerTest SARS-CoV-2+Flu A + Flu B (SARS-CoV-2, strip A) showed some cross reaction with SARS and very low with MERS.

For Strip B. Influenza A

Adenovirus	Coronavirus strain OC43	Legionella pneumophila	Salmonella typhi
Astrovirus	Entamoeba histolytica	Listeria monocytogenes	Salmonella typhimurium
Bovine Haemoglobin	Escherichia coli O157	MERS	Shigella boydii
Campylobacter jejuni	Giardia	Norovirus GI	Shigella dysenteriae
Clostridium difficile antigen GDH	Helicobacter pylori	Norovirus GII	Shigella flexneri
Clostridium difficile Toxin A	Human Calprotectin	Pig Haemoglobin	Shigella sonnei
Clostridium difficile Toxin B	Human Haemoglobin	Rotavirus	Streptococcus pneumococcal
Coronavirus SARS-1 (SARS)	Human Lactoferrin	Respiratory Syncytial Virus	Streptococcus pyogenes
Coronavirus SARS-2 (SARS-CoV-2)	Human Transferrin	Salmonella enteritidis	Yersinia Enterocolitica O:3
Coronavirus strain 229E	Influenza type B	Salmonella paratyphi A	Yersinia Enterocolitica O:9
Coronavirus strain NL63			

For Strip C. Influenza B

Adenovirus	Coronavirus strain OC43	Legionella pneumophila	Salmonella typhi
Astrovirus	Entamoeba histolytica	Listeria monocytogenes	Salmonella typhimurium
Bovine Haemoglobin	Escherichia coli O157	MERS	Shigella boydii
Campylobacter jejuni	Giardia	Norovirus GI	Shigella dysenteriae
Clostridium difficile antigen GDH	Helicobacter pylori	Norovirus GII	Shigella flexneri
Clostridium difficile Toxin A	Human Calprotectin	Pig Haemoglobin	Shigella sonnei
Clostridium difficile Toxin B	Human Haemoglobin	Rotavirus	Streptococcus pneumococcal
Coronavirus SARS-1 (SARS)	Human Lactoferrin	Respiratory Syncytial Virus	Streptococcus pyogenes
Coronavirus SARS-2 (SARS-CoV-2)	Human Transferrin	Salmonella enteritidis	Yersinia Enterocolitica O:3
Coronavirus strain 229E	Influenza type A	Salmonella paratyphi A	Yersinia Enterocolitica O:9
Coronavirus strain NL63			

CERTEST SARS-CoV-2 + Flu A + Flu B

One Step test to detect SARS-CoV-2, Influenza A and Influenza B antigens in combo card format/
Prueba combo de un solo paso para detección de SARS-CoV-2, Influenza A e Influenza B antígenos en formato cassette

Interferences

An evaluation was performed to determine the possible interferences of CerTest SARS-CoV-2 + Flu A + Flu B; no interferences against the substances tested were detected:

Potencial interfering substance	Concentration	Result	Potencial interfering substance	Concentration	Result
Metronidazole	3.0 mg/mL	Negative	Loperamide hydrochloride (<i>Fortasec</i>)	$1.5 \cdot 10^{-2}$ mg/mL	Negative
Ampicillin	3.0 mg/mL	Negative	Heparin (<i>Hibor</i>)	350.0 IU/mL	Negative
Oseltamivir	$3.0 \cdot 10^{-2}$ mg/mL	Negative	Almagato (<i>Almax</i>)	3.0 mg/mL	Negative
Amantadine	0.3 mg/mL	Negative	Fosfamycin (<i>Monurol</i>)	$3.0 \cdot 10^{-3}$ mg/mL	Negative
Ribavirin	3.0 mg/mL	Negative	Acetylcysteine (<i>Fluimucil</i>)	$3.0 \cdot 10^{-2}$ mg/mL	Negative
Codeine (<i>Toseina</i>)	0.2 mg/mL	Negative	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL	Negative
Benzocaine (<i>Angileptol</i>)	$3.0 \cdot 10^{-2}$ mg/mL	Negative	Levofloxacin	$3.0 \cdot 10^{-2}$ mg/mL	Negative
Cloperastine (<i>Flutox</i>)	0.3 mg/mL	Negative	Ciprofloxacin	0.3 mg/mL	Negative
Carbocysteine (<i>Iniston mucolítico</i>)	$3.0 \cdot 10^{-2}$ mg/mL	Negative	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL	Negative
Loratadine	0.3 mg/mL	Negative	Phenoxymethylpenicillin potassium	3.0 mg/mL	Negative
Dexchloropheniramine (<i>Polaramine</i>)	0.3 mg/mL	Negative	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL	Negative
Ebastine (<i>Ebastel</i>)	3.0 mg/mL	Negative	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL	Negative
Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	Negative	Lysine Carbocysteinate (<i>Pectox</i>)	$3.0 \cdot 10^{-2}$ mg/mL	Negative
Ibuprofen (<i>Espidifen</i>)	0.3 mg/mL	Negative	Hydroxyzine dihydrochloride	0.3 mg/mL	Negative
Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	Negative	Lorazepam	$3.0 \cdot 10^{-3}$ mg/mL	Negative
Metamizole (<i>Nolotil</i>)	5.0 mg/mL	Negative	Amoxicillin	3.0 mg/mL	Negative
Prednisone	0.3 mg/mL	Negative	Mercaptopurine	0.3 mg/mL	Negative
Omeprazole	$2.0 \cdot 10^{-3}$ mg/mL	Negative			

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1. Palestino, G.; García-Silva, I.; González-Ortega, O. and Rosales-Mendoza, S. **Can nanotechnology help in the fight against COVID-19?** *Expert Review of Anti-infective Therapy*. <https://doi.org/10.1080/14787210.2020.1776115>. Jun. 2020.
2. Moeis MR.; Rahayu AP.; Ihsani N. and Pertiwi W. **Coronavirus-SARS-CoV-2: Biology and Problems in rRT-PCR Detection**. *Borneo Journal of Pharmacy*. Vol. 3; special Issue 1:136-145. Jun. 2020. <https://doi.org/10.33084/bjop.v3iSpecial-1.1429>
3. Arnedo-Pena, A.; Sabater-Vidal, S.; Meseguer-Ferrer, N.; Pac-Sa, M. R.; Mañes-Flor, P.; Gascó-Laborda, J. C.; Larrea, R. M.; Tirado-Balaguer, M. D.; Romeu-García, M. A.; Gil-Fortuño, M.; Safont-Adsua, L.; Blasco, A.; Gomila-Sar, B.; Moreno-Muñoz, M. R. and Bellido-Blasco, J. **COVID-19 secondary attack rate and risk factors in household contacts in Castellon (Spain): Preliminary report**. *Revista Enfermedades Emergentes*. 2020;19(2):64-70.
4. Younes, N.; Al-Sadeq, D. W.; AL-Jighefee, H.; Younes, S.; Al-Jamal, O.; Daas, H. I.; Yassine, H. M. and Nasrallah, G. K. **Challenges in Laboratory Diagnosis of the Novel Coronavirus SARS-CoV-2**. *Viruses (Multidisciplinary Digital Publishing Institute)*. 12 (6): 582. <https://doi.org/10.3390/v12060582>. May. 2020.
5. Rincón, A.; Moreso, F.; López-Herradón, A.; Fernández-Robres, M.A.; Cidraque, I.; Nin, J.; Méndez, O.; López, M.; Pájaro, C.; Satorra, A.; Stuard, S. and Ramos, R. **The keys to control a coronavirus disease 2019 outbreak in a haemodialysis unit**. *Clinical Kidney Journal*. 2020: 1–8.
6. Fragkou, P. C.; Papaevangelou, V.; Antoniadou, A.; Kavvatha, D.; Ploussi, A.; Pantazis N.; Sirmipilantze, T.; Psarrakis, C.; Pournaras, S. A.; Tsiodras, S. and Kelekis, A. **Preliminary Data of a Quantitative Point of Care Test for SARS-CoV-2 Antibodies From Greece**. *In vivo*. 34: xxx-xxx (2020).
7. Ezhilan, M.; Suresh, I. and Nesakumar, N. **SARS-CoV, MERS-CoV and SARS-CoV-2: A Diagnostic Challenge**. *Measurement*. 168 (2021): 108335. <https://doi.org/10.1016/j.measurement.2020.108335>. Aug. 2020.
8. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. *Journal of Clinical Microbiology*, Mar. 2000, p. 1161–1165 Vol. 38, No. 3
9. Adriana Weinberg and Miranda L. Walker. Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza Virus A and B. *Clinical And Diagnostic Laboratory Immunology*, Mar. 2005, p. 367–370 Vol. 12, No. 3
10. Andreea C. Cazacu, Gail J. Demmler et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. *Journal of Clinical Microbiology*, Aug. 2004, p. 3661–3664 Vol. 42, No. 8.
11. Kenneth E. Irmen and James J. Kelleher. Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital. *Clinical and Diagnostic Laboratory Immunology*, May 2000, p. 396–403 Vol. 7, No. 3.
12. Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays. Interim guidance. 11 September 2020. World Health Organization.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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F-368 rev00

CERTEST SARS-CoV-2 POSITIVE CONTROL
CERTEST Influenza A POSITIVE CONTROL
CERTEST Influenza B POSITIVE CONTROL

INTENDED USE

The CerTest SARS-CoV-2 Positive Control is an external quality control for CerTest SARS-CoV-2 + Flu A + Flu B combo card test. **Only for the strip A (SARS-CoV-2).**

The CerTest Influenza A Positive Control is an external quality control for CerTest SARS-CoV-2 + Flu A + Flu B combo card test. **Only for the strip B (FLU A).**

The CerTest Influenza B Positive Control is an external quality control for CerTest SARS-CoV-2 + Flu A + Flu B combo card test. **Only for the strip C (FLU B).**

Use those positive control swabs to check that the extraction reagents and the test are working properly.

COMPOSITION

CerTest SARS-CoV-2 Positive Control: inactivated SARS-CoV-2 antigens dried in swab containing a preservative.

CerTest Influenza A Positive Control: inactivated type A Influenza antigens dried in swab containing a preservative.

CerTest Influenza B Positive Control: inactivated type B Influenza antigens dried in swab containing a preservative.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The reagent is stable through the expiration date printed on the sealed pouch. The reagent must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- The swab must be carried out within 2 hours of opening the sealed bag.
- Do not use the swab if the aluminium pouch has been damaged during storage.
- The swab should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.

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
- Disposable gloves
- Timer

QUALITY CONTROL TESTING PROCEDURE


Allow tests and positive control swabs to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Remove the CerTest SARS-CoV-2 + Flu A + Flu B combo card test from its sealed bag just before using (1).
2. Open the cap (2). This reagent is the same for both positive control swabs and for the samples. Use a separate pipette and testing tube for each control.


CerTest SARS-CoV-2 Positive Control swab procedure

3. Add 15 drops (approx. 500µL) (3) of reagent and immediately introduce the **SARS-CoV-2 Positive Control swab** into the testing tube.
4. Mix de solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (4). Discard the swab.
5. Dispense 3 drops (approx. 120µL) from the testing tube in the **circular window marked with the letter A (5).**
6. **Read the result at 10 minutes.** Do not read the test result after more than 10 minutes.

Wait 10 minutes


CerTest Influenza A Positive Control swab procedure


7. Take the reagent (the same used in the CerTest SARS-CoV-2 Positive Control procedure and samples procedure), add 15 drops (approx. 500µL) (6) of diluent into a different testing tube and immediately introduce the **Influenza A Positive Control swab** into the testing tube
8. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the tube as the swab is withdrawn (7). Discard the swab.
9. Dispense 3 drops (approx. 120µL) from the testing tube in the **circular window marked with the letter B (8).**
10. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

(6) Add 15 drops of the Reagent SARS-CoV-2 + Flu A + Flu B

(7) Introduce the **CerTest Influenza A Positive Control swab** into the testing tube rotating 1 minute and extract the liquid

Wait 10 minutes


(8) Add 3 drops in the circular window **Influenza A (strip B)**

CerTest Influenza B Positive Control swab procedure


11. Take the reagent (the same used in the CerTest SARS-CoV-2 and Influenza A Positive Control procedure and samples procedure), add 15 drops (approx. 500µL) (9) of diluent into a different testing tube and immediately introduce the **Influenza B Positive Control swab** into the testing tube
12. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the tube as the swab is withdrawn (10). Discard the swab.
13. Dispense 3 drops (approx. 120µL) from the testing tube in the **circular window marked with the letter C (11).**
14. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

(9) Add 15 drops of the Reagent SARS-CoV-2 + Flu A + Flu B

(10) Introduce the **CerTest Influenza B Positive Control swab** into the testing tube rotating 1 minute and extract the liquid

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VALID RESULT

Strip A: SARS-CoV-2 POSITIVE. A GREEN line (control line (C)) and a RED line (SARS-CoV-2 test line (T)) in the result window result appear at the same time. That means the test is working properly.

Strip B: Influenza A POSITIVE. A GREEN line (control line (C)) and a RED line (Influenza A test line (T)) in the result window appear at the same time. That means the test is working properly.

Strip C: Influenza B POSITIVE. A GREEN line (control line (C)) and a RED line (Influenza B test line (T)) in the result window appear at the same time. That means the test is working properly.

CerTest SARS-CoV-2 Positive Control is a control for strip A in CerTest SARS-CoV-2+Flu A + Flu B combo card test.


CerTest Influenza A Positive Control is a control for strip B in CerTest SARS-CoV-2 + Flu A + Flu B combo card test.

CerTest Influenza B Positive Control is a control for strip C in CerTest SARS-CoV-2 + Flu A + Flu B combo card test.


INVALID RESULTS

Any other result is considered to be invalid. Total absence of the control colored lines (GREEN) regardless the appearance or not of the test lines (RED), or only the control (C) coloured lines (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.


1. POSITIVE


2. INVALID

3. INVALID


4. INVALID


A: Green-Red→Positive SARS-CoV-2
B: Green-Red→Positive FLU A
C: Green-Red→Positive FLU B


5. INVALID


6. INVALID

7. INVALID

8. INVALID

	A SARS-CoV-2	B FLU A	C FLU B	Interpretation of the results
1	+	+	+	SARS-CoV-2 (A) POSITIVE , Influenza A (B) POSITIVE and Influenza B (C) POSITIVE : A GREEN line (control line (C)) and a RED line (test line (T)) in the three strips, in the results window, appear at the same time due to presence of antigens in the SARS-CoV-2 Positive Control swab in strip A, presence of antigens in Influenza A Positive Control swab in strip B and presence of antigens in Influenza B Positive Control swab in strip C.
2,3,4,5,6,7	Any other result			Invalid result, we recommend repeating the assay with another test and positive control swabs.

SYMBOLS FOR IVD PRODUCTS AND REAGENTS

IVD

In vitro diagnostic device

LOT

Batch code

Consult instructions for use

REF

Catalogue number

Keep dry

Σ

Contains sufficient for <n> tests

Temperature limitation

Manufacturer

Use by

DIL

Diluent

Control

Positive Control

Do not re-use

CE

CE marking

CerTest Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840, San Mateo de Gállego, Zaragoza
BIOTEC S.L. (SPAIN) www.cerTEST.es

CERTEST SARS-CoV-2 POSITIVE CONTROL
CERTEST Influenza A POSITIVE CONTROL
CERTEST Influenza B POSITIVE CONTROL

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