VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

Cryptosporidium, Giardia & E. histolytica

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ENGLISH

1. Intended use

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit is designed for specific detection and differentiation of *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in human stool samples from patients with signs and symptoms of gastrointestinal infection. This test is intended for use as an aid in the diagnosis of the *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in combination with clinical and epidemiological risk factors. DNA is extracted from stool specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica*.

2. Summary and Explanation

Enteric protozoa continue to be the most commonly encountered parasitic diseases affecting millions of people each year and causing significant morbidity and mortality worldwide. Among them, *Cryptosporidium* spp., *Entamoeba histolytica*, and *Giardia lamblia* are the major etiological agents.

Cryptosporidium can be acquired through person-to-person, contact with companion and farm animals, ingestion of contaminated food, drinking water and recreational water. *Cryptosporidium* can be present without symptoms or cause diarrhea and abdominal discomfort with weight loss and malabsorption. *Giardia lamblia* infection occurs by the faecal-oral route transmission, either through direct contact or by ingestion of contaminated food or water. Diarrhoea is the major symptom of the disease but the clinical manifestations vary; from non-symptomatic infections to severe diarrhea, abdominal cramps, bloating and flatulence often accompanied by nausea and weight loss. *Entamoeba histolytica* disease is spread primarily by food or water contaminated with cysts but may also be transmitted from person to person by faecal-oral route. Clinical features of amoebiasis range from asymptomatic colonization (90% of the cases) to amoebic dysentery (stomach pain, nausea or vomiting and severe diarrhea with bloody and slimy stools, fever) and invasive extraintestinal disease, most commonly in the form of liver abscesses.

Microscopy has been considered to be the gold standard method for diagnosis of these parasites, however is time-consuming, not very sensitive and unable to distinguish the invasive pathogenic *E. histolytica* from the commensal parasite *E. dispar*. Real-time PCR assay is less labor-intensive and has higher sensitivity and specificity, making it an attractive alternative.

3. Principle of the procedure

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit is designed for the diagnosis of the *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in human stool samples. After DNA isolation, the identification of *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* is performed by the amplification of a conserved region of the 18S rRNA gene for *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* respectively, using specific primers and a fluorescent-labelled probe.

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. *Cryptosporidium* DNA targets are amplified and detected in Cy5 channel, *Giardia lamblia* DNA targets are amplified and detected in FAM channel, *Entamoeba histolytica* DNA targets are amplified and detected in ROX channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1 and Table 2:

Reference	Reagent/Material	Description	Colour	Amount
VS-KGE1SL/ VS-KGE1SH	<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-KGE1C	<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit with Ref. VS-KGE106L, VS-KGE106H, VS-KGE112L and VS-KGE112H.

Reference	Reagent/Material	Description	Color	Amount
VS-KGE1PL/ VS-KGE1PH	<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-KGE1C	<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit with Ref VS-KGE113L and VS-KGE113H.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler) (to check compatibility see Annex I).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes.
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- For professional *in vitro* diagnostic use.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.

- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different kits and/or lots.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

8. Test procedure

8.1. Sample preparation

Stool samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

8.1.1. DNA extraction

For DNA extraction from human stool samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- QIAamp DNA Mini kit (QIAGEN).
- QIAamp DNA Stool Mini Kit (QIAGEN).
- Invisorb® Spin Universal Kit (Strattec).
- RIDA® Xtract (r-Biopharm).
- Maxwell®RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

8.2. Lyophilized positive control

Cryptosporidium, *Giardia* & *E. histolytica* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Cryptosporidium*, *Giardia* & *E. histolytica* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Cryptosporidium*, *Giardia* & *E. histolytica* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. Centrifuge briefly.

Load the plate or the strips in the thermocycler.

- 3) Set up your thermocycler.

Program your thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the Cy5 (*Cryptosporidium*), FAM (*Giardia lamblia*), ROX (*Entamoeba histolytica*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>Cryptosporidium</i> (Cy5)	<i>Giardia lamblia</i> (FAM)	<i>E. histolytica</i> (ROX)	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	<i>Cryptosporidium</i> , <i>Giardia lamblia</i> and <i>E. histolytica</i> Positive
-	-	-	+	-	+	<i>Cryptosporidium</i> , <i>Giardia lamblia</i> and <i>E. histolytica</i> Negative
+	-	-	+/-	-	+	<i>Cryptosporidium</i> Positive, <i>Giardia lamblia</i> and <i>E. histolytica</i> Negative
+	+	-	+/-	-	+	<i>Cryptosporidium</i> and <i>Giardia lamblia</i> Positive, and <i>E. histolytica</i> Negative
+	-	+	+/-	-	+	<i>Cryptosporidium</i> and <i>E. histolytica</i> Positive, and <i>Giardia lamblia</i> Negative
-	+	-	+/-	-	+	<i>Giardia lamblia</i> Positive, <i>Cryptosporidium</i> and <i>E. histolytica</i> Negative
-	+	+	+/-	-	+	<i>Giardia lamblia</i> and <i>E. histolytica</i> Positive, <i>Cryptosporidium</i> Negative
-	-	+	+/-	-	+	<i>E. histolytica</i> Positive, <i>Cryptosporidium</i> and <i>Giardia lamblia</i> Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

Table 4. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with human faecal samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit was testing using 172 faecal samples from symptomatic patients. These results were compared with those obtained by two commercial Real Time PCR Kit (RIDA®GENE Parasitic Stool Panel (R-biopharm) or FTD Stool parasites (fast-track DIAGNOSTICS)).

The results were as follows:

VIASURE <i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Real Time PCR Detection Kit	RIDA®GENE Parasitic Stool Panel II (r-Biopharm) + FTD Stool parasites (fast-track DIAGNOSTICS)			
		+	-	Total
	+	38	3*#	41
	-	1#	130	131
	Total	39	133	172

Table 5. Comparative results of *Cryptosporidium*

These 4 samples have been evaluated by an additional in-house Real Time PCR (Hadfield et al., 2011), confirming our results.

* The low amount of template DNA detected in these samples is below the detection limit of the method used and so it may produce random positive results.

VIASURE <i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Real Time PCR Detection Kit	RIDA®GENE Parasitic Stool Panel II (r-Biopharm) + FTD Stool parasites (fast-track DIAGNOSTICS)			
		+	-	Total
	+	58	2*	60
	-	0	112	112
	Total	58	114	172

Table 6. Comparative results of *Giardia lamblia*

* The low amount of template DNA detected in these samples is below the detection limit of the method used and so it may produce random positive results.

VIASURE <i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Real Time PCR Detection Kit	RIDA®GENE Parasitic Stool Panel II (r-Biopharm) + FTD Stool parasites (fast-track DIAGNOSTICS)			
		+	-	Total
	+	8 [‡]	1 ^{‡,*}	9
	-	0	163	163
	Total	8	164	172

Table 7. Comparative results of *Entamoeba histolytica*

* The low amount of template DNA detected in these samples is below the detection limit of the method used and so it may produce random positive results.

† Positive results were obtained for seven clinical stool specimens and 2 samples which were spiked with culture of *Entamoeba histolytica* prior to DNA extraction.

The results show a high sensitivity and specificity to detect *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* using VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit has a detection limit of ≥ 50 DNA copies per reaction for *Cryptosporidium* and ≥ 10 DNA copies per reaction for *Giardia lamblia* and *Entamoeba histolytica* (Figure 2, 3 and 4).

Figure 2. Dilution series of *Cryptosporidium* (10^7 - 10^2 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Cy5 channel).


Figure 3. Dilution series of *Giardia lamblia* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (FAM channel).



Figure 4. Dilution series of *Entamoeba histolytica* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (ROX channel)

12.3. Analytical specificity

The specificity of the *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* assays were confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine.

No cross-reactivity of *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* was detected against any of the following microorganisms tested.

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter lari</i>	-	<i>Enteropathogenic Escherichia coli</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter fetus</i>	-	<i>Klebsiella oxytoca</i>	-
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter coli</i>	-	<i>Listeria monocytogenes</i>	-
<i>Helicobacter heilmannii</i>	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-	<i>Candida albicans</i>	-
<i>Shigella flexneri</i>	-	<i>Campylobacter upsaliensis</i>	-	<i>Arcobacter butzleri</i>	-
<i>Shigella dysenteriae</i>	-	<i>Proteus vulgaris</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Salmonella typhi</i>	-	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	-	<i>Enterococcus faecalis</i>	-
<i>Salmonella paratyphi</i> A	-	<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3	-
<i>Salmonella paratyphi</i> B	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Salmonella typhimurium</i>	-	<i>Serratia liquefaciens</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella bongori</i>	-	<i>Vibrio parahaemolyticus</i>	-	Adenovirus serotypes 40/41	-
<i>Salmonella enteritidis</i>	-	<i>Clostridium difficile</i>	-	Rotavirus A	-
<i>Salmonella enterica</i> subsp. <i>enterica</i>	-	<i>Clostridium perfringens</i>	-	Norovirus Genotypes I and II	-
<i>Salmonella pullorum</i>	-	<i>Enterotoxigenic Escherichia coli</i>	-	Astrovirus Genotype I-VIII	-
<i>Salmonella gallinarum</i>	-				

Table 8. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for *Cryptosporidium* was evaluated against *Cryptosporidium parvum* showing positive results.

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for *Giardia lamblia* was evaluated against *Giardia lamblia* showing positive results.

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for *Entamoeba histolytica* was evaluated against *Entamoeba histolytica* strain DS4-868 showing positive result.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Fast Real-Time PCR System
Applied Biosystems	7500 Fast Dx Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
Bio-Rad	CFX96™ Real-Time PCR Detection System
Roche	LightCycler®480 Real-Time PCR System [#]
Roche	LightCycler®96 Real-Time PCR System [#]
Roche	Cobas z480 Analyzer [#]
Agilent Technologies	AriaMx Real-Time PCR System
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®*

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Bio-Rad	CFX96™ Deep Well Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Eppendorf	Mastercycler™ep realplex
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
Abbott	Abbott m2000 RealTime System
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DPrime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®*
VIASURE	VIASURE 48 Real Time PCR System †
VIASURE	VIASURE 96 Real Time PCR System †

Table A1/A2. Compatible low and high profile Real Time PCR systems.

* The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

Shell Frame grid plate which fits in these Roche qPCR System is necessary.

† See Annex 3 to check optical measurement exposure setting.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". Set exposition values as follow:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -150, HEX channel – 3000, ROX channel – 2000 and Cy5 channel - 1500.
- DTLite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 150, HEX channel - 500, ROX channel – 500 and Cy5 channel – 100.

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13. Bibliography/Bibliografía

1. J J.R. Limor. Detection and differentiation of *Cryptosporidium* parasites that are pathogenic for humans by real-time PCR. *Journal of Clinical Microbiology* 2002; 40(7): 2335-2338.
2. S.E. Stroup et al. Real-time PCR detection and speciation of *Cryptosporidium* infection using Scorpion probes. *Journal of Medical Microbiology* 2006; 55: 1217-1222
3. S.J. Hadfield et al. Detection and differentiation of *Cryptosporidium* spp. in human clinical samples by use of real-time PCR. *Journal of Clinical Microbiology* 2011; 49(3): 918-924
4. D.V. den Bossche et al. Comparison of four rapid diagnostic tests, ELISA, microscopy and PCR for the detection of *Giardia lamblia*, *Cryptosporidium* spp. and *Entamoeba histolytica* in feces. *Journal of microbiological methods* 2015; 110: 78-84

5. R.D. Adam. Biology of *Giardia lamblia*. *Clinical Microbiology Reviews* 2001; 14(3): 447-475.
6. J.J. Verweij *et al.* Real-time PCR for the detection of *Giardia lamblia*. *Molecular and Cellular Probes* 2003; 17(5): 223-225
7. R. Fotedar *et al.* PCR Detection of *Entamoeba histolytica*, *Entamoeba dispar*, and *Entamoeba moshkovskii* in Stool Samples from Sydney, Australia. *Journal of Clinical Microbiology* 2007; 45 (3): 1035-1037.
8. R. Haque *et al.* Multiplex real-time PCR assay for detection of *Entamoeba histolytica*, *Giardia intestinalis*, and *Cryptosporidium* spp. *The American Journal of Tropical Medicine and Hygiene* 2007; 76(4): 713-717.
9. Y. Lau *et al.* Real-time PCR assay in differentiating *Entamoeba histolytica*, *Entamoeba dispar*, and *Entamoeba moshkovskii* infections in Orang Asli settlements in Malaysia. *Parasites & Vectors* 2013; 6(1): 250.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cyclers, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid) and Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommended to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit ha sido validado en los siguientes equipos: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cyclers, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid) y Roche Molecular Diagnostics Cobas z480 Analyzer. Cuando se utiliza el equipo Applied Biosystems 7500 Fast con tiras, se recomienda colocar el soporte adecuado para reducir el riesgo de aplastar el tubo (Ref. PN 4388506).

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev00

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC