VIASURE

Real Time PCR Detection Kits


by CerTest
BIOTEC

Tick Borne Diseases

Handbook for the following references/

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VIASURE <i>Tick Borne Diseases</i> Real Time PCR Detection Kit 9 x 8-well strips, low profile	VS-TBD106L
VIASURE <i>Tick Borne Diseases</i> Real Time PCR Detection Kit 9 x 8-well strips, high profile	VS-TBD106H
VIASURE <i>Tick Borne Diseases</i> Real Time PCR Detection Kit 18 x 8-well strips, low profile	VS-TBD112L
VIASURE <i>Tick Borne Diseases</i> Real Time PCR Detection Kit 18 x 8-well strips, high profile	VS-TBD112H
VIASURE <i>Tick Borne Diseases</i> Real Time PCR Detection Kit 27 x 4-well strips, Rotor-Gene®	VS-TBD136

ENGLISH

1. Intended use

VIASURE Tick Borne Diseases Real Time PCR Detection Kit is designed for the specific identification and differentiation of viral RNA or genomic DNA specific for Tick Borne Encephalitis Virus (TBEV), *Rickettsia* spp., *Babesia microti*, *Babesia divergens*, *Ehrlichia chafeensis*, *Ehrlichia muris*, *Borrelia burgdorferi* sensu lato (s.l.), *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophilum* and/or *Coxiella burnetii* in blood, serum, tissue samples and microbiological culture from ticks, biopsy skin, cerebrospinal fluid (CSF) and synovial fluid from patients with signs and symptoms of Tick Borne diseases. This test is intended for use as an aid in the diagnosis of Tick Borne diseases in combination with clinical and epidemiological risk factors. RNA/DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific TBEV, *Rickettsia* spp., *Babesia microti*, *Babesia divergens*, *Ehrlichia chafeensis*, *Ehrlichia muris*, *Borrelia burgdorferi* sensu lato s.l., *Borrelia miyamotoi* and *Borrelia hermsii*, *Anaplasma phagocitophilum* and *Coxiella burnetii*.

2. Summary and Explanation

Tick Borne diseases comprise a group of infections transmitted to humans by the bite of ticks infected with bacteria, viruses, or parasites. Tick Borne diseases affecting humans include: Lyme disease, relapsing fever, babesiosis, anaplasmosis, ehrlichiosis, Q fever, Tick Borne encephalitis and spotted fever among others.

Lyme disease (or Lyme borreliosis) is the most common Tick Borne disease worldwide. It is caused by the spirochaete *Borrelia burgdorferi* and it is transmitted to humans through the bite of infected blacklegged ticks of the genus *Ixodes*. Typical symptoms of Lyme disease include fever, headache, fatigue and a characteristic skin rash called erythema migrans. If left untreated, infection can spread to joints, the heart, and the nervous system.

Borrelia miyamotoi and *B. hermsii* are species of spiral-shaped bacteria that is closely related to the bacteria that cause tick-borne relapsing fever (TBRF). First identified in 1995 in ticks from Japan, *B. miyamotoi* has also been detected in different *Ixodes* ticks species, whereas spirochete *Borrelia hermsii* is transmitted by its argasid tick vector, *Ornithodoros hermsi*. Patients with this infection were most likely to have fever, chills, and headache. Other common symptoms included body and joint pain and fatigue.

Anaplasmosis is caused by the bacterium *Anaplasma phagocytophilum* and is transmitted to humans through ticks of the genus *Ixodes*. Anaplasmosis shows high clinical variability, with symptoms like headache, fever, chills, malaise, muscle pain, nausea, cough, confusion and rash. If not treated correctly, anaplasmosis can evolve to severe clinical manifestations and even death in <1% of cases.

Q fever is a zoonosis caused by the bacteria *Coxiella burnetii*. Cattle, sheep, and goats are the primary reservoirs of *Coxiella burnetii*, and transmission to humans occurs primarily through inhalation of aerosols from contaminated soil or animal waste, though it can be transmitted through tick bites. Only about 50% of the infected people show clinical symptoms, which are flu-like symptoms such as headache, fever, chills, fatigue, muscle aches, nausea, cough, chest pain and weight loss. In severe cases people may develop pneumonia or hepatitis.

Babesiosis is produced by many species of protozoa of the genus *Babesia*, mainly *Babesia microti* and *Babesia divergens*. Babesiosis is also transmitted by ticks of the genus *Ixodes*, and it is frequently found as a coinfection with Lyme disease. Babesiosis usually shows none or flu-like mild symptoms, such as headache, fever, chills, body aches, loss of appetite, nausea or fatigue. However, since *Babesia* parasites infect and destroy red blood cells, babesiosis can lead to hemolytic anemia.

Ehrlichiosis is due to different species of bacteria of the genus *Ehrlichia*. Causative agents in humans are *Ehrlichia chaffeensis* and *Ehrlichia muris*. Ehrlichiosis is transmitted by the lone star tick (*Amblyomma americanum*). Typical symptoms of ehrlichiosis include headache, fever, chills, malaise, muscle pain, nausea, conjunctival infection, confusion and rash. Severe disease may present lethargy, myalgia, reduction of sodium levels and platelets and elevated liver enzymes, with fatal outcome in 3% of cases.

Tick Borne encephalitis is produced by the Tick Borne encephalitis virus (TBEV) of the family *Flaviviridae*. It is transmitted to humans through ticks of the genus *Ixodes*. The onset of the disease includes nonspecific symptoms such as fever, malaise, anorexia, muscle aches, headache, nausea and vomiting. In 20-30% of patients, a second phase is observed that involves the central nervous system with symptoms of meningitis, encephalitis, or meningoencephalitis, leaving neuropsychiatric sequelae in 10-20% of patients.

Spotted fever is caused by bacteria of the genus *Rickettsia* and is widely distributed by different geographical areas, being able to transmit depending on it by different ticks such as: wood tick, *Dermacentor Andersoni* (in the Rocky Mountain states, USA), dog tick, *Dermacentor variabilis* (other areas of the USA), ticks *Amblyomma cajennense* (in South America) and *Rhipicephalus sanguineus* (in Mexico). Spotted fever is characterized by darks scabs at the site of the tick bite (eschar), accompanied by nonspecific symptoms such as headache, fever, rash and muscle ache.

Since most Tick Borne diseases show similar symptoms, diagnosis can be problematic. Real-time PCR assays have been shown to be a sensitive and specific diagnostic tool for the detection of the causative agent.

3. Principle of the procedure

VIASURE Tick Borne Diseases Real Time PCR Detection Kit is designed for the diagnosis of TBEV, *Rickettsia* spp., *Babesia microti*, *Babesia divergens*, *Ehrlichia chaffeensis*, *Ehrlichia muris*, *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* in clinical samples. The detection of TBEV is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of the 3'UTR sequence of TBEV using specific primers and a fluorescent-labelled probe. After DNA isolation, the identification of rest of pathogens is performed by the amplification of a conserved region of the 23S rRNA gene (*Rickettsia* spp.), *CCT-eta* gene (*Babesia microti*), *hsp70* gene (*Babesia divergens*), *GroEl* gene (*Ehrlichia chaffeensis* and *Ehrlichia muris*), 23S rRNA gene (*Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and/or *Borrelia hermsii*), *msp2* gene (*Anaplasma phagocitophylum*) and *IS1111* gene (*Coxiella burnetii*), using specific primers and a fluorescent-labelled probe.

VIASURE Tick Borne Diseases Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE Tick Borne Diseases Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. Each kit includes three kinds of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophilum* and/or *Coxiella burnetii* (*Borrelia*, *Anaplasma* & *Coxiella* 4/8-well strips). *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and *Borrelia hermsii* DNA targets are amplified and detected in FAM channel, *Anaplasma phagocitophilum* DNA targets are amplified and detected in ROX channel, and *Coxiella burnetii* DNA targets are amplified and detected in Cy5 channel. The first strip (*Borrelia*, *Anaplasma* & *Coxiella* 4/8-well strips) contains an internal control (IC) which is amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2). The second strip contains the multiplex reaction mix for the detection of *Rickettsia* spp., *Babesia microti* and/or *Babesia divergens*, and *Ehrlichia chafeensis* and/or *Ehrlichia muris* (*Rickettsia*, *Babesia* & *Ehrlichia* 4/8-well strips). *Rickettsia* spp. DNA targets are amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex. 2), *Babesia microti* and *Babesia divergens* DNA targets are amplified and detected in ROX channel and *Ehrlichia chafeensis* and *Ehrlichia muris* DNA targets are amplified and detected in FAM channel. The third strip contains the monoplex reaction mix for the detection of TBEV (TBEV 4/8-well strips). TBEV RNA targets are amplified and detected in FAM channel.

4. Reagents provided

VIASURE Tick Borne Diseases Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1 and 2. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>Borrelia</i> , <i>Anaplasma</i> & <i>Coxiella</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	3/6 x 8-well strip
<i>Rickettsia</i> , <i>Babesia</i> & <i>Ehrlichia</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	White	3/6 x 8-well strip
TBEV 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	White	3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1/2 vials x 1.8 mL
Tick Borne Diseases Positive Control	Non-infectious synthetic lyophilized cDNA/DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE Tick Borne Diseases Real Time PCR Detection Kit with Ref. VS-TBD106L, VS-TBD106H, VS-TBD112L and VS-TBD112H.

Reagent/Material	Description	Colour	Amount
<i>Borrelia</i> , <i>Anaplasma</i> & <i>Coxiella</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	9x 4-well strip
<i>Rickettsia</i> , <i>Babesia</i> & <i>Ehrlichia</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	Transparent	9x 4-well strip
TBEV 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs in stabilized format	Transparent	9x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Tick Borne Diseases Positive Control	Non-infectious synthetic lyophilized cDNA/DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	27 x 4-cap strip

Table 2. Reagents and materials provided in VIASURE Tick Borne Diseases Real Time PCR Detection Kit with Ref. VS-TBD136. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Tick Borne Diseases* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA/DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Tick Borne Diseases* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-TBD136). Remove any excess air in the pouches prior to sealing.

- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For reference VS-TBD136 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Make sure to use a well for determining TBEV, another well for determining *Rickettsia* spp., *Babesia microti*/*Babesia divergens* and *Ehrlichia chafeensis*/*Ehrlichia muris* and another well for determining *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophilum* and/or *Coxiella burnetii*. Be careful not to mix them throughout the process.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Sample preparation

The proper collection and transport of clinical specimens is critical for the isolation, identification, and characterization of pathogens. Specimens (ticks, blood, serum, tissue samples, biopsy skin, cerebrospinal fluid (CSF) and synovial fluid) should be collected properly in a clean area and processed as soon as possible to avoid loss of viability of the etiological agents for suitable microbiological culture and/or to prevent nucleic acid degradation, as well as, to guarantee the quality of the test. We recommend to use fresh samples or immediately frozen.

For longer storage, the samples should be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenize sample as thoroughly as possible prior to preparation. Freezing and thawing cycles before isolating nucleic acids should be avoided.

Given that the number of *B. burgdorferi* spirochetes in infected tissues or body fluids of patients is very low, appropriate procedures for sample collection and transport and preparation of DNA from clinical samples are critical for yielding reliable and consistent PCR results.

Blood samples might be collected using vacutainer system (tubes stabilized with EDTA) and be stored at 4°C for up to one week. Blood components (serum and plasma samples) can be obtained after blood centrifugation and should be stored at 4 °C for up to one week or at -20°C indefinitely. Nucleic acids isolation from whole blood, serum or plasma could be performed using 100 µL and/or following the recommendations of the manufacturer.

Tissue samples should be collected in a clean container and stored immediately at -20°C or -80°C until use.

Ticks could be collected directly from patients using a mechanical method (forceps). For bacterial growth, tissues from each tick could be pooled and cultivated in primary cell line.

For biopsy skin, the peripheral border of the cutaneous lesion should be identified and under sterile conditions 4-mm-diameter punch biopsy specimen could be obtained from the peripheral aspect of the lesion. The specimens should be stored immediately at -20°C or -80°C. For microorganism isolation, each skin biopsy specimen could be placed in a polystyrene tube containing 6 ml of isolation/broth medium without antibiotics and incubated at corresponding temperature during the corresponding time. If the tissue is kept in BSK medium for over 24h, some spirochetes will have migrated from the skin biopsy to the culture medium. In this case, DNA should be prepared from both the skin biopsy and the medium.

Cerebrospinal fluid specimen (CSF) should be processed in a microbiology laboratory within 1 hour after collection or inoculated into Trans-Isolate (T-I) medium or similar for transport to the laboratory if processing within 1 hour is not feasible. These samples should be maintained at room temperature prior testing. Refrigeration is not recommended. If samples can be processed within 1 hour, centrifuge the CSF for 15 min at 1000 x g. Take the sediment and seed it in primary media. If not, incubate the T-I medium over-night.

Synovial fluids (SF) samples from patients with Lyme arthritis could be collected without additives and stored overnight initially at 4° and then transferred to -20°C or -80°C until use.

After specimen microbiology cultivation (ticks, biopsy skin, cerebrospinal fluid (CSF) and synovial fluid), acid nucleic could be extracted using 1 mL of the culture after being pelleted by centrifugation (2 min 8000 x g). To perform an extraction of bacterial DNA from gram positive Bacteria addition of lysozyme is needed.

This section summarized a brief description of sample preparation. To perform more suitable sample collection, transport and storage; we suggest following the manufacturer's recommendations appearing in the instructions for use of the extraction kit used.

8.1.1. RNA/DNA extraction

For RNA/DNA extraction from serum, blood, tissue and microbiological culture samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available RNA/DNA extraction kit and follow the manufacturer's instructions.

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

8.2. Lyophilized positive control

Tick Borne Diseases Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Tick Borne Diseases* Positive Control (red vial) by adding 300 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA/DNA extracted from each sample, reconstituted *Tick Borne Diseases* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*TBEV*, *Borrelia burgdorferi* s.l./*Borrelia miyamotoi*/*Borrelia hermsii* and *Ehrlichia chafeensis*/*Ehrlichia muris*), HEX, JOE or VIC (*Rickettsia* spp and Internal Control (IC)), ROX (*Babesia microti*/*Babesia divergens* and *Anaplasma phagocitophilum*), and Cy5 channels (*Coxiella burnetii*). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005PTM Real Time PCR System check

that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in negative control well and the presence of signal for *Tick Borne Diseases* Positive Control well. The analysis of the samples is done by the software itself of the used real time PCR equipment according to manufacturer's instructions.

Interpretation of results for *Borrelia*, *Anaplasma* & *Coxiella* 4/8-well strips:

- A sample is considered positive for ***Borrelia burgdorferi s.l./Borrelia miyamotoi/ Borrelia hermsii*** if there is an amplification signal in FAM channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for ***Anaplasma phagocitophylum*** if there is an amplification signal in ROX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for ***Coxiella burnetii*** if there is an amplification signal in Cy5 channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.
- In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

Interpretation of results for *Rickettsia*, *Babesia* & *Ehrlichia* 4/8-well strips:

- A sample is considered positive for ***Rickettsia spp.*** if there is an amplification signal in HEX/VIC/JOE channel and the Ct value obtained is less than 40.
- A sample is considered positive for ***Babesia microti/Babesia divergens*** if there is an amplification signal in ROX channel and the Ct value obtained is less than 40.
- A sample is considered positive for ***Ehrlichia chafeensis/Ehrlichia muris*** if there is an amplification signal in FAM channel and the Ct value obtained is less than 40.

- A sample is considered negative, if the sample shows no amplification signal in the detection system, but the internal control of *Borrelia*, *Anaplasma* & *Coxiella* assay is positive.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.

Interpretation of results for TBEV 4/8-well strips:

- A sample is considered positive for **TBEV** if there is an amplification signal in FAM channel and the Ct value obtained is less than 40.
- A sample is considered negative, if the sample shows no amplification signal in the detection system, but the internal control of *Borrelia*, *Anaplasma* & *Coxiella* assay is positive.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Borrelia*, *Anaplasma* & *Coxiella*).


Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Rickettsia*, *Babesia* & *Ehrlichia*).



Figure 3. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Monoplex reaction mix TBEV).

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control (if available) signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with RNA/DNA extracted from blood, serum, tissue samples and microbiological culture from ticks, biopsy skin, cerebrospinal fluid (CSF) and synovial fluid.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by contamination by Tick Borne diseases, either samples containing high concentrations of target RNA/DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE Tick Borne Diseases Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well included in the strip 1 (*Borrelia*, *Anaplasma* & *Coxiella* 4/8-well strips) confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE Tick Borne Diseases Real Time PCR Detection Kit for *Borrelia*, *Anaplasma* & *Coxiella* 8-well strip was tested using 95 DNA samples extracted from microbiological culture from ticks, biopsy skin, cerebrospinal fluid (CSF) and synovial fluid. A total of 17 well characterized *Borrelia* strains comprising 9 different *Borrelia burgdorferi* sensu lato genospecies, notably *B. japonica* (n=1), *B. burgdorferi* sensu stricto (n=2; B31 and PBr strains), *B. bavariensis* (n=1; PBi strain), *B. garinii* (n=5; PBr, PHei, PWudII, PRef and PLa strains), *B. bissettii* (n=1; PGeb strain), *B. afzelii* (n=2; PKo and PVPM strains), *B. lusitanae* (n=1; Poti B2strain), *B. spielmanii* (n=1; PSiglI strain), *B. valaisiana* (n=1; VS116 strain)) were included. Additionally, two relapsing fever control strains *B. hermsii* (n=1) and *B. miyamotoi* (n=1) were included as well as the potentially cross-reactive spirochaetes *Leishmania* spp. (n=2) and *Treponema* spp. (n=1) were tested. The VIASURE real-time PCR successfully detected all tested *Borrelia* sensu lato genospecies and the relapsing fever group strains *B. hermsii* and *B. miyamotoi*. No cross reactivity was observed with DNA from *Leptospira* and *Treponema* species for VIASURE assay.

VIASURE Tick Borne Diseases Real Time PCR Detection Kit for *Borrelia*, *Anaplasma* & *Coxiella* 8-well strip was evaluated with 3 INSTAND *Coxiella burnetii* & *Bacillus anthracis* panels from 2017 and 2018, as well as, 17 additional tissue samples. The results were compared with the final EQA program reports or with those obtained by a commercial qPCR assay (EXOone *Coxiella burnetii* (EXOPOL)). All *Coxiella burnetii* positive samples (6/12) from 3 INSTAND programs were detected and 15/17 tissue samples showed a positive result in the identification of *Coxiella burnetii*.

The clinical performance of VIASURE Tick Borne Diseases Real Time PCR Detection Kit for *Rickettsia*, *Babesia* & *Ehrlichia* 8-well strips and for TBEV 8-well strips was evaluated using 90 samples from 16 different QCMD panels (Tropical diseases and borreliosis panels) and 2 clinical specimens (serum and blood). VIASURE assay for TBEV 8-well strips found 8/90 positive samples for Tick Borne Encephalitis Virus.

In conclusion, the results show a high sensitivity and specificity to detect TBEV, *Rickettsia* spp., *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophilum* and *Coxiella burnetii* using VIASURE Tick Borne Diseases Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE Tick Borne Diseases Real Time PCR Detection Kit has a detection limit of ≥ 10 RNA/DNA copies per reaction (Figure 4, 5, 6, 7, 8, 9 and 10).

Figure 4. Dilution series of *Borrelia burgdorferi*/*Borrelia miyamotoi*/ *B. hermsii* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Borrelia*, *Anaplasma* & *Coxiella*, channel FAM)


Figure 5. Dilution series of *Anaplasma phagocitophilum* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Borrelia*, *Anaplasma* & *Coxiella*, channel ROX).


Figure 6. Dilution series of *Coxiella burnetii* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Borrelia*, *Anaplasma* & *Coxiella*, channel Cy5).


Figure 7. Dilution series of *Rickettsia* spp (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Rickettsia*, *Babesia* & *Ehrlichia*, channel HEX).


Figure 8. Dilution series of *Babesia microti*/*Babesia divergens* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Rickettsia*, *Babesia* & *Ehrlichia*, channel ROX).


Figure 9. Dilution series of *Ehrlichia chafeensis*/*Ehrlichia muris* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *Rickettsia*, *Babesia* & *Ehrlichia*, channel FAM).



Figure 10. Dilution series of TBEV (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Monoplex reaction mix TBEV, channel FAM).

12.3. Analytical specificity

The specificity of the Tick Borne Diseases assay was confirmed by testing a panel consisting of different microorganisms representing the most common Tick Borne pathogens. No cross-reactivity was detected against any of the following microorganisms tested:

Cross-reactivity testing			
<i>Anaplasma marginale</i>	-	<i>Borrelia garinii</i>	-/+
<i>Bartonella henselae</i> strain Houston-1	-	<i>Borrelia japonica</i>	-/+
<i>Borrelia hermsii</i>	-/+	<i>Borrelia miyamotoi</i>	-/+
<i>Borrelia lusitanae</i>	-/+	<i>Borrelia spielmanii</i>	-/+
<i>Borrelia valaisiana</i>	-/+	<i>Coxiella burnetii</i> strain Nine Mile Q	-/+
<i>Borrelia azfeli</i> strain P-Ko/1984	-/+	<i>Theileria annulata</i>	-
<i>Borrelia bavariensis</i>	-/+	<i>Rickettsia conorii</i> strain Moroccan	-/+
<i>Borrelia bisetti</i>	-/+	<i>Leptospira</i>	-
<i>Borrelia burgdorferi</i> sensu stricto-strain IRS	-/+	<i>Treponema phagedenis</i>	-
<i>Borrelia burgdorferi</i> sensu stricto strain B31	-/+	Tick Borne Encephalitis Virus (TBEV) strain Neudorfl	-/+

Table 4. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE Tick Borne Diseases Real Time PCR Detection Kit was evaluated against TBEV strain Neudorfl, *Rickettsia conorii* strain Moroccan, synthetic sequence of *Babesia microti*, synthetic sequence of *Babesia divergens*, synthetic sequence of *Ehrlichia chafeensis*, synthetic sequence of *Ehrlichia muris*, *Borrelia azfeli* (P-Ko/1984 and PVPM strains), *Borrelia bavariensis* (PBi strain), *Borrelia bisettii* (PGeb strain), *Borrelia bisettiae*, *Borrelia burgdorferi* sensu stricto (B31, IRS and PBre strains), *Borrelia garinii* (PHei, PWudII, PRef and PLa strains) and *Borrelia garinii* OspA Typ3 strains, *B. japonica*, *B. lusitanae* (Poti B2 strain), *B. spielmanii* (PSiglI strain), *Borrelia valaisiana* (VS116 strain), *B. hermsii* and *B. miyamotoi*, *Anaplasma phagocitophilum* and *Coxiella burnetii* strain Nine Mile Q showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía

1. V. Parikh et al. Infections of the nervous system. *International Journal of Critical Illness and Injury Science* 2012; 2(2): 82-97.
2. J. Brites-Neto et al. Tick Borne infections in human and animal population worldwide. *Veterinary World* 2015, 8(3):301-315.
3. J.J. Sormunen et al. Tick-borne bacterial pathogens in southwestern Finland. *Parasit Vectors*. 2016; 9: 168.
4. M.P. Nelder et al. Human pathogens associated with the blacklegged tick *Ixodes scapularis*: a systematic review. *Parasites & Vectors* 2016; 9:265
5. L. Michelet et al. High-throughput screening of tick-borne pathogens in Europe. *Frontiers in Cellular and Infection Microbiology* 2014; 4:103.

6. M.A. Diuk-Wasser *et al.* Coinfection by the Tick Borne pathogens *Babesia microti* and *Borrelia burgdorferi*: ecological, epidemiological and clinical consequences. *Trends in Parasitology* 2016; 32(1): 30–42.
7. M. Schwaiger *et al.* Development of a quantitative real-time RT-PCR assay with internal control for the laboratory detection of tick borne encephalitis virus (TBEV) RNA. *Journal of Clinical Virology* 2003; 27(2): 136-145.
8. C.Y. Kato *et al.* Assessment of Real-Time PCR Assay for Detection of *Rickettsia* spp. and *Rickettsia rickettsii* in Banked Clinical Samples. *Journal of Clinical Microbiology* 2013; 51(1): 314-317.
9. J. Liu *et al.* Development of a TaqMan Array Card for Acute-Febrile-Illness, Outbreak Investigation and Surveillance of Emerging Pathogens, Including Ebola Virus. *Journal of Clinical Microbiology* 2016; 54(1): 49-58.
10. Centers for Disease Control and Prevention ([https://www.cdc.gov/niosh/topics/Tick Borne/](https://www.cdc.gov/niosh/topics/Tick%20Borne/))
11. World Health Organization (<http://www.who.int/en/>).

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: May 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.cerfest.es

VIASURE online

F-362 rev01

VIASURE

CE

IVD

Real Time PCR Detection Kits

CerTest
BIOTEC