

### **Hoja de especificaciones del producto**

#### **TS/17-B**

Hisopo con mango flexible de polipropileno con punta de viscosa pura, empackado individualmente en bolsa aparte abre fácil

|                           |                                                                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Código del producto       | <b>TS/17-B</b>                                                                                                                                                   |
| Descripción               | Hisopo con mango flexible de polipropileno con punta de viscosa pura, empackado individualmente en bolsa aparte abre fácil                                       |
| Intención de uso          | Hisopo para cultivo; Muestreo de pacientes y uso general en el laboratorio.                                                                                      |
| Componentes               | Tapa de poliestireno punta de viscosa / mango de polipropileno (diámetro de la punta - 5,15 mm) bolsa aparte abre fácil pegamento no tóxico Etiqueta y película. |
| Contenido                 | No aplica                                                                                                                                                        |
| Método de esterilización  | Esterilizado por irradiación gamma a dosis de 25-46 kGy                                                                                                          |
| Vida útil                 | 24 meses desde la fecha de fabricación                                                                                                                           |
| Condiciones de almacenaje | 4-25 ° C; en zona seca y ventilada; proteger de la luz, el calor y la humedad en un recipiente original bien cerrado.                                            |
| Embalaje                  | 1000 en una caja                                                                                                                                                 |
| otros                     | Detalles de las pruebas de control de calidad en el producto final: i) Validación de la esterilidad por irradiación.                                             |

Esta especificación es sólo para orientación.

Usted siempre debe consultar a un profesional con experiencia en su aplicación prevista en caso de duda.

Edición #: 3 Fecha: 26/07/2013

ProductSpec\_TS17-B



## **Product Specification Sheet**

### **TS/17-B**

Swab with a flexible polypropylene shaft tipped with pure viscose, single in peel pouch

|                         |                                                                                                                      |
|-------------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Code            | TS/17-B                                                                                                              |
| Description             | Swab with a flexible polypropylene shaft tipped with pure viscose, single in peel pouch                              |
| Intended Use            | Culture swab; Patient sampling & general Laboratory use                                                              |
| Components              | Polystyrene cap viscose tip/polypropylene shaft swab (tip diameter – 5.15 mm) Peel pouch Non-toxic glue Label & film |
| Contents                | Not applicable.                                                                                                      |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-46kGy                                                                |
| Shelf Life              | 24 months from date of manufacture                                                                                   |
| Storage Conditions      | 4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container      |
| Pack Size               | 1000 in a box                                                                                                        |
| Other                   | Quality control testing details on final product: i) Validation of irradiation sterility.                            |

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 26.07.2013

ProductSpec\_TS17-B

