



# Microbiology

## MYCOFAST® REVOLUTION 2

Urogenital mycoplasma diagnosis  
**DIRECT METHOD**



### FULFILLING YOUR NEEDS IN THE MOST DEMANDING SITUATIONS

*Detection, Enumeration, Identification and Antimicrobial Susceptibility Testing*

- ◆ Complete diagnosis in an All-in-1 product
- ◆ Compatible for use with liquid-Amies collection and transport media
- ◆ Optimized management of co-infections
- ◆ Antibiotics compliant with worldwide recommendations
- ◆ The right pathologic threshold tested for each specimen type
- ◆ Fast results for clinicians
- ◆ Visual, colorimetric reading



**Labcare  
de Colombia**  
[www.labcarecolombia.com](http://www.labcarecolombia.com)



**ELITechGroup**  
EMPOWERING IVD

# MYCOFAST® REVOLUTION 2

## BENEFITS

### Complete diagnosis in 24h\*

- ◆ Detection, Enumeration
- ◆ Identification
- ◆ Antimicrobial Susceptibility Testing

### Cutting-edge performance

#### Sensitivity

- ◆ Detection and enumeration of *U. urealyticum* at  $10^3$ ,  $10^4$  and  $\geq 10^5$  CCU/mL
- ◆ Detection and enumeration of *M. hominis* at  $\geq 10^4$  CCU/mL
- ◆ The right pathologic threshold tested for each specimen type

### Specificity

- ◆ Selective growth media for each species tested separately

#### Reliability

- ◆ 5 antibiotics (LVX/MXF/E/CM/TE) tested in compliance with CLSI guidelines (Clinical & Laboratory Standards Institute)
- ◆ Doxycyclin added to the panel

#### Fast results

- ◆ Integrated growth activator for *M. hominis*


### Ease of use

- ◆ Ready-to-use reagents
- ◆ Each well pre-filled with growth medium
- ◆ Visual colorimetric reading
- ◆ No instrument required

### Flexibility and Safety

- ◆ Optimized transport medium: Specimens preserved up to 72 hours before analysis
- ◆ Clip-on closing system for the tray

## METHODOLOGY



\* Reading of all samples in 24h.  
Only read liquid samples having negative results within 24h.

## PRODUCT RANGE

REF.	NAME	QTY	PURPOSE
00080	MYCOFAST® Revolution 2	25 tests	Complete kit
00081	MYCOFAST® Revolution 2 AMIES	25 tests	Complete kit
00061	UMMt RevolutionN (optional)	50 vials	Additional transport medium
00083	UMMt AMIES RevolutionN (optional)	50 vials	Additional transport medium (for Amies specimens)
00082	COMPLEMENT MYCOFAST® Revolution 2 (optional)	25 tests	Additional test trays
00064	MYCOPLASMA Stabilizer (optional)	100 tests	Additive for sample storage > 56 hours
00900	MYCOPLASMA Control (optional)	12 tests	U.u. and M.h. control strains

OPTIONAL



#### WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.

#### Headquarters

Australia  
BELUX  
Brazil  
France  
Italy

#### +33 1 41 45 07 10

1800 815 098  
+32 92 820 531  
+55 27 3025 1415  
+33 4 83 36 10 82  
+39 011 67 61 913

New Zealand  
Serbia  
The Netherlands  
UK  
United States

0800 555 611  
+381 11 2467119  
+31 313 430 500  
+44 1442 869320  
+1 435 0752 6011



ELITECH Microbio (France)



ELITECH Group

EMPOWERING IVD



Labcare  
de Colombia  
www.labcarecolombia.com