



## HQ-Chex



HQ-Chex is a whole blood glucose and hemoglobin quality control assayed for select HemoCue systems. The control resembles the performance of patient samples and is available in three clinically significant levels. Its long shelf life reduces the number of shipments, associated costs and time spent performing lot-to-lot correlation. HQ-Chex is packaged in squeeze dropper vials for accurate sample dispensing, eliminating the need to pipette.

### Features

- + Whole blood control resembles patient sample
- + Available in three clinically significant levels
- + Packaged in plastic squeeze dropper vials
- + Compatible with HemoCue Glucose 201 Analyzer, Hb201+ Analyzer and B-Hemoglobin Photometer
- + 30-day open-vial stability
- + 180-day closed-vial stability

### Benefits

- + Helps clinical laboratories meet quality requirements
- + Established assay ranges provide less room for error
- + Long shelf life reduces the number of shipments and associated costs
- + Fewer lot-to-lot correlations saves time and resources
- + Access to *STATS*®, a free interlaboratory quality control program with peer group data comparison and online submission

### Ordering information

Level	Description	Catalog#
Level 1, 2, 3	3 x 2.5 mL	232749
Level 1	6 x 2.5 mL	232753
Level 2	6 x 2.5 mL	232754
Level 3	6 x 2.5 mL	232755

Level	Description	Catalog#
Level 1, 2, 3	6 x 2.5 mL	232756
Level 1, Level 3	6 x 2.5 mL	232757