



Demand, Ready-to-Use WasteWatR™

Item No. 743

One 500 mL whole-volume bottle is ready to analyze.



SPECIFICATIONS

Product Family	Environmental
Product Category	Water Pollution
Product Sub Category	Inorganics
Item Type	CRM
Accreditations	ISO 17034
Methods	EPA 405.1, EPA 410.1, EPA 410.2, EPA 410.3, EPA 410.4, EPA 415.1, EPA 415.2, EPA 9060, SM 5210 B, SM 5310 B, SM 5310 C, SM 5310 D, SM 5520 B
Shelf Life from date of manufacture (months)	2
Minimum Shelf Life on Ship Date (Months)	2
Container Type	500mL Amber glass bottle

Formulation	Liquid
Components/Pkg	1
Volume	500 mL
Weight	1.61 lbs
Application	Environmental