



NEW! **Locking** Microcentrifuge Tubes

Available In 0.5mL, 1.5mL, 2.0mL & 5mL MIDI Tube

New Globe Scientific Locking Microcentrifuge Tubes are engineered with a secure locking system that provides the security of a screw cap microtube with the simplicity of a snap cap tube.



Features at a Glance:

Secure – Maintains sample integrity with locking cap which prevents accidental cap opening and evaporation during centrifugation

Versatile – Single use tubes for laboratory procedures including centrifugation, boiling, incubation, freezing and sample storage

Quality – Tubes are made from optimally clear USP Class VI medical grade polypropylene

Certified – Manufactured in a class 100,000 clean room and are certified to be free from RNase and DNase, Human DNA, Pyrogens, PCR Inhibitors, and ATP

Ease of Use – Tubes feature graduations, writing area, pierceable cap, are autoclavable, and are packaged in zip, re-closable bags

Item #	Description	Unit
111552LK	0.5mL, polypropylene (PP), natural,	500
111562LK	1.5mL, polypropylene (PP), natural	500
111572LK	2.0mL, polypropylene (PP), natural	500
111582LK	5.0mL, polypropylene (PP), natural	100