

Certificate of Analysis

**For Research Use Only
Not for Human or Animal Consumption**

KPL 10X Tris-Glycine-SDS Buffer

PART CODE: 5460-0016 (51-11-02)
LOT NUMBER: 10683744
SIZE 5 L

DESCRIPTION Gel electrophoresis buffer. To prepare 1X running buffer for Tris-Glycine-SDS gel electrophoresis, mix 100 mL KPL 10X Tris-Glycine-SDS buffer and 900 mL reagent quality water. Do not adjust pH.

FORM/STORAGE Liquid. Store at room temperature.

CONTENTS Supplied as a 5L container(s).

STABILITY Product is stable for one year from date of receipt when stored at room temperature.

**DATE OF
MANUFACTURE:** 04 OCT 2023

Approval:

A handwritten signature in black ink, appearing to read 'E. Barker'.

17 OCT 2023

QA Verified By

Date