

Certificate of Analysis

Product Information

Product Name:	Certified Drug Free Urine
Product Number:	88121-CDF(L)
Date of Manufacture:	xx/xx/xx
Lot Number:	xxxx
Expiration Date:	xx/xx
Intended Use:	In vitro diagnostic
Homogeneity:	Sufficiently homogeneous for analytical testing

Precautions

Matrix:	Urine
Form:	Liquid
Anticoagulant:	N/A
Preservative:	Sodium Azide (0.01%)
pH:	4.5-8
Virology:	Non-reactive for HBsAg, HCV, and HIV
Microbial Growth:	<100 colony-forming units per milliliter

Storage and Stability

Storage Condition:	2-8°C (35-46°F)
Open Vial Stability:	Up to Expiration Date
Vial Volume:	1 gallon

ELISA Raw Material Screen

UTAK Laboratories, Urine Pool #xxx

Analyte	Qualitative Value	Cutoff Value	Unit
d-Amphetamine	Negative	10	ng/mL
d-Methamphetamine	Negative	5	ng/mL
Benzoylecgonine (BE)	Negative	8	ng/mL
Cotinine	Negative	10	ng/mL
Morphine	Negative	5	ng/mL
Oxazepam	Negative	10	ng/mL
(I)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	Negative	8	ng/mL

